DOMINIQUE DUTSCHER SAS

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**DOMINIQUE DUTSCHER SAS**

**INSTRUCTIONS FOR USE  
MULTI-TESTS FIELD EQUIPMENT**

DOMINIQUE DUISCHER SAS

# Summary

1.

■ Introduction p.23

2.

■ Carrying case External description p.24

3.

■ Description of the equipment p.25

4.

■ Installation of the Multi-Tests Field Equipment p.26

5.

■ Timer p.26 to 27

6.

■ Operation of the filtration manifold p.27 to 31

The waste bottle.....	p.27
Filtration stations.....	p.28
Usage of magnetized funnels and membrane supports...	p.29
Usage of disposable filtration units.....	p.29
Filtration.....	p.30
Revelation.....	p.31
Rinsing and purging the filtration manifold.....	p.31

7.

■ Incubation p.32  
Incubation process

8.

■ Colorimeter p.44 to 48

9.

■ Temperature recorder p.48 to 52

10.

■ Care & maintenance p.53

11.

■ Malfunctions p.54 to 56

## • IMPORTANT

This equipment is intended only for environmental analysis of Hot Domestic Waters, Drinking Waters and / or Cooling Tower Waters, and shall be used exclusively with Legio EZ reagents kits.

For technical assistance, please contact the technical support.  
Contact information on the back of this leaflet.

# Introduction

This case contains all the equipments to carry out the required analysis with **Legio EZ** reagents kits.

## RECOMMENDATIONS FOR AN OPTIMAL USE :

- You are required to carefully read the attached user manual before carrying out the analysis.

➔ **If you have any query about how to use your device, please contact the technical support.** ←

- The equipment is safe in normal conditions of use.
- It is the responsibility of the user to use this equipment in a well-lit area.
- The equipment shall be placed on a stable support at such a height that the plate of the equipment is at the level of the elbows.
  - Check that the power cord is not damaged before connecting it.
  - The power cord shall be already connected to the equipment before being connected into a 110/230V outlet. Similarly, the power cord shall be already disconnected from the 110/230V outlet before being disconnected from the equipment.
- The noise level emitted is less than 70dBA.
- This equipment shall be maintained or repaired by the renowned and skilled persons within the meaning of the local and currently valid legislation.

# Carrying case external description

➔ This carrying case is intended to protect equipment during transport.

*Designed in a high-quality plastic resin, the use of this material provides the carrying case increased strength and extreme impact resistance.*

## **Soft grip handle -**

This handle has been designed to promote the handling of the carrying case. The handle is coated with a soft material for greater comfort. Thanks to a «hard point» at the end, the handle is retained firmly in the stowed position to avoid any damage during other modes of transport.

## **Trolley with wheels -**

The carrying case is equipped with a trolley, foldable / retractable by simply pressing on its lock. The presence of ball bearings secures a frictionless rolling capacity that helps movement.

## **Easy Open Latches -**

Easy opening by pressing a button on the latches. Just press and release the latches at the same time.

## **Lockable Hasps -**

Twin mating hasp holes on the rim of the case allow attachment of padlocks. We recommend using «TSA» padlocks for air travel.

### **• IMPORTANT**

**No liquid should be present in the case during transport, whatever it is. You are required to carefully follow the instructions described in the chapter «Filtration manifold».**

# Description of the equipment

## LOCATION OF THE DIFFERENT COMPONENTS :

Colorimeter

Storage box for accessories including :

- Power cord
- Timer
- 10ml syringe
- Instructions for use Legio EZ Analysis Protocol
- Instructions for use for Multi-tests Field Equipment
- Instructions for use for Temperature recorder
- Depending on the Equipment chosen: a plasticized color chart, or a colorimeter with its instructions.


Un storyboard (visual memo)

Pressure gauge

Antibacterial filter

Waste bottle

6 holes for filtration equipped with 6 caps, 6 membrane supports and 6 magnetized funnels

Incubator (containing a Petri dish holder, equipped with a temperature recorder)

Blue pilot light POWER

USB port

3. Red pilot light "Hot Domestic Waters, Drinking Waters"

4. 3-positions power switch I/O/ II for the incubator

5. Red pilot light "Cooling Tower Waters"


7 8

7. Power switch ON / OFF for the filtration

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- Before use, the equipment shall be placed on a stable support at such a height that the plate of the equipment is at the level of the elbows.
- Take the power cord inside the storage box, and connect it first on the rear of the equipment.
- Then, connect the other side of the power cord into an outlet 110/230V – 50/60Hz.
- The blue pilot light **POWER** lights up confirming the equipment is powered.
- The USB port, located next to the blue pilot light is now powered allowing the connection of accessories such as Tablet or Smartphone. The maximum current is 3A.

## Timer

- ➔ The timer is delivered from factory with a battery already installed. To activate the battery, remove the transparent tab on the back of the timer. The display lights up.

### PROGRAMMATION

- Press M and S simultaneously to reset.
- Press M until the desired number of minutes is displayed.
- Press S until the desired number of seconds is displayed.
- If the desired number is exceeded, reset to zero.

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# Operation of the filtration manifold

## THE WASTE BOTTLE

- The waste bottle has a maximum capacity of 1 000 ml. **It is imperative to empty it completely at the end of each round of analyzes.**
- To empty it, first hold the bottle body in place while unscrewing the blue ring from the cap. Once this blue ring rotates freely, remove the cap without disconnecting the tubes and remove the bottle from its original place.
- Empty its content according to the best practices of your institution.
- Once the bottle is empty, replace it in its original place and screw again its cap. Screw the blue ring until the end-stop.
- The output of this bottle is equipped with an antibacterial filter. The filtration manifold must never be used without this filter.

### • IMPORTANT


The output of this bottle is equipped with an antibacterial filter. The filtration manifold must never be used without this filter.

## FILTRATION STATIONS

The equipment is provided with 6 filtration stations numbered from 1 to 6.


The **Legio EZ** analysis protocols require the use of reusable « magnetized funnels and membrane supports», supplied with the equipment, and afterwards the use of «disposable filtration units» supplied with the reagents kits.

Each filtration station is composed as follows:


To transport the equipment or in the case where one or more filtration stations are not used during an analysis, make sure to respect the following configuration :

**Cap screwed on the hole - Membrane support screwed on the plug - Magnetic funnel correctly placed on the membrane support.**

## USAGE OF MAGNETIZED FUNNELS AND MEMBRANE SUPPORTS

- Unscrew and remove 1 to 6 caps according to the number of analyzes desired.
- Always use the leftmost holes first, starting with No. 1.
- Unscrew the membrane supports from the caps you just removed and screw them with their funnels attached on each hole, **tightening without excess**.
- Remove the magnetized funnels by pulling them upwards.
- Place the membranes well centered and well oriented as indicated below. The tongue of the membrane should form an «L».
- Then place each funnel back on the membrane support.
- Proceed in reverse to remove the membranes.

### ATTENTION

After use, the membranes and funnels must be thoroughly **RINSED** with clean water and **DRIED**.


## USAGE OF THE DISPOSABLE FILTRATION UNITS

- Unscrew and remove the membrane supports with their funnels, screw them onto their respective caps and set them aside.
- Position the disposable filtration units by placing them in the center of the holes, then turning them a quarter of a turn while pressuring downward, until they are at the bottom end-stop.

## FILTRATION

- After having poured the liquid to be filtered into the funnel of each filtration unit, start the filtration by setting the FILTRATION power switch to “I”.
- A pressure gauge indicates the vacuum level.
- Stop the pump by pressing “0”, as presented in the analysis protocols instructions for use.
- The suction circuit is automatically returned to ambient pressure in a few seconds.
- Proceed the same way for the next filtrations.

### • IMPORTANT

Make sure the pressure gauge needle got back to “0” before pouring any other liquid into the funnels.


## REVELATION

At the end of the analysis, the protocol requires the installation of a collector in the filtration hole. This collector is delivered with each Legio EZ reagents kit.

Collector

- Take the collector and remove its cap.
- Then remove the filtration unit from its hole, place the collector inside the hole and replace the filtration unit in its original position.
- During the next filtration, the revealing liquid will be gathered directly in the collector. Once this operation is complete, after stopping the filtration, remove the filtration unit.

➡ **RESULTS** : please refer to the analysis protocol instructions for use. ⬅

## RINSING AND PURGING THE FILTRATION MANIFOLD

The filtration manifold must be rinsed and purged just after the end of the analysis.

- To do this, remove all filtration units and caps, and start the filtration.
- Then pour around 25 ml of drinking water into each hole. Once the holes are empty, keep the filtration running, and screw the caps one by one by screwing them without excess. Start with the 1, then the 2, etc.... until the 6.
- Once all caps are screwed, stop the filtration.
- Then, empty completely the waste bottle, according to the best practices of your institution. To empty the waste bottle, please refer to the instructions detailed in the “Waste bottle” part.

### • IMPORTANT

The equipment must not be transported with any liquids.

# Incubation

## INCUBATION PROCESS

The incubator allows a maximum of 6 Petri dish (delivered with each Legio EZ reagents kit).

- To do this, open the incubator and remove the Petri dish holder.
- Place the Petri dish(es) into the Petri dish holder upside down, black face up, beginning with the lowest position: dish N°1 at the bottom position, dish N°2 at the penultimate position etc...
- Then, turn the Temperature recorder on : press its button until both LEDs flash 3 times fast (see the recorder manual for further details).
- Replace the Petri dish holder into the incubator as shown on the following illustration and close it.
- Start the incubation, positioning the incubator switch according to the type of water you are testing (“I” for Domestic Hot Waters / Drinking Waters or “II” for Cooling Towers Waters).

### • IMPORTANT

Since the Domestic Hot Waters / Drinking Waters or Cooling Tower Waters incubation programs are specific to each type of water, it is not possible to simultaneously test Domestic Hot Waters / Drinking Waters and Cooling Tower Waters. All your samples must be of the same type at the risk of compromising the results of the analyzes.

The red pilot light indicates the start of the incubation, which will last 47 hours.

The red pilot light switches off as soon as the incubation cycle is over. Then, position the switch back to “0” and open the incubator to get back the Petri dishes and continue the analysis protocol.

## WARNING

The incubation cycle starts when the power switch is positioned on one of the two programs.

- Positioning the switch to “0” will stop the cycle.
- Positioning the switch to of the two programs (“I” or “II”) will relaunch a new incubation cycle of 47 hours.
- Any power failure will stop the incubation cycle. A new 47-hours incubation cycle will start as soon as the power returns.

**IT IS THEREFORE IMPORTANT TO CONNECT THE EQUIPMENT TO A RELIABLE CIRCUIT PROTECTED FROM ACCIDENTAL POWER FAILURES.**


- You are required to ensure the equipment is not likely to be unplugged by a third party.

## • IMPORTANT

Opening the incubator during the analysis cycle may compromise the results. Please make sure the incubation program went well by getting the temperature data from the recorder, and by comparing them to the corresponding temperature profile (further details in the Temperature recorder instructions for use).


# Colorimeter

## SECURITY

Although rugged, the colorimeter may be damaged if a substantial spill of liquid comes into contact with the electronic components.

Cuvette

### • IMPORTANT

Avoid spills, use the caps provided with the cuvettes.


- Wipe gently, without scratching, the outside of the cuvette before placing it in the colorimeter.
- In the event of a slight flow, the liquid is collected in the drain pan located at the bottom of the cuvette holder and can be easily cleaned with a cotton swab or a sheet of absorbent paper previously folded.

### WARNING


- The colorimeter readings may be affected by temperature changes. Avoid exposing the device to the sun's rays.
- Avoid working in an environment subject to significant temperature fluctuations.

## BUTTONS


- The power button turns the device on or off.
- When used stand-alone, the device turns off after 2 minutes to save battery power.

## CAL

The **CAL** button calibrates the device at 0.000 (absorbance).

## SCREEN

The screen displays absorbance values from 0.000 to 1.999.

### Non-numeric readings

- The display can not display negative values or values above 1.999 (absorbance).
- The following non-numeric displays are possible :

729

Absorbance is negative, i.e. less than the absorbance of the solution used for calibration. **Recalibrate** the device using water or the least absorbent solution.


H<sub>1</sub>

➔ Absorbance is greater than 1.999.

bAt

➔ The batteries are weak and need to be replaced. This event is normally very rare.

## POWER


- The device is powered 6V by four AA batteries, inserted by removing the compartment screw.
- When the batteries need to be replaced, the screen displays : **bAt** at startup.
- Battery life with AA 1000 mAh batteries is estimated at 100 hours of continuous use.

**WARNING**

**The replacement of the batteries must occur very rarely.**

## CUVETTES

- The samples to be measured are placed in standard 4.5 ml plastic cuvettes. The light path is 10 mm long. Satisfactory results can be obtained with volumes of 2ml (half-filled cuvette).
- Some cuvettes have four transparent faces, but many models have two transparent faces and two frosted or ribbed faces.
- For reading correctly, position the cuvette in the colorimeter so that the transparent faces are oriented left and right, and not front and back.

## USE OF COLORIMETER

### Calibration

└ Calibration is performed by adjusting the absorbance to «zero» with clear, limpid water.

- Press the power button 
- Place a cuvette containing water in the colorimeter.
- Make sure that the transparent faces are oriented right and left.
- Press **CAL**. The display shows **CAL** for a few seconds, then 0.000 (A).

### Taking readings

- └ Replace the « blank » with a cuvette containing the liquid to be measured, close the lid and read the value on the screen.
- In stand-alone mode, the colorimeter turns off automatically after 2 minutes to conserve battery life, while keeping calibration in mind.

### After use

└ To turn off the device, press the power button.

Check the inside of the cuvette holder and wipe off any traces of moisture with a sheet of pre-folded paper towel.

The colorimeter was subjected to a series of tests that established its compliance with the following standards :

- EN 60950-1: 2006+A11:2009+A1:2010+A12:2011
- EN 301 489-01 V1.9.2: 2011
- EN 301 489-17 V2.2.1: 2012
- EN 300 328 V1.8.1: 2012

The test certificates can be consulted on the website [www.mystrica.com](http://www.mystrica.com).

FCC ID: 2AC8AMYS2  
Copyright 2016 Mystrica Ltd

# Temperature recorder

MX100 Bluetooth Recorder Quick

Mounting Loop

Alarm LED

Status LED

Start Button

Your equipment has been delivered with an MX100 temperature recorder to ensure that the incubation cycle has gone perfectly during the entire incubation period :

No power failure, no inadvertent opening of the incubator . This recorder is attached to the Petri dishes holder.

We already configured this temperature recorder according to the parameters below so its implementation get facilitated :

- Name of the recorder: Serial number of your equipment (SN: 8102XXX)
- Recording start: " When pressing the start button "
- Recording stop : " Automatic after 3 days "
- Recording interval : " every 2 minutes"
- Memory Saturation : " Never " (loop recording)
- Alarms: " Off "
- Led Display: " On "
- Bluetooth always off: " Off ".

➔ **Thus configured, your recorder device will have a duration of more than 2 years.** ←

## RECORDING START

Press the start button the time indicated below until the 2 LEDs (green and red) flash 3 times quickly :

- 4 seconds to start a new recording.
- 10 seconds to reset and restart a recording in progress.

## RECORDING STOP

THE RECORDING WILL STOP AUTOMATICALLY AFTER 3 DAYS.

## VISUALIZATION OF THE INCUBATION CYCLE

- 1 To perform the first visualization, you will need to download the HOBOMobile® application on your tablet or smartphone (Apple Store for IOS - Play Store for Android)  
If the application is already loaded, skip to step 2.
- 2 Open HOBOMobile®. Activate Bluetooth® on your device if requested by the application.
- 3 Press the "HOBO" icon at the bottom of the screen.
- 4 Then select your recorder to connect to it. If it does not appear, or if you detect several recorders, briefly press the recorder button to make it appear at the top of the list. Make sure you are in an accessible signal range (30m in free field).
- 5 Once connected, press " Play " to collect the data. After the message " Successful Reading ", press " OK ".
- 6 At the bottom of the screen, select " Data file ". After the message " Close the connection ? Press " Yes ".
- 7 Temperature curves are displayed on the screen, from the most recent to the oldest.
- 8 The incubation was successful if your curve corresponds to the following typical profile :

➔ For the " Domestic Hot Waters/ Drinking Waters " incubation program:

A rise in temperature to about **37°C** then a stabilization, for a **total duration of 47 hours**.

➔ For the " Cooling Towers Waters " incubation program:

A rise in temperature to about **52°C** then down to **37°C** and stabilization for a **total duration of 47 hours**.

Any other form of profile will indicate that the incubation has been disturbed by external elements (power failure, opening of the incubator during the incubation, ...). If so, the analysis results are not interpretable.

## CAUTION

**Do not open, incinerate, heat above 85°C (185°F), or recharge lithium batteries. Batteries may explode if the recorder is exposed to extreme heat, or conditions that could damage or destroy the battery envelope. Do not place the recorder or batteries in a fire. Do not expose the batteries to water. Recycle batteries according to local recommendations.**

DOMINIQUE DUTSCHER SAS

# Care and maintenance

In normal conditions of use, this equipment requires very few maintenance. You are required to keep the equipment clean. The filtration circuit must be rinsed, purged at each analysis as indicated in the chapter “Rinsing and purging the filtration manifold”.

- If any liquid is spilled on the equipment, wipe it up immediately with a damp sponge and let it dry.
- You are required to replace the antibacterial filter every 50 tests or every 36 months (first term expired).

## ATTENTION

This filter (Ref: 00 088) has a direction of connection. The “IN” side must be connected to the waste bottle side.

- If the equipment is not used for more than one week, leave the carrying case open with the bottle open and the incubator open for a few hours to avoid stagnant moisture..
- In the event of a malfunction, refer to the table below.

Any intervention on the equipment, other than those listed in this manual are prohibited, **END THE GUARANTEE OF THE EQUIPMENT AND COMPROMISE THE RESULTS OF ANALYSIS.**

# Malfunctions

## ➔ POWER UP : THE PILOT LIGHT POWER DOES NOT SWITCH ON

POSSIBLES CAUSES	SOLUTIONS
The power cord is not correctly connected to the carrying case	Check that the power cord is correctly connected on the rear of the carrying case.
The power cord is not correctly connected into the 110/230V – 50/60Hz outlet	Check that the power cord is correctly connected into the 110/230V – 50/60Hz outlet
No power into the outlet	If possible, check the outlet by connecting another equipment.
Damaged power cord	Unplug the power cord from the outlet, check the integrity of the power cord, and replace the power cord if it is damaged.
Equipment Failure	Contact technical support.

## ➔ FILTRATION : THE PUMP DOES NOT START (NO SOUND)

POSSIBLES CAUSES	SOLUTIONS
The carrying case is not powered	Check that the pilot light POWER switches on. Otherwise, refer to the table above.
The switch is not correctly pressed	Press « 0 », then « I » correctly.
Equipment Failure	Contact technical support.

## FILTRATION : THE PUMP STARTS, BUT THE FILTRATION IS INSUFFICIENT OR NON-EXISTENT

(The pressure gauge indicates a depression less than 0.4 bars)

POSSIBLES CAUSES	SOLUTIONS
Air leak	Check that the caps are correctly screwed on the unused holes.
Air leak	Check that the membrane supports are correctly screwed, and that the funnels and the membranes are correctly positioned.
Air leak	Check that each filtration unit is correctly placed on its respective hole.
Air leak	Check that each filtration unit has its cap correctly positioned on its side.
Air leak	Check that the cap of the waste bottle is correctly screwed.
Air leak	Check that the two tubes are correctly connected to the filter and are not damaged.
Internal air leak	Contact technical support.

## FILTRATION : THE PUMP STARTS, BUT THE FILTRATION IS INSUFFICIENT OR NON-EXISTENT

(The pressure gauge indicates a depression higher than 0.4 bars)

POSSIBLES CAUSES	SOLUTIONS
Antibacterial filter plugged	Replace the filter by a new one.
Clogged membrane(s)	<p><b>The event occurs at the beginning of the analysis</b> The water sample is too heavily loaded with impurities, pre-filter your sample on a 10 µm filter.</p> <p><b>The event occurs after incubation</b> The sample is too charged with microorganisms, contact technical support.</p>

## INCUBATION : THE RED PILOT LIGHT OF THE INCUBATOR DOES NOT LIGHT UP

POSSIBLES CAUSES	SOLUTIONS
The carrying case is not powered	Check that the pilot light POWER is on. Otherwise, refer to the table above.
The incubation switch is not correctly pressed	Press « 0 », then « I » or « II » according to the desired incubation program.
Equipment Failure	Contact technical support.

## INCUBATION : INCUBATION IS STILL RUNNING AFTER 47 HOURS

POSSIBLES CAUSES	SOLUTIONS
Unfinished incubation (Incubation time is about 47h)	Wait one hour before moving on to the next possible cause.
There was a power failure that interrupted the incubation and started it back from the beginning	The analysis result can be compromised. Check the profile of the temperature curve (temperature recorder). Ensure to connect the case to a reliable circuit.
The power switch was manipulated during the incubation	The incubation cycle has restarted from the beginning. The result of the analysis can be compromised. Check the profile of the temperature curve (temperature recorder). Ensure to protect the access (closed room, or lock of the carrying case).

# DECLARATION OF CONFORMITY

We hereby : **C4Hydro SAS**  
**Grand Luminy Technopole**  
**Case 922 – 163 Avenue de Luminy**  
**13009 Marseille, FRANCE**

Declare that the hereunder specified equipment :

Product designation : **MULTI-TESTS FIELD EQUIPMENT**

Model : **Field Equipment**

Product type : **MT**

Serial # : **810XXXX**

Complies with the following European Directives :

**Directive 2006/42/CE (Machinery)**

**Directive 2014/35/UE (Low Voltage)**

**Directive 2014/30/UE (Electromagnetic Compatibility)**

And conforms with the following harmonized standards :

**EN 12100 (Machinery)**

**EN 60204-1 (LVD)**

**EN61000-6-1:2007; EN61000-6-3:2007 (EMC)**

Person authorized to compile the technical file:

**M. Eric Agostini**  
**Grand Luminy Technopole**  
**Case 922 – 163 Avenue de Luminy**  
**13009 Marseille, France**

Place of issue : **Marseille**

Date of issue : **January 8<sup>th</sup>, 2018**

Signatory: **M. Eric Agostini**

Quality of signatory : **General Manager**

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