DOMINIQUE DUTSCHER SAS

Legio

 EZ-Test™

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Legio

 EZ-Test™

**INSTRUCTIONS FOR USE
UNI-TEST FIELD EQUIPMENT**

Summary

1.

■ Introduction p.18

2.

■ Description of the equipment p.18

3.

■ Timer p.19

4.

■ Filtration p.19 to 20

5.

■ Révélation p.21

6.

■ Incubateur p.22

Installation of the Petri dish..... p.22

Start the incubator..... p.23

Temperature Adjustment..... p.23

Stop the incubator..... p.24

7.

■ Colorimeter p.31 to 35

8.

■ Malfunctions p.35 to 37

• IMPORTANT

This equipment is intended only for environmental analysis of Hot Domestic Waters and Drinking Waters.

For technical assistance, please contact the technical support.

Introduction


This case contains all the equipments to carry out the required analysis with **Legio EZ-Test** reagents kits dedicated to Hot Domestic Waters and Drinking Waters.

Recommendations for an optimal use :

You are required to carefully read the attached user manual before carrying out the analysis. If in doubt, contact the technical support. The equipment is safe in normal conditions of use.

Description of the equipment

THE CARRYING CASE 1 CONTAINS :

Timer

The timer is delivered from factory with a battery already installed. To activate the battery, remove the transparent tab on the back of the timer. The display lights up.

PROGRAMMING :


- Press M and S simultaneously to reset.
- Press M until the desired number of minutes displayed.
- Press S until the desired number of seconds displayed.
- If the desired number is exceeded, reset to zero.
- The timer is now set to the desired time. At the selected time, start the countdown by pressing START/STOP.
- The timer sounds when the time is up.
- Press START/STOP to stop the alarm.


Filtration

FILTRATION BEFORE INCUBATION :

The filtration before incubation is carried out using the vacuum pump, the bottle, the membrane support, the magnetized funnel, and the membrane delivered with the **Legio EZ-Test** reagents kit.

- Remove the cap from the bottle, take the membrane support and position it on the bottle by pushing it to the end-stop.
- Remove the cap from the bottle, take the membrane support and position it on the bottle by pushing it to the end-stop.


- Position the assembled equipment on a stable surface. Connect the power supply connector to the pump connector.
- Connect the power supply to a 110/230V outlet.

At the beginning of the analysis, you will be asked to position the free membrane (supplied with the disposable filtration unit) in the center of the membrane support (see detailed information in the **Legio EZ-Test** analysis protocol instructions for use).

Once the membrane is in place, carefully position the magnetized funnel on the membrane support.

• **IMPORTANT**

The filtration must be carried out on a plan clean and horizontal surface. The vacuum pump shall be started only indoors.

- To start the filtration, start the pump by pressing ON.
- Press OFF to stop the filtration of the vacuum pump.

ATTENTION

After use, the membrane support, the magnetic funnel and the bottle must be disconnected and RINSED thoroughly with clear water and DRIED.

FILTRATION BEFORE INCUBATION :

In the second part of the analysis, the filtration is carried out using the vacuum pump, the bottle, and the disposable filtration unit supplied with the Legio EZ reagents kit.

- Firmly screw the bottle under the disposable filtration unit.
- Connect the suction tube to the disposable filtration unit by removing its cap.

- You are required to check the filtration unit is stable on the work plan.
- To start the filtration, start the pump by pressing ON.
- Press OFF to stop the filtration of the vacuum pump.

Revelation

- The analysis protocol requires the installation of a collector delivered with each **Legio EZ** reagents kit.
- Take the collector and remove its cap.
- Unscrew the disposable filtration unit from the bottle and place the collector inside the bottle.
- Then screw the disposable filtration unit on the bottle.
- During the next filtration, the revealing liquid will be gathered directly in the collector.
- Once this operation is complete, after stopping the filtration, remove the disposable filtration unit.
- At this point, the disposable filtration unit can be discarded.
- Get the collector using the forceps and reseal it.

• IMPORTANT


When you will discard the content of the bottle at the end of the analysis protocol, you are required to carefully check the internal ring is in its original position. Otherwise, you won't be able to install the collector for the next analysis.

Incubator

The incubator allows to hold the Petri dish (delivered with each **Legio EZ-Test** reagents kit) at a temperature of 37°C during 47 hours.

IT IS COMPOSED OF :

- Display
- Two pilot lights : WORK and SET
- Three buttons : SET, + and -

• IMPORTANT

The incubator shall be started only :


- Indoor
- On an horizontal, stable and dry surface

The incubator must not operate for more than 1 minute without its cover.

If necessary, stop the incubator and restart it only when its cover is replaced.

INSTALLATION OF THE PETRI DISH

- Remove its cover of the incubator.
- Hold the Petri dish, **TURN IT UPSIDE DOWN** and slide it into the holder, its cover facing downwards as shown in the diagram below.
- Replace the cover on the incubator.

START THE INCUBATOR

- Connect the power cord to the power connector of the incubator.
- You are required to check the power cord fits perfectly.
- Connect the other side of the power cord into an outlet 110/230V.
- At the connection, the incubator fan switches on, the display shows the actual temperature under the cover.
- The set-up temperature should be 37°C. The incubator is delivered already set at 37°C by default. There is therefore no adjustment to be made.
- When starting the process, the WORK pilot light (left pilot light) switches on indicating that the heating is active.
- The WORK pilot light switches off as soon as the temperature of 37°C is reached, and switches back on again as soon as the temperature has dropped below 37°C.
- This alternative mode of operation is normal.

TEMPERATURE ADJUSTMENT

- If the temperature has been accidentally de-adjusted, it can be adjusted by pressing SET.
- The SET pilot light (right pilot light) switches on indicating the setting mode is active: the display screen shows the setting temperature.
- Adjust the temperature by pressing + and -. Once the temperature is set, press SET or wait 5 seconds to exit the setting mode.
- The SET pilot light switches off when you exit setting mode. Even disconnected, the incubator keeps the last temperature setting in memory.
- The temperature is regulated only by heating. There is no cooling system.

• IMPORTANT

The ambient temperature must be permanently higher than 15°C and lower than 37°C.

STOP THE INCUBATOR

- Stop the incubator by disconnecting the power cord from the outlet 110/230V.
- There is no need to disconnect the power cord from the incubator power connector if you do not have to carry the incubator.


• IMPORTANT

The incubator does not stop automatically.
You must be present 47 hours after having started your incubation to stop the cycle manually.

DOMINIQUE DUTSCHER SAS


Colorimeter

SECURITY

Although rugged, the colorimeter may be damaged if a substantial spill of liquid comes into contact with the electronic components.

• IMPORTANT

Avoid spills, use the caps provided with the cuvettes.


- Wipe gently, without scratching, the outside of the cuvette before placing it in the colorimeter.
- In the event of a slight flow, the liquid is collected in the drain pan located at the bottom of the cuvette holder and can be easily cleaned with a cotton swab or a sheet of absorbent paper previously folded.

WARNING

- The colorimeter readings may be affected by temperature changes. Avoid exposing the device to the sun's rays.
- Avoid working in an environment subject to significant temperature fluctuations.

BUTTONS


- The power button turns the device on or off.
- When used stand-alone, the device turns off after 2 minutes to save battery power.

CAL

The **CAL** button calibrates the device at 0.000 (absorbance).

SCREEN

The screen displays absorbance values from 0.000 to 1.999.

Non-numeric readings

- The display can not display negative values or values above 1.999 (absorbance).
- The following non-numeric displays are possible :

729

Absorbance is negative, i.e. less than the absorbance of the solution used for calibration. **Recalibrate** the device using water or the least absorbent solution.


H₁

➔ Absorbance is greater than 1.999.

bAt

➔ The batteries are weak and need to be replaced. This event is normally very rare.

POWER


- The device is powered 6V by four AA batteries, inserted by removing the compartment screw.
- When the batteries need to be replaced, the screen displays : **bAt** at startup.
- Battery life with AA 1000 mAh batteries is estimated at 100 hours of continuous use.

WARNING

The replacement of the batteries must occur very rarely.

CUVETTES

- The samples to be measured are placed in standard 4.5 ml plastic cuvettes. The light path is 10 mm long. Satisfactory results can be obtained with volumes of 2ml (half-filled cuvette).
- Some cuvettes have four transparent faces, but many models have two transparent faces and two frosted or ribbed faces.
- For reading correctly, position the cuvette in the colorimeter so that the transparent faces are oriented left and right, and not front and back.

USE OF COLORIMETER

Calibration

└ Calibration is performed by adjusting the absorbance to «zero» with clear, limpid water.

- Press the power button 
- Place a cuvette containing water in the colorimeter.

- Make sure that the transparent faces are oriented right and left.
- Press **CAL**. The display shows **CAL** for a few seconds, then 0.000 (A).

Taking readings

└ Replace the « blank » with a cuvette containing the liquid to be measured, close the lid and read the value on the screen.

- In stand-alone mode, the colorimeter turns off automatically after 2 minutes to conserve battery life, while keeping calibration in mind.

After use

└ To turn off the device, press the power button.

Check the inside of the cuvette holder and wipe off any traces of moisture with a sheet of pre-folded paper towel.

The colorimeter was subjected to a series of tests that established its compliance with the following standards :

EN 60950-1: 2006+A11:2009+A1:2010+A12:2011
EN 301 489-01 V1.9.2: 2011
EN 301 489-17 V2.2.1: 2012
EN 300 328 V1.8.1: 2012

The test certificates can be consulted on the website www.mystrica.com.

FCC ID: 2AC8AMYS2
Copyright 2016 Mystrica Ltd

Malfunctions

FILTRATION : THE PUMP DOES NOT START (NO SOUND)

POSSIBLE CAUSES	SOLUTIONS
The pump is not connected	Check that the pump is correctly connected to the power supply unit, and check the power supply unit is correctly connected to an outlet 110/230V. Check the outlet is correctly powered.
The button is not correctly pressed	Press « OFF » correctly, then « ON ».
Equipment failure	Contact technical support.

FILTRATION : THE PUMP STARTS, BUT THE FILTRATION IS INSUFFICIENT OR NONEXISTENT

POSSIBLE CAUSES	SOLUTIONS
Air leak	Check that the membrane support is pressed firmly into the bottle. If necessary, remove it and put it back on correctly. Check that the waste bottle is correctly screwed under the disposable filtration unit. If needed, unscrew it completely and screw it back again.
Air leak	Check that the connector is correctly connected to the filtration unit (or the membrane support).
Air leak	Check that the tube connecting the pump to the disposable filtration unit/membrane support is not damaged. Change it if necessary.
Internal air leak	Pinch firmly the tube connecting the pump to the disposable filtration unit/membrane support. If the sound of the pump does not change, contact technical support.
Clogged membrane(s)	<p>The event occurs at the beginning of the analysis : The water sample is too heavily loaded with impurities, pre-filter your sample on a 10 µm filter.</p> <p>The event occurs after incubation : The sample is too charged with microorganisms, contact the technical support.</p>

INCUBATION : THE DISPLAY OF THE INCUBATOR DOES NOT LIGHT UP

POSSIBLE CAUSES	SOLUTIONS
The incubator is not powered	Check that the power cord is correctly connected to the incubator on the one hand, and to an outlet 110/230V on the other hand. Check that the outlet is correctly powered.
Equipment failure	Contact technical support.

INCUBATION : THE DISPLAY OF THE INCUBATOR LIGHTS UP BUT THE INCUBATOR FAN DOES NOT SWITCH ON

POSSIBLE CAUSES	SOLUTIONS
Equipment failure	Contact technical support.

INCUBATION : THE DISPLAY LIGHTS UP BUT THE INCUBATOR IS NOT HEATING

POSSIBLE CAUSES	SOLUTIONS
Equipment failure	Contact technical support.

INCUBATION : CONDENSATION APPEARS INSIDE THE COVER OF THE INCUBATOR

CAUSES POSSIBLES	SOLUTIONS
Normal phenomenon	Continue the incubation.

DECLARATION OF CONFORMITY

We hereby : **C4Hydro SAS**
Grand Luminy Technopole
Case 922 – 163 Avenue de Luminy
13009 Marseille, FRANCE

Declare that the hereunder specified equipment :

Product designation : **UNI-TEST FIELD EQUIPMENT**

Model : **Field Equipment**

Product type : **UT**

Serial # : **810XXX**

Complies with the following European Directives :

Directive 2006/95/CE (Low Voltage)

Directive 2014/30/UE (Electromagnetic Compatibility)

And conforms with the following harmonized standards :

EN 60204-1 (LVD)

EN61000-6-1:2007; EN61000-6-3:2007 (EMC)

Person authorized to compile the technical file:

M. Eric Agostini
Grand Luminy Technopole
Case 922 – 163 Avenue de Luminy
13009 Marseille, France

Place of issue : **Marseille**

Date of issue : **January 8th, 2018**

Signatory: **M. Eric Agostini**

Quality of signatory : **General Manager**

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