DOMINIQUE DUTSCHER SAS

Legio

PROTCOLE D'ANALYSE LEGIO EZ-COUNT

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Legio

 EZ-Count™

**INSTRUCTIONS FOR USE
ANALYSIS PROTOCOL**

DOMINIQUE DEBISCHER SAS

Summary

1.

Introduction p.28 to 29

- Technology used
- Recommendations for an optimal use

2.

Contents of the reagents kit Legio EZ-Count p.30

3.

Preparations before analyses p.31 to 32

4.

Step by step procedure p.33

Step A.....	p.33 à 34
Step B.....	p.35
Step C.....	p.36
Step D.....	p.36
Step E.....	p.37
Step F.....	p.38
Step G.....	p.39
Step H.....	p.40
Step I.....	p.41
Step J.....	p.41
Step K.....	p.42
Step L.....	p.43
Step M.....	p.43

5.

Results p.44 to 45

Introduction

Legio EZ-Count * is a simple and fast kit which detects only bacteria *Legionella pneumophila* capable of being cultivated (as well as regulatory procedure NF T90-431) in Cooling Towers Waters samples..

Legio EZ-Count identifies in 48 hours Cooling Towers Waters contain less than 1 000 CFU/L, between 1 000 and 100 000 CFU/L or more that 100 000 CFU/L of *L. pneumophila*.

Technology used

This solution uses a patented innovative technology** which combines:

- (i) culturing on selective agar GVPC of the sample to be analyzed,
- (ii) assimilation only by the *L. pneumophila* of a modified sugar,
- (iii) revelation of this modified sugar present after the assimilation in the part O-antigenic of lipopolysaccharide of *L. pneumophila*, through a reaction based on « click » chemistry for realizing colorimetric and enzymic reaction.

Recommendations for an optimal use

You are required to carefully read the protocol before carrying out the analysis. If in doubt, contact the service support.

- The protocol is safe in normal condition of use. As the water sample may potentially contain *L. pneumophila*, you are required to take the usual precautions throughout the analysis.
- Avoid any contact with skin, eyes and splashing.

- Keep the reagents kits away from children.
- Store the **Legio EZ-Count** reagents kit between +4°C and +8°C (ideally at +4°C) upon reception. A punctual exposition at ambient temperature (especially during transport, ...) is not detrimental.
- Do not expose to temperatures lower than 0°C.
- The **Legio EZ-Count** reagents kit is exclusively intended to be used with the following equipment:
 - Multi-tests Field Equipment
 - Filtration Manifold
- You are required to carefully read the instructions for use of the equipment you choose to use before beginning the analysis.
- The compliance of the protocol, of the indicated times and of the expiry date appearing on the box guarantee the performance of the kit.
- Do not expose the reagents to the sun.
- The analysis must be conducted quickly after taking your water sample.
- The water sample must not be kept cool before performing the analysis.
- At the conclusion of the analysis, remove the waste in accordance with the best practices from your establishment.

***Legio EZ-Test is a registered trademark of the company Diamidex.*

***Identification of living Legionella pneumophila using species-specific metabolic lipopolysaccharide labeling. Mas Pons J, Dumont A, Sautejeau G, Fugier E, Baron A, Dukan S, Vauzeilles B. Angew Chem Int Ed Engl. 2014 Jan 27;53(5):1275-8.*

• **IMPORTANT**


Use this kit only for environmental analyses of Drinking Waters and Hot Domestic Waters.

For technical assistance, please contact the technical support.
Contact information on the back of this leaflet

Contents of the reagents kit Legio EZ-Count

Box containing :

Each **Legio EZ-Count** reagent kit allows to carry out 1 analysis, for 1 Cooling Towar water sample.

To carry out new analyzes, please use new **Legio EZ-Count** reagents kits, each kit containing all the consumables mentioned above.


• IMPORTANT

From the receipt, store the Legio EZ-Count reagents kits in the refrigerator (+4°C and +8°C). The disposable filtration unit can be stored separately in a clean area, away from heat and moisture.

Disposable filtration unit

The Disposable Filtration Unit is a non-reusable consumable, used for the 2nd part of the protocol.

This consumable is supplied with each **Legio EZ-Count** reagents kit.

Preparations before analysis

IMPORTANT

The analysis should be carried out on a horizontal and dry surface.

The analysis protocol consists of 3 parts:

PART 1 : PREPARATION OF THE SAMPLE.

Step A to Step E : About 10 minutes of handling.

Incubation :

Step F : This part will require 47 hours of incubation. Be sure to take all the precautions indicated in the instructions for the equipment used, so that no interruption of the incubation cycle occurs, otherwise the results of the analysis can be compromised. Make sure you can **continue the analysis at the end of the 47 hours of incubation.**

PART 2 : REVELATION AND READING OF THE RESULT.

Step G to Step M : This part will require 50 minutes of handling.

• IMPORTANT


Please plan these steps carefully, taking into account their duration, in order to guarantee the results of the analysis.

[FILTER]

or

NB : In the rest of this protocol, the [FILTER] symbol will be associated with the following manipulations:

- Turn on the pump
- Wait until the liquid is completely sucked
- **Important** : Even if the liquid seems to have been completely drawn in, the pump must remain running for at least 40 seconds.
- Turn off the pump
- Once the pump is off, wait 10 more seconds before proceeding to the next step.

Indeed, you must be certain that the aspiration is complete before pouring another reagent. If the aspiration has not been completed, the next reagent may be aspirated before it can take effect.

The ways of turning on / off the pump are specific to the equipment chosen for the analyzes. Please refer to the corresponding instructions for use.

WARNING

In the case of a sample of water heavily loaded with impurities, we recommend prefiltering your sample with a 10µm filter.

If a filtration problem is encountered, refer to the «malfunctions» chapter of the chosen equipment instructions for use.


Step by step procedure

STEP A


- Take the membrane support and the magnetized funnel provided with the equipment chosen to perform the analysis.
- Position the membrane support on the equipment.
- Remove the **Legio EZ-Count** reagents kit from the refrigerator.
- Remove the disposable filtration unit from its bag while holding in place its cover.

NB : This part will not require the use of the disposable filtration unit, only the membrane located inside. The disposable filter unit will be retained for Part 2 of the protocol.

- Using the forceps, grab the free membrane of the disposable filtration unit by its lateral tab.
- ➔ **Do not grab it in any other place as the forceps could damage the membrane and compromise the results of the analysis.**

- Position the membrane on the membrane support so that:
 - It forms an « L » as shown on the diagram below.
 - It is perfectly centered on the membrane support, in the space provided for this purpose.


• IMPORTANT


The membrane must be centered and placed permanently on the same face throughout the analysis, otherwise the results of the analysis will be compromised.

- Position the magnetized funnel on the membrane support. Make sure that the magnets of the funnel and the membrane support engage properly.

- Put the cover back on the disposable filtration unit and store it in its bag.

- Pour the water sample into the collector until it is full (the volume is 25 ml, the water sample must be at room temperature).


- If you perform multiple analyzes simultaneously, be sure to use a different collector for each one of your samples, so they won't contaminate each-other and compromise the results of the analyzes.
- Then replace the collector in the **Legio EZ-Count** reagents kit, it will be used again at the end of the analysis protocol.

[FILTER]

STEP B


- Once the filtration is complete and the vacuum pump is stopped, open vial N°1 and slowly pour all of its contents into the funnel.
- Close surely and throw away vial N°1 in the plastic bag.
- Start the timer for 5 minutes.

During these 5 minutes, carry out step C.

STEP C

- Open the vial N°2 (pod) by tilting its cap on the side until it breaks.
- Open the vial N°3.
- By pressing on the vial N°2, pour 6 drops into the vial N°3.
- Close surely vial N°3.

➔ Keep vials N°2 and N°3 nearby. ←

- After the 5 minutes :

[FILTER]


STEP D

- Slowly pour the whole vial N°2 into the center of the funnel.

- Throw away vial N°2.

[FILTER]

STEP E

a/

- Remove the Petri dish from its pouch. A notch on the edge of the bag allows easy opening.
- Remove its cover from the Petri dish.

b/

- Do not touch the inside of the Petri dish to avoid contaminating it.
- Pour the whole vial N°3 into the center of the Petri dish.
- Do not spread the liquid.


• IMPORTANT

Do not move the Petri dish, the deposited liquid must remain in the center.

- Remove the magnetized funnel from the membrane support.


c/

- Using the forceps, grab by its lateral tab the free membrane located onto the membrane support.
- Carefully place it as shown in the diagrams opposite.
- Place the edge of the membrane at the edge of the drop previously deposited at the center of the Petri dish, so that the liquid is in contact with the entire surface of the membrane.
- Make sure there is no bubble underneath the membrane.
- With the tab still forming an « L ».

d/

- Replace the cover of the Petri dish. Keep the box perfectly horizontal.
- Apply the adhesive tape all around the box making 2 turns.
- Keep the forceps for further analysis.
- Turn the Petri dish upside down so that its cover is underneath (black face up).

STEP F

INCUBATE THE PETRI DISH ACCORDING TO THE EQUIPMENT USED:

EQUIPEMENT	INCUBATION
Multi-tests Field Equipment	Cooling Tower Mode
Laboratory Incubator	45 minutes at 52°C then 47 hours at 37°C

Put back the **Legio EZ-Count** reagents kit in the refrigerator.

CLEAN THE MEMBRANE SUPPORT AND THE MAGNETIZED FUNNEL AS DESCRIBED IN THE INSTRUCTIONS SUPPLIED WITH YOUR EQUIPMENT.

• **IMPORTANT**


Note the day and time of the beginning of the incubation and come back after 47 hours to carry out the 2nd part of the analysis.


47H

STEP G


- After 47 hours, remove the **Legio EZ-Count** reagents kit from the refrigerator.
- Open the vial N°6.
- Take precisely 5ml of the vial W by using the syringe provided with the equipments (avoid the presence of air in the syringe)

- **Carefully transfer these 5 ml into the vial N°6.**
- Close surely the vial N°6.


- Tap the bottle on the worktop, so that the liquid passes well under the powder.
- Reposition vial N°6 at its numbered location in the consumable box and let stand until Step K.


STEP H

- **TAKE THE DISPOSABLE FILTRATION UNIT.**
- Remove the Petri dish from the incubator.
- Place the disposable filtration unit on the equipment you are using to conduct the analysis.
- Remove its cover from the disposable filtration unit.
- Remove the adhesive tape and the cover from the Petri dish.
- Using the forceps, grab the free membrane by its lateral tab and place it into the center of the funnel.

a/

• IMPORTANT

At this point, and in most cases, the membrane will have few (or no) elements visible to the eye. If you notice many spots on the membrane, this may indicate that your water sample was heavily loaded with microorganisms and prevent the next filtration from working. If the following filtration does not work, please refer to the examples at the end of this leaflet.

• FOR REMINDER

Respect the diagram below: **the tongue forming an «L»**. The membrane must be perfectly **centered at the bottom of the funnel**.

- Throw away the Petri dish and the forceps in the plastic bag.
- Take **10 ml from the vial W** with the syringe.
- Empty the contents of the syringe into the center of the funnel.
- Gently rub the entire surface of the membrane with the spreader **for 30 seconds**.
- Throw away spreader in the plastic bag.

b/

c/

- Take a second time 10 ml from the vial W with the syringe.

- Empty the contents of the syringe into the center of the filtration unit.

[FILTER]

STEP I

a/

- The pump being turned off, make sure the pressure gauge needle is at «0» (if your equipment has a pressure gauge).
- Open the vial N°4 and pour all its content into the center of the funnel.
- Respect a waiting time of 5 minutes.
- After the 5 minutes : **[FILTER]** 

b/

- Take 10 ml from the vial W with the syringe.
- Empty the contents of the syringe into the center of the funnel.

[FILTER] 

STEP J

a/

- With the pump still off, open the N°5 vial and pour all its content into the center of the funnel.
- Start the timer for 1 minute.
- Close surely vial N°5 and throw it away in the plastic bag.

After 1 minute :

[FILTER] 

b/

- Take 10 ml from the vial W with the syringe.
- Empty the contents of the syringe into the center of funnel.

[FILTER] 

STEP K

a/

- Take vial N°6 and shake until completely dissolved (the solution should be homogeneous).

➔ To do this, keep the cap in place, and shake it gently upside down several times in a row.

- With the pump still off, pour the contents from vial N°6 into the center of the funnel.
- Close surely the filtration unit with its cover.
- Respect a waiting time of 30 minutes.
After these 30 minutes :

[FILTER] 

30 min


b/

- Take 10 ml from the vial W with the syringe.
- Empty carefully the contents of the syringe into the center of funnel.

[FILTER] 

- Take a second time 10 ml from the vial W with the syringe.
- Empty the contents of the syringe into the center of funnel.

[FILTER] 


(bis)


STEP L

- Pour the contents of the remaining vial N°2 (pod) into the center of the funnel.

[FILTER] 


STEP M

- Place the collector under the disposable filtration unit according to the type of equipment you are using for conducting the analysis.
- With the pump still off, pour the contents of the vial N°7 into the center of the funnel.
- Replace its cover on the disposable filtration unit.
- Respect a waiting time of 5 minutes.
- Close surely vial N°7 and throw it away into the plastic bag.
- After these 5 minutes:
[FILTER] 
- Respect a waiting time of 5 minutes again.

RESULTS

- Once the 5 minutes have passed, get the collector back using the forceps and transfer a portion of the liquid of the collector into the cuvette using the pipette provided. The cuvette must be filled by more than half to allow reading.
- Place the cuvette in the colorimeter, and read the absorbance.

• IMPORTANT

The cuvette must be set in a certain direction so that the absorbance reading can be performed correctly. Please refer to the operating instructions of the colorimeter supplied with your equipment.

- If the number displayed on the colorimeter is **less than 0.24** then the result is **NEGATIVE** (the water sample tested contains less than 1000 CFU/L of *Legionella pneumophila*).
- If the number displayed is **between 0.24 and 0.6** then the water sample tested contains **between 1 000 CFU/L and 100 000 CFU/L** of *Legionella pneumophila*.
- If the number displayed is **greater than 0.6** then the water sample tested contains **at least 100 000 CFU/L** of *Legionella pneumophila*.


END OF ANALYSIS, MAINTENANCE OF EQUIPMENT

- Once the analysis is complete, dispose of the waste (plastic bag containing the empty bottles, disposable filtration unit ... and liquid waste) according to the best practices of your institution.
- Clean the equipment in accordance with the maintenance recommendations detailed in the instructions for use supplied with any one of them.
- The syringe must be rinsed with clean water and dried before being stored in your equipment.


• IMPORTANT

No liquid should remain inside the equipment, which must be clean and dry before being transported and/or stored.

The filtration step after incubation did not work, the membrane presents a large number of spots. Case examples:

Presence of an interfering flora too important to interpret the concentration of *Legionella pneumophila* in your sample. In this type of case, most methods of detection and enumeration of *Legionella pneumophila* will give uninterpretable results.

WARNING

Your sample most likely contains a very large amount of *Legionella pneumophila* (> 100 000 CFU/L).

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

Ref : 00-388-A

Protocole d'Analyse Legio EZ-Count | Analysis Protocol Legio EZ-Count

© Diamidex 2019