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Fax +49.551.308.3289

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

Specifications subject to change without notice.

Printed and copyrighted by  
Sartorius Stedim Biotech GmbH  
W402.08 · G

Publication No.: SM-6047-p05085

## Directions for Use of Gelatin Filters

### 1. Packing Specifications

The sterile filters are sealed in units of five each in a polyethylene bag. Each filter is separately inserted in a pocket-like wrapper of special paper for added protection and convenience. One box contains 50 filters. They have been presterilized by gamma radiation.

#### Order numbers:

12602-037-ALK (37 mm diameter)  
12602-047-ALK (47 mm diameter)  
12602-050-ALK (50 mm diameter)  
12602-080-ALK (80 mm diameter)

### 2. Storage

Gelatin filters can be stored at temperatures ranging from 2°C to 35°C. Avoid exposing the filters to moisture and solvent vapors or to the vapors of other chemicals.

### 3. Labeling

The label on the front side of the box specifies the type, number per box, diameter, pore size and the lot number. Please state the lot number in any letter of inquiry you send to us.

### 4. How to Use

Water-soluble gelatin filters, 12602, have been specially designed for detection and analysis of airborne microbes. The air flow rate for these filters is approximately  $2.7 \pm 0.5$  l/min/cm<sup>2</sup> at a differential pressure of 0.05 bar (50 mbar). To remove a filter, cut open the polyethylene bag on the side where the fold of the pocket is located. Now you can remove each pocket by hand without risking contamination of the other filters. As you remove the filter, please avoid contaminating it by touching it with your fingers or allowing it to come into contact with foreign objects, for example. Handle it carefully, because gelatin filters are somewhat brittle by nature. To use, place each filter in a dry, sterile filter holder, such as the open-faced filter holder supplied for the Sartorius Stedim Biotech MD8 airscan Air Sampler (16746, 16747 or 16748) or MD8 AirPort (16757), by carefully letting the filter slide out of the pocket onto the filter support of the aluminum filter holders. Make sure to load the filter under sterile conditions (for example, in a laminar flow bench). If you use forceps to aid positioning of the filter, do not grasp the filter too tightly, and avoid bending it. After sampling, detach the top part of the aluminum filter holder from the base by turning it counterclockwise. Place the prepared, agar-filled petri dish base over the aluminum filter holder. The gelatin filter adheres to the agar surface so that the petri dish along with the gelatin filter can be lifted off. Then, immediately cover the petri dish base with the lid. The gelatin filter dissolves due to the moisture in the agar culture medium.

### 5. Incubation and Evaluation for Airborne Microbe Determination | Detection

After transferring the filter to the culture medium, incubate it in an incubator with the lid of the petri dish facing up. To prevent liquid from collecting on the agar surface, always use predried agar that is not freshly prepared. Choose the time, temperature and type of culture medium according to your microbiological application or test procedure; for instance, Standard, Caso or Plate Count Agar are suitable for determining the colony count (total CFU count); Sabouraud, Malt Extract or Wort Agar for detecting yeasts and molds; and Blood Agar for detecting pathogenic microbes causing hemolysis. By comparing the number of colonies in relation to the volume of air originally sampled, you will obtain the quantity of colony-forming units per cubic meter of air (CFU/m<sup>3</sup>).

## 6. Additional Application Notes

a. Dissolving the gelatin filters after sampling airborne microbes  
Gelatin filters can also be dissolved in sterile liquid warmed to 35–40°C, such as physiological saline or 0.1% peptone water. Stirring the solution using a sterile magnetic stirrer accelerates dissolving of the filter. Afterwards, process the solution according to Koch's pour plate method or the membrane filtration technique. During the dissolving and stirring process, the colony-forming units are separated into individual microbes so you will obtain a higher CFU count with this indirect method than that by directly placing the exposed filter on a culture medium. It is strongly recommended to dissolve the gelatin filter (indirect method) whenever you take samples in an area sprayed with disinfectants (or if antibiotic airborne particles are present), for instance. The reason: all traces of disinfectants, which could inhibit the growth of microbes on the culture medium, can be removed as follows: dissolve the gelatin filter, then filter this liquid through a 0.45 µm pore size membrane filter, for example, 11406; afterwards, add sterile water to the filter holder to rinse the filter.

b. Dissolving gelatin filters after sampling virus and phage aerosols  
To determine the airborne virus concentration, the gelatin filter must be dissolved in order to cultivate the vi-ruses and phages. For detailed information, please see the Application Notes "Collecting Airborne Viruses and Phages using the Sartorius Stedim Biotech Gelatin Membrane Filter Method" (publication no. SLF4028-e) as well as the two special reprints from the BioTec journal concerning the same subject (publication nos. FM-8030-e "Sampling Virus Aerosols using the Gelatin Membrane Filter" and FM-8034-e).

Further details will be provided at your request.

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen, Germany

Phone +49.551.308.0

Fax +49.551.308.3289

[www.sartorius-stedim.com](http://www.sartorius-stedim.com)

Specifications subject to change without notice.

Printed and copyrighted by  
Sartorius Stedim Biotech GmbH  
W402.08 · G

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