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DOMINIQUE DUTSCHER SAS

Table of contents	Page
Delivery extent	5
Warranty	5
Anaerobic jars	7
Cleaning	7
Operation	9
Maintenance	11
Ordering information	11
Accessories	12
Special equipment for connection to Anoxomat	15
Spare-parts	17

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DOMINIQUE DUTSCHER SAS

Delivery extent

- Anaerobic jar with or without valve/stopcock and manometer
- Operating instructions
- Warranty Card with reply card

Warranty

The manufacturer grants a one-year guarantee on the Anaerobic jars.

To establish your right for warranty coverage, send in the reply card accompanying the Warranty Card to the manufacturer within 14 days.

Device faults and malfunctions resulting from improper use, especially due to over pressure or autoclaving of the lid are not covered by the warranty.

Only original accessories and spare-parts supplied by the manufacturer may be used.

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Anaerobic jars

The Anaerobic jars made of stainless steel are characterised by a break-proof construction. The Anaerobic jars made of transparent polycarbonate (PC) allow for monitoring the growth of colonies optimally during incubation.

The lids are made of black plastic or transparent polycarbonate (PC).

The two corner valves with tube clip (5 mm inner-Ø) mounted to the lid allow for evacuating the jars up to -1 bar and the following filling with the requested gas for cultivation of anaerobic or microaerophilic microorganisms.

The manometer mounted to the lid permits control of the vacuum or over pressure inside the Anaerobic jar at any time. The over pressure may not exceed 0.2 bar as otherwise the tongs will deform, which may cause leakage of the lid.

The Anaerobic jar “eco” has been especially designed for operation with GasPack-Kits for chemical production of the gas mixture. It provides a ventilation screw at the lid for pressure compensation if a vacuum is reached inside the jar during incubation.

All Anaerobic jars may use GasPack-Kits for reaching anaerobic conditions by chemical reduction of oxygen. We recommend our racks, which provide holders for gas pack and catalyst.

Cleaning

The Anaerobic jars may be cleaned with commercial cleaning solutions or alcohol as well as other suitable disinfectants. The Anaerobic jars “crystal” may only be wiped with alcohol as longer residence times may expand the material and cause cracks.

Lids made of plastic with manometer are not suitable for temperatures exceeding 60 °C and may not be autoclaved (high-temperature manometers available on request).

Lids made of plastic without manometer may be autoclaved (Caution: Lids may not be installed on stainless steel container – danger of implosion!).

Lids made of Polycarbonat (PC) may not be autoclaved.

The stainless steel racks may be autoclaved.

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Operation

Upon inoculation, the Petri dishes or tubes are placed inside the optionally available rack.

For Anaerobic jar with valve/stopcock and manometer, methods 1-3 are suitable for reaching anaerobic conditions. For Anaerobic jar models "eco" as well as jars with lids not providing valves/stopcocks and manometer, only method 3 is suitable for reaching anaerobic conditions.

Method 1:

Evacuate and Fill with gas

- Insert the rack in the Anaerobic jar.
- If need is, place an anaerobiosis indicator in the holder at the lid*.
- Close the lid. Make sure the tongs are inserted in the ring groove at the lid's surface and clamp tightly.
- Connect the vacuum pump (ideal with <math><50\text{ mbar}</math> residual pressure) with a thick-walled hose to the tube clip of one corner valve at the lid.
- Connect a gas bottle containing the requested gas mixture to the second corner valve. The pressure reduction valve is to be set to max. 0.2 bar over pressure at the outlet side.
- Open the valve connected with the vacuum pump and extract the gas from the Anaerobic jar.
- Close the valve connected with the vacuum pump and open the valve connected with the gas bottle.
- Fill the Anaerobic jar with the requested gas mixture up to max 0.2 bar over pressure.
- Depending on the sensitivity of the microorganisms, the process of evacuation and filling with gas may be repeated once or twice.
- After the last filling with gas, close both valves and remove the hoses.
- Check the pressure at the manometer and if there is any over pressure left, open the valves to drain it.

Method 2:

Flow with gas

- Insert the rack in the Anaerobic jar.
- If need is, place an anaerobiosis indicator in the holder at the lid*.
- Close the lid. Make sure the tongs are inserted in the ring groove at the lid's surface and clamp tightly.
- Connect a gas bottle containing the requested gas mixture to one of corner valves. The pressure reduction valve is to be set to max. 0.2 bar over pressure at the outlet side. The other valve should be connected with a tube, which allows to transfer the gas (i.e. flue).
- Open both valves.
- Open the gas supply and flow the Anaerobic jar with the requested gas mixture.
- Flow the Anaerobic jar with the requested gas mixture up to max 0.2 bar over pressure.
- Depending on the sensitivity of the microorganisms, the Anaerobic jars may be flowed from 5 to 15 minutes.
- After flowing with gas, close both valves and remove the hoses.
- Check the pressure at the manometer and if there is any over pressure left, open the valves to drain it.

*Only with transparent lid (models "eco" and "crystal").

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Method 3:

Operation with gas generating (GasPack-kits)

- Follow the instructions of the manufacturer enclosed to the gas packs.
- If need is, place an anaerobiosis indicator in the holder at the lid*.
- If need is, regenerate the catalyst overnight at 150 °C and insert it in the top coil of the rack (only requested for gas packs of other manufacturers using hydrogen as reducer).
- Hold a sufficient number of gas packs available
(Anaerobic jar “large“ = 2 packs, all other Anaerobic jars = 1 pack).
- Load the rack with the samples.
- Open the plastic packing of the GasPack, take out the pack and insert it in the bottom coil of the rack.
- Insert the rack in the Anaerobic jar.
- Immediately close the lid. Make sure the tongs are inserted in the ring groove at the lid’s surface and clamp tightly.
- Inside the Anaerobic jar, the oxygen is now reduced or the requested gas mixture produced.
- If during the gas production an over pressure of 0.2 bar is reached, open the valves to drain the over pressure (Caution: The gas may be inflammable!).

The Anaerobic jars may now be placed inside the incubator for incubation.

When using the GasPack-kit “anaerobic”, the anaerobiosis indicator should be discoloured when opening directly after incubation. After 5-10 seconds, the indicator turns blue. In case the indicators is already blue when opening the jars, the reduction might have been too low or there has been a leakage of the Anaerobic jar (see “Maintenance”).

Maintenance

In case the sealing rings of the Anaerobic jar’s lid are flattened from intensive use or if cracks are caused from ageing of the material, which leads to leakage of the Anaerobic jars, the sealing rings may easily be replace. Please contact us.

Ordering information

The offered 5 models of Anaerobic jars differ in size, material and the methods for reaching anaerobic conditions inside the jar. The most important features are presented in the following chart.

On request, the lids are available as custom-made versions, e.g. made of other material, without valves and manometer or with high-temperature manometer.

*Only with transparent lid (models "eco" and "crystal").

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Special equipment for connection to Anoxomat

All anaerobic jar types available with or without valve or manometer. Equipped with special coupling for connection to Anoxomat.

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schuett-biotec GmbH