Revisions


Rev from	Rev to	ECO #
2010/06	2011/03	5703-11

Notes:


1. BD Cat. Number Various
2. Blank (Sheet) Size: Length: 17" Width: 27.5"
 Number of Pages: 20 Number of Sheets: 1
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 4.25" x 2.75"
3. Style (see illustrations below): # 7


#1


#2


#3


#4

#5

#6

#7

4. See Specification Control Number VS-80-1041-2PR for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled by BD Caribe, LTD

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8010412		Category and Description Package Insert, GasPak EZ Gas Generating Container Systems	Sheet: 1 of 21 Scale: N/A	A

BD GasPak™ EZ Gas Generating Container Systems

English: pages 1 – 3 Italiano: pagine 9 – 11
Français : pages 3 – 6 Español: páginas 11 – 14
Deutsch: Seiten 6 – 9 Japanese: pages 15 – 17

8010412
2011/03

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instruccie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The GasPak™ EZ Gas Generating Container Systems are multi-use systems that produce atmospheres suitable to support the primary isolation and cultivation of anaerobic, microaerophilic, or capnophilic bacteria by use of gas generating sachets inside multi-use incubation containers.

SUMMARY AND EXPLANATION

In 1965, Brewer and Allgeier introduced a disposable hydrogen generator envelope, which was later modified to include carbon dioxide generation and an internal catalyst.^{1,2} The GasPak EZ Gas Generating Container Systems offer bio-performance and reliability without the need for a catalyst or an activation procedure. The GasPak EZ Gas Generating Sachet contains all the ingredients needed to create a specific atmosphere for specimen incubation. The sachet and specimens are placed in the clear, stackable, plastic GasPak EZ Incubation Container. The container is sealed and incubated.

PRINCIPLES OF THE PROCEDURE

The GasPak EZ Gas Generating Sachet consists of a reagent sachet containing inorganic carbonate, activated carbon, ascorbic acid and water. When the sachet is removed from the outer wrapper, the sachet becomes activated by exposure to air. The activated reagent sachet and specimens are placed in the GasPak EZ Incubation Container and the container is sealed. The sachet rapidly reduces the oxygen concentration within the container. At the same time, inorganic carbonate produces carbon dioxide.

For the cultivation of anaerobic bacteria, the GasPak EZ Anaerobe Container System Sachets produce an anaerobic atmosphere within 2.5 h with less than 1.0% oxygen, and greater than or equal to 13% carbon dioxide within 24 h.

For the cultivation of microaerophilic bacteria, the GasPak EZ Campy Container System Sachets produce a microaerophilic atmosphere within 2 h with approximately 6 – 16% oxygen, and 2 – 10% carbon dioxide within 24 h.

For the cultivation of carbon dioxide requiring organisms, the GasPak EZ CO₂ Container System Sachets produce an atmosphere within 24 h with greater than or equal to 2.5% carbon dioxide.

REAGENTS

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Observe established precautions against microbiological hazards throughout all procedures. Prior to discarding, decontaminate specimen containers and other contaminated material by autoclaving.

After use, the sachet may remain warm. Discard sachet only after reaction is complete (when cool).

DO NOT use sachet if outer foil is damaged or open in any manner. Examine carefully before use.

Avoid direct sunlight and excessive temperatures.

DO NOT stack sachets when using more than one.

DO NOT open container until incubation is complete.

GasPak EZ Container System Sachets are intended for use in GasPak EZ Incubation Containers, and BBL GasPak 100 System or BBL GasPak 150 Large Anaerobic System. Their use in other systems may yield misleading results.

GasPak EZ Incubation Containers are designed for GasPak EZ Container System Sachets. Do not use with other gas generating sachets or envelopes (e.g., BBL GasPak/GasPak Plus Envelopes).

Storage Instructions: On receipt, store sachets at 2 – 25°C in a dry environment.

GasPak EZ Container System Sachets are ready to use. The expiration date is for unopened, intact and properly stored sachets. Do not open until ready to use.

Product Deterioration: Do not use reagent sachets if outer package has been damaged or opened. Do not use containers that show cracks, chips, damaged gaskets or other irregularities. Containers to be discarded shall be treated as biohazardous waste.

Active Ingredients: Ascorbic acid, activated carbon and water

Warning: Activated carbon is irritating to eyes, respiratory system and skin.

PROCEDURE

Materials Provided: GasPak EZ Container System. Each system requires the use of a container and sachet(s) (See "Availability").

Container	# of Sachets	# of BBL Stacker Petri dishes
GasPak EZ Small Incubation Container (Cat. No. 260002)	1	10 – 12
GasPak EZ Standard Incubation Container (Cat. No. 260671)	2	15 – 18
GasPak EZ Large Incubation Container (Cat. No. 260672)	3	30 – 33*
BBL GasPak 100 System (Cat. No. 260626)	1	12*
BBL GasPak 150 Large Anaerobic System (Cat. No. 260628)	3	36*

*Or a combination of Petri dishes and tubed media.

Materials Required But Not Provided: GasPak EZ Incubation Container Rack (See "Availability"), BBL GasPak tube holder, ancillary culture media, reagents, indicators, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure:

1. Place the desired plates in the rack inside the GasPak EZ Incubation Container. (If tubes are being used in the BBL GasPak Jar, use the BBL GasPak tube holder, Cat. No. 260630.)
2. Use the above table to determine the number of sachet(s) needed per container. Remove the GasPak EZ Container System Sachet(s) from the carton. Remove the outer foil packaging.
3. Place the activated sachets in the GasPak EZ Incubation Container or BBL GasPak Jar with the plates. The sachet(s) should be placed along side the plates between the plate rack and the outside of the container. If using the GasPak EZ Anaerobe Sachet for incubation, add an anaerobic indicator (Cat. No 271051) to the container or jar at this time. If using GasPak EZ w/Indicator, there is no need to add an additional O₂ indicator. The O₂ indicator will be white when the sachet is removed from the foil packet, gradually becoming blue as the sachet is exposed to an oxygenated atmosphere; it will then revert to white when the environment is reduced.
4. Close the GasPak EZ Container by seating the lid on the container bottom. Snap all four latches completely closed. If using the BBL GasPak Jar, refer to the package insert for the jar.
5. Incubate the GasPak EZ Container System at a temperature appropriate for the organism being cultured, but not to exceed 45°C. The GasPak EZ Incubation Containers may be stacked to conserve incubator space.
6. After incubation, open the container, remove plates and dispose of the GasPak EZ Container System Sachets in the appropriate manner. Ensure sachets are cool to the touch before discarding.

Maintenance Recommendations:

Proper care of the GasPak EZ Incubation Container consists of the following:

1. Avoid contact with abrasives.
2. Avoid contact with solvents and strong detergents.
3. Rinse and dry thoroughly after cleaning with a mild detergent/disinfectant.
4. Never autoclave container or lid.

User Quality Control:

Each GasPak EZ Gas Generating Container System should be tested periodically for its ability to provide adequate conditions for the growth of appropriate bacteria.

System	QC organism
GasPak EZ Anaerobe Container System	<i>Bacteroides fragilis</i> ATCC™ 25285
GasPak EZ Campy Container System	<i>Campylobacter jejuni</i> ATCC 29428
GasPak EZ CO ₂ Container System	<i>Neisseria gonorrhoeae</i> ATCC 19424

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Anaerobic conditions are achieved within 2.5 h, with greater than or equal to 13% carbon dioxide within 24 h at 35°C. Visible condensate should occur within 30 min of activation. Blood-containing agar plates appear reduced within 2 – 4 h at 35°C. The integrated anaerobic indicator should appear reduced (white) within 9 h at 35°C.

Microaerophilic conditions are achieved within 2 h with an oxygen concentration of approximately 6 – 16% and a carbon dioxide concentration of 2 – 10% at 35°C within 24 h.

Carbon dioxide enriched conditions are achieved within 24 h with a carbon dioxide concentration of greater than or equal to 2.5% at 35°C.

PERFORMANCE CHARACTERISTICS

GasPak EZ Anaerobic Container System Sachets (260678)

Prior to release, representative samples of each lot of GasPak EZ Anaerobic Container System Sachets are tested for performance characteristics.

Each sample of the GasPak EZ Anaerobic Container System Sachets is activated and placed in the incubation container along with an anaerobic indicator. The container is sealed according to label directions and observed for the time required for visible condensation to form. The container is then incubated at 35 ± 2°C. After approximately 2.5 h, a 22-gauge 1.5" needle is inserted into the incubation container and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of oxygen present in the incubation container. The percent of oxygen for each container is less than or equal to 1% and the mean percent for all samples tested for oxygen is less than 0.7%. After approximately 24 h incubation at 35 ± 2°C, another sample is taken from the container and analyzed for the amount of carbon dioxide. The percent of carbon dioxide is dependent upon container size but is greater than or equal to 13%. Within 9 h, the anaerobic indicators are white (reduced). The average time for visible condensation to appear inside the incubation containers is within 20 min after activation.

GasPak EZ Anaerobic System Sachets with Indicator (260001)

Prior to release, representative samples of each lot of **GasPak EZ Anaerobic System Sachets with Indicator** are tested for performance characteristics.

Each sample of the **GasPak EZ Anaerobic System Sachets with Indicator** is activated and placed in the incubation container. The container is sealed according to label directions and observed for the time required for visible condensation to form. The container is then incubated at $35 \pm 2^\circ\text{C}$. After approximately 2.5 h, a 22-gauge 1.5" needle attached to a 3 mL syringe is inserted into the incubation container and a sample of the gas is removed, placed in a gas chromatograph and analyzed for the amount of oxygen present in the incubation container. The percent of oxygen for each container is less than or equal to 1% and the mean percent for all samples tested for oxygen is less than 0.7%. After approximately 24 h incubation at $35 \pm 2^\circ\text{C}$, another sample is taken from the container and analyzed for the amount of carbon dioxide. The percent of carbon dioxide is dependent upon container size but is greater than or equal to 13%. Reduction of the Anaerobic Indicator pill attached to the sachet must occur within 24 h of incubation at $35 \pm 2^\circ\text{C}$. The average time for visible condensation to appear inside the jars is within 20 min after activation.

GasPak EZ Campy Container System Sachets (260680)

Prior to release, representative samples of each lot of **GasPak EZ Campy Container System Sachets** are tested for performance characteristics.

Each sample of the **GasPak EZ Campy Container System Sachets** is activated and placed in the incubation container. The container is sealed according to label directions and observed for the time required for visible condensation to form. The container is then incubated at $35 \pm 2^\circ\text{C}$. After approximately 2 h, a 22-gauge 1.5" needle is inserted into the incubation container and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of oxygen present in the incubation container. The percent of oxygen for each container is 6 to 16%. After approximately 24 h incubation at $35 \pm 2^\circ\text{C}$, another sample is taken from the container and analyzed for the amount of carbon dioxide. The percent of carbon dioxide for each container is 2 to 10%. The average time for visible condensation to appear inside the incubation containers is within 20 min after activation.

GasPak EZ CO₂ Container System Sachets (260679)

Prior to release, representative samples of each lot of **GasPak EZ CO₂ Container System Sachets** are tested for performance characteristics.

Each sample of the **GasPak EZ CO₂ Container System Sachets** is activated and placed in the incubation container along with 2 filled Petri dishes. The container is sealed according to label directions and incubated at $35 \pm 2^\circ\text{C}$. After approximately 24 h, a 22-gauge 1.5" needle is inserted into the incubation container and a sample of the gas is removed, placed in a gas chromatograph, and analyzed for the amount of carbon dioxide present in the incubation container. The percent of carbon dioxide for each container is greater than or equal to 3%.

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DOMINIQUE DUTSCHER SAS
Benex Limited