




Remove both caps of the system. Using the tip of inoculating needle, pick up a well isolated colony from a culture medium.

DOMINIO

DOMINIO
LITSHER SAS



*Inoculate **EnteroPluri-Test** turning and withdrawing the needle throughout the sectors of the system. Reinsert the needle with a turning movement until the breakage notch; break the inoculating needle folding it in correspondence with the notch.*

DOMINIQUE DUTSHER SAS



The portion of the needle remaining inside the system keeps anaerobic conditions necessary for reactions of the sectors **Glucose/Gas, Lysine** and **Ornithine**.



Use the broken portion of the needle to punch the plastic film on the holes of the sectors **Adonitol, Lactose, Arabinose, Sorbitol, VP, Dulcitol/PA, Urea, Citrate.**

DOMINIQUE DOUILLIER SAS



Screw again both caps and incubate **Enteropluri-Test** at 36 ± 1 °C for 18-24 hours, putting it on its flat surface or vertically.



*Inject with a syringe 3 or 4 drops of Kovac's Reagent in the sector **H₂S/Indole**.*



*Inject with a syringe 3 drops of alpha-naphtol and 2 drops of potassium hydroxide in the sector **VP**.*



Form the 5-digit code following the instructions and identify the bacterium using the **Codebook**.