INSTRUCTIONS FOR USE

CE 0123

DOMINIQUE DITSCHER SAS

Swabs with Stuart Gel, Amies Gel transport medium Instructions for use

INTENDED USE

Swabs with transport medium are intended for the collection, transport and preservation of clinical samples for bacteriological testing. These swabs are sterile and ready-for-use devices.

NUMBER	TRANSPORT MEDIUM	APPLICATOR SWAB	INTENDED USE/SAMPLING SITE *
26070	Stuart Agar Gel with Charcoal	Regular single plastic applicator	Mouth, Throat, Vagina, Wounds
26066	Stuart Agar Gel without charcoal	Regular single plastic applicator	Mouth, Throat, Vagina, Wounds
26068	Amies Agar Gel without charcoal	Regular single plastic applicator	Mouth, Throat, Vagina, Wounds
26071	Amies Agar Gel with charcoal	Regular single plastic applicator	Mouth, Throat, Vagina, Wounds

* This table is just a suggested reference. Please consult your GLP procedures to choose the most appropriate device for the specific sampling site.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This is made possible thanks to the use of swabs with transport medium. Each device consists of a sterile peelpouch containing an applicator swab used to collect the sample, and a tube containing the transport medium in which the swab applicator is placed after sampling.

The applicator swab is available with various transport media: Amies gel and Stuart gel. These transport media are non-nutritious, buffered with phosphate, and provide a reduced environment, due to their formulation with sodium thioglycolate^{1,2,3,4,5}.

The medium is designed to maintain the viability of organisms during transit to the laboratory. The viability of resistant bacteria such as *Neisseria gonorrhoeae* can be prolonged due to the presence of charcoal in the Amies Gel with charcoal medium. Once a swab sample is collected, it should be placed inside the tube with the medium and taken to the laboratory as soon as possible, then cultured on a suitable isolation medium (Blood Agar, Laked Blood Agar, MacConkey, etc.) For specific recommendations with regards to the collection of samples for microbiological analyses and primary isolation techniques, please refer to the internal instructions of the laboratory.

The single devices are packaged in peelpouches made of paper film.

REAGENTS

The nominal formula for each medium is as follows:

Stuart Agar Gel Transport Medium Without Charcoal

Sodium Glycerophosphate
Calcium Chloride
Sodium Thioglycollate
Agar
Distilled Water

Stuart Agar Gel Transport Medium with Charcoal

Sodium Glycerophosphate
Calcium Chloride
Sodium Thioglycollate
Agar
Charcoal
Distilled water

Amies Agar Gel Transport Medium Without Charcoal

Sodium Chloride
Potassium Chloride
Calcium Chloride
Magnesium Chloride
Monopotassium Phosphate
Disodium Phosphate
Sodium Thioglycollate
Agar
Distilled Water

Amies Agar Gel Transport Medium With Charcoal

Sodium Chloride
Potassium Chloride
Calcium Chloride
Magnesium Chloride
Monopotassium Phosphate
Disodium Phosphate
Sodium Thioglycollate
Charcoal
Agar
Distilled Water

TECHNICAL NOTES

The formula of the transport medium of the devices contains thioglycolic acid, an important component for product performance and viability of the microorganisms.

Thioglycolic acid has a natural sulphurous odour. Do not be alarmed should you perceive this odour upon opening the product, since it is perfectly normal and causes no hazards. The yellowing effect of the tube or of the fibre tip of the applicator swab is a well-known phenomenon tied to the product sterilisation treatment and has no contraindications in terms of product quality or performance.

PRECAUTIONS

1. This is a single-use, disposable product; reuse may cause a risk of infection and/or inaccurate results.
2. For in vitro use.
3. The device is certified as a Class IIa device under the classification terms of the European Medical Device Directive EC 93/42. In particular, the applicator swab is qualified for short transient contact with the patient in order to collect a specimen. This short contact is made with the patient's external surfaces or internal surfaces via normal body orifices such as the nose, throat or vagina or via surgical wounds.
4. When collecting swab samples from patients, care should be taken to not use excessive force or pressure, as this may result in breakage of the swab shaft.
5. The fibre attachment to the applicator stick was designed for instantaneous sample collection; prolonged contact must be avoided as it may result in the detachment of fibre.
6. Instructions for use must be followed carefully. The manufacturer shall not be held responsible for any unauthorised or unqualified use of the product.
7. When the swab sample is cultured in the laboratory, the shaft must be cut; it is recommended to use sterile scissors for a safe and clean break.
8. It must be assumed that all specimens contain infectious microorganisms; therefore, all specimens should be handled with the appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.
9. Swab sample processing should be performed inside a protective safety cabinet or protective hood. Protective laboratory clothing and eyeglasses should be worn at all times when processing culture swab samples.
10. The product is supplied sterile, and cannot be re-sterilised.
11. Product use with specific molecular platforms, antigen kits (Elisa or Rapid Antigen Kit) or other uses not mentioned in this document must be validated and qualified by the kit manufacturer or by the final user.

STORAGE AND STABILITY

Store the devices at a temperature between 5 and 25°C. Do not freeze or overheat. Do not use after the expiration date, which is clearly printed on the outer box, on each pack of swabs, on each individual sterile swab pouch and on the specimen transport tube label. If the product is stored at temperatures other than the indicated ones, it will not maintain the stated product specifications and performance.

PRODUCT DETERIORATION

The contents of unopened or undamaged units are guaranteed sterile. Do not use if they show evidence of damage, dehydration or contamination. Do not use if expiration date has passed.

MATERIALS SUPPLIED


Each cardboard box contains 100 sterile pieces. Each individual swab pouch contains an applicator and a plastic tube containing transport medium.

MATERIALS REQUIRED BUT NOT SUPPLIED

Material needed for the isolation, differentiation and cultivation of aerobic and anaerobic bacteria.

INSTRUCTIONS FOR USE

The instructions for use are printed on each unit along with descriptive diagrams, and are summarized as follows:

1. Peel open the sterile pouch of the device where it says "Peel Here".
2. Remove cap from transport tube.
3. Remove applicator swab and collect specimen. During specimen collection, the applicator tip should touch only the area where the infection is suspected to be, to minimize potential contamination.
4. Place applicator swab in transport tube and replace cap firmly to completely seal.
5. Record patient's name and information on tube label.
6. Send specimen to the laboratory for immediate analysis.

Precaution - When collecting swab samples from patients, care should be taken to not use excessive force or pressure which may result in breakage of the swab shaft.

QUALITY CONTROL

All raw materials, swab components and batches of finished products are subjected quality control. As part of these test procedures, a panel of control bodies is used to check compliance with product specifications. Sterility and quality control certificates indicate some microorganisms with which the product is tested during quality control.

RESULTS

The survival of bacteria in a transport medium depends on many factors. These include the type of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium.

The devices guarantee the viability of many microorganisms for 24 hours.

For resistant bacteria, such as *Neisseria gonorrhoeae*, *Streptococcus pneumoniae* and anaerobes, swab specimens should be plated directly onto culture medium or transported immediately to the laboratory and cultured within 24 hours. Published studies have demonstrated that swabs with Amies and Stuart transport medium are able to maintain the viability of clinically important aerobic and anaerobic bacteria for 24 hours⁶⁻¹⁶.

LIMITATIONS

Swabs with Stuart Gel and Amies gel are intended exclusively for the collection and transport of bacteriological samples.

Samples containing viruses or chlamydia should be collected and transported using alternative specific transport systems.

The device has not been validated for environmental sampling and sterility tests.

PERFORMANCE CHARACTERISTICS

Survey studies were performed using the swab with Amies gel and Stuart gel culture medium with various bothersome aerobic and anaerobic organisms. Swabs were dosed with inoculum and inserted into the transport tube containing medium. The tubes were stored at room temperature prior to subculturing onto appropriate media.

Aerobic organisms evaluated were *Escherichia coli* (NCTC 9001 and ATCC 25922), *Haemophilus influenzae* (ATCC 10211), *Neisseria gonorrhoeae* (ATCC 43069), *Neisseria meningitidis* (NCTC 10025 and ATCC 13090), *Pseudomonas aeruginosa* (NCTC 9332 and ATCC 27853), *Staphylococcus aureus* (NCTC 5532 and ATCC 25923) and *Streptococcus pyogenes* (ATCC 19615). Anaerobic organisms evaluated were *Bacteroides fragilis* (ATCC 25285), *Bacteroides levii* (ATCC 29147), *Bacteroides thetaiotaomicron* (ATCC 29741), *Bacteroides vulgatus* (ATCC 8482), *Clostridium difficile* (ATCC 9689), *Clostridium perfringens* (ATCC 13124), *Clostridium sporogenes* (ATCC 3584), *Clostridium tertium* (ATCC 19405), *Fusobacterium necrophorum* (ATCC 25286), *Fusobacterium nucleatum* (ATCC 25586), *Peptostreptococcus anaerobius* (ATCC 27337), *Peptostreptococcus magnus* (ATCC 29328), *Porphyromonas gingivalis* (ATCC 33277), *Prevotella melaninogenica* (ATCC 25845), and *Propionibacterium acnes* (ATCC 6919). All organisms tested remained viable for more than 24 hours when maintained at room temperature.

Please refer to symbol table at the end of the instructions for use

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BIBLIOGRAPHY

1. Stuart R.D. The diagnosis and control of gonorrhoea by bacteriological cultures. *Glasgow M. J.* 27: 131-142, May 1946.
2. Moffett M. Young J.L. and Stuart R.D. Centralised gonococcus culture for dispersed clinics. *Brit. M. J.* 2: 421-424. Aug. 28, 1948.
3. Stuart R.D., Toshach S.R. and Patsula T.M. The problem of transport of specimens for culture of gonococci. *Cand. J. Health*, 45: 73-83, February 1954.
4. Stuart R. D. Transport Medium for specimens in Public Health Bacteriology. *Public Health Reports* 74: No. 5, 431-438, May 1959.
5. Amies C.R. A modified formula for the preparation of Stuart's medium. *Canadian Journal of Public Health*, July 1967. Vol. 58, 296 - 300.
6. Balows A., Hausler, Jr. W. J., Herrmann K.L., Isenberg H. D., Shadomy H.J. *Manual of Clinical Microbiology*. Fifth Edition. American Society for Microbiology, Washington DC, 1991.
7. Isenberg H. D. (Editor in Chief). *Clinical Microbiology Procedures Handbook*. American Society for Microbiology, Washington DC, 1992.
8. Van Horn K., Toth C. and Wegienek J. Viability of aerobic microorganisms in four swab systems. Poster Session 249/C Abstract C-436. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
9. Zavala M. K., Citron D. M., Goldstein E. J. C. Evaluation of a novel specimen transport system for anaerobic bacteria. *Clinical Infectious Diseases*, 1997; Vol. 25 (Suppl 2): S132 - 3.
10. Perry J.L. Assessment of swab transport systems for aerobic and anaerobic organism recovery. *Journal of Clinical Microbiology*, May 1997. Vol. 35, 1269 - 1271.
11. Van Horn K., Toth C. and Warren K. Comparison of Copan Amies agar swab and BBL Port-a-cul swab for recovery of anaerobic bacteria. Poster Session 249/C Abstract C-437. 98th General Meeting of American Society for Microbiology, Atlanta, GA, May 1998.
12. Arbique J.C, Forward K. R., and Le Blanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. *Diagnostic Microbiology and Infectious Disease*, 2000 Vol. 26, 163-168.
13. Van Horn K., Warren K., and Toth C. Comparison of four swab systems for the recovery of anaerobic bacteria. Abstract C-32. 99th General Meeting of American Society for Microbiology, Chicago, IL, May 1999.
14. Olsen C. C., Schwebke J. R., Benjamin W. H., Beverly A., and Waites K. B. Comparison of direct inoculation and Copan transport systems for the isolation of *Neisseria gonorrhoeae* from endocervical specimens. *Journal of Clinical Microbiology* Vol. 37, No. 11, Nov. 1999, p. 3583 - 3585.
15. Thompson D. S., French S. A. Comparison of commercial Amies transport systems with in-house Amies medium for recovery of *Neisseria gonorrhoeae*. *Journal of Clinical Microbiology* Vol. 37, No. 9. Sept. 1999, p. 3020-3021.
16. Hetchler C., Brown C., and Galbraith J. Comparison of 3 Amies Gel transport systems for the recovery of 12 clinically significant organisms. Abstract C-152. 100th General Meeting of American Society for Microbiology, Los Angeles, CA, May 2000.
17. Summanen P., Baron E.J., Citron D. M., Strong C.A., Wexler H.M. and Finegold S. M. *Wadsworth Anaerobic Bacteriology Manual*. Fifth Edition. Star Publishing Company, Belmont, CA, 1993.
18. Ronald M. *Atlas, Handbook of Microbiological Media*. 4th Edition, CRC Press Taylor and Francis Group 2010

INDEX OF SYMBOLS:

Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos	Meaning / Significato / Signification / Bedeutung / Sens / Significado
	Manufacturer / Fabbricante / Fabricante/ Hersteller / Fabricant / Fabricante
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Número de identificación del organismo notificado / Kennnummer der benannten Stelle / Identification de l'organisme notifié / Número de identificação do organismo notificado
STERILE R	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não reutilizar
REF	Numero di catalogo / Catalogue number / Número de catálogo / Bestellnummer / Numéro de catalogue / Número de catálogo
	Temperature limits / Limiti di temperatura / Límites de temperatura / Temperaturbegrenzung / Limites de température / Limites de temperatura
	Use before / Utilizzare entro / Fecha de caducidad / Verwendbar bis / À utiliser avant / Utilizar antes de
	Consult the instructions for use / Consultare le istruzioni per l'uso / Consultar las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Rasgar e abrir
LOT	Batch code (lot) / Codice del lotto (partita) / Código de lote (lote) / Chargenbezeichnung (Los) / Numéro de lot (lot) / Código do lote (remessa)
	Contents sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Inhalt ausreichend für <n> Prüfungen / Contenu suffisant pour <n> tests / Conteúdo suficiente para <n> testes
	Do not use if the packaging is damaged / Non utilizzare in caso di confezionamento danneggiato / No utilizar si el envase está dañado / Nicht verwenden, wenn Verpackung beschädigt / Ne pas utiliser si l'emballage est endommagé / Não utilizar se a embalagem estiver danificada

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia Italy
Tel +39 030 2687211
Fax +39 030 2687250

Distributed by:
FL MEDICAL s.r.l. Unipersonale
Via Enrico Mattei, 20
35038 Torreglia (PD) - Italy
Tel: +39 049 5211475 / +39 049 5212835
Fax: +39 049 5212566

E-mail: info@flmedical.com
Website: www.flmedical.com

Copyright © 2019 Copan Italia S.p.A.
All rights reserved