**EBA 280**

**EBA 280 S**

Ⓓ	Bedienungsanleitung.....	9
Ⓔ	Operating Manual.....	47
Ⓕ	Mode d'emploi .....	85
Ⓘ	Istruzioni per l'uso .....	124

Andreas Hettich GmbH & Co. KG  
Föhrenstraße 12, D-78532 Tuttlingen / Germany  
Phone +49 (0)7461 / 705-0  
Fax +49 (0)7461 / 705-1125  
info@hettichlab.com, service@hettichlab.com  
www.hettichlab.com

**Hettich**

© 2014 by Andreas Hettich GmbH & Co. KG

All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.

Änderungen vorbehalten! , Modifications reserved! , Sous réserve de modifications ! , Con riserva di modifiche!

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We hereby declare under our sole responsibility that the designated device and its accessories, which are listed in the technical documentation for this device and whose conformity has been assessed together with the device, conform to the Directive 98/79/EC on in vitro diagnostic medical devices.

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The conformity evaluation process was performed in accordance with appendix III of Directive 98/79/EC.

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**The following additional European directives and ordinances have been applied:**

- Machinery Directive 2006/42/EU
- EMC directive 2014/30/EU
- Low voltage directive 2014/35/EU
- RoHS II Directive 2011/65/EU (without involvement of a notified body)
- Ordinance (EC) No. 1907/2006 (REACH) (without involvement of a notified body)

Applied standards:

According to the list of applied standards, which is part of the product file.

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Kanada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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---

## Standards and regulations which apply to this device

The device is a high-end technical product. It is subject to extensive testing and certification procedures according to the following standards and regulations in their respectively valid version:

### Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Excess-voltage category II)
- IEC 61010-2-010 „Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of Materials" (only valid for centrifuges with heating)
- IEC 61010-2-011 „Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-011: Particular requirements for refrigerating equipment" (only valid for centrifuges with cooling)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"
- IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment"

### Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

### Risk management:

- DIN EN ISO 14971 "Application of risk management to medical devices"

### Restriction of Hazardous Substances (RoHS II):

- EN 50581 "Technical documentation for assessing electric and electronic devices with regard to the restriction of hazardous substances"

### European directives applied for conformity assessment procedures:

- In vitro diagnostic device directive 98/79/EG  
EC conformity assessment procedure according to annex III "EC DECLARATION OF CONFORMITY" – self-declaration by the manufacturer
- Directive 2011/65/EU for the restriction of use of certain hazardous substances in electric and electronic devices. Carrying out the EC conformity assessment process is the sole responsibility of the manufacturer, without the involvement of a notified body.

### Applied medical device regulations outside Europe:

- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

### Certified quality management system according to

- ISO 9001 "Quality management systems – Requirements"
- ISO13485 "Medical devices - Quality management systems - Requirements for regulatory purposes"

### Environmental management system according to

- ISO 14001 "Environmental management systems - Requirements with guidance for use"

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### Restrictions relatives aux substances dangereuses (RoHS II):

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- **USA** : QSR, 21CFR 820 "CFR Title 21 - Food and Drugs : TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada**: CMDR, SOR/98-282 "Medical Devices Regulations"

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## Contents

1	How to use these operating instructions .....	49
2	Symbol meanings .....	49
3	Intended use .....	50
4	Remaining risks .....	50
5	Technical specifications .....	51
6	Notes on safety .....	52
7	Transport and storage .....	53
7.1	Transport .....	53
7.2	Storage .....	53
8	Scope of delivery .....	53
9	Removing the transport securing device .....	54
10	Initial operation .....	54
11	Opening and closing the lid .....	55
11.1	Opening the lid .....	55
11.2	Closing the lid .....	55
12	Inserting and removing the rotor .....	56
13	Loading the rotor .....	57
14	Operating and display elements .....	58
14.1	Displayed symbols .....	58
14.2	Control panel keys .....	58
14.3	Setting options .....	59
15	Programming .....	60
15.1	Inputting or changing programs .....	60
15.2	Calling up programs .....	62
15.3	Saving programs .....	62
16	Centrifugation .....	63
16.1	Centrifugation with preset time .....	63
16.2	Continuous run .....	64
16.3	Short-term centrifugation .....	65
17	Emergency Stop .....	65
18	Rotor detection .....	66
19	Cycle counter .....	67
20	Settings and queries .....	68
20.1	Querying system information .....	69
20.2	Acoustic signal .....	70
20.3	Optical signal after ending the centrifugation run .....	71
20.4	Automatic unlocking of the lid after the centrifugation run .....	72
20.5	Backlighting of the display .....	73
20.6	Displayed centrifugation parameters after switching on the device .....	74
20.7	Activating/deactivating the "Dual time" function .....	75
20.8	Querying the hours of operation and the number of centrifugation runs .....	76
20.9	Resetting the cycle counter to zero .....	77
21	Relative centrifugal force (RCF) .....	78

22	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm <sup>3</sup> .....	78
23	Emergency unlocking .....	78
24	Maintenance and servicing .....	79
24.1	Centrifuge (housing, lid and centrifuging chamber).....	79
24.1.1	Surface cleaning and care .....	79
24.1.2	Surface disinfection .....	79
24.1.3	Removal of radioactive contaminants .....	79
24.2	Rotors and Attachments .....	80
24.2.1	Cleaning and care.....	80
24.2.2	Disinfection .....	80
24.2.3	Removal of radioactive contaminants .....	80
24.2.4	Rotors and accessories with limited service lives .....	80
24.3	Autoclaving .....	81
24.4	Centrifuge containers.....	81
25	Faults .....	82
26	Change mains input fuses.....	84
27	Returning Devices .....	84
28	Disposal .....	84
29	Anhang / Appendix .....	162
29.1	Rotoren und Zubehör / Rotors and accessories.....	162

## 1 How to use these operating instructions

- Before using the centrifuge, read the operating instructions and observe them.
- The operating instructions are a part of the device. They must always be kept readily available.
- If the device is set up at a different location, the operating instructions must be provided with it.

## 2 Symbol meanings

Symbol on the device:  
Attention, general hazard area.

Symbol on the device:  
Observe operating instructions.  
This symbol indicates that the user must observe the operating instructions provided.

Symbol in this document:  
Attention, general hazard area.  
This symbol refers to safety relevant warnings and indicates possibly dangerous situations.  
The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the device and in this document:  
Beware of biohazard.

Symbol in this document:  
This symbol refers to important circumstances.

Symbol on the device:  
The symbol is on the housing, near the hinge.  
Direction of rotation of the rotor.

Symbol on the device:  
The symbol is located on the front of the housing, below the borehole.  
The direction of rotation of the Allen key to unlock the lid in an emergency.

Symbol on the device and in this document:  
Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).  
Applies in the countries of the European Union, as well as in Norway and Switzerland.

Symbol on the shipping carton label:  
This way up.

Symbol on the shipping carton label:  
The shipping packaging must be transported and handled within the indicated humidity range (10% - 80%).

Symbol on the shipping carton label:  
The shipping packaging must be stored, transported and handled within the temperature range shown (-20°C - +60°C).

Symbol on the shipping carton label:  
The shipping packaging must be kept away from rain and kept in a dry environment.

Symbol on the shipping carton label:  
Fragile, handle with care.

### 3 Intended use

This device is a medical product (laboratory centrifuge) within the context of the IVD Directive 98/79/EC.

The centrifuge was designed only for the separation of substances or mixtures with a density of no more than 1.2 kg/dm<sup>3</sup>, particularly of samples for preparation for human medical in-vitro diagnostics.

This especially includes the determination of potassium in blood serum. The running behaviour allows for gentle separation, thereby preventing the accumulation of additional erythrocytic potassium in the serum. This prevents the test results from being falsified.

The centrifuge is only meant for this purpose.

Another use or one which goes beyond this, is considered to be non-intended. The company Andreas Hettich GmbH & Co. KG is not liable for damage resulting from this.

Observing all information in the operating instructions and complying with the measures described therein is also a part of the intended use.

If the centrifuge is installed in another device or in a system, the manufacturer of the overall system is responsible for its safety.

### 4 Remaining risks

The device is built according to the state-of-the-art and the recognized safety regulations. If used and handled improperly, there could be life-threatening danger to the user or third parties, or the device could be impaired or there could be other property damage. The device is only to be used for its intended purpose and only when it is in safe working condition.

Malfunctions which could affect safety must be corrected immediately.

## 5 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	EBA 280		EBA 280 S	
Type	1101	1101-01	1102	1102-01
Mains voltage ( $\pm 10\%$ )	200 - 240 V 1~	100 - 127 V 1~	200 - 240 V 1~	100 - 127 V 1~
Mains frequency	50 - 60 Hz	50 - 60 Hz	50 - 60 Hz	50 - 60 Hz
Connected load	185 VA	185 VA	330 VA	330 VA
Current consumption	0.85 A	1.75 A	1.6 A	3.0 A
Max. capacity	6 x 50 ml			
Allowed density	1.2 kg/dm <sup>3</sup>			
Speed (RPM)	6000			
Force (RCF)	4146		5071	
Kinetic energy	2700 Nm			
Obligatory inspection (BGR 500)	no			
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 2°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.			
– Set-up site				
– Altitude				
– Ambient temperature				
– Humidity				
– Excess-voltage category (IEC 60364-4-443)	II			
– Pollution degree	2			
Device protection class	I			
Not suitable for use in explosion-endangered areas.				
EMC	EN / IEC 61326-1, Class B		FCC Class B	
– Emitted interference, Interference immunity	EN / IEC 61326-1, Klasse B		FCC Class B	
Noise level (dependent on rotor)	$\leq 51$ dB(A)		$\leq 56$ dB(A)	
Dimensions				
– Width	326 mm			
– Depth	400 mm			
– Height	242 mm			
Weight	approx. 11 kg			

## 6 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge, make sure you check that the rotor is engaged on the motor shaft.
- When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- For centrifuges without temperature control, when the room temperature is increased and/or if the device is frequently used, the centrifuging chamber could be heated up. Therefore, it can't be ruled out that the sample material might be changed due to the temperature.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). Before centrifuge vessels are used which are not listed in the chapter "Appendix, Rotors and accessories", the user must make sure they can be used by asking the manufacturer.
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
  - inflammable or explosive materials
  - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.  
For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").  
Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.  
No biosafety systems are available for this centrifuge.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:  
EN / IEC 61010-1 EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
  - the centrifuge is operated in accordance with the operating instructions,
  - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations.
- Meeting the country-specific requirements concerning occupational safety with regard to the use of laboratory centrifuges at the workplaces provided for this purpose by the user is the responsibility of the user.

## 7 Transport and storage

### 7.1 Transport

 Before transporting the device, the transport securing device must be installed.

When the device and accessories are transported, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 10% to 80%, non-condensing

### 7.2 Storage

 The device and the accessories may only be stored in closed and dry rooms.

When the device and accessories are stored, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 10% to 80%, non-condensing

## 8 Scope of delivery

The following accessories are delivered with the centrifuge:

- 1 connecting cable
- 2 fuses
- 1 hex head wrench
- 1 operating manual
- 1 information sheet, transport securing device
- 1 information sheet, emergency unlocking

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

## 9 Removing the transport securing device


It is imperative that the transport securing device be removed.

Keep the transport securing device in a safe place since it must be installed again before transporting the device.

The device may only be transported with the transport securing device installed.

To protect the device from damage during transport, the motor is fixed in place.  
This transport securing device must be removed when the device is put into operation.

- Remove the two screws (a) and spacer sleeves (b).

The transport securing device is installed in the opposite order.

## 10 Initial operation

- Remove the transport securing device from the housing floor; see the "Removing the transport securing device" chapter.
- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

**When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**

- Ventilation openings may not be blocked.  
A distance of 300 mm must be maintained from the ventilation slots and openings of the centrifuge.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the power cord to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Switch on the mains switch (switch position "I"). The LEDs will briefly light up.  
The following displays appear one after the other:
  1. The centrifuge model
  2. The type number and program version
  3. The rotor code (R) and the maximum speed of the rotor (maxRPM) of the last rotor detected by the rotor detector
  4. The centrifugation parameters of the last used program or program 1

If the lid is closed, the message "Open the lid" is displayed.  
In this case, open the lid to display the centrifugation parameters.

## 11 Opening and closing the lid

### 11.1 Opening the lid

 The lid can only be opened if the centrifuge is switched on and the rotor is stationary. If this shouldn't be possible, see the "Emergency unlock" chapter.

 When the cycle counter is activated, after a centrifugation run, while opening the lid, the remaining number of running cycles (centrifugation runs) is briefly displayed.

Example:

- Press the key.  
The lid is unlocked by the motor.  
 : Lid unlocked.  
The left LED in the key will go out.

Example:

### 11.2 Closing the lid

 Do not reach with your fingers between the lid and housing.  
Do not slam the lid closed.

- Put the lid on and lightly press the front edge of the lid down.  
The lid is locked by the motor.  
 : Lid locked.


The left LED in the  key is illuminated.

Example:

## 12 Inserting and removing the rotor

### Rotor with unlocking knob

#### Inserting the rotor:


Dirt particles between the motor shaft and rotor prevent the rotor from having a perfect seat and cause it to run unsteadily. If necessary, clean the motor shaft (B) and the rotor hub before inserting the rotor.

- Lift up the rotor by the unlocking knob (A) and place horizontally on the motor shaft (B). The rotor will automatically engage on the motor shaft.
- Hold the rotor firmly at the left and right and pull up lightly to check whether it has properly engaged on the motor shaft.

#### Removing the rotor:

- Pull up the unlocking knob (A) on the rotor, hold onto it and lift the rotor from the motor shaft.

### Rotor without unlocking knob

#### Inserting the rotor:

Dirt particles between the motor shaft and rotor prevent the rotor from having a perfect seat and cause it to run unsteadily. If necessary, clean the motor shaft (B) and the rotor hub before inserting the rotor.

- Lift up the rotor by the handle (C), place it horizontally on the motor shaft (B) and press down as far as it can go.

#### Removing the rotor:


- Hold the rotor firmly by the handle (C) and lift it off of the motor shaft.

### 13 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Hold the rotor firmly at the left and right and pull up lightly to check whether it has properly engaged on the motor shaft.
- With swing-out rotors all rotor positions must be lined with **identical** hangers.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".  
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded


**Not permitted!**

Rotor is not evenly loaded


- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- The centrifuge containers may only be filled outside of the centrifuge.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.

In the case of angle rotors, the centrifuging vessels may only be filled so far that no fluid can be expelled from them while the centrifuge is running.

- When loading the angle rotors, no liquid may get into the angle rotors or centrifuging chamber.
- When loading the hangers and swing-out rotors and when swivelling out the hangers while the centrifuge is running, no liquid may enter the hangers or the centrifuging chamber.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

## 14 Operating and display elements

### 14.1 Displayed symbols

Lid unlocked.

Lid locked.

Rotation display. The rotation display lights up, rotating counter clockwise, as long as the rotor is turning.

### 14.2 Control panel keys

- Calling up programs
- Save program.

- Switching between the RPM display (RPM) and RCF display (>RCF<). RCF values are displayed in parentheses >(<.<. The LED above the key will light up if the RCF display is activated.

RPM Speed  
RCF : Relative centrifugal acceleration

- Selection key for activating the individual parameters. Every time the key is pressed, the next parameter is activated.
- Call up the "MACHINE MENU" (keep the key pressed for 8 seconds).
- In the "Machine Menu", select the menus "-> Info", "-> Settings" and "-> Time & Cycles".
- Scroll forward in the menus.

- Change the value of a previously selected parameter. If the key is kept pressed, the value changes with increasing speed.

- Start the centrifugation run. The LED in the key will be illuminated during the centrifugation run as long as the rotor is turning.
- Short-term centrifugation. Centrifugation is run as long as the key is kept pressed. The LED in the key will be illuminated during the centrifugation run as long as the rotor is turning.
- Save inputs and changes.
- Call up the menus "-> Info", "-> Settings" and "-> Time & Cycles".

- Finish the centrifugation run. The rotor runs down with a pre-selected brake stage. The right LED in the button is illuminated until the rotor is still. After the rotor stops, the left LED in the key will light up. Pressing the key twice triggers the Emergency Stop.
- Unlock the lid. The left LED in the key will go out.

### 14.3 Setting options

**PROG RCL** **Program location of the program to be called**  
The program locations 1 to 9 and # can be called.

**t/min** **Runtime.** Adjustable from 1 - 99 min, in steps of 1 minute.

**t/sec** **Runtime.** Adjustable from 1 - 59 s, in steps of 1 second.

**Continuous run "--:--".** Set the parameters **t/min** and **t/sec** to zero.

**Begins at START** **Counting the runtime.**  
The runtime begins counting immediately after the centrifugation run begins.

**Begins at SPEED** The runtime only begins counting after reaching the set speed.

 This setting is only possible when the "Dual time" function is activated. See the "Activating/deactivating the 'Dual time' function" chapter. This function is activated ex works.

**RPM** **Speed**  
A number value from 200 rpm to the maximum speed of the rotor can be set.  
Settable in increments of 10.  
For the maximum rotor speed, see the chapter "Appendix, Rotors and accessories".

**RAD/mm** **Centrifuging radius**  
Adjustable from 10 mm to 330 mm, in steps of 1 mm.  
For the centrifuging radius, see the chapter "Appendix, Rotors and accessories".

 It is only possible to input the centrifuging radius if the RCF display (➤RCF◀) is activated.

**RCF** **Relative centrifugal acceleration**  
A number value can be set which results in a speed between 200 rpm and the maximum rotor speed.  
Settable in increments of 1.

 It is only possible to input the relative centrifugal acceleration (RCF) if the RCF display (➤RCF◀) is activated.  
The relative centrifugal acceleration (RCF) depends on the centrifuging radius (RAD). The relative centrifugal acceleration RCF) and the set centrifuging radius (RAD) are displayed.

 **Starting steps 1 - 9.** Step 9 = shortest starting time, Step 1 = longest starting time.

 **Brake steps 0 - 9.** Step 9 = shortest run-down time, Step 1 = long run-down time, Step 0 = longest run-down time (brakeless run-down).

**PROG STO** **Program location where the program is stored.**  
9 different programs can be stored (program locations 1 to 9).

 No programs should be saved at program location #.  
The program location # serves as intermediate memory for changed centrifugation parameters.  
The centrifugation parameters at program location # are overwritten every time a centrifugation run is carried out with changed centrifugation parameters which were not saved to a program location.

## 15 Programming

 No programs should be saved at program location #.

The program location # serves as intermediate memory for changed centrifugation parameters.

If centrifugation parameters are changed, but are not saved to a program location afterwards, the display shows a dash "-" instead of the program location number. After starting the centrifugation run, these centrifugation parameters are then automatically saved at program location #.

The centrifugation parameters at program location # are overwritten every time a centrifugation run is carried out with changed centrifugation parameters which were not saved to a program location.

### 15.1 Inputting or changing programs

 In order to set the continuous run, the parameters **t/min** and **t/sec** must be set to zero. In the time display (t/min:s), "--:--" appears.

Example:

The relative centrifugal force (RCF) depends on the centrifuging radius (RAD). A correct calculation is only possible if the right centrifuging radius is set. During the input of the RCF, the set centrifuging radius is displayed.

Example:

It is only possible to display and input the relative centrifugal force (RCF) and the centrifuging radius (RAD) if the RCF display (>RCF<) is activated.

If no key is pressed for 8 seconds after selection or during parameter input, the previous values are shown on the display. The parameters must then be input again.

Entering parameters can be cancelled at any time by pressing the  key. In this case, the settings are not stored.

 • Press the key.

**PROG RCL** : Program location of the displayed centrifugation parameters.

↓

 • With the keys, set the program location whose centrifugation parameters are to be displayed.

↓

 • Press the key.

The centrifugation parameters are displayed.

↓

Continued on next page

**Example:**  
RPM display (RPM)

**Example:**  
RCF display (>RCF<)

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Keep pressing the key until the desired centrifugation parameter is displayed. Example: **t/min** : Runtime, minutes.

- Set the desired value with the keys.

- Input all the desired centrifugation parameters in succession.

- Press the key until the following is displayed:  
**PROG STO** : Program location where the centrifugation parameters are stored.

- Set the desired program location with the keys.

- Press the key. The centrifugation parameters are stored at the program location.

**Example:  
RPM display (RPM)**

**Example:  
RCF display (>RCF<)**

DOMINIQUE DUISSEL

### 15.2 Calling up programs

- 
  - Press the key.  
**PROG RCL** : Program location of the displayed centrifugation parameters.
- ↓
- 
  - With the keys, set the program location whose centrifugation parameters are to be displayed.
- ↓
- 
  - Press the key.  
The centrifugation parameters are displayed.

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

### 15.3 Saving programs

- 
  - Press the key twice.  
**PROG STO** : Program location where the centrifugation parameters are stored.
- ↓
- 
  - With the keys, set the program location where the centrifugation parameters are to be stored.
- ↓
- 
  - Press the key.  
The centrifugation parameters are stored at the program location.

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

## 16 Centrifugation

 In accordance with EN / IEC 61010-2-020, no persons, hazardous substances or objects may be within a safety range 300 mm around the centrifuge during a centrifugation run.


 After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.

Example:

If the permissible weight difference within the rotor load is exceeded, the centrifugation run is aborted during the start-up and the following display appears:

A centrifugation run can be aborted at any time by pressing the  key.

During a centrifugation run, centrifugation parameters can be selected and changed. The changed values are stored at program location #.

With the  key, you can toggle between the RPM display (RPM) and RCF display (>RCF<) at any time. If the RCF display (>RCF<) is worked with, the centrifuging radius must be input.

The following is displayed:

Further operation of the centrifuge is then only possible after opening the lid.

Operation errors and malfunctions will be shown (see the chapter on "Malfunctions").

- Plug in the mains switch. Switch position I.
- Load the rotor and close the centrifuge lid.

### 16.1 Centrifugation with preset time

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Input the desired centrifugation parameters (see the "Programming" chapter).

- Press the key. The centrifugation run is started. During the centrifugation run, the rotor speed or the resulting RCF value and remaining time are displayed.

Continued on next page

Example:  
RPM display (RPM)


Example:  
RCF display (>RCF<)


- After the time elapses or if the centrifugation run is aborted by pressing the  key, the run-out occurs with the set brake stage. The brake stage is displayed. Example  $\sim 7$

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

## 16.2 Continuous run

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Input the desired centrifugation parameters. Set the parameters **t/min** and **t/sec** to zero (see the "Programming" chapter).

- Press the key. The centrifugation run is started. During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.

- Press the key. The centrifugation run is ended. The run-out occurs with the set brake stage. The brake stage is displayed. Example  $\sim 7$

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

### 16.3 Short-term centrifugation

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Input the desired centrifugation parameters (see the "Programming" chapter).

- Press the key and keep it pressed. During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.

- Release the key again. The centrifugation run is ended. The run-out occurs with the set brake stage. The brake stage is displayed. Example  $\sim 7$ .

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

### 17 Emergency Stop

- Press the key twice. During the Emergency Stop, the run-out occurs at brake stage 9 (shortest run-out time). Brake stage  $\sim 9$  is displayed.

Example:  
RPM display (RPM)

Example:  
RCF display (>RCF<)

## 18 Rotor detection

The centrifuge is equipped with automatic rotor detection.

Rotor detection is carried out after starting each centrifugation run.

If the rotor code of a rotor is read in for the first time by the rotor detection, the centrifugation run then stops automatically after rotor detection. The rotor code (R) and the maximum speed (maxRPM) of the rotor permitted for this centrifuge are displayed.

Example:

After the rotor comes to a standstill (the left LED in the  key is illuminated), press any key. Afterwards, the maximum permissible number of run cycles is displayed.

Example:

Now, the maximum permissible number of run cycles of the rotor/hangers must be set (see the "Cycle counter" chapter).

When using a rotor whose rotor code has already been read in by the rotor detection, after starting the centrifugation run, the maximum speed (ROTOR MAX) of the rotor permissible for this centrifuge is briefly displayed.

Example:

 If the set speed is greater than the maximum speed of the rotor permissible for this centrifuge, the centrifugation run is aborted after rotor detection and the following displays alternately appear:

Example:

In this case, the speed on the display must be checked and adapted to the maximum rotor speed.

## 19 Cycle counter

 The period of use of certain rotors and hangers is limited to a maximum number of run cycles (centrifugation runs). The rotors and hangers in question are labeled with the maximum permissible number of run cycles.

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs) of the various rotor codes (see also the "Rotor detection" chapter).

After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.

Example:

PROG | RCF< | RPM | t/min:s  
**RemCycles= 16703**

If the maximum permissible number of running cycles has been exceeded, the following is displayed every time the centrifugation run is started, and the centrifugation run must be restarted.

PROG | RCF< | RPM | t/min:s  
**Cycles passed**

 If this display appears:

PROG | RCF< | RPM | t/min:s  
**Cycles passed**

the rotor or hanger in question must be immediately exchanged for a new rotor/hanger for safety reasons.

After the rotor/hanger has been exchanged, the cycle counter of the respective rotor code must be reset to zero again (see "Resetting the cycle counter to zero" chapter).

If the rotor code of a rotor is read in for the first time by the rotor detection, the centrifugation run then stops automatically after rotor detection. The rotor code (R) and the maximum speed (maxRPM) of the rotor permitted for this centrifuge are displayed.

Afterwards, the maximum permissible number of run cycles of the rotor/hanger must be set.

 In the case of swing-out rotors, the maximum permissible number of run cycles of the hangers must always be set.

The maximum permissible number of run cycles of the rotor/hanger must be set as follows.

### Example:

Display after rotor detection.

**R** : Rotor code of the rotor.

**maxRPM** : The maximum speed of the rotor permitted for this centrifuge.

PROG | RCF< | RPM | t/min:s  
**R 3 maxRPM 5000**

Press any key.

- Wait until the left LED in the  key is illuminated.
- Press any key.  
The maximum permissible number of run cycles will be displayed.

PROG | RCF< | RPM | t/min:s  
**Cyc lim=< 50000**

Continued on next page

## Example:

- Check whether the displayed maximum number of run cycles matches the specifications on the rotor/hangers.

If not, the maximum permissible number of run cycles specified on the rotor/hangers must be set.

- ☞ For rotors and hangers whose period of use is not limited to a maximum number of run cycles (centrifugation runs), zero must be set.


- Set the correct value with the keys.

- Press the key. The maximum permissible number of run cycles is saved.

## 20 Settings and queries

### Query:

- System information (section 20.1)

### Setting:

- Acoustical signal (section 20.2)
- Automatic unlocking of the lid after the centrifugation run (section 20.4)
- Automatic switch-off of the background lighting (section 20.5)
- Optical signal after ending the centrifugation run (section 20.3)
- Displayed centrifugation parameters after switching on the device (section 20.6)
- Activating/deactivating the "Dual time" function (section 20.7)

### Query:

- Hours of operation (section 20.8)
- Number of carried-out centrifugation runs (section 20.8)

### Setting:


- Cycle counter (section 20.9)


## 20.1 Querying system information

The following system information can be queried:

- Centrifuge model
- The maximum speeds of the various rotor codes,
- Centrifuge program version
- Frequency converter type
- Frequency converter program version

With the rotor at a standstill, the query can proceed as follows:

 The procedure can be aborted at any time by pressing the  key.

	<ul style="list-style-type: none"> <li>• Press and hold the button for eight seconds.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key. Centrifuge model.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• Press the key. <b>R</b> : Rotor code of the rotor. <b>*</b> : The rotor code of the currently used rotor is marked with an asterisk. <b>RPMmax</b> : The maximum speed of the rotor permitted for this centrifuge.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• The maximum speeds of the various rotor codes can be displayed with the keys.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• Press the key. Centrifuge program version.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• Press the key. Frequency converter type.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• Press the key. Frequency converter program version.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• Press the key twice to exit the "-&gt; Info" menu or press three times to exit the "* MACHINE MENU *".</li> </ul>	<p>Example:</p>  

## 20.2 Acoustic signal

The acoustic signal sounds:

- after a malfunction occurs in a 2-second interval.
- after ending the centrifugation run and the rotor is at a standstill in a 30-second interval.

The acoustic signal is ended by pressing any key.

If the rotor is at a standstill, the acoustic signal can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

	<ul style="list-style-type: none"> <li>• Press and hold the button for eight seconds.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key until the following is displayed:</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key.  <b>End beep</b> : Acoustic signal after ending the centrifugation run.  <b>off</b> : Signal deactivated.  <b>on</b> : Signal activated.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• With the keys, set <b>off</b> or <b>on</b>.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• Press the key.  <b>Error beep</b> : Acoustic signal after a malfunction occurs.  <b>off</b> : Signal deactivated.  <b>on</b> : Signal activated.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• With the keys, set <b>off</b> or <b>on</b>.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• Press the key.  <b>Beep volume</b> : Volume of the acoustic signal.  <b>min</b> : soft  <b>mid</b> : medium  <b>max</b> : loud</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• With the keys, set <b>min</b>, <b>mid</b>, or <b>max</b>.</li> </ul>	  

Continued on next page

- Press the key to save the setting.

- Press the key once to exit the "-> Settings" menu or press twice to exit the "\* MACHINE MENU \*".

Example:

### 20.3 Optical signal after ending the centrifugation run

The backlighting of the display flashes after the centrifugation run to visually signalize that the centrifugation run has finished.

The optical signal can be switched on or off when the rotor is at a standstill:

The procedure can be aborted at any time by pressing the key.

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until one of the following is displayed:  
**End blinking** : Backlighting of display flashes after the centrifugation run.  
**off** : Backlighting does not flash.  
**on** : Backlighting flashes.

- With the keys, set **off** or **on**.

- Press the key to save the setting.

- Press the key once to exit the "-> Settings" menu or press twice to exit the "\* MACHINE MENU \*".

Example:

### 20.4 Automatic unlocking of the lid after the centrifugation run

It can be set whether the lid should automatically unlock or not after the centrifugation run.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

	<ul style="list-style-type: none"> <li>Press and hold the button for eight seconds.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>Press the key until the following is displayed:</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>Press the key.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key until one of the following is displayed:  <b>Lid AutoOpen</b> : Automatic unlocking of the lid after the centrifugation run.  <b>off</b> : Lid does not unlock automatically.  <b>on</b> : Lid unlocks automatically.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>With the keys, set <b>off</b> or <b>on</b>.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key to save the setting.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key once to exit the "-&gt; Settings" menu or press twice to exit the "* MACHINE MENU *".</li> </ul>	
		<p>Example:</p> 

## 20.5 Backlighting of the display

To save energy, it can be set that, after a centrifugation run, the backlighting of the display switches off after 2 minutes.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.


	<ul style="list-style-type: none"> <li>Press and hold the button for eight seconds.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key until the following is displayed:</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key until one of the following is displayed:  <b>Power save</b> : Automatic shutdown of the backlighting.  <b>off</b> : Automatic shutdown deactivated.  <b>on</b> : Automatic shutdown activated.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>With the keys, set <b>off</b> or <b>on</b>.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key to save the setting.</li> </ul>		
↓			
↓			
	<ul style="list-style-type: none"> <li>Press the key once to exit the "-&gt; Settings" menu or press twice to exit the "* MACHINE MENU *".</li> </ul>		Example: 

### 20.6 Displayed centrifugation parameters after switching on the device

It can be set whether the centrifugation parameters of the last used program or of program 1 should be displayed after switching on the device.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

	<ul style="list-style-type: none"> <li>Press and hold the button for eight seconds.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>Press the key until the following is displayed:</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>Press the key.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key until one of the following is displayed:  <b>Start Pr</b> : Displayed program after switching on.  <b>First</b> : Program 1.  <b>Last</b> : Last used program.</li> </ul>	 
↓		
 	<ul style="list-style-type: none"> <li>Set <b>First</b> or <b>Last</b> with the keys.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key to save the setting.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>Press the key once to exit the "-&gt; Settings" menu or press twice to exit the "* MACHINE MENU *".</li> </ul>	

Example:

## 20.7 Activating/deactivating the "Dual time" function

In order to be able to set the runtime counter (parameter **Begins at START** and parameter **Begins at SPEED**), the "Dual time" function must be activated. This function is activated ex works.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

	<ul style="list-style-type: none"> <li>Press and hold the button for eight seconds.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key until the following is displayed:</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key until one of the following is displayed:  <b>Dual time</b> : Counting the runtime.  <b>on</b> : Function activated.  <b>off</b> : Function deactivated.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Set with the <b>off</b> or <b>on</b> keys.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key to save the setting.</li> </ul>		
↓			
	<ul style="list-style-type: none"> <li>Press the key once to exit the "-&gt; Settings" menu or press twice to exit the "* MACHINE MENU *".</li> </ul>		<p>Example:</p> 

## 20.8 Querying the hours of operation and the number of centrifugation runs

The hours of operation are divided into internal and external hours of operation.

Internal hours of operation: Total time the device was switched on.

External hours of operation: Total time of the previous centrifugation runs.

With the rotor at a standstill, the query can proceed as follows:

 The procedure can be aborted at any time by pressing the  key.

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.  
**TimeExt** : External hours of operation.

Example:

- Press the key.  
**TimeInt** : Internal hours of operation.

Example:

- Press the key.  
**Starts** : Number of all centrifugation runs.

Example:

- Press the key twice to exit the "-> Time & Cycles" menu or press three times to exit the "\* MACHINE MENU \*".

Example:

## 20.9 Resetting the cycle counter to zero

After the rotor/hanger has been exchanged, the cycle counter of the respective rotor code must be reset to zero again.

 The cycle counter may only be reset to zero if the rotor/hanger has been exchanged for a new rotor/hanger first.

With the rotor at a standstill, the cycle counter can be reset as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

	<ul style="list-style-type: none"> <li>• Press and hold the button for eight seconds.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key until the following is displayed:</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key.</li> </ul>	<p>Example:</p> 
↓		
	<ul style="list-style-type: none"> <li>• Press the key until the following is displayed: <b>Cyc sum</b> : Number of completed run cycles.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key. The number of completed run cycles is reset to zero.</li> </ul>	
↓		
	<ul style="list-style-type: none"> <li>• Press the key to save the setting.</li> </ul>	 
↓		
	<ul style="list-style-type: none"> <li>• Press the key twice to exit the "→ Time &amp; Cycles" menu or press three times to exit the "* MACHINE MENU *".</li> </ul>	<p>Example:</p>  

## 21 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left( \frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge container. For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

## 22 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm<sup>3</sup>

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm<sup>3</sup>.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm<sup>3</sup>

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

## 23 Emergency unlocking


In the event of a power failure, the lid cannot be unlocked with the motor. Emergency unlocking must be done by hand.

To unlock in an emergency, disconnect the centrifuge from the mains.  
Open the lid only when the rotor is at a standstill.

**CAUTION!** Damage to the lock during emergency unlocking by turning the hexagon Allen key in clockwise direction (to the right).  
The Allen key may only be turned counter clockwise (to the left); see figure.

- Switch off the mains switch (switch setting "0").
- Look through the window in the lid to make sure that the rotor is at a standstill.
- Insert the Allen key horizontally in the bore (A) and turn carefully counter clockwise (to the left) until the lid opens.
- Pull the Allen key back out of the bore.

## 24 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
  - have a pH between 5 - 8
  - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

### 24.1 Centrifuge (housing, lid and centrifuging chamber)

#### 24.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents:  
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

#### 24.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants:  
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

#### 24.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:  
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

## 24.2 Rotors and Attachments

### 24.2.1 Cleaning and care

- In order to avoid corrosion and changes in materials, the rotors and accessories have to be cleaned regularly with soap or with a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:  
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- The rotors and accessories have to be checked weekly for wear and corrosion.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

### 24.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:  
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

### 24.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:  
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

### 24.2.4 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

- usable until end of month/year: 10/2024 ("usable until end of month/year: 10/2024")
- "max. cycles: 40000" (max. run cycles: 40000)

For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs) of the various rotor codes. For a description, see the "Cycle counter" chapter.

### 24.3 Autoclaving

The following accessory can be autoclaved at 121°C / 250°F (20 min):

- Swing-out rotors without hangers
- Reductions

No statement can be made about the degree of sterility.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

### 24.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts, glass splinters and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

## 25 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

 <b>Perform a MAINS RESET:</b> <ul style="list-style-type: none"> <li>- Switch off the mains switch (switch position "0").</li> <li>- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").</li> </ul>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Message / fault		Cause	Remedy
No display		No voltage Mains input fuses defective.	<ul style="list-style-type: none"> <li>- Check distribution voltage.</li> <li>- Check mains power input fuse, refer to Chapter "Change mains input fuse".</li> <li>- Mains switch ON.</li> </ul>
IMBALANCE		The rotor is unevenly loaded.	<ul style="list-style-type: none"> <li>- Open the lid after the rotor is at a standstill.</li> <li>- Check the loading of the rotor, see chapter "Loading the rotor".</li> <li>- Repeat the centrifugation run.</li> </ul>
RPM > ROTOR MAX		Speed in the selected program greater than the maximum speed of the rotor.	<ul style="list-style-type: none"> <li>- Check the set speed.</li> <li>- Reduce the set speed</li> </ul>
MAINS INTERRUPT		Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> <li>- Open the lid after the rotor is at a standstill.</li> <li>- Push  button.</li> <li>- Repeat the centrifugation run if necessary.</li> </ul>
Rxxx WRONG ROTOR (xxx = Anzeige des Rotorcodes)		The used rotor is not approved for this device.	<ul style="list-style-type: none"> <li>- Use a rotor which is approved for this device.</li> </ul>
KEYBOARD ERROR		Error / defect electronics	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
TACHO ERROR	1, 2, 96.1	Failure of speed impulses during operation. No rotor inserted. Error / defect electronics	<ul style="list-style-type: none"> <li>- After the displayed waiting time expires (150 seconds), carry out a MAINS RESET.</li> <li>- Check whether the rotor is inserted.</li> </ul>
LID ERROR	4.1 – 4.127	Error in lid locking or lid closure.	<ul style="list-style-type: none"> <li>- Perform a MAINS RESET when the rotor has been stationary.</li> </ul>
OVER SPEED	5	Rotation too fast	
ROTOR ERROR	10.1 – 10.6	Incorrect rotor coding	
VERSION ERROR	12	Incorrect centrifuge model recognized. Error / defect electronics	
UNDER SPEED	13	Rotation too slow	


Message / fault		Cause	Remedy
CTRL ERROR	22 – 25.4	Error / defect electronics	– Perform a MAINS RESET when the rotor has been stationary.
CRC ERROR	27, 27.1		
COM ERROR	31 – 36		
FC ERROR	60, 61.1 – 61.21, 61.64 – 61.142		
TACHO ERR	61.22	Speed measurement error	– The device must not be switched off as long as the "wait" message is being displayed. After the message "wait ...!" is no longer displayed, carry out a MAINS RESET.
FC ERROR	61.23		
FC ERROR	61.153	Error / defect electronics	<ul style="list-style-type: none"> <li>– Perform a MAINS RESET.</li> <li>– Check the loading of the rotor, see chapter "Loading the rotor".</li> <li>– Repeat the centrifugation run.</li> </ul>
VERS. ERR	61.154		

## 26 Change mains input fuses

Switch off the mains switch and separate the device from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the device to the mains supply.

Model	Type	Fuse	Order no.
EBA 280	1101	T 3,15 AH/250V	E997
EBA 280 S	1102	T 3,15 AH/250V	E997
EBA 280	1101-01	T 6,3 AH/250V	2266
EBA 280 S	1102-01	T 6,3 AH/250V	2266

## 27 Returning Devices

Before returning the device, a transport securing device has to be installed.

If the device or its accessories are returned to Andreas Hettich GmbH & Co. KG, in order to provide protection for people, the environment and materials, it has to be decontaminated and cleaned before being shipped.

We reserve the right to refuse contaminated devices or accessories.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

## 28 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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1) Sample temp. rise during maximum speed and 1 hour running time  
 2) Observe the tube manufacturer's instructions







1142		1127-A							
<b>Ausschwingrotor 12-fach / Swing out rotor 12-times</b>   60°									
		0553							
									
Kapazität / capacity	ml	5	1,6 - 5	2,6 - 3,4	2,7 - 3				
Maße / dimensions $\varnothing \times L$	mm	13 x 75	13 x 75	13 x 65	11 x 66				
Anzahl p. Rotor / number p. rotor		12							
Drehzahl / speed	RPM	5000							
RZB / RCF	<sup>2)</sup>	2963							
Radius / radius	mm	106							
 9 (97%)	sec	10							
 9	sec	12							
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	10							

1148		1131-A				1132-A			
<b>Ausschwingrotor 8-fach / Swing out rotor 8-times</b>   90°									
		0553				0501		2079	
									
Kapazität / capacity	ml	5	1,6 - 5	2,6 - 3,4	2,7 - 3	6	4 - 5,5	10	4 - 7
Maße / dimensions $\varnothing \times L$	mm	13 x 75	13 x 75	13 x 65	11 x 66	12 x 82	15 x 75	17 x 70	16 x 75
Anzahl p. Rotor / number p. rotor		8							
Drehzahl / speed	RPM	5000							
RZB / RCF	<sup>2)</sup>	2991							
Radius / radius	mm	107							
 9 (97%)	sec	8							
 9	sec	10							
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	9							

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


- 1) Sample temp. rise during maximum speed and 1 hour running time  
 2) Observe the tube manufacturer's instructions




1137	1634		1633		1635				
<b>Winkelrotor 6-fach / Angle rotor 6-times</b>  36°	—								
	0521	0546	0519	0545					
Kapazität / capacity	ml	50	50	25	30	9-10	10	15	4 - 7
Maße / dimensions Ø x L	mm	34 x 100	29 x 107	24 x 100	26 x 95	16 x 92	15 x 102	17 x 100	13 x 100
Anzahl p. Rotor / number p. rotor		6							
Drehzahl / speed	RPM	6000							
RZB / RCF	<sup>2)</sup>	4025	3904	3703	3703	3783	3783	3783	3783
Radius / radius	mm	100	97	92	92	94	94	94	94
 9 (97%)	sec	20							
 9	sec	17							
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	7							

1137	1632	1641	1631	1635 + 1054-A			
<b>Winkelrotor 6-fach / Angle rotor 6-times</b>  36°							
	0578	0513	0509				
Kapazität / capacity	ml	7	50	15	1,6 – 5,0		
Maße / dimensions Ø x L	mm	12 x 100	29 x 115	17 x 120	13 x 75		
Anzahl p. Rotor / number p. rotor		18	3	6	6		
Drehzahl / speed	RPM	6000					
RZB / RCF	<sup>2)</sup>	3944	3824	3824	2978		
Radius / radius	mm	98	95	95	74		
 9 (97%)	sec	20					
 9	sec	17					
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	7					

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

- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions

1139		1054-A											
<b>Winkelrotor 12-fach / Angle rotor 12-times</b>  35°													
		0553	--									0518	
				 3)									
Kapazität / capacity	ml	5		12	1,6 – 5,0	1,1 – 1,4	2,6 – 3,4	2,7 – 3	4,5 – 5	4,9	4	15	
Maße / dimensions Ø x L	mm	12 x 75	13 x 75	17 x 102	13 x 75	8 x 66	13 x 65	11 x 66	11 x 92	13 x 90	12 x 60	17 x 100	
Anzahl p. Rotor / number p. rotor		12											
Drehzahl / speed	RPM	6000											
RZB / RCF	<sup>2)</sup>	3300	4146	3300	3300	3300	4146	3260	4146				
Radius / radius	mm	82	103	82	82	82	103	81	103				
 9 (97%)	sec	16											
 9	sec	16											
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	6											

1139					1058				6305		1063-6	
<b>Winkelrotor 12-fach / Angle rotor 12-times</b>  35°												
								0509	0507			2078
												
Kapazität / capacity	ml	7,5 – 8,5	9 – 10	10	4 - 7	8	8,5 - 10	15	4	0,5	1,5	2,0
Maße / dimensions Ø x L	mm	15 x 92	16 x 92	15 x 102	13 x 100	16 x 125	16 x 100	17 x 120	17 x 100	10 x 88	10,7 x 36	11 x 38
Anzahl p. Rotor / number p. rotor		12			12	6	12	6	12	12	12	
Drehzahl / speed	RPM	6000										
RZB / RCF	<sup>2)</sup>	4146			4146	4146	4146	4146	3502	2777	2737	
Radius / radius	mm	103			103	103	103	103	87	69	68	
 9 (97%)	sec	16										
 9	sec	16										
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	6										

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- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions.
- 3) without lid

1133							
<b>Winkelrotor 12-fach / (mit Dekantierhilfe) Angle rotor 12-times (with decanter aid)</b>  $\sphericalangle$ 35°		0553	0501	0578			
							
Kapazität / capacity	ml	5	6	7			
Maße / dimensions $\varnothing$ x L	mm	12 x 75	12 x 82	12 x 100			
Anzahl p. Rotor / number p. rotor		12					
Drehzahl / speed	RPM	5000					
RZB / RCF	<sup>2)</sup>	2879					
Radius / radius	mm	103					
 9 (97%)	sec	8					
 9	sec	10					
Probenerwärmung/Sample temp. rise	K <sup>1)</sup>	18					

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- 1) Sample temp. rise during maximum speed and 1 hour running time
- 2) Observe the tube manufacturer's instructions.