

Cefixime Tellurite Supplement (CT) ISO

Cat. 6064

A selective supplement for the isolation of *E. coli* O157:H7.

Practical information

Applications	Categories
Selective isolation	<i>Escherichia coli</i> O157

Industry: Clinical / Food

Regulations: ISO 16654



Principles and uses

Cefixime Tellurite Supplement is added to standard Sorbitol MacConkey agar media, which is an internationally accepted culture media for the selective isolation of *E. coli* O157:H7.

Macconkey Agar with Sorbitol (CT-SMAC) (Cat. 1099) is based on the formula developed by Rappaport & Hening. This medium is recommended for the research of *E. coli* O157:H7 in clinical and food testing. The composition is similar to MacConkey Agar (Cat. 1052) but the lactose has been substituted with sorbitol for differentiating enteropathogenic *E. coli* serotypes. These strains are typically sorbitol-negative. On standard MacConkey Agar containing lactose, this strain cannot be distinguished from other lactose fermenting *E. coli*.

E. coli O157:H7 does not ferment Sorbitol and therefore produces colorless colonies. As most of the other *E. coli* do ferment it, their colonies are pink.

The addition of cefixime and tellurite significantly reduces the number of sorbitol nonfermenters that need to be screened during the attempted isolation of *E. coli* O157:H7. Cefixime inhibits *Proteus* spp. and tellurite inhibits non-O157 *E. coli* and other organisms, thus improving the selectivity of the medium.

Formula per vial

Potassium tellurite (mg)	1,25	Cefixime (mg)	0,025
--------------------------	------	---------------	-------

Preparation

Aseptically reconstitute 1 vial with 5 ml of sterile distilled water. Mix gently until complete dissolution and aseptically add to 500 ml of MacConkey Agar with Sorbitol (CT-SMAC) (Cat. 1099), autoclaved and cooled to 45-50 °C. Mix well and distribute into sterile containers.

When Cefixime Tellurite Supplement (CT) is required to be added to another media, like *E. coli* O157:H7 Chromogenic Agar Base (Cat. 1588), refer to the specific instructions of the medium for the quantity of Cefixime Tellurite Supplement (CT) that should be added.

Instructions for use

- » For clinical diagnosis, the type of sample is feces.
- Inoculate on the surface. Parallel striae with the handle or hyssop.
- Incubate in aerobic conditions at 35±2 °C for 18-24 hours.
- Reading and interpretation of the results.
- » For other uses not covered by the CE marking:
Detection of *Escherichia coli* O157 according to ISO 16654:
 - Prepare the initial suspension adding the test portion to Trypticasein Soy Broth Modified with Novobiacin (Cat. 1292) prewarmed to 41,5 °C to obtain a ratio of 1/10.
 - Separate and concentrate the microorganisms by immunogenic particles coated with antibodies to *E. coli* O157.
 - Incubate for 6 hours then a further 12 hours to 18 hours at 41,5 °C.
 - Subculture the immunomagnetic particles with the bacteria adhered on MacConkey Agar with Sorbitol (Cat. 1099) and a second selective isolation agar of choice by the laboratory. The optimal incubation temperature for *E. coli* O157 is 37±1 °C for 18-24 hours.
 - Confirm by indol production in Tryptophan Culture Broth (Cat. 1237) and agglutination with the serum anti *E. coli* O157.

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
w/o rests	White lyophilized	N/A	Transparent	N/A

Microbiological test

According to ISO 11133:

Incubation conditions: (37±1 °C / 21±3 h).

Inoculation conditions: Productivity qualitative (10³ - 10⁴ CFU) / Selectivity (10⁴- 10⁶ CFU).

Microrganisms	Specification	Characteristic reaction
Escherichia coli ATCC 25922	Inhibición parcial (1)	Growth of some pink colonies
Staphylococcus aureus ATCC 25923	Total inhibition (0)	
Escherichia coli O157:H7 ATCC 700728	Good growth (2)	Transparent colonies with a pale yellowish-brown appearance and a diameter ~1 mm

Storage

Temp. Min.:2 °C

Temp. Max.:8 °C

Bibliography

Doyle, M.P. and J.L. Schoeni. 1987. Applied Environmental Microbiology 53:2394-2396.

J. G Wells et al, 1991. Isolation of Escherichia coli serotype O157:H7 and other Shiga-like-toxin-producing E. coli from dairy cattle.