

KPC Chromogenic Medium

Cat. 2063

Chromogenic medium for detection of Gram-negative with reduced susceptibility to most of the carbapenem agents.

Practical information

| Applications | Categories |
|--------------|-----------------------------------|
| Detection | Bacteria resistant to carbapenems |

Industry: Clinical



Principles and uses

Klebsiella pneumoniae carbapenemase (KPC)-producing bacteria are a group of emerging highly drug-resistant Gram-negative bacilli causing infections associated with significant morbidity and mortality. Carbapenem antibiotics are generally not effective against KPC-producing organisms.

Although *K. pneumoniae* remains the most prevalent bacterial species carrying KPCs, the enzyme has been identified in several other gram-negative bacilli. Infections caused by bacteria-producing carbapenemases are becoming an increasingly significant problem worldwide because they are often not detected by routine susceptibility screening and possess an exceptional potential for dissemination. Infections caused by these organisms present clinicians with serious treatment challenges, due to limited antibiotic options.

Peptones and growth factors provide nitrogen, vitamins, minerals and amino acids essential for growth. Chromogenic mixture allows the identification of Gram-negative bacteria with a reduced susceptibility to the carbapenem agents. The supplement inhibits the growth of all the KPC non-producing bacteria.

Characteristics of KPC colonies:

- *Escherichia coli*: colonias pink.
- *Enterobacter aerogenes*: dark blue.
- *Klebsiella pneumoniae*: dark blue.

Formula in g/L

| | | | |
|----------------------|----|---------------------|----|
| Bacteriological agar | 16 | Chromogenic mixture | 3 |
| Peptone | 14 | Growth factors | 15 |

Preparation

Suspend 48,0 grams of medium in one liter of distilled water. Mix well and dissolve by heating with frequent agitation. Boil for one minute until complete dissolution. Sterilize in autoclave at 121 °C for 15 minutes. Cool to 50 °C and aseptically add two vials of KPC supplement (Cat. 6044) previously reconstituted with 5 ml of sterile distilled water. Mix well and dispense into plates.

Instructions for use

For clinical diagnosis, the type of sample is urine, lung aspirations and rectal samples.

- Inoculate on the surface making parallel striae with the handle or hyssop.
- Incubate in aerobic conditions at 35±2 °C for 18-24 hours.
- Reading and interpretation of results.

Note: It's important to notice that, as it happens in other chromogenic media, bacteria with atypical KPC enzyme may produce anomalous reactions in this medium.

Quality control

| Solubility | Appearance | Color of the dehydrated medium | Color of the prepared medium | Final pH (25°C) |
|------------|------------|--------------------------------|------------------------------|-----------------|
|------------|------------|--------------------------------|------------------------------|-----------------|

Microbiological test

Incubation conditions: (35±2 °C / 18-24 h).

| Microorganisms | Specification | Characteristic reaction |
|----------------------------------|--------------------|-------------------------|
| Klebsiella pneumoniae ATCC 13883 | Total inhibition | |
| Klebsiella BAA 1705 | Good growth | Blue colonies |
| Enterococcus faecalis ATCC 19433 | Partial inhibition | Light blue colonies |
| Escherichia coli ATCC 2469 | Good growth | Pink colonies |
| Escherichia coli ATCC 25922 | Total inhibition | |
| Staphylococcus aureus ATCC 25923 | Total inhibition | |
| Proteus mirabilis ATCC 25933 | Total inhibition | |

Storage

Temp. Min.:2 °C

Temp. Max.:8 °C

Bibliography

- Ryan S. Arnold, MD, Kerri A. Thom, MD, MS, [...], and Daniel J. Morgan, MD Emergence of Klebsiella pneumoniae Carbapenemase (KPC)-Producing Bacteria.
- Paterson DL, Bonomo RA. Extended-spectrum beta-lactamases: a clinical update. *Clinical Microbiology Reviews*. 2005;18:657–686. [PMC free article] [PubMed]
- Osterblad M, Kirveskari J, Koskela S, et al. First isolations of KPC-2-carrying ST258 Klebsiella pneumoniae strains in Finland, June and August 2009. *Euro Surveill*. 2009;14(40):19349. [PubMed]
- Martín-Gil J, Villa FM, Ramos-Sánchez MC, Martín-Gil FJ. "Studies on beta-lactam antibiotics - Differential thermal-analysis of Cephalosporins". *J. Thermal Anal Cal*, 1984, 29 (6): 1351-1357.
- Sebastian Droguett Perez (Dr.2)(2004)*Rossi S. (Ed.) (2004). *Australian Medicines Handbook 2004*. Adelaide: Australian Medicines Handbook. ISBN 0-9578521-4-2. (* Rossi S (Ed.) (2004).