

Reference: 0974

Technical Data Sheet Product: COLUMBIA CNA + 5% SHEEP BLOOD AGAR

Specification

Selective medium from isolation of gram-positive cocci of clinical samples.

Presentation

Shelf Life Storage **Packaging Details 20 Prepared Plates** 90 mm 2-14 ºC 3 months 1 box with 2 packs of 10 plates/pack. Single cellophane.

with: 21 ± 2 ml

Composition

Composition (g/l):	
Casein pancreatic digest	10.000
Meat peptic digest	5.000
Heart Pancreatic digest	.3.000
Yeast Extract	5.000
Sodium chloride	5.000
Starch	1.000
Agar	15.000
Colistin sulphate	.0.010
Nalidixic ac	0.015
Defibrinated Sheep Blood	.50 ml

Description / Technique

Description:

This medium is a modification of the Blood Columbia Agar with the selective antimicrobial agents Colistin sulfate and Nalidixic Acid added. These antibiotics inhibit the growth of Enterobacteriaceae and Pseudomonas.

Techniques recommended use:

Inoculate the samples directly on the surface of agar, streaking to obtain isolated colonies. Some stab inoculations should also be carried out to deposit Beta-haemolytic streptococci deep in the medium as this subsurface growth allows manifestation of both oxygen-stable and oxygen-labile streptolysin activity, giving clear haemolytic reactions.

The plates are incubated in (aerobic, anaerobic or 5-10% CO, enriched atmosphere) according to laboratory protocol, for each sample type. After incubation for 18 to 24 /48 hours at 37°C the plates are examined for growth and, subsequently, for haemolytic reactions:

- Alpha-haemolysis (a) is the reduction of haemoglobin to methaemoglobin in the medium surrounding the colony, producing a green halo.
- Beta-haemolysis (b) is the total lysis of the blood erythrocytes producing a clear zone around the colony.
- Gamma-haemolysis (g) is indicated by no haemolysis: No change in the environment.
- Alpha-prime-haemolysis (a) presents as a zone of complete lysis next to the colony surrounded by an area of partial lysis.

The haemolytic effect of streptococci depends on many factors. Ruoff (1995) noted that incubation in atmospheres enriched in (5-10%) CO2 optimized the action of beta-haemolytic streptococci and some strains of streptococci, (Lancefield group D) behave differently depending on the animal origin of the blood used in the medium: In Blood Agar with horse, human, or rabbit blood, beta-haemolytic action is manifested and with sheep blood alpha-haemolytic action is best observed.

Different animal blood source, greater incubation times, humidity or larger percentage of carbon dioxide in atmosphere,... may be required depending on the sample, on the specifications of the laboratory, the expected isolations to be found.

Quality control

Physical/Chemical control

Red pH: 7.2 ± 0.2 at 25°C Color:

Microbiological control

Spiral Spreading: Practical range 100 ± 20 CFU. min. 50 CFU (productivity) / 104-106 CFU (selectivity). Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020

Microaerofilic incubation at 35 ± 2 °C for 18-24h

Microorganism

Streptococcus pyogenes ATCC® 19615 Stph. aureus ATCC® 25923, WDCM 00034 Streptococcus pneumoniae ATCC® 49619 Proteus mirabilis ATCC® 12453

Sterility Control

Incubation 48 hours at 30-35 °C and 48 hours at 20-25 °C: NO GROWTH. Check at 7 days after incubation in same conditions.

Growth

Good (Beta-hemolysis) Good (Beta-hemolysis) Good (Alfa-hemolysis) Inhibited

Page 1 / 2 Revision date: 13/12/21



Reference: 0974

CE IVD

Technical Data Sheet

Product: COLUMBIA CNA + 5% SHEEP BLOOD AGAR

Bibliography

- · ATLAS, R.M. & L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press. London.
- · BARON, E.J., L.R. PETERSON & S.M. FINEGOLD (1994) Bailey & Scott's Diagnostic Microbiology. 9th ed. Mosby-Year Book Inc. St Louis, MO. USA.
- · CASMAN, E. (1947) A non-infusion blood agar base for neiseriae, pneumococci and streptococci. Am. J. Clin. Path. 17:281-289.
- · ELLNER, P.D., C.J. STOESSEL, E. DRAKEFORD, & F. VASI (1966) A new culture medium for medical bacteriology. Amer.J.Clin.Path 45:502-504.
- · ESTEVEZ, E.G. (1984) Bacteriological Plate Media: review of mechanisms of action. Lab. Med. 15:258-262.
- · ISENBERG H.D. (1992). Clinical Microbiology Procedures Handbook. Vol 1 ASM Washington DC, USA.
- . ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- · RUOFF, K.L. (1995) Streptococcus p. 299-305. En Manual of Clinical Microbiology 6th ed. Por Murray, Baron, Pfaller, Tenover y Yolker (editors) ASM. Washington DC. USA.

Page 2 / 2 Revision date: 13/12/21