## Violet Red Bile Lactose + Neutralizing

Selective medium for environmental and personnel hygiene monitoring with inactivation of disinfectants.

Instructions For Use

ENGLISH

### DESCRIPTION

Violet Red Bile Lactose (VRBL) + Neutralizing is a prepared plated medium used for detection and enumeration of coliform bacteria in surfaces testing and air microbial control by dynamic sampling technique.

The composition of the base culture medium complies with the specifications given by the American Public Health Association and ISO 4832 with the addition of neutralizing agents, which are included in the medium to inactivate residual disinfectants allowing detection of microorganisms surviving after treatment of surface and material with antiseptics.

### TYPICAL FORMULA\*

	(g/litre)
Enzymatic Digest of Animal Tissues	7.0
Yeast Extract	3.0
Lactose	10.0
Sodium Chloride	5.0
Bile Salts	1.5
Neutral Red	0.03
Crystal Violet	0.002
Agar	14.0
Histidine	1.0
Lecithin	0.7
Polysorbate 80	5.0
Sodium Thiosulfate	0.5

Final pH 7.4 ± 0.2 at 25°C

\*Adjusted and/or supplemented as required to meet performance specifications.

### METHOD PRINCIPLE

Enzymatic digest of animal tissues provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Lactose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium. Bile salts and Crystal violet are selective agents effective against Gram-positive cocci. Neutral red is the pH indicator. Agar is the solidifying agent. Histidine inactivates aldehydes. Lecithin neutralizes quaternary ammonium compounds. Histidine, lecithin, polysorbate 80 (Tween 80) and sodium thiosulfate are neutralizers which have demonstrated to be effective for neutralisation of disinfectant residues containing the following active agents: aldehydes, quaternary ammonium compounds, biguanides, phenolic and related compounds, mercurial derivatives, halogens and oxidizing compounds.

### TEST PROCEDURE

**For Active Air Sampling**, place the plate, lid up, in the air sampler, remove the lid of the plate and sample a specific volume of air according to laboratory procedure.

**For Surface Testing**, gently press the agar surface on the test surface for about 10 seconds with a steady pressure. Do not move laterally the plate. Residues of the medium should be subsequently removed from the area tested.

After sampling close the plate with the lid and incubate in ambient air at 30°C or 37°C, depending on the organisms under study, for 24 ± 2 hours. Other temperatures and times of incubation may be chosen to suit specific requirements of target organisms but such conditions should be validated under the specific protocol used for environmental monitoring testing.

### INTERPRETING RESULTS

Observe for the formation of purplish-red colonies with a diameter of at least 0.5 mm and record the number of CFU per plate.

Colonies should be further isolated and identified with appropriate biochemical tests and/or serological procedures.

### APPEARANCE

Slightly opalescent, reddish-purple.

**STORAGE**

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

Avoid quick temperature shifts to prevent condensation.

**SHELF LIFE**

9 months.

**QUALITY CONTROL**

To check the performance of the medium the following microbial reference strains can be used.

Control strain		Inoculum	Incubation	Criteria	Specification
<i>Escherichia coli</i>	WDCM 00012	50-100 CFU	24 ± 2 h / 30 ± 1°C	P <sub>R</sub> ≥ 0.5	Purplish-red colonies with or without precipitation halo
<i>Escherichia coli</i> + 50 µl Aerodesin 2000	WDCM 00012				
<i>Enterococcus faecalis</i>	WDCM 00009 or WDCM 00087	10 <sup>4</sup> -10 <sup>6</sup> CFU		Total inhibition	—
<i>Pseudomonas aeruginosa</i>	WDCM 00025	10 <sup>3</sup> -10 <sup>4</sup> CFU		—	Colourless to beige colonies

A productivity ratio (P<sub>R</sub>) of 0.5 is equivalent to a recovery rate of 50%.

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**BIBLIOGRAPHY**

See the references at the end of this document.

**TABLE OF SYMBOLS**

See the table of symbols at the end of this document.

**The product is available in the configurations listed below.** There may be additional product ref. numbers as well. For an updated listing of available products, visit [liofilchem.com](http://liofilchem.com)

Product Description	Format	Packaging	Ref.
Violet Red Bile Lactose + Neutralizing	Contact Plate	20 plates	15361

This document is available from the online Support Center:

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[liofilchem.com/ifu-sds](http://liofilchem.com/ifu-sds)

**BIBLIOGRAPHY / BIBLIOGRAFIA / BIBLIOGRAPHIE**

1. EN ISO 11133:2014+Amd1:2018. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. ISO 4832:2006. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – Colony count technique.
3. Davidson, Roth, and Gambrel-Lenarz (2004) In Wehr and Frank (ed.) Standard methods for the microbiological examination of dairy products, 17<sup>th</sup> ed. American Public Health Association, Washington, D.C.
4. Kornacki and Johnson (2001) In Downes and Ito (ed.) Compendium of methods for the microbiological examination of foods, 4<sup>th</sup> ed. American Public Health Association, Washington D.C.

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