Millipore®

SigmaAldrich.com

User Guide

Millex®-SV Filter Unit (Sterile)

REF **SLSV025LS**

English
Français
Italiano
Deutsch
Español
Português
Nederlands
Dansk
Svenska
Latviski
Magyar
Česky
Polski
繁體中文
Norsk
Türkçe
Slovenčina

CE 0459

℞only

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

PR06219 03/25

Date of issue: 2025-03

Symbol	Definition	Symbol	Definition
STERILE EO	Sterilized using ethylene oxide		Use-by date
	Non-pyrogenic		Date of manufacture
	Consult instructions for use		Manufacturer
	Do not use if package is damaged and consult instructions for use		Caution
	Do not resterilize	CE 0459	CE conformity marking
	Do not re-use	℞only	Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.
REF	Catalogue number	CH REP	Swiss authorised representative
LOT	Batch code	UDI	Unique device identifier
	Temperature limit		

Made in Ireland

Merck Millipore Ltd.
Tullagreen, Carrigtwohill,
Co. Cork, Ireland
T45 KD29
an affiliate of Merck KGaA,
Darmstadt, Germany

CH REP
MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

Brazilian Registration Holder (BRH)
Merck S/A
Mirante de Cajamar Business Center, South Access Road,
Km 30/Anhanguera Highway,
KM 29.5-Km 29 plus 503 M, Warehouse Module B4 – Warehouse 01,
Polvilho – Cajamar/SP/Brazil
SP/ Brazil
ZIP CODE: 07790-330
CNPJ no 33.069.212/0008-50

The initial M, Millipore, Millex, Durapore and Sigma Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

English

Millex®-SV Filter Unit (Sterile)

SLSV025LS Qty: 50 pack

• 25 mm • Sterile • Single use only • Non-pyrogenic


Intended Purpose
The Millex®-SV Filter Unit (Sterile) is a medical device intended for use as a bidirectionally supported syringe filter to clarify low volume solutions in direct patient care and pharmacy admixture applications. Device is intended for medical and/or laboratory professional use only.

Introduction
This filter is a low-binding, bidirectionally supported syringe-operated filter for the filtration of aqueous, proteinaceous solutions where minimal sample loss is desired. Clinical benefits include the removal of particles, precipitates, and undissolved powders larger than the membrane's rated pore size from solutions for clinical use. This sterile, single-use filter is non-pyrogenic and non-toxic. It consists of a 5.0 µm Durapore® membrane filter sealed in a polyvinyl chloride (PVC) housing. Typical applications include particle or aggregate removal from protein solutions such as monoclonal antibodies, extracted cellular or tissue-based solutions, hormones, and drugs.

Chemical Compatibility
This filter is compatible with most aqueous solutions. Compatibility of solution should be tested before use. Visit SigmaAldrich.com/MedMillex for more information.

WARNINGS
- Filter solutions with a temperature limit of 45 °C (113 °F).
- Syringe filters are for manual use only; do not use on automated systems.
- Not suitable for filtering high viscosity solutions or blood.
- Do not filter solutions containing 5 mg or less of active drug materials unless binding studies have been performed.
- To ensure sterility, do not use this product if the package is damaged.
- Use aseptic technique to maintain sterility when removing product from primary package.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter and/or personal injury.
- Do not use for direct neuraxial applications which may come in contact with cerebrospinal fluid.
- Do not use the same syringe filter to filter solutions in both directions.
- Sudden loss of pressure could indicate failure of the filter.
- Single use only; do not re-use or resterilize.
- Make sure to wet the filter membrane thoroughly before injecting the solution; improperly wetted filters can become airlocked.
- Not suitable for use with medical pressure infusion equipment or gravity fed intravenous (IV) administration sets.
- Discard appropriately after single use. See "Disposal" below.

Directions for Use

1 Wearing gloves using aseptic technique: fill syringe with solution to be filtered.
2 Remove cover from package.
3 Attach suitable syringe to filter inlet and remove assembly from package.
4 Hold syringe with filter pointing up and top off by pushing a few drops through to wet filter. ⚠ Excess solution may be hazardous and should be disposed of with care.
5 Depress syringe plunger to filter solution into suitable container or compatible attachment on outlet.

Specifications and Performance Characteristics

Materials	
Membrane	Hydrophilic Durapore® polyvinylidene fluoride (PVDF), type SV
Pore size	5.0 µm
Housing	Polyvinyl chloride (PVC)
Dimensions	
Inlet to outlet	25 mm (0.98 in.)
Diameter	29 mm (1.14 in.)
Filtration area	4 cm² (0.62 in²)
Temperature limit	45 °C (113 °F)
Pressure limit at 21 °C	5.2 bar (75 psi) inlet and differential
Filtration volume	1–100 mL
Hold-up volume	≤0.1 mL after air purge
Sterilization method	Ethylene oxide gas
Connections	Female Luer-Lok™ inlet, male Luer-slip outlet
Flow rate at 2.1 bar (30 psi), 21 °C	≥600 mL/min

Notice
We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Contact Information
For the location of the office nearest you, go to SigmaAldrich.com/offices.

Disposal
Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance
Visit the tech service page on our web site at SigmaAldrich.com/techservice.

Any serious incident involving this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty
The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

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