DOMINIQUE DUTSCHER SAS

AGITADOR ROTATIVO «MOVIL ROD»

TUBE SHAKER «MOVIL ROD»

7001723

J.P. SELECTA s.a. Ctra. NII Km 585.1 Abrera 08630 (Barcelona) España Tel 34 937 700 877 Fax 34 937 702 362
e-mail: selecta @jpselecta.es - website: http://www.jpselecta.es

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2 Installation

Unit description:

Shaker with 7 rollers with double motion action: rotation and swing, this provide a suitable homogenize effect.

Suitable for use in liquid suspension, biology and haematology

Application examples:

- Suitable for mixing blood and anticoagulant.
- Mixing of substances as dry powder or freezed.
- Liquid mixing.

Location:

- Place on flat and levelled surface. Suitable for unit weight.
- Let free 20cm, at least, around the unit.
- Near to power socked (1m).

Connection to power supply:

- Check that power socked have suitable voltage, according unit plate.
- Power socked have earth connector. Unit must be grounded.
- Unit could be dangerous is not properly grounded.

3 Technical specification

Rolling and vibrating speed:	From 10...80 rpm.
Vibrating range:	22mm±1mm
Power consumption:	12W
Fuse:	2A (Ø5 x 20mm)
Voltage:	230V AC ±10% (50/60Hz)
Weight:	3.6Kg
Dimensions:	394 x 266 x 98 mm

4 Operation

- Place the tubes to stirr.
- Press the main switch
- Adjust the speed.

Near once a year, lubricate ,slightly, the bearings.

5 Warranty

The warranty covers the defective components or faulty assembling of the unit for 1 year since the purchasing date. Keep purchasing bill as the only valid document for the purchasing date.

The warranty does not covers the damages caused for a unit malfunction.

Return goods are not accepted without previous J.P. SELECTA,s.a. approval.

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