Operating Instructions

CO_2 Incubator

MCO-173AIC
MCO-173AICUV
MCO-173AICUVH

DOMINIQUE DUTSCHER SAS

Please read the operating instructions carefully before using this product and keep the operating instructions for future use.

See page 92 for the model number.

CONTENTS

1. BEFORE USING

INTRODUCTION 4
INTENDED USE AND PRECAUTIONS 4
SAFETY PRECAUTIONS 5
LABELS ON INCUBATOR 11
INCUBATOR COMPONENTS 12
- Main body 12
- Accessories 16

INSTALLATION 17
- Installation site 17
- Installing the incubator 18
- Connecting to the power inlet connector 19
- When stacking an incubator 20

CLEANING 21
- Removing the components 21
- Cleaning the chamber and components 23
- Cleaning the control panel 24
- Installing the cleaned components 24

PREPARATION BEFORE USING 26
- Connecting a CO_2 gas cylinder 26
- Installing a remote alarm device 29
- Filling the humidifying pan 30

START-UP PROCEDURE 31
FOR OPTIMAL CULTIVATION 32
- Precautions for cultures 32
- Preventing contamination 32

2. TOUCH SCREEN OPERATION

USING THE TOUCH SCREEN 33
- Basic operation 33
- Entering alphanumeric values and symbols 33
- Other operations 35

HOME SCREEN 36
MENU BAR BUTTONS 39

3. SETTINGS AND OPERATION

SETTING FOR BASIC OPERATION 42
- Setting temperature, CO_2 density, and high temperature limit 42
- Configuring alarm settings 43
- Alarm and buzzer stop button 45

DISINFECTION BY UV-LED 46
- Setting the UV-LED 46
- Using UV-LED 47

SCREEN LOCK 48
- Setting screen lock 48
- Operation when screen lock is ON 49
- Canceling screen lock setting 50

DOMINIQUE DUTSCHER SAS

AUTOMATIC DOOR LOCK 51
Setting user ID and password 51
Changing the password 53
Deleting user ID and password 53
Setting automatic door lock 54
Unlocking the outer door 55
Unlocking and locking the outer door during a power failure 56
Canceling the automatic door lock 57
SETTING DATE, TIME, AND TIME ZONE 58
SETTING BRIGHTNESS AND SLEEP 60
ANALOG OUTPUT 61
Setting the data to be output 61
Connecting to external device 61
Analog data output specification 62
SETTING UNIT NAME 63
REFERRING TO THE OPERATING INSTRUCTIONS USING MOBILE DEVICE 64
EVENT LOG 65
Viewing the event log 65
Exporting event log data 67
OPERATION LOG CHART 69
Displaying operation log chart 69
Exporting the operation log chart data displayed on the screen 71
EXPORTING LOG DATA 73
GAS AUTO CHANGER (OPTION) 75
Connecting CO_2 gas cylinders 75
Automatic CO_2 gas supply line switching 76
Manual CO_2 gas supply line switching 77

4. MAINTENANCE & TROUBLESHOOTING

ROUTINE MAINTENANCE 78
Cleaning the exterior, interior, and inner attachments 78
H_2O_2 decontamination 79
DURING/AFTER POWER FAILURE 84
Operation during power failure 84
Operation after recovery from power failure 84
ALARMS AND SELF-DIAGNOSIS 85
Alarm messages 85
Notification messages 86
TROUBLESHOOTING 87
SPECIFICATIONS 90
PERFORMANCE 92
SAFETY ENVIRONMENTAL CONDITIONS 93
EMC SPECIFICATION 94
SAFETY CHECK SHEET 97

DISTRIBUE DUTSCHER SAS

INTRODUCTION

■ Read the operating instructions carefully before using the product and follow the instructions for safe operation.

■ PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.

■ Keep the operating instructions in a suitable place so that they can be referred to as necessary.

■ The operating instructions are subject to change without notice for improvement of performance or function.

■ We accept no liability for any damages, including indirect damages such as losses due to differences in instructions, which may arise from discrepancies between the operating instructions on our website and those included with the product you purchased.

■ Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.

■ No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

IMPORTANT NOTICE

PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for:

- any loss or damage to the contents of the product
- any indirect damage caused by data damage or loss

INTENDED USE AND PRECAUTIONS

■ For medical use:
<Intended Use>
This equipment is designed for culture of cells, tissues, organs, and embryos.

<Intended User>
Intended for healthcare professionals.

It is subject to only Austria, Belgium, Cyprus, France, Germany, Ireland, Italy, Liechtenstein, Luxembourg, Malta, the Netherlands, Poland, Spain, Switzerland, the UK, and Malaysia.

■ For laboratory use:
This equipment is designed for cell and tissue culture for laboratory use.

■ The adapted culture condition depends on the kind of sample. It is necessary to determine the culture temperature, CO_2 density, and culture period that is required.

■ For culture of embryos, low O_2 density is desirable. Using an O_2/CO_2 incubator is recommended.

■ For IVF/ART use, special attention should be paid to traceability, since incidents may be disclosed several months or years later, at the baby's birth, or during its later life. It is therefore recommended that the following data is maintained: product serial number, incubation term, and incubation parameters. (Refer to details on MEDDEV 2.2/4.)

DOMINIQUE DUTSCHER SAS

SAFETY PRECAUTIONS

Be sure to observe the operating instructions as they contain important safety advice. If the product is used in a manner not specified by the operating instructions, the protection provided by the product may be impaired.

For correct and safe use of the product, follow the precautions and procedures in these operating instructions carefully. Failure to do so could result in injury or damage to the product.

Precautions are illustrated in the following way:

⚠ WARNING

Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

⚠ CAUTION

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

The following symbols are used in this document and some of them are attached to the unit.

	Actions are prohibited.
	Actions are mandatory.
	Caution must be taken.
	This symbol indicates the possibility of an electric shock. High-voltage electrical components are placed under the covers. Only a qualified engineer or service personnel should be allowed to open these covers.
	This symbol indicates a risk of burns from contact with hot surfaces.
	This symbol indicates a risk of an electric shock by touching the product with wet hands.
	This symbol indicates an earth terminal. Connect the earth terminal to the ground to prevent an electric shock.
	This symbol indicates a risk of an electric shock by an electric leakage caused by a wet product.
	This symbol indicates a risk of injuries such as an electric shock caused by disassembling the product.
	This symbol indicates that the user must disconnect the mains plug for the purposes of maintenance, in the case of malfunction or when left unattended.

DOMINIQUE DUTSCHER SAS

SAFETY PRECAUTIONS

	This symbol indicates the “ON” for the power switch.
○	This symbol indicates the “OFF” for the power switch.

In the event of a serious incident (accident) occurring in relation to this product, contact PHC Europe B.V. and your national competent (regulatory) authority.

⚠WARNING

As with any equipment that uses CO_2 gas, there is a likelihood of oxygen depletion in the vicinity of the equipment. It is important that you assess the work site to ensure there is suitable and sufficient ventilation. If restricted ventilation is suspected, then other methods of ensuring a safe environment must be considered. These may include atmosphere monitoring and warning devices.

DOMINIQUE DUTSCHER SAS

⚠WARNING

Installation

- **Do not use the unit outdoors.**
 Exposure to rain may cause an electric leakage and/or an electric shock.
- **Do not install the unit in a location where flammable or volatile substances are present.**
 Installing the unit in such a location may cause explosions and/or a fire.
- **Do not install the unit in a location where there are high levels of moisture or where it may be splashed with water.**
 This may cause the insulation to deteriorate, leading to an electric leakage and/or an electric shock.
- **Do not install the unit in a location where corrosive gases such as acids are present.**
 Deterioration of the insulation due to corrosion of the electric components may cause an electric leakage or an electric shock.
- **Do not install the unit in a location where there are water or steam pipes above it.**
 This may cause the insulation to deteriorate, leading to an electric leakage and/or an electric shock.
- **Do not leave the plastic bags used for packing in a place where they can be reached by children.**
 This may result in unexpected accidents such as suffocation.

- **Only qualified engineers or service personnel should install the unit.**
 Installation by unqualified personnel may cause water leakage, an electric shock, or a fire.
- **Install the unit in a location capable of bearing the total combined weight (product + optional accessories + stored items). After installing the unit, be absolutely sure to take precautions to prevent the unit from falling over.**
 If the unit is installed in a location which is not strong enough or if the proper precautions are not taken, the unit may fall over and cause injuries.
- **Connect the unit to a dedicated power source as indicated on the rating label attached to the unit.**
 Use of any other voltage or frequency other than that on the rating label may cause a fire or an electric shock. Also, using an extension cable or power strip may cause a fire resulting from abnormal heating.
- **When handling harmful samples (for example, those which consist of toxic, pathogenic or radioactive substances), install the unit inside a designated isolation facility.**
 If the unit is installed in a location other than an isolation facility, there may be detrimental effects on both people and the natural environment.
- **Install the unit in a well-ventilated (airy) location.**
 As with any equipment that uses CO_2 gas, there is a likelihood of oxygen depletion in the vicinity of the equipment. It is important that you assess the work site to ensure there is suitable and sufficient ventilation. If lack of ventilation is suspected, then other methods of ensuring a safe environment must be considered. These may include atmosphere monitoring systems and warning devices with alarms.
- **This unit must be plugged into a dedicated circuit protected by branch circuit breaker.**
 Otherwise it may cause an electric shock or a fire.
- **When stacking incubators, use the hooks on the back panel of the incubator and secure it to the wall by passing a strong rope or chain through the hooks.**
 The incubator may fall over and cause injuries.

- **Be absolutely sure to earth (ground) the unit to prevent an electric shock.**
 Failure to earth the product may cause an electric shock. If necessary, ask a qualified contractor to do this work.
- **Do not connect the earth wire to a gas pipe, water pipe, or lightning rod when earthing the unit.**
 Earthing the unit improperly may cause an electric shock.

DOMINIQUE DUTSCHER SAS

SAFETY PRECAUTIONS

⚠WARNING

Power supply plug and cord

- **Never damage the power supply cord or power supply plug (by breaking, adapting, placing near a source of heat, bending with force, twisting, pulling, adding weight, or binding).**
 A damaged power supply cord or power supply plug may cause an electric shock, a short circuit, or a fire. Contact our sales representative or agent for repairing the cord and plug.
- **Do not place this unit or other equipment in a position where it is difficult to plug in or disconnect the power supply plug.**
 Failure to disconnect the power supply plug may cause a fire in the event of a problem or malfunction.
- **Do not pull the power supply cord, the power supply plug, or the product side plug (inlet part) hard.**
 If the cord breaks, it may cause electric shock or fire.
- **Do not use the power supply cord that came with this product for other electrical equipment.**
 It may cause a fire or an electric shock.

- **Do not plug or unplug the power supply plug with wet hands.**
 This may cause an electric shock.

- **Remove dust from the power supply plug periodically.**
 Dust on the power supply plug may lead to an insulation failure due to moisture and thus cause a fire. Disconnect the power supply plug and wipe it with a dry cloth.
- **Make sure the power supply plug is pushed fully in.**
 Faulty insertion of the power supply plug may cause an electric shock or a fire due to a generation of heat. Never use a damaged power supply plug or loose power outlet.
- **Grip the power supply plug when disconnecting the power supply cord from the outlet.**
 Pulling the power supply cord may cause an electric shock or a short circuit.
- **Always use the detachable power supply cord that came with the product.**
 Other power supply cord may cause electric shock or fire.

- **Disconnect the power supply plug before moving the unit.**
 Take care not to damage the power supply cord. A damaged power supply cord may cause an electric shock or a fire**.**
- **Disconnect the power supply plug when the unit is not in use for long periods.**
 Keeping the unit connected may cause an electric shock, an electric leakage, or a fire due to the deterioration of insulation.
- **Before proceeding with maintenance or checking the unit, turn off the power switch (if the switch is provided), and disconnect the power supply plug.**
 Performing the work while power is still being supplied to the product or while the power supply plug is still connected may cause an electric shock and/or injuries.

When something is wrong with the unit

- **Never disassemble, repair, or modify the unit yourself.**
 A high-voltage area is located inside the unit. Any work carried out by unauthorized personnel may result in an electric shock. Contact our sales representative or agent for maintenance or repair.

- **Turn off the power switch and disconnect the power supply plug if something is wrong with the unit.**
 If the unit keeps running under such conditions, there may be a risk of an electric shock or a fire. Contact our sales representative or agent immediately for maintenance or repair.

DOMINIQUE DUTSCHER SAS

⚠WARNING

When using the unit

- **Never put containers with liquid on top of the unit.**
 This may cause an electric shock or short circuit if the liquid is spilled.
- **Never insert metal objects such as pins and wires into any vent, gap, or outlet on the unit.**
 This may cause an electric shock or injury by accidental contact with moving parts.
- **When opening the incubator door, do not directly inhale the air in the chamber.**
 It may be harmful to your health, causing gas poisoning or lack of oxygen.
- **Do not look directly at UV light.**
 UV light is harmful to the eyes.

- **Never splash water directly onto the unit.**
 This may cause an electric shock or short circuit.

- **Never store volatile or flammable substances in this unit except in a sealed container.**
 Such substances may cause an explosion or a fire if they leak.

- **Do not give strong shock or vibration when moving or using the unit.**
 The piping may be damaged, causing gas poisoning or lack of oxygen due to gas leakage.

When using gas cylinders

- **When connecting a gas cylinder to the incubator, check the gas type.**
 Using an incorrect gas type may result in explosion or fire or may cause gas poisoning or oxygen deprivation in the event of gas leakage.
- **Check that the gas cylinder is connected securely and that the gas will not leak.**
 Gas leakage may result in explosion or fire or may cause gas poisoning or oxygen deprivation.
- **When using CO_2 gas, ventilate the room.**
 A poorly ventilated confined space may cause gas poisoning or oxygen deprivation. When ventilating, exhaust air should be directed outdoors whenever possible. Avoid letting exhaust air stagnate indoors.
- **Be sure to use the specified pressure.**
 Using an incorrect pressure may result in explosion or fire or may cause gas poisoning or oxygen deprivation due to gas leakage.

When performing H_2O_2 decontamination

- **Use the reagent specified by our company for H_2O_2 decontamination.**
 Using other reagents may result in explosion or damage to the incubator.
- **When performing H_2O_2 decontamination, close the inner and outer doors securely.**
 Failure to do so may be harmful to health due to leakage of H_2O_2 gas.
- **During H_2O_2 decontamination, cover the access port (from the inside and back side of the chamber) with the silicon caps that came with the product.**
 Failure to do so may be harmful to health due to leakage of H_2O_2 gas.
- **Do not unlock and open the outer door using the accessory key during H_2O_2 decontamination or during H_2O_2 gas decomposition.**
 Doing so may cause harm to health from H_2O_2 gas leakage.
- **Wear safety goggles and rubber gloves when handling the H_2O_2 reagent.**
 Direct contact with the H_2O_2 reagent may result in inflammation of the skin.

DOMINIQUE DUTSCHER SAS

SAFETY PRECAUTIONS

⚠WARNING

When storing and disposing of the unit

- **If the unit is to be stored unused in an unsupervised area for a long period, ensure that children do not have access to the unit and that doors cannot be closed completely.**
 There may be a risk of child entrapment.
- **Disconnect the power supply plug when the unit is not used for long periods.**
 Keeping the connection may cause electric shock, current leakage, or fire due to the deterioration of insulation.
- **Use a qualified contractor to carry out disassembly/disposal of the unit and do not leave CO_2 incubator in a location that can be accessed by third parties.**
 This may result in unexpected accidents (e.g. the incubator may be used for unintended purposes).

⚠CAUTION

- **Do not climb on top of the unit or put any objects on the unit.**
 Falling from the unit may cause injury; falling objects may cause damage to the unit.
- **Do not damage the glass or give it a shock.**
 The inner door of the CO_2 incubator is made of tempered glass, but it can be broken and cause injury if it is used incorrectly.
- **Do not stack three or more incubators.**
 Falling incubators may cause injuries.

- **Turn the leveling feet to level the unit.**
 If the unit is not level, the unit may inadvertently move out of position when its door is opened or closed, causing injury.
- **Never store corrosive substances such as acids or alkalis in this unit except in a sealed container.**
 These may be harmful to your health and may cause corrosion of internal components or electrical parts.
- **Wear rubber gloves and a mask when cleaning the chamber.**
 You may get injured by the corners of interior parts. Also, touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.
- **When stacking the product, the floor should be strong enough to bear the weight.**
 If the products fall, it may cause injury.


- **When moving the unit, be sure to take precautions to prevent it from falling over.**
 Moving the unit with too much force may cause it to fall over, possibly resulting in damage or injury.
- **Set up or remove the shelf with both hands.**
 If the shelf falls down, it may cause injury.

- **Be careful not to burn yourself.**
 During operation, the inside of the outer door becomes hot. Touching the hot parts may result in burns.

DOMINIQUE DUTSCHER SAS

LABELS ON INCUBATOR

Users are advised to read carefully the warnings and cautions contained on stickers at key locations on the interior and exterior of the incubator.


Possible Danger	**Warning/Caution Type** Location of Danger	Warning/Caution Label	Description of Danger
Burns	**Hot Surface** Door & Chamber		During operation, avoid touching the inside of the door and chamber, which can reach high temperatures and may cause burns.
Personal injury	**Hazardous UV Light** Interior		The UV light is hazardous. Do not look directly at the UV light.
Personal injury	**Gas Poisoning or Oxygen Deprivation** Environment		When using CO_2 gas, ensure there is adequate ventilation. Using CO_2 gas in a small room without adequate ventilation may cause gas poisoning or oxygen deprivation. In addition, when opening the incubator doors, do not directly inhale the air in the chamber.
Personal injury	**Gas Poisoning or Oxygen Deprivation** Interior		Excessive pressure may cause gas supply lines inside the incubator to come loose, which may result in gas poisoning or oxygen deprivation due to the escape of gas.
Personal injury	**Electric shock** Rear cover		Never remove this cover. Only service personnel should remove the cover to avoid the electric shock.

DOMINIQUE DUTSCHER SAS

INCUBATOR COMPONENTS

Main body

Front and left side

1. **USB port:**
 A USB flash drive can be inserted into either of these ports when exporting log data. If USB flash drives are inserted into both USB ports, the one on the control panel will take precedence.


 Notes:
 ▪USB flash drives with a capacity of 32 GB or less that employ the FAT16/FAT32 file system are supported. USB flash drives that require passwords cannot be used. We do not guarantee the correct operation of all USB flash drives even if these conditions are satisfied. If issues occur, consider using a different type or brand of USB flash drive.
 ▪Do not insert devices other than USB flash drives into the USB port.

2. **Outer door:**
 Ensure the outer door is fully closed. If the outer door is left ajar, the incubator may not perform as intended. Always open and close the door gently. Forceful handling may lead to spillage of the culture medium, incomplete closure, or damage to the gasket. The door's opening direction is reversible. Contact our sales representative or agent to switch the door hinge from left to right or vice versa.

3. **Analog output terminals:**
 Used to obtain chamber temperature and CO_2 density data as current output. For details, refer to pages 61 and 62.


DOMINIQUE DUTSCHER SAS

Back

4. Access port (back side):
Used when routing cables, such as the sensor cable of a measurement instrument, into the chamber. Place the silicon cap on the port when it is not in use.

Right side


5. Switch cover:
Prevents shutting down of the unit in case of accidental pressing of the main power switch or prevents the gas tube from accidentally being pulled out. The cover is fixed by screws at its right and left sides.

6. Sample air outlet:
Used for obtaining sample air from the chamber. When this outlet is not used, cover the outlet with the sample air outlet cap. If not covered, condensation may form, causing leakage of liquid or corrosion.

7. Power switch:
Main switch for the incubator (ON-"I," OFF-"O").

8. CO_2 gas connection port A:
The gas tube is connected here. Refer to pages 26-28 for gas cylinder connection.

9. Position for power cord cover plate (accessory):
Prevents the power supply cord from being disconnected. Fix the cover plate after the power supply cord is connected to the power inlet connector.


10. CO_2 gas connection port B:
When the optional gas auto-changer (MCO-21GCP) is installed, both ports A and B are available. Refer to pages 75-77 for optional gas auto changer.

DOMINIQUE DUTSCHER SAS


INCUBATOR COMPONENTS

Inside

11. Inner door latch

You can open or close the inner door using this inner door latch. When closing the inner door, make sure to engage the latch. Otherwise, the incubator may not achieve its intended performance.

12. Inner door

Close the inner door properly. If the inner door is partially open, the incubator may not achieve its intended performance. Although the inner door is made of tempered glass, do not subject the glass to excessive impact. The glass may break and cause injuries.
You can replace the inner door with the optional small door (MCO-173ID). For purchase and installation, contact our sales representative or agent.

13. Tray runners

Trays are inserted to fit into the integrated shelf runners.

14. Humidifying pan cover

Prevents the UV light from leaking. Always use it. Using without the cover may have an impact on the chamber temperature distribution and humidity recovery.


15. Door switch

Detects the opening/closing status of the outer door. When the door is open, “Door:Open” blinks in orange on the Home screen, and audible alarm sounds after set door open alarm sound delay time (1 to 30 minutes) (pages 43-44).

16. Electric lock*


When the electric lock is equipped (either as a standard feature or added as an optional feature), the outer door is automatically locked by the electric lock. During a power failure, if you need to unlock the outer door manually using the key, insert the key into the key slot located at the bottom center of the outer door.
*: Available for 173AICUVH, or 173AICUV/173AIC with the optional electric lock (MCO-170EL).

Note:
The accessory key is used only for emergencies.
Do not use the key to lock the outer door in ordinary use.


DOMINIQUE DUTSCHER SAS

Inside (viewed from front)

[With inner attachments]

[Without inner attachments]

17. Access port (inside):

Used when routing cables such as the sensor cable of a measurement instrument into the chamber. Place the silicon cap when the port is not used.

18. CO_2 density measurement port:

The CO_2 sensor inside measures CO_2 density. Make sure not to spray ethanol for disinfection into the port when cleaning the chamber.

19. Humidity control bar:

Reduces condensation in the chamber caused by the external environment and frequent door openings.

20. UV-LED*1

The UV-LED does not generate ozone. For details about UV-LED, refer to pages 46-47. For replacement, contact our sales representative or agent.

21. H_2O_2 decontamination connector terminal*2

When performing H_2O_2 decontamination, the connector from the H_2O_2 generator is connected here (page 80).







*1: In the case of MCO-173AICUVH/173AICUV, or MCO-173AIC with the optional UV system set (MCO-LUVS).
*2: In the case of MCO-173AICUVH, or MCO-173AICUV/173AIC with the optional H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL).

DOMINIQUE DUTSCHER SAS

INCUBATOR COMPONENTS

Accessories

Check that following accessories are included with this product. If any accessories are missing, contact our sales representative or agent.

Name	Qty	Appearance	Usage
Power supply cord	2		Detachable power supply cords (UK/EU Type) for the incubator.
Power cord cover plate	1		Prevents the power supply cord from being accidentally pulled out of the unit.
Tray	4		Trays are set in the chamber and used for putting plates etc. on them.
Gas tube	1		Connects a gas cylinder and the incubator.
Humidifying pan	1		Used to maintain chamber humidity by filling this pan with sterile distilled or pure water.
Conversion joint	1		Connects the tube from a gas cylinder and the tube from the incubator.
Key*	2		For manually unlocking the outer door locked by the electric lock during a power failure.

* In the case of MCO-173AICUVH.

DOMINIQUE DUTSCHER SAS

INSTALLATION

Installation site

To run the incubator properly, install the incubator in a location which meets all the conditions described below.

■ Appropriate air environment
Install the incubator in an appropriate air environment.

■ A well-ventilated (airy) location
To ensure proper ventilation around this product, leave at least a 10 cm gap on the left, right, back, and top (if stacked, 10 cm above the upper unit). Blocking the ventilation may degrade the product performance or cause the product to malfunction. Also, at the front of the unit, allow 77 cm of space for opening and closing the door.

■ A location not exposed to direct sunlight
Do not install the incubator in a location where it will be exposed to direct sunlight. If the incubator is operated in direct sunlight, the product performance will be adversely affected.

■ A location away from heat sources
Do not install the incubator near significant heat sources, such as heaters, boilers, ovens, or autoclaves. Heat will adversely affect the performance of the incubator.

■ Ambient temperature at least 5°C lower than set temperature
The control temperature of the incubator is at least 5°C higher than the ambient temperature. For example, if the chamber is controlled at 37°C, the ambient temperature must be 32°C or lower. Do not allow the ambient temperature to be too high.

■ Firm and level floor
Select a site with a strong and level floor. If the floor is uneven or the installation is not level, the incubator will be unstable, and this may cause accident or injury. Also, this may adversely affect cultures. To avoid vibration and noise, ensure that the incubator is stable. An unstable surface may result in vibration or noise.

■ Away from products that cause vibration
Do not install the incubator near products that cause vibration. Vibration may adversely affect cultures.

■ Low humidity
Select a site with a relative humidity of 80% RH or lower. Using the incubator in high humidity may result in current leakage or electric shock.

■ No inflammable or corrosive gas
Never install the incubator in a location where it will be exposed to inflammable or corrosive gas. Doing so may result in explosion or fire. In addition, insulation may deteriorate due to corrosion of protective casing, resulting in current leakage or electric shock.

■ No falling objects
Do not install the incubator in a location where there is the possibility of objects falling from above. Doing so may result in damage or an accident.

■ Quiet place
Do not install the incubator in a high noise level location. If the incubator status becomes abnormal, users may not notice the audible alarm, causing damage to the cultures.

Note:
Keep electric products which emit electromagnetic waves away from the unit. Interference from the electromagnetic waves may cause the unit to malfunction.

DOMINIQUE DUTSCHER SAS

INSTALLATION

To prevent contamination:

■ **Avoid locations with high temperatures or humidity** where the air may contain more microorganisms.

■ **Do not place the incubator near doors, air conditioners, fans, etc,** where fast-moving air currents and high-traffic may increase the risk of microorganisms entering the chamber.

■ **Keep the incubator out of direct airflows from air conditioners or fans.** Airflow from an air conditioner may cause condensation or lead to possible contamination.

■ **If possible, use a cleanroom.**
To achieve a better culture, it is recommended that a cleanroom be used if one is available.

Installing the incubator

1. Remove the packing tape.


Remove all the tapes securing the doors and the internal components. Take all accessories out of the chamber. Open the doors to ventilate the chamber for a while. If the outer panels are dirty, use a cloth to wipe them with a diluted neutral detergent (undiluted detergent can damage the plastic components. For the dilution, refer to the instructions on the detergent).
Wipe off the residual detergent with a piece of gauze moistened with water and then wipe off any moisture with a dry cloth.


Note: Remove the cable tie around the power supply cord to prevent corrosion of the cord coating.

2. Set the humidifying pan and trays.


First, set the humidifying pan cover and humidifying pan.

Set the humidifying pan cover by lifting its front side a little and engaging the rear side tabs into the pins on both sides of the duct.


Lift the front side of the humidifying pan cover slightly and insert the humidifying pan under the humidifying pan cover.


Then, place the four trays at the desired levels, ensuring that the bent-down part of each tray is at the front.

DOMINIQUE DUTSCHER SAS


Note:
When the optional small door installation kit (MCO-173ID) is installed, we recommend placing the three trays to align with the lower edge of the small doors to make the most efficient use of the interior space.

3. Adjust the leveling feet.

Check that the incubator is level by placing a spirit level on a tray. Adjust the leveling feet by turning them to level the incubator.

Note: Incubating in a sloping condition may adversely affect the cultures.

4. Earth the incubator.

Earth the incubator during installation to prevent electric shock. If there is no earth wire at the location, consult with qualified service personnel.


When an earth must be installed:

If an earthed 3-pole outlet is not available, then an earth must be installed. Consult with qualified service personnel.

Connecting to the power inlet connector


Connect the included power cord to the unit by following the steps below.

1. Remove the switch cover.

2. Insert the unit-side plug of the included power cord into the power inlet connector on the lower right side of the unit.

3. Screw the power cord cover plate to prevent the cord from being accidentally disconnected.


DOMINIQUE DUTSCHER SAS

INSTALLATION

When stacking an incubator

⚠ WARNING	● **When stacking incubators, use the hooks on the back panel of the incubator and secure it to the wall by passing a strong rope or chain through the hooks.** If not, the incubator may fall over and cause injuries.

The optional double stacking bracket MCO-170PS or the stacking plate MCO-170SB or MCO-230SB is required to stack an incubator on this product.
Contact our sales representative or agent to purchase these options and fix them.

Notes:

- When stacking an incubator, fix the upper incubator to the wall by passing a strong rope or chain through the hooks on the back side of the unit to prevent falling over.
- When stacking our other CO_2 incubator or O_2/CO_2 incubator with this product, the brackets/plates to be used vary depending on the combination. Refer to Table 1 on page 91.

DOMINIQUE DUTSCHER SAS

CLEANING

Before using the incubator for the first time, or setting the incubator again after moving the incubator, always thoroughly clean the chamber, inner attachments, and humidifying pan (accessory) to remove dirt (tape residue, oil, etc.). Cleaning these items is essential to ensure the best performance of the incubator. Follow the steps described in the following sections to properly clean the incubator.

CAUTION	● **Wear rubber gloves and a mask when cleaning the chamber.** You may get injured by the corners of interior parts. Also, touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.

Removing the components


Before cleaning the chamber, remove the inner attachments in the chamber. Be careful not to damage the UV-LED and the humidity control bar in the duct.


Removing the inner attachments

1. Turn off the power to the incubator.


2. Open the outer door and inner door, then pull out the parts in the following order.

Pull out all the trays. → Lift the humidifying pan cover and pull out the humidifying pan. → Lift the humidifying pan cover a little to disengage the tabs of the humidifying pan cover from the pins on the duct at the back of the chamber and pull it out.

3. Remove the connector cap.*


*: In the case of MCO-173AICUVH, or MCO-173AICUV/173AIC with the optional H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL).


DOMINIQUE DUTSCHER SAS

CLEANING


4. Lift the fan cover to detach its tabs from the duct and pull out the fan cover.

5. Lift the duct and disengage it from the pins on the rear side and pull it out.


6. Pull out the fan.


7. Remove the silicon caps on the interior side and exterior side.

[Exterior side]
Loosen the screws at the access port cover and remove the cover by lifting it. Then, remove the silicon cap.

[Interior side]
Remove the silicon cap at the rear top left corner of the chamber.


DOMINIQUE DUTSCHER SAS

Cleaning the chamber and components

Turn off the incubator before cleaning.


1. Put on rubber gloves, and then disinfect the surface of the rubber gloves with ethanol for disinfection.

2. Thoroughly wipe clean the inner walls of the chamber, the removed inner attachments, the temperature sensor, the humidity control bar, and the UV-LED using gauze moistened with a proper amount* of ethanol for disinfection. Make sure to take particular care in cleaning corners and joints of the product.

3. Thoroughly wipe clean the inner door, inner door gasket (contact points between the unit and the gasket, which is back side of the gasket), handle, and inner door glass using gauze moistened with a proper amount* of ethanol for disinfection. Especially, take special care of the back side of the gasket since an unclean gasket may lead to contamination.


Note:
If the incubator starts operation with the inner door gasket not in place, the humidified air in the chamber will leak and will cause condensation between the unit and outer door. After wiping, confirm that the inner door gasket is securely in place and does not have any creases.

* The amount that cannot form droplets on the surface.

Notes:
- The most effective procedure to prevent contamination is to clean each part with a cloth or sponge moistened with a neutral detergent diluted to 5% or less, and then wipe down each part with gauze moistened with distilled water or pure water. You can also use autoclave sterilization (121°C for 20 minutes) for the silicon caps for access ports (2 pcs) and the fan.
- If the result of cultivation is undesirable due to dirt (tape residue, oil, etc.) on the removed inner attachments, clean the removed inner attachments following the steps below.
 1) Wash with a neutral detergent.
 2) Rinse well with distilled water or pure water.
 3) Wipe with gauze, sterile non-woven cloth, or sterile paper.
- After cleaning, ensure the interior is completely dry and free of any ethanol odor before turning the power back on at a later time. Turning the power on while the interior is still wet may result in damage to the CO_2 sensor.

Do not use acidic, alkaline or chlorine solvents or disinfectants. Use of these solutions may cause discoloration, corrosion or rust.
Be careful not to spray solutions such as ethanol for disinfection to the ①fan motor shaft bearing, ④inner sample gas access port, ⑤H_2O_2 decontamination connector terminal*1, ⑥CO_2 density measurement port, ⑦temperature sensor, and ⑧CO_2 gas injection port. Also, do not wash the ②humidity control bar and ③ UV-LED*2 using detergent. Use of these solutions may cause damage.

*1: In the case of MCO-173AICUVH, or MCO-173AICUV/173AIC with the optional H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL).
*2: In the case of MCO-173AICUVH/173AICUV, or MCO-173AIC with the optional UV system set (MCO-LUVS).

DOMINIQUE DUTSCHER SAS

CLEANING

Cleaning the control panel

When sterilizing and cleaning the control panel, follow the precautions below.

- Do not spray liquid on the control panel directly. The panel may become damaged.
- When sterilizing and cleaning, wipe the surface using a piece of gauze moistened with a proper amount of ethanol for disinfection (the amount that cannot form droplets).


Installing the cleaned components

Install all the inner attachments by reversing the steps used to remove them.
If the inner attachments are not inserted properly, the intended air velocity cannot be achieved, and it may cause culture failure or insufficient H_2O_2 decontamination.


Installing the inner attachments

1. Put the silicone caps on the interior side and exterior side.

[Rear side]


[Chamber side]



2. Install the fan.

When installing the fan, position the center hole of the fan over the projection of the motor shaft and push it fully in. Manually turn the fan lightly to check that it does not touch the rear wall. Then, pull the fan lightly to check that it does not come off.


3. Install the duct by fitting the 4 pins at the rear side of the chamber (shown by ➡ in the figure below) into the 4 holes of the duct.


DOMINIQUE DUTSCHER SAS

4. Attach the fan cover. Be sure to insert the two hooks on the upper part of the fan cover into the holes on the top of the duct.

5. Attach the connector cap.*


*: In the case of MCO-173AICUVH, or MCO-173AICUV/173AIC with the optional H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL).


6. Set the humidifying pan cover and humidifying pan.

Set the humidifying pan cover by lifting its front side a little and engaging the rear side tabs into the pins on both sides of the duct.

Set the humidifying pan under the humidifying pan cover in a longitudinal direction with its shorter side touching the chamber's rear wall. If the humidifying pan is not set correctly, re-condensed water droplets will fall outside of the humidifying pan. This can cause condensate to accumulate on the bottom surface.

7. Insert the tray so that the bent-down part of the tray is at the front.

DOMINIQUE DUTSCHER SAS

PREPARATION BEFORE USING

Before running the incubator, you need to connect a CO_2 gas cylinder to the incubator, install a remote alarm device (if necessary), and fill sterile distilled water or pure water into the humidifying pan.


Connecting a CO_2 gas cylinder

Note: Before connecting a CO_2 gas cylinder, read “SAFETY PRECAUTIONS” on pages 5 to 10.

1. Prepare a CO_2 gas cylinder and attach the optional gas regulator MCO-010R to it.


Notes:

- Use a liquefied-CO_2 gas cylinder (at least 99.5 % pure). The siphon (dip tube) type cannot be used.
- When MCO-010R is not available, attach a two-stage gas pressure regulator rated at 25 MPa(G) (250 kgf/cm^2(G), 3,600 psi(G)) for the primary side, and 0.25 MPa(G) (2.5 kgf/cm^2(G), 36 psi(G)) for the secondary side.

2. Connect one end of the gas tube that came with this product to the CO_2 gas connection port A. When the CO_2 gas auto-changer MCO-21GCP (option) is installed, refer to page 75 for the connection.

Lower right side of the unit

About tube fitting

This product employs tube fittings for the connection ports. When connecting and removing a tube, follow the instruction below.

1) Connection

Insert the tube into the gas connection port all the way to the end so that the tube is secured by the inner locking stopper and sealed with elastic sleeve around it.

* In order to make sure that the tube is securely connected, pull the tube after connection. If the tube comes off, pull the release ring, and then reconnect the tube and check if it does not become detached.

2) Removal

While pushing the release ring, pull the tube to remove it.


DOMINIQUE DUTSCHER SAS

3. Connect the other end of the gas tube to the gas regulator (MCO-010R). There are following two ways to connect to the gas regulator.


Notes:
- Use the gas tube that came with this product for the connection between the unit and the conversion joint.
- Check that the connected tubes are not bent.
- Do not connect multiple incubators to one CO_2 gas cylinder by branching the gas tube. Doing so may not obtain the intended performance.

When connecting to the hose joint of the gas regulator

General view


1) Insert the other end of the gas tube into the release ring of the conversion joint that came with this product.

2) Connect the hose joint of the gas pressure regulator and the other side of the conversion joint using the connection tube (A) between them and tighten the tube ends using the tube bands (B).

Parts to be prepared by the customer (Reference)	A. Connection tube: φ6 inner diameter/Tygon tube B. Tube bands


DOMINIQUE DUTSCHER SAS


PREPARATION BEFORE USING

When connecting using a tube fitting

General view


1) Remove the hose joint from the gas regulator.

2) Attach the prepared conversion bush (B) to the gas regulator using the packing or seal tape (A) and then, attach the tube fitting (C) to the conversion bush.

Parts to be prepared by the customer (Reference)	A. Packing: 6A (ID13 mm x OD 18 mm) or seal tape B. Conversion bush: bush A M12 x 1, R1/8 C. Tube fitting: female straight, adaptive tube OD 6 mm, Rc1/8

3) Insert the other end of the gas tube described in step 2 to the release ring of the tube fitting (C).

Note: If you use seal tape, wind it around the thread part only. Make sure that the seal tape does not stick out of the thread.

4. After connecting the gas tube, make sure that no gas is leaking (e.g. by using a gas leak detector spray).

5. Set the secondary pressure for the CO_2 gas to 0.03 MPa(G) to 0.1 MPa(G) (0.3 kgf/cm^2(G) to 1 kgf/cm^2(G), 4.4 psi(G) to 14.5 psi(G)) for gas injection. Recommended pressure: 0.03 MPa(G) (0.3 kgf/cm^2(G), 4.4 psi(G)).

Note: Do not set the secondary pressure too high. As the pressure increases, the CO_2 gas density will fluctuate widely. Also, excessive pressure may cause gas supply lines inside the incubator to come loose, which may result in gas poisoning or oxygen deprivation due to leaking gas. If gas lines come loose, the incubator must be repaired.

6. When there is no CO_2 gas left and the CO_2 gas empty alarm is activated, replace the empty gas cylinder with a new one.
When the optional gas auto-changer MCO-21GCP is installed, it automatically switches the empty CO_2 gas supply line to the other cylinder. Refer to page 76.

Notes:
- The gas tubes connected to the incubator will degrade over time. If any deterioration or abnormalities are found during inspection, replace the tubes immediately.
- Close the valve of the CO_2 gas cylinder when the CO_2 gas is not in use.


DOMINIQUE DUTSCHER SAS

Installing a remote alarm device

CAUTION	**Before connecting an alarm device to the remote alarm terminals, turn off and unplug both the unit and the device.** Connecting while the power switch is on may cause electric shock.

The alarm condition of the incubator and the door open/close status can be transmitted to a remote location by connecting a commercially available external alarm device to the remote alarm terminals. The use of a remote alarm device is especially recommended if the incubator is placed in an unattended environment to ensure that the alarm condition is promptly notified to the administrator. For installation of the alarm device, contact our sales representative or agent.

The remote alarm terminals (shown in the right figure) are located on the lower left side of the incubator. Alarm (N.O., COM., and N.C. terminals) and door status (DOOR STATUS terminals) signals are output from these terminals as non-voltage contact outputs. Contact capacity is DC 30 V, 2 A.
The following tables show the contact output conditions of the terminals.

Contact output from alarm terminals

"Remote Alarm" setting (pages 43-44)	Connection terminal	Normal condition	Abnormal condition (Including in the event of power outage and when power plug is pulled out.)	
				When Buzzer Stop button is tapped
ON: Not linked with Buzzer Stop button	COM.-N.C.	Close	Open	Open (indicates abnormal condition)*
	COM.-N.O.	Open	Close	Close (indicates abnormal condition)*
OFF: Linked with Buzzer Stop button	COM.-N.C.	Close	Open	Close (goes back to normal condition)
	COM.-N.O.	Open	Close	Open (goes back to normal condition)

*When the alarm is "A1:CO2 Gas Empty", the incubator goes back to normal condition.

Contact output from the DOOR STATUS terminals

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Notes:
- For the types of alarms that can be transmitted to the remote alarm device, refer to pages 85 to 86.
- The signal from DOOR STATUS terminals are not linked with the "N4:Outer Door Open notification", which is indicated by the incubator itself using the incubator software. Therefore, you cannot use the Buzzer Stop button on the control panel to stop the door opening alarm indicated by the remote alarm device.
- Use a commercially available alarm device when using the remote alarm function. Creating a new device or modifying an existing one yourself may result in electric shock or malfunction.
- When connecting an alarm device to the remote alarm terminals, turn off the power switch and unplug the incubator. Failing to do so may result in electric shock.
- Use shielded wire for connection. A maximum length of 30 meters of cable is recommended to connect to external devices. Cables longer than 30 meters are not covered under warranty.

Use a twisted shield wire for the connection.
Type: UL 2343, UL 2448, UL 2464, UL 2552, UL 2623.
Length: 30 m max.


FOMINOUE DUTSCHER SAS

PREPARATION BEFORE USING

Filling the humidifying pan

To fill the humidifying pan with sterile distilled water or pure water, or to change the water in the humidifying pan, follow the steps below.


1. If possible, lift the humidifying pan cover and pull out the humidifying pan

2. Dispose of the remaining water in the pan. Then, rinse it thoroughly with distilled water or pure water and wipe it using ethanol for disinfection.

3. Wipe the moisture on the bottom of the chamber.

4. Put back the humidifying pan under the humidifying pan cover and pour sterile distilled water or pure water (approx. 1.5 L, preheated to the incubator's temperature set point).

5. Set the humidifying pan in a longitudinal direction with its shorter side touching the chamber's rear wall. Then, close the inner door and the outer door.

Notes:

- If you run the incubator without sterile distilled water or pure water in the humidifying pan, the chamber temperature may temporarily increase, resulting in a higher temperature than the set temperature.
- Use sterile distilled water or pure water and change the water once a week.
- Before pouring sterile distilled water or pure water, preheat it to the chamber temperature. Adding cool water will lower the temperature and humidity in the chamber.
- Refill the water when the water in the pan is running low.
- Do not mix any reagent in the water for humidifying pan. Doing so may affect the cultures.
- After cleaning the chamber using ethanol for disinfection, change the water in the humidifying pan. The ethanol for disinfection that got into the pan may influence the correct control of CO_2 density.
- Do not use a humidification water bag (commercially available). Using such a water bag may not obtain the product's intended performance.
- Ensure that the humidifying pan is set firmly against the back wall inside the chamber. The humidity control bar inside the duct is always kept at a low temperature to recondense the water vapor inside the chamber. If the humidifying pan is not pushed all the way to the rear wall of the chamber, the recondensed water droplets will fall outside of the humidifying pan, causing condensation to accumulate. This could lead to contamination or mold growth. The accumulated water may evaporate and leave white marks; however, this is not considered abnormal. To clean, wipe thoroughly with gauze or non-woven fabric moistened with ethanol for disinfection, or if the stains are severe, polish them off with a cream cleanser.
- When the incubator is not used, dispose of the water in the humidifying pan, eliminate moisture in the chamber, and close the doors once the interior is dry. The presence of humidity may lead to failure of the incubator.

DOMINIQUE DUTSCHER SAS

START-UP PROCEDURE

Follow the procedure below to start test run or start operation of this product. Before starting, check that the incubator is turned off and the power plug is disconnected.

Cleaning



1. Clean the chamber by following the steps (removing components, cleaning the chamber, and installing the components) on pages 21-25.

Preparation & Start-up

2. Prepare the incubator for cultivation (pages 26-29).

3. Insert the power plug into the dedicated outlet.

Note:
Two removable power supply cords are provided. Select the appropriate code based on your region of residence.


For UK	For EU countries
	

4. Turn on the power switch for the unit.

5. Fix the switch cover.

Setting after start-up

6. Set the date, time, and time zone.
When you turn on the incubator for the first time, the following Notice message is displayed on the Home screen. Tap the OK button to close the window and then, set current date, time, and time zone by referring to pages 58 and 59.

Note: If you do not tap the Apply button after selecting each value, the setting is not saved, and this Notice window will be displayed every time when the Home screen is displayed.

Configuration

7. Configure the settings for cultivation (pages 42).
Set the chamber temperature, CO_2 density in the chamber, and high temperature limit for overheat alarm.

Note:
When running the incubator for the first time, or after a long unused period, set the temperature to desired temperature, and CO_2 density to 0%. After filling the humidifying pan (page 30), run the incubator with the setting for about 4 hours or more until the temperature and CO_2 sensor become stable and the message “N1: Gas Sensor initializing.” (page 86) in the message display area disappears. Then, set the desired CO_2 density, and high temperature limit for overheat alarm.

8. Configure other settings such as alarm, UV-LED, screen lock, and key lock etc., if necessary.

Cultivation

9. Fill the humidifying pan (page 30).

10. Put cultures into the chamber.

DOMINIQUE DUTSCHER SAS

FOR OPTIMAL CULTIVATION

Precautions for cultures

- **Leave space between culture containers.**
 Always leave space for ventilation between culture containers (petri dishes, flasks, etc.). Inadequate spacing may result in uneven temperature distribution and CO_2 gas density.
- **Do not place harmful materials in the chamber.**
 Never place samples that release acidic, alkaline, or corrosive gas into the chamber. Doing so may cause damage resulting from discoloration or corrosion of the interior, internal components, temperature sensor, and CO_2 sensor.
- **Close the inner door.**
 Make sure to engage the inner door latch properly to close the inner door. Failure to close the inner door will adversely affect the performance of the incubator even if the outer door is closed.
- **Open and close the doors gently.**
 Closing the doors forcefully may cause spillage of the culture medium, incomplete door closing, or damage to the gasket.
- **Be careful when closing the outer door.**
 Use the handle when closing the outer door. Not doing so may cause injury by getting fingers caught in the door. Do not hang on to the outer door. Doing so may result in injury from the outer door coming loose or the incubator falling over or may cause current leakage or electric shock.
- **Be careful of the inside of the outer door.**
 There is a heater in the outer door. The inner side of the outer door may become hot.
- **Avoid using excessive force on the inner door.**
 Do not hit the glass, poke it with sharp objects, or apply excessive force. Doing so may result in injury from breaking the glass.
- **Check the cause of any alarm.**
 If an audible alarm sounds while the incubator is in use, check the cause of the alarm immediately. For details on what may cause an audible alarm to sound, refer to pages 85-86. If the cause of an audible alarm is uncertain, contact our sales representative or agent.
- **Vibration of a shaker.**
 If cultures requiring a shaker and those not requiring a shaker are carried out together, the vibrations from the shaker may transfer and affect the growth of the cells.

Preventing contamination

- **Use clean containers.**
 Dirty containers are the greatest cause of contamination for cultures. Be careful not to get containers or trays dirty when taking them in and out.
- **Keep the chamber clean.**
 Wipe off any fingerprints. If water spills from the humidifying pan, or if the doors are left open for a long time, condensation may form on the inner side of the inner door. Wipe off condensation with a piece of dry sterile gauze. In particular, if the culture medium is spilled, clean and disinfect the chamber immediately (pages 23-24). Also, clean the chamber and inner attachments at least once a month.
- **Use sterile distilled water or pure water (do not use ultrapure water) for the humidifying pan.**
 Use sterile distilled water or pure water that has been purified through methods such as deionization or desalination for the humidifying pan and replace the water regularly (once a week).
 Note that if ultrapure water is used or if the humidifying water is not replaced for longer than the above-mentioned periods, red rust-like particles may appear in the humidifying pan. If these red rust-like particles appear, promptly discard the water, thoroughly clean and rinse the pan, and then refill it with fresh sterile distilled water or pure water.
- **Take care when putting items that make the temperature high into the chamber.**
 Putting items such as shaker, which make the temperature high, into the chamber may cause condensation and lead to possible contamination.




DOMINIQUE DUTSCHER SAS

2. TOUCH SCREEN OPERATION

USING THE TOUCH SCREEN

For controlling the incubator, an LCD touch screen is provided at the front of the incubator.

Basic operation


Tap Lightly touch the screen with your finger and immediately release the finger from the screen.	**Slide** With lightly touching the screen with your finger, move to the destination and then release the finger from the screen.	**Long tap** Lightly touch the screen for a few seconds and then release the finger from the screen.
		

Notes:
- Use your finger or a capacitive stylus pen to protect the screen.
- Do not use sharp pointed hard items such as the tip of a ballpoint pen for touching the screen.

Entering alphanumeric values and symbols

Entering numerical values:


1. Tap the numeric input field.
►The numeric input window is displayed.

2. Tap the numeric keys or Up/Down keys to enter a numerical value.


Details of keys in the numeric input window
- Numeric keys (0 to 9): Input numbers directly.
- Up/Down keys (▲/▼): Increase or decrease the numerical value displayed in the numeric input field.
- Clear button: Deletes the numerical value displayed on the numeric input box.

Note:
∞, +, and - keys are unavailable.

3. Tap the OK button.

Notes:
- In some situations, the Up/Down keys are not displayed.
- When the numeric input window is displayed, the Top and Back buttons are unavailable.


DOMINIQUE DUTSCHER SAS

USING THE TOUCH SCREEN


Entering alphanumeric characters and symbols:

1. Tap the alphanumeric input field.
►The alphanumeric input window is displayed.

2. Tap the alphabetic keys, numeric keys, or symbol keys to enter alphanumeric characters and symbols.
You can switch the uppercase letters/numbers window to lower case letters/symbols window by tapping the A⇌a button. The screen switches in the order of the uppercase letters/numbers window → lowercase letters/symbols window.

[Uppercase letters, numbers]

[Lowercase letters, symbols]


Details of keys in the alphanumeric input window
• Alphabetic keys (A to Z, a to z, Space): Input alphabetic characters or spaces.
• Numeric keys (0 to 9): Input numbers.
• Symbol keys: Input symbols.
• Left/Right buttons (◄/►): Move the cursor to the left/right.
• Delete button: Deletes an alphanumeric character or symbol immediately to the left of the cursor.

3. Tap the OK button.

Note:
While the alphanumeric input window is open, the Home and Back buttons are unavailable.

Entering date and time:


1. Tap the date input field or time input field.
►The input window is displayed.


2. Tap the ▲/▼ buttons to enter the values.

Note:
The time is expressed using a 24-hour clock.


[Date] Year/Month/Day

[Date] Day/Month/Year

[Time] Hour/Minute/Second

3. Tap the OK button.

Note:
While the input window is open, the Home and Back buttons are unavailable.

DOMINIQUE DUTSCHER S.A.S.

Other operations

1. Common button operation on the touch screen:

Button	Name/Details
OK	OK button Confirm the selected setting and return to the previous screen (the setting has not been saved yet).
Cancel	Cancel button Cancel the selected setting and return to the previous screen.
Apply	Apply button Apply and save the setting into the unit and return to the previous screen.
Abort	Abort button Abort ongoing operation.
Home *	Home button Return to the Home screen.
Back *	Back button Return to the previous screen.
▼ ▲ ◀ ▶	Arrow buttons for selecting numbers, a time zone, or moving the cursor in dialog windows Tapping the button increases or decreases the value one step at a time or moves one position at a time. Long tapping the button continuously increases or decreases the value or makes the cursor move or scroll continuously.

* If you tap these buttons on the screens for selecting/setting values before tapping the Apply button, the values are not saved into the unit.

2. Auto return function

When there is no touch screen operation for about 90 seconds on screens other than the Home screen, the display automatically returns to the Home screen. Therefore, when you are in the middle of setting some values and do not touch the screen for more than about 90 seconds, the screen returns to the Home screen without saving the data.

When the sleep function is set to ON (page 60) and there is no touch screen operation for about 90 seconds without any alarms or errors, the display automatically returns to the Home screen and enters the sleep state.


DOMINIQUE DUTSCHER SAS

HOME SCREEN

The following screen is displayed on the LCD display after turning on the power switch for the incubator. This screen is known as “Home screen” and shows basic information about the status of the incubator. It will take about 10 seconds before the Home screen is displayed after turning on the power switch.

Note:

During warm-up, the message “N1: Gas Sensor initializing.” appears on the message display area. Also, “--.-“ is displayed on the current CO_2 density display area.

Temperature section


1. **Current temperature in the chamber**
2. **Temperature set point**

CO_2 section

3. **Current CO_2 density:**
 When CO_2 density is set to “0.0,” nothing is displayed here.
4. **CO_2 density set point**
5. **CO_2 gas supply line indicator A and B*:**
 When the optional CO_2 gas auto changer MCO-21GCP is installed, the indicator of the CO_2 gas supply line (A or B) that is currently supplying CO_2 gas is turned on (pages 76 and 77).

* The indicator becomes active only when the optional gas auto changer (MCO-21GCP) is installed.

DOMINIQUE DUTSCHER SAS

6. Overheat alarm indicator:

When overheat status occurs, “Over Heat” blinks in dark orange.

7. UV-LED status indicator*:

UV-LED ON: “UV : On” is displayed in orange.
UV-LED OFF: “UV : Off” is displayed.

* In the case of MCO-173AICUVH/173AICUV or MCO-173AIC with the optional UV system set (MCO-LUVS).

8. Outer door (open/close) status indicator:

Open: “Door : Open” blinks in orange.
Closed: “Door : Closed” is displayed.

9. 2D barcode

When the incubator is in an alarm or notification state, a 2D barcode is displayed here. By scanning the code with a device such as a smartphone or tablet, you can view the countermeasure for the alarm or notification displayed in the message display area.

10. Message selection button: When there are multiple messages relating to Alarm or Notification status, tapping these buttons can switch to the previous or next message.

11 Status indicator

The status of the incubator is indicated by the following indicator ring color.

Color	Status	Details
Green	Normal	No abnormality.
Orange	Notification	An error not related to the incubation and decontamination performance has occurred. Refer to page 86 for measures to address the error. If the status remains unsolved, stop the incubator operation and contact our sales representative or agent.
Red	Alarm	An error that affects the incubation and decontamination performance has occurred. Refer to pages 85 to 86 for measures to address the error. If the status remains unsolved, or the countermeasure is not described, stop the incubator operation and contact our sales representative or agent.

DOMINIQUE DUTSCHER SAS

HOME SCREEN


12. Buzzer stop button

When the audible alarm sounds (intermittent tone), “Buzzer Stop” is displayed in the center of the Buzzer Stop button. Tap this button to stop the sound. However, when the alarm sound is a continuous tone (Alarm: A38: Door Lock Failure or A97: Over Heat), “Buzzer Stop” indication is not displayed, and you cannot stop the sound.
When the ring back function is set to ON and the alarm status still continues after the ring back set time elapsed, the audible alarm sounds again (page 44).

13. Message display area

Alarm or Notification information is displayed.


DOMINIQUE DUTSCHER SAS

MENU BAR BUTTONS

Using the buttons in the menu bar, you can switch the screens for unlocking the door, viewing data of the incubator, exporting log data, cleaning, or settings.


1. Home button

Tapping this button on a screen other than the Home screen returns to the Home screen.


2. Setting button

Tapping this button displays the Setting screen where various settings such as temperature set point, alarm, screen lock can be configured.


The Setting screens consist of two pages. You can navigate between pages by:

- tapping the ◁ (previous page) or ▷ (next page) icon on the left and right sides of the screen or
- tapping one of the page indicator ● ○ for the desired page (current page is indicated in orange).

First page

Second page

Buttons in the Settings screen	Details	Page
Stand-by	Temperature set point, CO_2 density, high temperature limt setting	42
Auto-Lock*1	Turning ON or OFF the auto-lock setting	54-57
Auto-Lock User*1	Adding or delating users for automatic door lock function	51-53
Screen Lock	Turning ON and OFF the screen lock setting, password setting	48-50
Alarm	Setting for alarm (temperature and CO_2 gas), door delay, ring back, and remote alarm	43, 44
UV*2	UV setting	46, 47
Date & Time	Date, time, and time zone setting	58, 59
Brightness/Sleep	Screen brightness and sleep mode setting	60

*1: In the case of MCO-173AICUVH, or MCO-173AICUV/173AIC with the optional electric lock (MCO-170EL).

*2: In the case of MCO-173AICUVH/173AICUV, or MCO-173AIC with with the optional UV system set (MCO-LUVS).

"In other cases, the buttons will not be displayed.

DOMINIQUE DUTSCHER SAS

MENU BAR BUTTONS


Buttons in the Settings screen	Details	Page
DAQ	This button is for service staff use only. Do not tap it.	-
System Info.	Information of the model No.and serial No., unit name setting, and software version information.	63
Analog Output	Turning ON or OFF the analog output setting	61,62
Operating Instructions	Displays 2D barcode for accessing online operation instructions for the product.	64

3. Cleaning button*

Tapping this button displays the Cleaning screen where you can start H_2O_2 decontamination.

* The cleaning function can be used in either of the following cases. If neither condition is met, the Cleaning button is grayed out and becomes unavailable.


- MCO-173AICUVH
- MCO-173AICUV/173AIC equipped with the optional accessories: H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL)

Buttons in the Cleaning screen	Details	Page
H_2O_2 Decontamination	Starts H_2O_2 decontamination.	79-83

4. Event log button


Tapping this button displays the Event Log screen where you can check event log data saved in the unit and can export the data to a USB flash drive.

DOMINIQUE DUTSCHER SAS


5. Chart button

Tapping this button displays the Chart screen where the chart of the operation log data can be shown. Also, on the Chart screen, selected log data displayed on the screen can be exported to a USB flash drive (pages 69-72).

6.Export button


Tapping this button displays the Export screen where you can export log data (pages 73-74).

7. Door unlock button

Tapping this button when the automatic door lock function* is ON displays a screen for entering a user ID and password for opening the outer door (pages 51-55).

* The automatic door lock function can be used in either of the following cases. If neither condition is met, the Door unlock button is grayed out and becomes unavailable.
 - MCO-173AICUVH
 - MCO-173AICUV/173AIC with the optional electric lock (MCO-170EL).

The button icon changes as follows depending on the status of the incubator.

Button icon	Description	Button control
	The automatic door lock function is set to ON and the outer door is locked. To unlock the outer door, tap this button.	Yes
	The automatic door lock function is set to OFF, or the automatic door lock function is set to ON and the outer door is unlocked.	No


DOMINIQUE DUTSCHER SAS


SETTING FOR BASIC OPERATION

Setting temperature, CO_2 density, and high temperature limit

The chamber temperature, CO_2 density, and high temperature limit for overheat alarm can be set by following the steps below.


1. Tap the Setting button in the menu bar.
►The Setting screen is displayed.

2. Tap the Stand-by button on the Settting screen.
►Stand-by screen is displayed.

3. Tap the numeric input field to enter each value or slide the ON/OFF buttons to configure the settings.

Note:
When running the incubator for the first time, or after a long unused period, fill the humidifying pan (page 30), then set the temperature to desired temperature and CO_2 density to 0%. Run the incubator with the setting for about 4 hours or more until the temperature becomes stable and the message “N1: Gas Sensor initializing.” (page 86) in the message display area disappears. Then, set the desired CO_2 density, and high temperature limit.

Details of entry

Item	Details	Settable range	Default value
Temperature	Temperature set point of the chamber. **Note:** Set a value of at least 5°C or higher than the ambient temperature.	0.0°C to 50.0°C	37.0°C
High Limit	When the chamber temperature exceeds the high temperature limit set here, “Over Heat” alarm activates to notify the user of the incubator’s over temperature status. This function works independently of the high/low temperature alarm. **Notes:** • This setting should be done after the chamber temperature reached the set temperature set point. • Set a value of at least 1°C higher than the temperature set point.	20.0°C to 53.0°C	53.0°C
CO_2 density	To control the CO_2 density in the chamber, turn ON the slide switch by sliding it to the right. Then tap the numeric input field to enter desired CO_2 density value (initially 5.0 is displayed). When you turn OFF the slide button, the value is deleted, and the CO_2 density control is turned off. Also, when you enter “0” in the numeric input field, the button is automatically turned to “OFF.”	ON/OFF 0.0% to 20.0%	OFF

4. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

5. Tap the Home button in the Menu bar.
►The screen returns to the Home screen.


DOWNLOADED FROM DUTSCHER.COM


Configuring alarm settings

On the Alarm #1 and Alarm #2 screens, various alarm settings can be configured.


1. Tap the Setting button in the menu bar.
►The Setting screen is displayed.

2. Tap the Alarm button.
►The Alarm #1 screen is displayed.

3. Tap the numeric input field to enter each value.

Details of entry for “Alarm” fields


Item	Details	Settable range	Default value
Temperature	When the chamber temperature gets higher or lower than the set temperature set point plus or minus the value set here, the temperature alarm is activated.	±1.0 to ±5.0°C	±1.0°C
CO_2	When the CO_2 density gets higher or lower than the set CO_2 density plus or minus the value set here, the CO_2 density alarm is activated.	±0.5 to ±5.0%	±1.0%

Details of entry for “Alarm Delay” fields

Item	Details	Settable range	Default value
Temperature	Alarm delay is a function that when the alarm status occurs, instead of activating audible alarm or remote alarm immediately, activates the alarm after the elapse of the set alarm delay time. **Note:** When the incubator recovers from the alarm status within the set alarm delay time, audible alarm or remote alarm does not activate after the elapse of the alarm delay time.	0 to 15 min	15 min
CO_2			

DOMINIQUE DUTSCHER SAS

4. Tap the right page indicator or ▷ button to move to the next page.
►The Alarm Setting #2 screen is displayed.

5. Tap the numeric input field to enter a value or slide the ON/OFF buttons to configure the settings.

Details of entry


Item	Details	Settable range	Default value
Door Delay	Door delay is a function that when the incubator is in N4: Outer Door Open notification (the outer door is open) status, instead of sounding immediately, the alarm will sound after the elapse of the set door delay time. **Note:** When the incubator recovers from the alarm status within the set door alarm delay time, the door alarm sound does not activate after the elapse of the door alarm delay time.	1 to 30 min	2 min
Ring Back	Ring back is a function that the audible alarm sounds again when an alarm status still continues after the ring back set time elapsed even though the alarm sound was stopped by tapping the Buzzer Stop button. The Ring Back can be turned ON or OFF by sliding the ON/OFF button. When turning it on, enter the ring back time in the right numeric input field. **Note:** For A1 (CO_2 gas empty) or the N4 (Outer Door Open), the alarm does not activate again (page 45).	ON/OFF 1 to 99 min	ON 30 min
Remote Alarm	Selects behavior of the remote alarm when the Buzzer Stop button is tapped to stop the audible alarm sound of the incubator. When "ON" (not linked with the Buzzer Stop button) is selected by sliding the button to the right, the alarm indication by the remote alarm device does not stop even if the Buzzer Stop button is tapped. ON: Not linked with the Buzzer Stop button (remote alarm continues). OFF: Linked with the Buzzer Stop button (remote alarm stops).	ON/OFF	ON

6. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

7. Tap the Home button in the Menu bar.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS

Alarm and buzzer stop button

The behaviors of the alarm indication and the ring back function when tapping the Buzzer Stop button are shown in the following tables.
Tapping the Buzzer Stop button while the audible alarm is sounding will stop the sound. However, tapping this button alone does not resolve the cause of the alarm, except for "A1: CO2 Gas Empty*". Refer to pages 85 to 86 for information about the cause of the alarm and how to resolve the issue.

* Tapping the Buzzer Stop button while the "A1: CO2 Gas Empty" alarm is active stops the alarm sound and clears the alarm.

Behaviors of alarm indication when tapping the Buzzer Stop button

<Setting of the incubator> <Alarm behavior>

Remote alarm setting	Ring back setting	Alarm sound from incubator		Remote alarm	
		When Buzzer Stop button is tapped	After the set ring back time has elapsed	When Buzzer Stop button is tapped	After the set ring back time has elapsed
ON: Not linked with Buzzer Stop button	ON	OFF (Alarm is not canceled)	ON	ON	ON (Alarm status is continuing)
	OFF		OFF		
OFF: Linked with Buzzer Stop button	ON		ON	OFF (Alarm is not canceled)	ON
	OFF		OFF		OFF

Note:
You cannot stop the audible alarm sound for the A97: Over Heat and A38: Door Lock Failure.

Behaviors of "A1: CO2 Gas Empty" alarm indication when tapping the Buzzer Stop button

<Setting of the incubator> <Alarm behavior>

Remote alarm setting	Ring back setting	Alarm sound from incubator		Remote alarm	
		When Buzzer Stop button is tapped	After the set ring back time has elapsed	When Buzzer Stop button is tapped	After the set ring back time has elapsed
ON: Not linked with Buzzer Stop button	ON	OFF (Alarm is canceled)	OFF (Alarm is already canceled)	Alarm is canceled	OFF
	OFF				
OFF: Linked with Buzzer Stop button	ON				
	OFF				

Note:
When the A1: CO2 Gas Empty alarm occurs, connect a new CO_2 gas cylinder properly and tap the Buzzer Stop button to stop the audible alarm. When the optional MCO-21GCP gas auto-changer is installed, and the gas supply is switched to the reserve gas cylinder, tap the Buzzer Stop button to stop the audible alarm, and replace the empty gas cylinder.

For the door opening notification "N4: Outer Door Open":
Tapping the Buzzer Stop button while the "N4: Outer Door Open" notification is active stops the alarm sound from the incubator. The "N4: Outer Door Open" notification is controlled by the incubator software and indicated by the incubator itself. However, the remote door status for the external remote alarm device is triggered by a signal from the DOOR STATUS terminals on the remote alarm terminals (page 29). Therefore, these behaviors are not linked. For this reason, you cannot stop the remote door opening indication on the external remote alarm device by tapping the Buzzer Stop button. In addition, for the remote door status indication, the door delay, remote alarm, ring back settings on the control panel are unavailable.

DOMINIQUE DUTSCHER SAS

DISINFECTION BY UV-LED

The water in the humidifying pan can be disinfected by UV-LED light. You can set if the UV-LED is turned on or not each time the outer door is closed. When set to ON, in order to reduce the heat impact on the chamber by the UV light, the UV-LED is automatically turned on or off during the automatically calculated time.


* The UV disinfection function can be used in either of the following cases. If neither condition is met, the UV button is not displayed on the Setting screen.
 - MCO-173AICUVH/173AIUV
 - MCO-173AIC with the optional UV system set (MCO-LUVS).

Setting the UV-LED

In this setting, you can select if the UV-LED is turned on or off after the outer door is closed and select the period for which the UV-LED is ON (Normal or Extended).
The UV light ON period is automatically calculated based on the set chamber temperature set point, selected “UV Light Time” (Normal or Extended), and the value automatically calculated and displayed at “UV Timer Ext.”


1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.

2. Tap the UV button.
►The UV screen is displayed.

3. Slide the ON/OFF buttons for “UV Light Mode,” and select Normal or Extended for “UV Light Time” if you turned ON the UV Light Mode.

Details of each item

Item	Details	Settable range	Default value
UV Light Mode	Slide the slide button to ON to set the UV-LED light to turn on after closing the outer door. When ON is selected, the UV-LED light ON period is automatically calculated and displayed in the right box (users cannot enter a value). The default ON period differs depending on the setting[*1].	ON/OFF	ON[*1]
UV Light Time	Select the duration of the UV-LED light ON period: Normal: For usual use. Extended: When longer disinfection time is needed. Based on this selection, actual ON period is calculated and displayed in the box at the UV Light Mode.	Normal/Extended	Normal
UV Total Time	Indicates the UV-LED’s total ON period up to now. If “1” is displayed in the box, it means that the UV-LED has been turned on for a total of 100 hours (1 x 100h). (Users cannot enter a value.)	-	0

DOMINIQUE DUTSCHER SAS

Details of each item

Item	Details	Settable range	Default value
UV Timer Ext.	Indicates how much UV-LED ON period is extended[*2]. The amount of UV irradiation decreases as the total UV-LED light ON period increases. To compensate for this, the UV-LED light ON period is automatically extended as the total UV-LED light ON period increases (this value is used by the incubator to calculate the time indicated at "UV Light Mode"). (Users cannot enter a value.)	-	0[*3]
UV ID No.	Left box: UV-LED ID number of the UV-LED currently used for the incubator is displayed. Right box: The UV-LED ID number of the latest UV-LED that can be used for the incubator is displayed.	-	1/1[*3]

*1: Default UV light ON period:

Temp. Set Point	UV Light Time	UV Timer Ext.	Default UV light ON period (min)
Temp. set point ≤ 40.0°C	Normal	0	30
	Extended		45
Temp. set point > 40.0°C	Normal		70
	Extended		105

*2: For example, the UV-LED light ON period is 30 minutes and the UV Timer Ext. is 10%, actual UV-LED ON period becomes 33 minutes.

*3: The default values differ depending on the types of UV-LED.

4. Tap the Apply button.

►The values are saved, and the screen returns to the Setting screen.

5. Tap the Home button in the menu bar.

►The screen returns to the Home screen.

Using UV-LED

WARNING

- **Do not look directly at UV light.**
 UV light is harmful to the eyes.

1. With all the inner attachment installed correctly, put in the cultivation samples on the trays in the chamber.

Notes:

- Be sure to install all the inner attachment correctly.
- Never turn the UV-LED on when the humidifying pan cover is removed.
- Always use the humidifying pan cover even when using the incubator with the UV-LED turned off. Using without humidifying pan cover may affect the chamber temperature distribution and humidity recovery.

2. Close the inner door and outer door.

►"UV : On" is displayed at the UV-LED status indicator on the Home screen, and the UV-LED turns on and off repeatedly for the time automatically calculated by the incubator.

Notes:

- If the outer door is opened while the UV-LED is lit, the light turns off. Then, when the door is closed, the UV-LED turns on and off repeatedly for the time automatically calculated by the incubator.
- If the outer door is opened and closed repeatedly, it may cause condensation in the chamber, and may affect the chamber temperature distribution because of the long-time heat generation by the UV-LED.
- If the outer door is not opened for more than 2 consecutive hours, the UV-LED automatically lights for 1/6 of the automatically preset ON time (e.g. for five minutes when the UV-LED light ON period is displayed as "30 min").
- If the UV-LED trouble occurs, a code number and message related to the trouble is displayed on the message display area. Refer to pages 85 and 86 for the details of the error and contact our sales representative or agent.

DOMINIQUE DUTSCHER SAS

SCREEN LOCK


Screen lock function prevents unauthorized access or changes to the incubator setting by prompting the user to enter a password when tapping a button in the menu bar.

Setting screen lock


The screen lock function can be turned ON as follows.

1. Tap the Setting button in the menu bar.
►The Setting screen is displayed.



2. Tap the Screen Lock button.
►The Screen Lock screen is displayed.

3. Configure the screen lock settings.

Details of entry

Item	Details	Settable range	Default value
Screen Lock	Slide the slide button to turn ON or OFF the screen lock function. To turn ON, slide the button to the right.	ON/OFF	OFF
Password	Enter a password to be used for unlocking the screen lock. You can use a 4 to 6 digit number for the password. Tap the Password input field to enter the password. **Note:** If the entered password is less than 4 digits, the following message is displayed. Tap the OK button and enter a 4 to 6 digit password. 	4 to 6 digit number	-
Confirm Password	Tap the Confirm Password input field to enter the same password as you entered in the Password input field and tap the OK button. **Note:** If the password does not match the password entered in the Password field, the following message is displayed. Tap the OK button and enter the same correct passwords in both input fields again. 	4 to 6 digit number	-

4. Tap the Apply buton.
►The values are saved, and the screen returns to the Setting screen.

5. Tap the Home button in the menu bar.
►The screen returns to the Home screen.


Note:
Be sure to manage the passwords properly.

DOMINIQUE DUTSCHER SAS


Operation when screen lock is ON

When the screen lock is set to ON, if you try to tap a button in the Menu bar, you are asked to enter a password.


1. On the Home screen, tap a button in the Menu bar.
►The Password input window is displayed.

2. Enter the screen lock password and tap the OK button.
►The screen corresponding to the tapped button is displayed.

Note:
When the password you entered is incorrect, the Notice message is displayed. Tap the OK button and return to Step 1 to enter the correct password.


DOMINIQUE DUTSCHER SAS

SCREEN LOCK


Canceling screen lock setting

You can turn off the screen lock setting by following the steps below.


1. Tap the Setting button in the menu bar.
►The Password input window is displayed.


2. Enter the password and tap the OK button.
►The Setting screen is displayed.

3. Tap the Screen Lock button.
►The Screen Lock screen is displayed.

4. Slide the Screen Lock slide button to the left.
►The screen lock is set to OFF.

5. Tap the Apply button.
►The setting is saved, and the screen returns to the Setting screen.

Note:
The registered password is deleted.

6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS

AUTOMATIC DOOR LOCK

When you set the automatic door lock function to ON, the outer door is automatically locked by the electric lock after it is closed, and a preset amount of time has elapsed. To unlock the outer door, you can select one of the following unlock modes:

• Quick mode: Tap the Unlock button in the menu bar.
• User-ID mode: After tapping the Unlock button in the menu bar, enter the User-ID and password.

The automatic door lock function can be used in either of the following cases. If neither condition is met, the Door unlock button in the menu bar is grayed out and becomes unavailable. Also, the Auto-Lock button and Auto-Lock User button are not displayed on the Setting screen.

- MCO-173AICUVH
- MCO-173AICUV/173AIC with the optional electric lock (MCO-170EL).


Note:
The outer door is also locked by the automatic door lock function when a power failure occurs, or the power cord is disconnected during use.

Setting user ID and password


To use the user-ID mode for unlocking the outer door, register the user ID and the password beforehand.

1. Tap the Setting button in the menu bar.
►The Setting screen is displayed.

2. Tap the Auto-Lock user button.
►The Auto-Lock User screen is displayed.

3. Enter a user ID and passwords.


Details of entry

Item	Details	Settable range	Default value
User ID	Tap the input field to enter a user ID to be used for unlocking the outer door.	1 to 8 digit alphanumeric characters	-
Password	Tap the input field to enter a password to be used for unlocking the outer door.	4 to 6 digit number	-
Confirm Password	Tap the input field to enter the same password as you entered into the “Password” input field and tap the OK button.	4 to 6 digit number	-


DOMINIQUE DUTSCHER SAS

Notes:


• If the entered user ID is the same as the user ID that is already registered, the following Confirm message is displayed. Tapping the Yes button overwrites the old password with the new one, and tapping the No button keeps the old password.


• If the entered password is less than 4 digits, the following message is displayed. Tap the OK button and enter a 4 to 6 digit password.


• If the password entered into the Confirm Password input field does not match the password entered in the Password field, the following Notice message is displayed. Tap the OK button and enter the same correct passwords in both input fields again.


• Up to 99 user IDs and passwords can be registered. If you try to register the 100th user ID, the following Notice message is displayed. In such a case, tap the OK button, and delete an unused user ID by referencing pages 53 and 54 before registering them.

4. Tap the Add button.
►The values are saved, and the following Information message is displayed.

Note:
If an error occurs, the following Notice message is displayed. Tap the OK button and register the user ID and password again.

5. Tap the OK button on the Notice message window.
►The Auto-Lock User screen is displayed.

6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.


Note:
To prevent unauthorized use of IDs and passwords, they should be properly managed by only qualified administrators.

DOMINIQUE DUTSCHER SAS

Changing the password


You can change the registered password for a user ID.

1. On the Auto-Lock User screen, enter the user ID already registered to the incubator. Then, enter a new password in both of the Password and Confirm Password input fields and tap the Add button.
►The following Confirm message is displayed.


2. Tap the Yes button.
►The new password is registered, and the Information message is displayed.

Note:
If you tap the No button, the password is not changed, and the screen returns to the Auto-Lock User screen.


3. Tap the OK button on the Information messsage window.
►The Auto-Lock User screen is displayed.

4. Tap the Home button on the menu bar.
►The screen returns to the Home screen.


Deleting user ID and password

1. On the Auto-Lock User screen, enter the user ID and password you want to delete into each input field. Then, tap the Delete button.
►The following Confirm message is displayed.


2. Tap the Yes button.
►The user ID and password are deleted, and the Information message is displayed.

Note:
If you tap the No button, the user ID and password are not deleted, and the screen returns to the Auto-Lock User screen.


DOMINIQUE DUTSCHER SAS

3. Tap the OK button.
►The Auto-Lock User screen is displayed.

Note:
If you delete all user IDs, and the User-ID mode has been set to “ON” on the Auto-Lock screen (pages 54-55), the following Information message notifying you that the User-ID mode is tuned off is displayed. Tap the OK button to return to the Auto-Lock User screen.


4. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

Setting automatic door lock


To turn on or off the automatic door lock function and to select the unlock mode, follow the steps below.

1.Tap the Setting button in the menu bar.
►The Setting screen is displayed.

2. Tap the Auto-Lock button.
►The Auto-Lock screen is displayed.

3. Configure each setting.

Details of entry


Item	Details	Settable range	Default value
Auto-Lock	Slide the slide button to turn ON or OFF the automatic door lock function. If you select to turn ON, tap the input field next to the slide button to enter the time period from the outer door is closed until the outer door is automatically locked by the electric lock. **Note:** Sliding to OFF while the outer door is locked unlocks the outer door.	ON/OFF 1 to 60 min	ON 1 min
User-ID	Select the unlock mode: ON: User ID mode You need to enter your user ID and password to unlock the outer door after tapping the Unlock button in the Menu bar (page 55). OFF: Quick mode You can unlock the door by just tapping the Unlock button in the Menu bar (page 55).	ON/OFF	OFF

DOMINIQUE DUTSCHER SAS

4. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

5. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

Notes:
• When no User-ID is registered, the following Notice message is displayed. Tap the OK button, and then register a User-ID and its password with reference to pages 51-52).


• In the User-ID mode, User-ID entered to unlock the outer door is saved as a part of log data that shows the open/close state of the outer door (pages 65-67).
• Even if the User-ID mode is turned OFF, registered User-IDs are not deleted.
• When all the registered User-IDs are deleted, the User-ID mode is turned OFF (pages 53-54).

Unlocking the outer door

To unlock the outer door, use either of the following depending on the set unlock mode.

Quick mode (“User-ID” setting is OFF)


Tap the Unlock button in the Menu bar.
►The outer door is unlocked.


User-ID mode (“User-ID” setting is ON)


1. Tap the Unlock button in the Menu bar.
►The User ID input window is displayed.

2. Enter your user ID and password and then tap the OK button.
►The outer door is unlocked.

Notes:
• The user ID entered to unlock the outer door is saved as a part of log data that shows the open/close state of the outer door (pages 65-67).
• If a wrong user ID or password is entered, one of the following Notice messages is displayed. Tap the OK button and enter a correct user ID or password.

• The outer door is automatically locked again after unlocking and closing the door, and a preset amount of time has elapsed.


DOMINIQUE DUTSCHER SAS

AUTOMATIC DOOR LOCK

Unlocking and locking the outer door during a power failure

In the event of a power failure, the outer door is locked by the electric lock. To unlock the outer door during a power failure, use the key that comes with the product.

Unlock

1. Insert the key into the key slot.

2. Turn the key 30° counterclockwise.
►The outer door is unlocked.

3. Open the outer door.


4. Turn the key to the position at which it is inserted and remove it.

Lock

1. With the outer door opened, check that the key is not inserted in the key slot (if you lock the outer door when the key is left in the slot, the electric lock may be damaged).

2. Set the latch part so that it comes in the central position.

3. Close the outer door.
►The outer door is locked.

Set the latch part so that it comes in the central position.
(Not the right or left side.)

Note:
The keys should be kept in a safe place and managed carefully. If the key is lost, contact our sales representative or agent.

DOMINIQUE DUTSCHER SAS

Canceling the automatic door lock


The automatic door lock function can be turned off by following steps below.

1. Tap the Setting button in the menu bar.
►The Setting screen is displayed.


2. Tap the Auto-Lock button.
►(User-ID mode) The User ID input window is displayed. Proceed to Step 3.
►(Quick mode) The Auto-Lock screen is displayed. Proceed to Step 4.

3. Enter your user ID and password, then tap the OK button.
►The Auto-Lock screen is displayed.

4. Slide the Auto-Lock slide button to the left.
►The Auto-Lock slide button changes to OFF.

5. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.


DOMINIQUE DUTSCHER SAS

SETTING DATE, TIME, AND TIME ZONE

You can set date, time, and time zone for the incubator.


1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.


2. Tap the ► button or the right page indicator.
►The second page of the Setting screen is displayed.

3. Tap the Date & Time button.
►The Date & Time screen is displayed.



4. Configure the settings for date and time.

Details of entry

Item	Details	Settable range	Default value
Format	Select the date format by tapping the radio button.	Year/Month/Day or Day/Month/Year	Day/Month/Year
Date	Enter the current date by tapping each input field and tapping the ▲/▼ buttons or numeric keys on the input window. Then, tap the OK button.	-	The date and time has been temporarily set at the factory. Set the correct date and time.
Time	Enter the current time by tapping each input field (hour: minute: second) and tapping the ▲/▼ buttons or numeric keys on the input window. Then, tap the OK button. The time is expressed using a 24-hour clock. **Note:** It is advisable to set the right time regularly since the error of about 1 minute may be observed within a month.	-	

DOMINIQUE DUTSCHER SAS

Item	Details	Settable range	Default value
Time Zone	Select the time zone of your country: 1. Tap the Time Zone input field. ►The UTC* Select window is displayed.  2. Scroll up or down by tapping the ▲/▼ buttons to locate your country's time zone and tap it. ►The selected time zone is highlighted in orange.  3. Tap the OK button. ►Selected time zone is confirmed, and the screen returns to the Date & Time screen.	Refer to the time zone table below.	The date and time has been temporarily set at the factory. Set the correct date and time.

*UTC: Coordinated Universal Time

List of time zones

Time zone	
UTC -12:00: International Date Line West	UTC +04:00: Gulf, Samara
UTC -11:00: Samoa	UTC +04:30: Kabul
UTC -10:00: Hawaii-Aleutian	UTC +05:00: Pakistan, Yekaterinburg
UTC -09:30: Marquesas islands	UTC +05:30: New Delhi, Mumbai, Colombo
UTC -09:00: Alaska	UTC +05:45: Kathmandu
UTC -08:00: Pacific	UTC +06:00: Bangladesh, Omsk
UTC -07:00: Mountain	UTC +06:30: Yangon
UTC -06:00: Central	UTC +07:00: Indochina, Novosibirsk
UTC -05:00: Eastern	UTC +08:00: China, Australian Western
UTC -04:30: Venezuela	UTC +08:45: Eucla
UTC -04:00: Atlantic	UTC +09:00: Japan, Korea
UTC -03:00: Brasilia, Greenland	UTC +09:30: Adelaide
UTC -02:00: Coordinated Universal Time-02	UTC +10:00: Australian Eastern, Vladivostok
UTC -01:00: Azores, Cape Verde	UTC +10:30: Lord Howe Island
UTC 00:00: Coordinated Universal Time, GMT	UTC +11:00: Solomon Islands, Magadan
UTC +01:00: Central European	UTC +12:00: New Zealand, Fiji
UTC +02:00: Eastern European	UTC +12:45: Chatham Island
UTC +03:00: Moscow, Eastern African	UTC +13:00: Coordinated Universal Time +13
UTC +03:30: Tehran	UTC +14:00: Christmas Island

5. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.


6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.


FURUNO FRANCE DUTSCHER SAS

SETTING BRIGHTNESS AND SLEEP

You can configure the brightness of the control panel and the sleep function.

1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.

2. Tap the ▷ button or the right page indicator.
►The second page of the Setting screen is displayed.

3. Tap the Brightness/Sleep button.
►The Brightness/Sleep screen is displayed.

4. Configure each setting on the Brightness/Sleep screen.

Details of entry

Item	Details	Settable range	Default value
Brightness (Active)	Brightness of the control panel in normal state. You can set the value using the slider bar or tapping the input field to enter the value.	40 to 100%	70%
Sleep	Sleep is a function that decreases the brightness of the control panel to save electricity when there is no touch screen operation. Slide the Sleep slide button to the right to turn ON the Sleep function. Then, enter the time period from the last touch screen operation to the time the Sleep state starts.	ON/OFF 1 to 5 min	ON 2 min
Brightness (Sleep)	Brightness of the control panel in sleep state. You can set the value using the slider bar or tapping the input field to enter the value.	0 to 40%	40%

5. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

6. Tap the Home button in the Menu bar.
►The screen returns to the Home screen.


DOMINIQUE DUTSCHER SAS


ANALOG OUTPUT

This product has a function to output the temperature and CO_2 density data as electric current data from the analog output terminals provided on the lower left side to an external device (recorder or data logger).

Setting the data to be output


1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.


2. Tap the ▷ button or the right page indicator.
►The second page of the Setting screen is displayed.

3. Tap the Analog Output button.
►The Analog Output screen is displayed.

4. Select ON or OFF for each item by sliding the slide button.

- Temperature: temperature in the chamber
- CO_2: CO_2 concentration in the chamber


5. Tap the Apply button.
►The values are saved, and the screen returns to the Setting screen.

6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

Connecting to external device

Connect the CO_2 and TEMP terminals to a commercially available recorder using commercially available cables. Be sure to connect the positive (+) and negative (-) terminals correctly.

Notes:
- When connecting the terminals, turn off the equipment and devices and unplug them. Otherwise, it may cause an electric shock.
- When connecting to the voltage terminal of the recorder, connect a commercially available shunt resistor (high precision) to the voltage terminal of the recorder.
 The formula for converting voltage reading to electric current value:
 Electric current (mA) = voltage reading (V) / shunt resistor value (Ω) x 1000
- Use shielded wire for connection. A maximum length of 30 meters of cable is recommended.

ANALOG OUTPUT

HUM and O_2 terminals are unavailable.

DOMINIQUE DUTSCHER SAS

ANALOG OUTPUT

Analog data output specification

Each control element is as follows.

	Range	Output
Temperature	25°C to 50°C	4 mA to 12 mA (0.32 mA/°C)
CO_2 density	0% to 20%	4 mA to 20 mA (0.8 mA/%)

Analog output characteristic (accuracy: ±3.0%FS)

Temp [°C]	Analog Output [mA]
25.0	4.0
27.5	4.8
30.0	5.6
32.5	6.4
35.0	7.2
37.5	8.0
40.0	8.8
42.5	9.6
45.0	10.4
47.5	11.2
50.0	12.0

CO_2 [%]	Analog Output [mA]
0	4.0
2	5.6
4	7.2
6	8.8
8	10.4
10	12.0
12	13.6
14	15.2
16	16.8
18	18.4
20	20.0

(Example)
Output of temperature 37°C: (37°C - 25°C) x 0.32 mA/°C + 4 mA = 7.84 mA
Output of CO_2 density 5%: 5% x 0.8 mA/% + 4 mA = 8 mA

DOMINIQUE DUTSCHER SAS

SETTING UNIT NAME

The unit name is used at the beginning of the file name when exporting the log data to a USB flash drive, and it is also recorded in the top part of the CSV file. By default, the serial number of the unit is set for the unit name.

(Example) When the unit name is “A00001”

File name: A00001_20251001_Data.csv


Top part of the CSV file: A00001

To change the unit name, follow the steps below.


1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.


2. Tap the ▷ button or the right page indicator.
►The second page of the Setting screen is displayed.

3. Tap the System Info. button.
►The System Info. screen is displayed.

4. Tap the input field for Unit Name.
►The alphanumeric input window is displayed.

5. Enter a unit name in alphanumeric characters or symbols (1 to 10 characters) and tap the OK button.
►The System Info. screen is displayed.

6. Tap the Apply button.
►The values are saved, and the screen returns to the Settings screen.

7. Tap the Home button in the sidebar menu.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS


REFERRING TO THE OPERATING INSTRUCTIONS USING MOBILE DEVICE


You can read the online operating instructions for this product on our website using mobile devices such as smartphones and tablets by scanning the 2D barcode displayed on the screen.

Notes:

- The content of the operating instructions on our website may be changed without prior notice due to product specification changes, regulatory changes, or other reasons. Therefore, the instructions provided on the website may differ from those included with the product you purchased.
- Please note that we accept no liability for any damages, including indirect damages such as losses due to differences in instructions, which may arise from discrepancies between the website instructions and the instructions included with your purchased product.


1. Tap the Setting button in the menu bar.
►The first page of the Setting screen is displayed.

2. Tap the ▷ button or the right page indicator.
►The second page of the Setting screen is displayed.

3. Tap the Operating Instructions button.
►The Operating Instructions screen is displayed

4. Scan the 2D barcode displayed on the screen by your mobile device.
►The operating instructions for the product are displayed on your device.

5. Tap the Home button.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS

EVENT LOG


This product is provided with a function to store event log data such as alarm messages.
- Log data type: Alarm messages, Notification messages, outer door status, settings, and system related events.
- Up to 40,000 logs can be stored. When the number of logs exceeds the maximum number, the oldest log data is overwritten with the new data.
- Log data is not recorded during a power failure.

Viewing the event log


You can view the event log data by selecting the duration and event types.

1. The event Log screen is displayed by tapping the Event Log button in the menu bar.
►The Event Log screen is displayed.

Event Log screen


No.	Item		Details
1	Event log indication		Indicates event log information stored in the incubator (date/time, details of an event). Tapping the ∨ mark on the right side of the Event column expands the indication area (2-line indication). Tapping the ∧ mark on the right side of the Event column collapses the indication area (1-line indication). The following information is displayed: The first line: event name The second line: details of the event (the second line is not displayed for system related event)
2	Period setting	Number of days	Tap the input field to enter the period of days you want to display the events that occurred on and before the date specified at 3 (base date) shown in the above figure. The settable range is 1 to 99 days.
3		Base date	Tap the input field to enter the base date for selecting the period of days you want to display events. Tapping ◄ and ► buttons can display the previous date and next date, respectively.
4	Filtering status indication		All: Indicates all event logs are displayed. Filter On: Indicates selected logs are displayed.
5	Export button		Tapping this button exports the selected event log data to the USB flash drive inserted into the USB port in CSV format.

DOMINIQUE DUTSCHER

EVENT LOG


No.	Item	Details
6	Filtering item selection pull-down	Tapping ∨ button displays the pull-down lists from where you can select the types of event log you wish to display. Select the types of event log by tapping the checkbox.
7	Scroll buttons	Tapping the ▲ button scrolls up one line. Dragging the scroll box up and down scrolls through the window. Tapping the ▼ button scrolls down one line.

2. Tap the base date input field.
►The date input window is displayed.


3. Tap the ▲/▼ buttons to enter the date and tap the OK button.
►The Event Log screen is displayed.
Note:
Tapping ◄ and ► buttons can display the previous date and next date, respectively.

4. Tap the Number of days input field.
►The numeric input window is displayed.

5. Enter the number of days you want to display the events that occurred on and before the base date specified in step 3.


6. Tap the filtering item selection pull-down and select the event types you want to display.


Items in the Filtering item selection list

Item	Details
All	All event logs
Alarm	Alarm message log
Notification	Notification message log
Door	Outer door open/close status log
Setting	Setting button access log
System	System related event log

- You can select multiple items.
- To deselect an item, tap its checkbox again.


7. Tap the filtering item selection pull-down ∧ again to collapse the pull down.
►The selected types of events occurred during the selected period of days are displayed.

DOMINIQUE DUTSCHER SAS

Notes:

- When all the event log items are selected, “All” is displayed at the filtering status indication. When specific event log items are selected, “Filter On” is displayed there.
- You can scroll up and down the screen to view the items on the next or previous page using the scroll buttons or scroll box.


8. After checking the logs, tap the Home button in the sidebar menu.

Exporting event log data


The event log data displayed on the Event Log screen can be exported to a USB flash drive inserted into the USB port in CSV format.

1. With the desired data displayed on the Event Log screen, insert a USB flash drive into the USB port.
For supported USB flash drives, refer to “1. USB port” on page 12.

2. Tap the Export button at the upper right corner of the Event Log screen.
►Export of the selected event log data starts.

3. Tap the OK button on the Information message indicating that the export has been completed.
►The Event Log screen is displayed.
A folder named “log” is created in the USB flash drive and the exported data file is saved in the folder in CSV format.
The file name consists of the unit name, date (YYYY/MM/DD or DD/Month/YYYY), and type (AllEvent or Event). The date format specified on the Display screen is used for the date (see page 58).

DOMINIQUE DUTSCHER SAS

EVENT LOG

Examples of the file names are as follows:

The data is exported with the unit name “A00001” and the period is last 4 days from October 4 in 2025.





Filtering item	YYYY/MM/DD	DD/Month/YYYY
All Filter On	A00001_20251001_20251004_Event.csv	A00001_01Oct2025_04Oct2025_Event.csv

The data is exported with the unit name “A00001” and the period is last 1 day from October 1 in 2025.

Filtering item	YYYY/MM/DD	DD/Month/YYYY
All Filter On	A00001_20251001_Event.csv	A00001_01Oct2025_Event.csv

Notes:

- If a unit name is not specified, the name of the exported file will be [unit’s serial number]_[date]_Data.csv.
- If there is an existing same name file, the old file is overwritten with a new file.
- The event log data stored in the incubator is not deleted even if the event log data is exported to the USB flash drive.
- If an error occurs, one of the following Error or Notice messages will appear.

● The USB flash drive is not inserted into the USB port. ● The USB flash drive is not inserted properly. ● The USB flash drive is not formatted in FAT16 or FAT32 format.  ►Tap the OK button and check if the USB flash drive is inserted into the USB port correctly and check the file system format.	● There is not enough free space in the USB flash drive.  ►Tap the OK button and check the memory capacity of the USB flash drive.
● The data during the specified period does not exist.  ►Tap the OK button and check the specified data.	● An Error has occurred during the export process.  ►Tap the OK button and retry.

4. Remove the USB flash drive from the USB port.

5. Tap the Home button in the Menu bar.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS

OPERATION LOG CHART

This product is equipped with a function to store operation log data (chamber temperature, CO_2 density, and outer door open/close status).


- Up to about 200 day's operation log data (fixed interval of every two minutes) can be stored. When log data reaches the maximum storage amount, the oldest log data is overwritten with the new data.
- Log data is not recorded during a power failure.

Displaying operation log chart

Operation log data (chamber temperature, CO_2 density, and outer door open/close status) stored in the incubator can be displayed on the screen in chart format. On the chart, up to two types of log can be displayed using the primary and secondary vertical axes.


1. Tap the Chart button in the menu bar.
►The Chart screen is displayed.

Note:
Chart with Actual Temp (°C) log data is displayed at this point .


2. Tap the list button.
►The Chart / Data Select screen is displayed.


3. Select each log for the primary axis and the secondary axis.

1) Primary vertical axis:
Tap the 1st radio button and select one type of log from the list.
►① appears on the left of the selected item and the item is highlighted in orange.

Actual Temp.(°C): Temperature log
Actual CO2 Level (%): CO_2 concentration log
Door Opening: Door open and close log


2) Secondary vertical axis:
Tap the 2nd radio button and select one type of log from the list.
►② appears on the left of the selected item and the item is highlighted in green.
Note:
You can choose not to select the secondary axis and display only the log type for the primary vertical axis.


DOMINIQUE DUTSCHER SAS

OPERATION LOG CHART


4. Tap the Date input field to enter the date of the operation log you want to display.
►The numeric input window is displayed.

5. Tap the Apply button.
►Operation log chart on the selected date and types is displayed.

Details of the Chart screen

No.	Item		Details
1	Period setting	Base date	Tapping the input field displays the data input window for entering the base date for selecting the period of days you want to see log data. Tapping ◄ and ► buttons can display the previous date and next date, respectively.
2		Period selection radio button	Select the period of days you want to see the log data on and before the date specified for base date (1) shown in the above figure. 1d: Only the base date 2d: The base date and the previous day 4d: Last four days including the base date
3	Export button		Tapping this button exports the operation log data displayed on the chart to a USB flash drive inserted in the USB port in CSV format.

DOMINIQUE DUTSCHER SAS

No.	Item	Details
4	Log name	The operation log names are displayed.
		Indicates name of the log selected for primary vertical axis. Tapping the list button () displays the Chart/Data Select screen where you can change the log type for the primary and secondary vertical axes.
		Indicates name of the log selected for secondary vertical axis. Tapping the log change button () can change the log for secondary vertical axis. Changes the log data currently selected for the secondary vertical axis to the one displayed one line up in the Card/Data Select screen. Changes the log data currently selected for the secondary vertical axis to the one displayed one line down in the Card/Data Select screen.
5	Chart display area	Chart is displayed in this area. The data series for the primary vertical axis is shown with an orange line, and the data series for the secondary vertical axis is shown with a green line.
6	Scroll buttons	Tapping the ▲ ▼ buttons can vertically scroll up and down the scale of the chart.
7	Scale of secondary vertical axis	Indicates scale of secondary vertical axis.
8	Time axis	Indicates time. When 2 Day or 4 Day duration is selected, date is also displayed under the time.
9	Scale up/down buttons	Tapping the + button enlarges the scale of the vertical axis. Tapping the - button reduces the scale of the vertical axis.
10	Scale of primary vertical axis	Indicates scale of primary vertical axis.
11	Alarm setting value indication	Selecting this checkbox displays pink lines in the chart that represent the upper and lower alarm activation values for the data series associated with the primary vertical axis, which are set on the Alarm #1 screen (page 43).


6. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

Exporting the operation log chart data displayed on the screen

The operation log data during the selected period on the Chart screen can be exported to a USB flash drive inserted into the USB port in CSV format. All operation log data is exported in a single file. You cannot export each data separately.


1. Insert a USB flash drive into the USB port.
For supported USB flash drives, refer to “1. USB port” on page 12.

2. Tap the Export button on the screen displaying the operation log data to be exported.
►Export of the operation log data starts.

DOMINIQUE DUTSCHER SAS

3. Tap the OK button when the Information message notifying the completion of the export is displayed.

►The screen returns to the Chart screen.
A folder named “log” is created in the USB flash drive and the exported data file is saved in the folder in CSV format.
The file name consists of the unit name, date (YYYY/MM/DD or DD/Month/YYYY), and data name (Data). The date format specified on the Date & Time screen is used for the date (see page 58).

e.g. When the log data for October 4, 2025, with 4d and the unit name A00001 selected is exported, the file name is as follows.

YYYY/MM/DD	DD/Month/YYYY
A00001_20251001_20251004_Data.csv	A00001_01Oct2025_04Oct2025_Data.csv

e.g. When the log data for October 1, 2025, with 1d and the unit name A00001 is exported, the file name is as follows.

Year/Month/Day	Day/Month/Year
A00001_20251001_Data.csv	A00001_01Oct2025_Data.csv

Notes:
- If a unit name is not specified, the name of the exported file will be [unit’s serial number]_[date]_Data.csv.
- If there is an existing same name file, the old file is overwritten with a new file.
- The operation log data stored in the incubator is not deleted even if the operation log data is exported to the USB flash drive.
- If an export error occurs, one of the Error or Notice messages shown in the table on page 68 appears.

4. Remove the USB flash drive from the USB port.

5. Tap the Home button in the menu bar.
►The screen returns to the Home screen.

DOMINIQUE DUTSCHER SAS


EXPORTING LOG DATA

Log (event log and operation log) data stored in the incubator can be exported to a USB flash drive inserted into the USB port in CSV format.


1. Insert a USB flash drive into the USB port.
For supported USB flash drives, refer to ”USB port” on page 12.

2. Tap the Import/Export button in the menu bar.
►The Import/Export screen is displayed.


3. Select types of the log data (Item) and date range (Range) of the data to be exported.


Details of setting:

<table>
<tr><th>Item</th><th>Type</th><th colspan="2">Details</th></tr>
<tr><td rowspan="2">Item</td><td>Status</td><td colspan="2">When this checkbox is selected, operation log (temperature, CO_2 density, and door open/close status log) data is exported.</td></tr>
<tr><td>Event</td><td colspan="2">When this checkbox is selected, event log data is exported. If “All” is selected, all event log data is exported. To select the event types you want to export, tap the pull-down button ∨ and choose the desired types from the pull-down list. You can select multiple event types. Once the desired types are selected, the indication on the pull-down menu changes to “Filter On.”
<table>
<tr><td>All</td><td>All event logs</td></tr>
<tr><td>Alarm</td><td>Alarm message log</td></tr>
<tr><td>Notification</td><td>Notification message log</td></tr>
<tr><td>Door</td><td>Outer door open/close status log</td></tr>
<tr><td>Setting</td><td>Setting button access log</td></tr>
<tr><td>System</td><td>System related event log</td></tr>
</table>

</td></tr>
<tr><td rowspan="2">Range</td><td>All</td><td>All data stored</td><td rowspan="2">Tap the All or Custom radio button.
If you tap the Custom radio button, specify the start date and end date by tapping the date input field.</td></tr>
<tr><td>Custom</td><td>Data during specified period</td></tr>
</table>

DOMINO DUE TO SCALE

EXPORTING LOG DATA

4. Tap the Export button at the upper part of the Import/Export screen.
►Export of the selected type of data starts.

5. Tap the OK button on the Information message indicating that the export has been completed.
►The Import/Export screen is displayed.
A folder named “log” is created in the USB flash drive and the exported data file is saved in the folder in CSV format.
The file name consists of the unit name, date (YYYY/MM/DD or DD/Month/YYYY), and type (AllData or Event). The date format specified on the Date & Time screen is used for the date (see page 58).

Examples of the file names are as follows:

The data is exported with the unit name “A00001” and ALL (October 1 in 2025 to January 30 in 2026) selected.

Item	Date format	File name
Status	YYYY/MM/DD	A00001_20251001_20260130_AllData.csv
	DD/Month/YYYY	A00001_01Oct2025_30Jan2026_AllData.csv
Event	YYYY/MM/DD	A00001_20251001_20260130_Event.csv
	DD/Month/YYYY	A00001_01Oct2025_30Jan2026_Event.csv

The data is exported with the unit name “A00001” and Custom (January 30 in 2026) selected.

Item	Date format	File name
Status	YYYY/MM/DD	A00001_20260130_AllData.csv
	DD/Month/YYYY	A00001_30Jan2026_AllData.csv
Event	YYYY/MM/DD	A00001_20260130_Event.csv
	DD/Month/YYYY	A00001_30Jan2026_Event.csv

Notes:
- If a unit name is not specified, the name of the exported file will be [unit’s serial number]_[date]_[AllData or Event].csv.
- If there is an existing same name file, the old file is overwritten with a new file.
- The log data stored in the incubator is not deleted even if the log data is exported to the USB flash drive.
- If an export error occurs, one of the Error or Notice messages on page 68 appears.

6. Remove the USB flash drive from the USB port.

7. Tap the Home button in the menu bar.
►The Home screen is displayed.

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GAS AUTO CHANGER (OPTION)


When the optional gas auto changer MCO-21GCP is installed, two connection ports for CO_2 gas pipes A and B are available, and you can connect two CO_2 gas cylinders. This enables the CO_2 gas supply line to be automatically switched when one of the CO_2 gas cylinders becomes empty.

Connecting CO_2 gas cylinders

- For purchase and installation of the gas auto changer MCO-21GCP, contact our sales representative or agent.
- Before connecting CO_2 gas cylinders, read “SAFETY PRECAUTIONS” on pages 5-10.
- For the tubes, use each gas tube that came with this product and optional gas regulator.
- Make sure that the connected tubes are not bent.
- Do not connect multiple incubators to a single gas cylinder. If the CO_2 gas is supplied to multiple incubators from a single gas cylinder, the intended performance may not be achieved.

1. Prepare two CO_2 gas cylinders (CO_2 gas cylinder A and B) and attach the optional gas regulators MCO-010R to both of the CO_2 gas cylinders.

Notes:
- Use a liquefied-CO_2 gas cylinder (at least 99.5 % pure). The siphon (dip tube) type cannot be used.
- When MCO-010R is not available, attach a two-stage gas pressure regulator rated at 25 MPa(G) (250 kgf/cm^2(G), 3,600 psi(G)) for the primary side, and 0.25 MPa(G) (2.5 kgf/cm^2(G), 36 psi(G)) for the secondary side.

2. Connect the CO_2 gas connection port A and the gas regulator for the CO_2 gas cylinder A using the gas tube and conversion joint that came with this product.

3. Connect the CO_2 gas connection port B and the gas regulator for the CO_2 gas cylinder B using the gas tube and conversion joint that came with the gas auto changer MCO-21GCP.
Note: For the connection of the product and the gas cylinder, refer to pages 26-28.

4. After connecting the gas tube, make sure that no gas is leaking (e.g. by using a gas leak detection spray).


5. For both CO_2 gas cylinders A and B, adjust the secondary side pressure of the gas regulator to 0.03MPa(G) – 0.1 MPa(G) (0.3 kgf/cm^2(G) – 1 kgf/cm^2(G), 4.4 psi(G) –14.5 psi(G)) while CO_2 gas is injecting.
Note: Do not set the secondary pressure too high. As the pressure increases, the CO_2 gas density will fluctuate widely. Also, excessive pressure may cause gas supply lines inside the incubator to come loose, which may result in gas poisoning or oxygen deprivation due to leaking gas. If gas lines come loose, the incubator must be repaired.

DOMINIQUE DUTSCHER SAS

GAS AUTO CHANGER (OPTION)

Automatic CO_2 gas supply line switching

When the optional gas auto changer MCO-21GCP is installed, the following indicator is displayed on the Home screen.

CO_2 gas supply line indicator A and B:

The icon of the CO_2 gas line currently supplying CO_2 gas to the incubator turns orange. When the current CO_2 gas cylinder becomes empty and the gas supply line automatically switches to the other one, the gas line indicator of the empty gas cylinder is displayed in reverse video and blinks. Also, you can switch the CO_2 gas supply line manually by sliding this button to A or B.

When the CO_2 density level in the chamber remains unchanged even though the CO_2 gas valve in the incubator is open, the incubator assumes the currently connected CO_2 gas cylinder to be empty and automatically switches the CO_2 gas supply line to the other one. The following table shows the automatic switching of the line.

Automatic switching from CO_2 gas cylinder A to B

	Status	CO_2 gas			CO_2 gas supply line indicator		
		Supply line	Cylinder A	Cylinder B		A	B
1	CO_2 gas is being supplied from cylinder A.	A	Remaining	Remaining		ON	OFF
2	CO_2 density in the chamber does not increase even if CO_2 gas valve is open.	A	Empty	Remaining		ON	OFF
3	If the above status continues for 2 or 3 minutes, the CO_2 gas supply line is automatically switched to cylinder B. At this time, the CO_2 gas empty alarm activates, and an audible alarm sounds with the indication of the message “A1: CO2 Gas Empty” on the message display area on the Home screen.	B	Empty	Remaining		Reverse video and blink	ON
4	If you tap the Buzzer Stop button, the audible alarm stops, indicator A turns off, and the message disappears.	B	Empty	Remaining		OFF	ON
5	Empty gas cylinder A is replaced with a new one.	B	Remaining	Remaining		OFF	ON


Notes:

- The switching of CO_2 gas supply line is determined by an increase of CO_2 density in the chamber. The switching of CO_2 gas cylinders may also occur in situations where the gas tube is clogged, the gas is leaking, the gas pressure has dropped down, or the valve on the gas cylinder is insufficiently open etc., even though the CO_2 gas cylinder being used is not empty. Check the remaining amount of gas in the gas cylinder before removing it.
- When the Buzzer Stop button is not tapped and the cylinder is not replaced with a new one when both CO_2 gas cylinder A and B gets empty, the switch operation between CO_2 gas supply lines A and B will be repeated. In this case, replace both CO_2 gas cylinders A and B with new ones immediately, and tap the Buzzer Stop button.

DOMINIQUE DUTSCHER SAS

Manual CO_2 gas supply line switching

The CO_2 gas supply lines A and B can be switched manually at any time by sliding the CO_2 gas supply line indicator button displayed on the Home screen.

The following table shows the behavior when the CO_2 gas supply line is switched by the CO_2 gas automatic switching function and then manually switched back to the empty CO_2 gas supply line before tapping the Buzzer Stop button.

Manual switching from cylinder B to empty cylinder A

Status		CO_2 gas			CO_2 gas supply line indicator		
		Supply line	Cylinder A	Cylinder B		A	B
1	CO_2 gas is automatically switched to the supply line B.	B	Empty	Remaining		Reverse video and blink	ON
2	Without tapping the Buzzer Stop button, the CO_2 gas supply line select button is long tapped.	A	Empty	Remaining		Blink	OFF

In both 1 and 2 cases, the CO_2 gas empty alarm activates, and an audible alarm sounds with the indication of the message “A1: CO2 Gas Empty” on the message display area on the Home screen.

DOMINIQUE DUTSCHER SAS

ROUTINE MAINTENANCE

Clean the exterior and interior of the incubator and inner attachments once a month. Regular cleaning keeps the incubator in good condition.

When the chamber is contaminated or when cleaning the chamber prior to starting a culture, you can perform H_2O_2 decontamination (pages 79-83).

Also, when the product will not be used for a long periods of time, dispose of the water in the humidifying pan. Then, eliminate the moisture in the chamber and dry it enough before closing the door. Not using the incubator for a long period of time leaving the chamber in high humidity status may cause the failure of the incubator.

Cleaning the exterior, interior, and inner attachments

CAUTION	**Wear rubber gloves and a mask when cleaning the chamber.** You may get injured by the corners of interior parts. Also, touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.

1. Take the humidifying pan out of the incubator, dispose of the water in it, and turn off the power to the incubator.

2. Take the inner attachments out of the chamber (pages 21-22).

3. Clean the exterior, inner attachments, and chamber (pages 23-24).
 If the outside panels are dirty, clean them by wiping them with a gauze moistened with a diluted neutral detergent (using an undiluted solution of detergent may cause the unit's plastic areas to crack. Follow the directions on the detergent for details of dilution). After that, be absolutely sure to wipe the surfaces using a cloth moistened with clean water to remove traces of the detergent and wipe the surfaces with a dry cloth. After that, wipe the surfaces with a cloth moistened with ethanol for disinfection.

4. Put back the inner attachments to the original positions (pages 24-25).

5. Open the outer and inner doors for a while to dry the chamber.

6. Wipe the control panel with a dry cloth.
 If the stain will not come out, wipe it using a cloth moistened with a proper amount of water (the amount that cannot form droplets) and then, wipe thoroughly with a dry cloth. If the water enters the control panel, it may cause the panel to malfunction.

Notes:
- If the stain is stubborn, consult our sales representative or agent.
- Do not use a brush, acid, thinner, laundry soap, powder detergent, or boiling water for cleaning. These may cause damage to painted surfaces or cause perishing of plastic and rubber components. Moreover, do not wipe plastic and rubber components with a volatile material.
- Do not clean the entire incubator unit using peracetic acid or hydrogen peroxide, or do not disinfect the unit by formalin fumigation.
- Never pour water onto or into the unit. Doing so can damage the electric insulation and cause failure.

DOMINIQUE DUTSCHER SAS

H_2O_2 decontamination

When the chamber is contaminated or when cleaning the chamber prior to starting a culture, you can perform H_2O_2 decontamination to decontaminate the chamber

* The H_2O_2 decontamination function can be used in either of the following cases. If neither condition is met, the Cleaning button in the menu bar is grayed out and becomes unavailable.
 - MCO-173AICUVH
 - MCO-173AICUV/173AIC with the optional H_2O_2 decon kit (MCO-HK) and electric lock (MCO-170EL).

To perform H_2O_2 decontamination, you have to purchase optional H_2O_2 generator (MCO-50HP) and H_2O_2 reagent (MCO-20H2O2).

Notes:
- Before performing decontamination, read the safety precautions on pages 5-10.
- The door is locked by the electric lock during the decontamination procedure. However, if the key that came with this product has been inserted into the key slot and left in the unlocked position, the door cannot be locked by the electric lock. When performing decontamination, be sure to take the key out of the key slot and store and manage the key in a safe place.
- If the optional small door (MCO-173ID) is attached to the incubator, you can perform H_2O_2 decontamination with each door closed.

Preparation for H_2O_2 decontamination

Before starting H_2O_2 decontamination, the following preparation is required.

Note:
Be sure to turn off the power to the incubator before the preparation.

1. Clean the chamber and inner attachments.
1) Take out the trays, humidifying pan, humidifying pan cover, fan cover, fan, duct, and connector cap from the chamber (pages 21-22).
2) Dispose of the water in the humidifying pan, and wipe the inner attachments removed from the chamber with a gauze containing ethanol for disinfection.
3) Wipe the inner wall of the chamber with a gauze containing ethanol for disinfection.
4) Attach the fan, duct, and fan cover (pages 24-25).

2. Set the H_2O_2 generator.
1) Pour a whole bottle of the H_2O_2 reagent (MCO-20H2O2) into the H_2O_2 generator (MCO-50HP) and then cover it with the lid. Do not use chemicals other than the H_2O_2 reagent specified by our company. Doing so may cause malfunction.

Observe the following precautions when using the optional H_2O_2 generator MCO-50HP and MCO-20H2O2 (H_2O_2 reagent).
● Cautions for handling
- Wear protective equipment such as protective glasses and rubber gloves when handling the H_2O_2 reagent.
- Do not recklessly use fire in an area where H_2O_2 reagent is handled.
- Do not leave H_2O_2 reagent in the bottle after or during use.
- Do not place combustible or inflammable objects near areas where the H_2O_2 reagent is handled.

● Cautions for storage
- Once the bottle of H_2O_2 reagent is opened, it must not be stored for future use.
- Store the H_2O_2 reagent in a cool, dark place or in a refrigerator.
- Make sure the bottle for H_2O_2 reagent is free of damage such as corrosion or cracking.
- Store the bottles for H_2O_2 reagent with the mouth facing upwards, and make sure that they cannot topple or fall.


● Cautions for disposal
- Dilute the remaining H_2O_2 reagent with tap water for disposal.


DOMINIQUE DUTSCHER SAS

2) Remove the connector cap on the H_2O_2 generator. Connect one end of the waterproof cable (included with the H_2O_2 generator) to the H_2O_2 generator and then, set the positioning fittings of the H_2O_2 generator into the slits of the duct.


Notes:

- Ensure that the H_2O_2 generator is properly set. An improper setting may cause inadequate decontamination.
- Keep the connector cap for later use.

3) Connect the other end of the waterproof cable to the H_2O_2 decontamination connector terminal at the rear central left part of the chamber.

3. Set the inner attachments for decontamination as follows.

1) Check that the duct, fan cover, and H_2O_2 generator are set properly.


2) Set the trays on the 2nd, 3rd, 4th, and 6th tray runners from the top.

Note:
Be sure to set the trays that came with the incubator. If the optional half trays (MCO-25ST) or trays for previous models are used, decontamination may not be sufficiently effective.

3) Set the humidifying pan and the humidifying pan cover at the position shown in the right figure.

Notes:

- On H_2O_2 decontamination, only the chamber and inner attachments can be decontaminated. Do not place anything other than the designated inner attachments into the chamber during decontamination.
- Position the humidifying pan cover on the back part of the tray so that it does not come into contact with the inner door.

4) Close the inner door and outer door.

Note:
When the optional small door (MCO-173ID) is attached to the incubator, check that each gastight split door is closed.


DOMINIQUE DUTSCHER SAS

Steps for performing H_2O_2 decontamination

After the preparation, turn on the incubator and follow the steps below to start H_2O_2 decontamination.


1. Tap the Cleaning button in the menu bar.
►The Cleaning screen is displayed.

2. Tap the H_2O_2 Decontamination button.
►The screen showing where to attach the connector cap is displayed.


3. Check that the inner attachments and H_2O_2 generator are properly set by looking at the description on the screen. Then, tap the Next button.
►The process of H_2O_2 decontamination Step 1 starts.

Note:
Perform H_2O_2 decontamination with the inner attachments arranged as specified by our company. Arranging them in a different way may result in insufficient decontamination.

Step 1 (Step indicator: )

The outer door is locked by the electric lock, and System check is performed automatically.
►After the system check finishes successfully, the process proceeds to Step 2.
If there is a problem, a message is displayed. In such a case, refer to the description of “System check error when starting H_2O_2 decontamination” for the cause of the problem on page 86.

About step indicator:
During the process, the step indicator is displayed at the upper right of the window. It indicates the progress of the process by showing the current step.

Gray: Completed step

Inactive: Incomplete step


Blue: Ongoing step

DOMINIQUE DUTSCHER SAS

Step 2 (Step indicator:)

4. Tap the OK button.

► H_2O_2 decontamination is performed (Step 3 to Step 5).

CAUTION **Do not unlock the outer door using the accessory key during H_2O_2 decontamination or during H_2O_2 gas resolve (decomposition).** Doing so may be harmful to health due to H_2O_2 gas leakage.


Step 3 (Step indicator:)

The chamber is warmed up to 48°C.

►When the temperature reaches 48°C, the process proceeds to Step 4.

Note:

If you tap the Abort button during Step 3, Step 3 is finished incompletely, and the "H_2O_2 Decontamination Successful" (Step 6) screen is displayed.


Step 4 (Step indicator:)

The H_2O_2 generator generates H_2O_2 gas to decontaminate the chamber.

►After H_2O_2 decontamination, the process proceeds to Step 5.

Note:

If you tap the Abort button to abort the decontamination during Step 4, Step 4 is finished incompletely and proceeds to H_2O_2 gas resolve (decomposition) step (Step 5).


Step 5 (Step indicator:)

The H_2O_2 gas in the chamber is resolved (decomposed).

►On completion of the gas resolve, "H_2O_2 Decontamination Successful" (Step 6) is displayed, the audible alram sounds, and the outer door is unlocked.

Note:

You cannot abort during Step 5.

Step 6 (Step indicator:)

5. Open the outer door and inner door to remove the waterproof cable connector from the H_2O_2 decontamination connector terminal and take out of the H_2O_2 generator.

►The message "Please remove H_2O_2 unit." disappears.

Note:

Wear protective glasses and rubber gloves to perform this step, and also, wipe the liquid including H_2O_2 left at the bottom of the chamber, where the H_2O_2 generator was positioned, with a dry cloth (e.g. unwoven cloth).


DOMINIQUE DUTSCHER SAS

6. Tap the OK button.


►The screen returns to the Home screen.

Notes:


●When the H_2O_2 generator is connected to the incubator ("Please remove H_2O_2 unit." is displayed on the Step 6 screen), the OK button is grayed out and unavailable.

●If an error occurs during the decontamination process, its code number and message are displayed on the message display area. Take an appropriate measure by referring to the table "Errors related to H_2O_2 decontamination" on page 86.


●If the decontamination process is canceled for any reasons, its code number and message are displayed on the message display area. Refer to 86 and retry the H_2O_2 decontamination process from the beginning.

●If a power failure occurs during the decontamination process, refer to page 84 for the operation after recovery from the power failure.

7. Ventilate the chamber sufficiently and put the inner attachments back to the original positions.

After performing H_2O_2 decontamination, properly cover the H_2O_2 decontamination connector terminal with the connector cap. Improper covering may cause the connector to rust.


Notes:

- After completing H_2O_2 decontamination, residual H_2O_2 liquid may remain at the bottom of the chamber and under the H_2O_2 generator duct. Be sure to wear protective goggles and rubber gloves, and wipe it off using non-woven fabric or similar materials.
- Wait at least one hour before starting cultivation until the incubation environment becomes stable.

Cleaning the H_2O_2 generator after decontamination

- Discard the remaining H_2O_2 reagent in the generator, wipe the exterior with a dry cloth, then dry thoroughly, and store it in a clean environment outside the device. At this time, do not touch the spray nozzle of the H_2O_2 generator or pour liquid into it. This may cause errors due to damage.
- Dispose of the discarded H_2O_2 reagent by diluting it with water. If dilution is insufficient, it may cause rust in the drainage destination.
- Please avoid cleaning the H_2O_2 generator with disinfectant ethanol.
- When handling the H_2O_2 generator, do not hold the positioning bracket. This may cause deformation or damage.

H_2O_2 generator

DOMINIQUE DUTSCHER SAS

DURING/AFTER POWER FAILURE

Operation during power failure

When a power failure occurs, the incubator behaves as follows.

- All of the incubator operations stop.
- Nothing is displayed on the control panel
- The outer door is locked by the electric lock.
- The operation log and event log are not recorded during a power failure.
- The clock does not stop.

Operation after recovery from power failure

During a power failure, the setting values are stored in the nonvolatile memory. Therefore, after the power comes back, the incubator resumes the operation using the setting values before the power failure. Reconfiguration for setting values is not required but check the settings such as operation settings and alarm settings.

Notes:
- After the power comes back, an alarm may be activated due to CO_2 density or increased or decreased temperature during the power failure.
- If the temperature at the CO_2 sensor has decreased during a power failure, after the power comes back, the CO_2 control does not start and indicates “N1: Gas Sensor initializing” message on the message display area until the CO_2 sensor becomes stable.

When a power failure occurs during the following operation, the operation resumes from the step described below.

Power failure when performing H_2O_2 decontamination

- During Step 1 to Step 3 (pages 81-82):
 The decontamination procedure is terminated, and the Step 6 screen is displayed. At this time, alarm A35 is displayed on the message display area, and the audible alarm sounds (page 86). The decontamination procedure has not completed. Try the steps again from the beginning.
- During Step 4 and Step 5 (page 82):
 The decontamination procedure restarts from Step 3 and when the chamber temperature reaches the temperature set point, the procedure goes to Step 5 to perform H_2O_2 gas resolve in Step 5. After that, the Step 6 screen is displayed with alarm A35 displayed on the message display area. The decontamination procedure has not been completed. Try the steps again from the beginning.
 Note:
 In this case, the Abort button on the Step 3 screen is grayed out and unavailable.
- During Step 6 (page 82):
 After the power comes back, the Step 6 screen is displayed. At this time, alarm A35 is displayed on the message display area, and the audible alarm sounds (page 86). In this case, the decontamination procedure has been completed.

DOMINIQUE DUTSCHER SAS

ALARMS AND SELF-DIAGNOSIS

When the incubator encounters an error, one of the following code numbers and alarm or notification messages will be displayed on the message display area of the control panel.
In that case, check the solution described below and also refer to the relevant section in the "TROUBLESHOOTING" for additional guidance. If the issue persists, or if an error occurs that is not covered in the description provided, stop using the product and contact our sales representative or agent.

Alarm messages

Indicates that a problem affecting cultivation and decontamination performance has occurred.

Message display area Code number: message	Type	Situation and solution	Audible alarm	Status indicator color	Remote alarm
A1: CO2 Gas Empty.	CO_2 gas cylinder empty	The CO_2 density does not increase when the CO_2 valve is opened. →See A1 section in TROUBLESHOOTING.	Intermittent tone	Red	ON
A5: Temp Sensor Open.	Chamber temp. sensor error	The chamber temperature sensor has been disconnected.			
A6: Temp Sensor Short.		The chamber temperature sensor has short-circuited.			
A7: CO2 Box Temp Sensor Open.	Temp. sensor for sensor box error	The temperature sensor for the sensor box has been disconnected.			
A8: CO2 Box Temp Sensor Short.		The temperature sensor for the sensor box has short-circuited.			
A9: AT Sensor Open.	Ambient temp. sensor error	The ambient temperature sensor has been disconnected.			
A10: AT Sensor Short.		The ambient temperature sensor has short-circuited.			
A11: CO2 Sensor Vref Abnormal.	CO_2 sensor error	The output voltage on the reference side of the CO_2 sensor is abnormal.			
A12: CO2 Sensor Vgas Abnormal.		The output voltage on the gas side of the CO_2 sensor is abnormal.			
A21: CO2 Sensor Connect Abnormal.		Communication error with CO_2 sensor has occurred.			
A13: Main Heater Abnormal.	Heater disconnection, circuit error	Main heater has been disconnected or the main heater circuit has short-circuited.			
A14: Bottom Heater Abnormal.		Bottom heater has been disconnected or the bottom heater circuit has short-circuited.			
A15: Door Heater Abnormal.		Door heater has been disconnected or the door heater circuit has short-circuited.			
A16: CO2 Box Heater Abnormal.		Overheat alarm has been activated, the sensor box heater has been disconnected, or the sensor box heater circuit has short-circuited.			
A17: Heater SSR Open.		The overheat alarm has been activated, or circuit for the main heater, bottom heater, door heater, or sensor heater has short-circuited.			
A18: UV-LED Abnormal.	UV-LED failure	UV-LED ON was detected even if the UV light was OFF.	Continuous tone		
A20: Door Lock Failure.	Door lock error	The outer door was opened when the automatic door lock function was ON. →See A20 section in TROUBLESHOOTING.	Intermittent tone		
A39: Fan Motor Error.	Fan error	Rotational speed of fan motor is abnormal.			
A65: Peltier Fan Error.	Peltier fan error	Rotational speed of peltier fan motor is abnormal.			
A56: Communication Failure.	Communication error*1	Communication between the control panel and control circuit board has been disconnected or unstable.			-----
A74: UV-LED Circuit Abnormal.	UV-LED circuit board error	The electric current value has exceeded or fallen below the defined value, or the internal heater has short-circuited.	Intermittent tone*2		ON*2
A75: UV-LED Short Circuit.		The UV-LED circuit board has short-circuited.			
A76: UV-LED Open Circuit.		The UV-LED circuit board has been disconnected.			
A77: UV-LED Parameter Error.		The drive parameter for the installed UV-LED is not supported by this product.	Intermittent tone		ON
A81: Analog Output Error.	Analog output error	Analog output terminals are not connected. →See A81 section in TROUBLESHOOTING.			
A82: Analog Output System Error.		Analog output circuit error has occurred.			

DOMINIQUE DUTSCHER SAS

ALARMS AND SELF-DIAGNOSIS

Message display area Code number: message	Type	Situation and solution	Audible alarm	Status indicator color	Remote alarm
A31: H2O2 Unit Not Connected.	System check error when starting H_2O_2 decontamination	The H_2O_2 generator is not connected. →See A31 section in TROUBLESHOOTING.	Intermittent tone (3 times)	Red	-----
A32: Low H2O2 Level.		The H_2O_2 reagent has not been put in the H_2O_2 generator, or the H_2O_2 level sensor has failed (or has been disconnected). →See A32 section in TROUBLESHOOTING.			
A33: Outer Door Open.		The outer door is not closed. →See A33 section in TROUBLESHOOTING.			
A34: H2O2 Volume.	Failure during H_2O_2 decontamination	The H_2O_2 reagent is left in the H_2O_2 generator, thus H_2O_2 has not been fully generated. →See A34 section in TROUBLESHOOTING.	Intermittent tone		ON
A35: Power Failure.		There was a power failure during the H_2O_2 decontamination procedure. →See A35 section in TROUBLESHOOTING.			
A38: Door Lock Failure.	Failure during H_2O_2 gas resolve	The outer door was opened. →See A38 section in TROUBLESHOOTING.	Continuous tone (Linked with opening of outer door)		ON (Linked with opening of outer door)
A97: Over Heat.	Overheat alarm*3	Temperature in the chamber has exceeded the value set for the high temperature limit for overheat alarm. →See A97 section in TROUBLESHOOTING.	Continuous tone		ON
A89: High Temp. or A90: Low Temp.	Temperature alarm*4	Temperature in the chamber has become higher or lower than the set temperature set point plus or minus the value (±1.0°C to ±5.0°C) set for the temperature alarm. →See A89 or A90 section in TROUBLESHOOTING.	Intermittent tone after preset alarm delay time (0 to 15 min)		Activates after preset alarm delay time (0 to 15 min)
A91: High CO2 Density. or A92: Low CO2 Density.	CO_2 density alarm*4	CO_2 density in the chamber has become higher or lower than the set CO_2 density plus or minus the value (±0.5% to ±5.0%) set for the CO_2 density alarm. →See A91 or A92 section in TROUBLESHOOTING.			
A106: Please Calibrate UV-LED.	UV-LED not calibrated	UV-LED is not calibrated.	Intermittent tone		ON

*1: When a communication error occurs, the control panel becomes unavailable except for the Event Log button.
*2: When the UV Light Mode is set to “OFF,” the audible alarm and remote alarm do not activate.
*3: A16: CO2 Box Heater Abnormal and A17: Heater SSR Open are generated when a certain period of time has elapsed after the overheat alarm activated.
*4: When the fan motor speed becomes slow due to malfunction or the end of its life, these alarms may activate because of uneven distribution of temperature, CO_2 density, or humidity in the chamber.

Notification messages

Indicates status of the incubator.

Message display area Code number: message	Type	Situation and solution	Audible alarm	Status indicator color	Remote alarm
N1: Gas Sensor initializing.	Gas control preparation	Waiting for the temperature to be stable for controlling the gas after the power to the incubator has been turned on, or after the chamber temperature set point has been changed. →See N1 section in TROUBLESHOOTING.	-----	Orange	-----
N3: USB Over Current.	USB overcurrent	Overcurrent was applied to the USB connector. →See N3 section in TROUBLESHOOTING.			
N4: Outer Door Open. “Door: Open” blinks in orange.	Door open	The outer door is open. →See N4 section in TROUBLESHOOTING.	Intermittent tone after preset alarm delay time (1 to 30 min)		-----*

* The remote door status for an external alarm device is triggered by a signal from the Door Status terminals. It is not linked to the N4 notification issued by the incubator software.

DOMINIQUE DUTSCHER SAS

TROUBLESHOOTING

If the incubator does not seem to be working properly, check the following solutions before calling for service.

Problems related to Alarm messages

Code number	Cause	Solution
A 1	CO_2 gas supply tube is disconnected or there is a gas leak.	Check the connection between the CO_2 gas cylinder and the product (refer to pages 26-28).
	CO_2 gas cylinder is empty.	Replace the CO_2 gas cylinder (when the primary pressure is 3.8 MPa(G) (38 kgf/cm²(G)) or lower, it indicates low gas remaining, so replace it promptly).
A20	The outer door, which is supposed to be locked with an electric lock, is open for some reason, such as being manually unlocked with the physical key that came with the product.	▪Close the outer door and check that the control panel displays "Door: Closed." ▪Normally, the outer door is locked with the electric lock, but it cannot be locked if the provided physical key is inserted into the keyhole. Ensure that the physical key is not inserted into the keyhole.
A31	The optional H_2O_2 generator (MCO-50HP) is not installed.	Optional H_2O_2 generator is required for H_2O_2 decontamination. For purchase, contact our sales representative or agent.
	The cable of the H_2O_2 generator is not connected properly.	Ensure that the waterproof cable of the H_2O_2 generator is correctly connected to the main unit and the H_2O_2 decontamination connector terminal.
A32	The amount of H_2O_2 decontamination reagent (MCO-20H2O2) is insufficient.	Use the entire contents of one bottle of the H_2O_2 decontamination reagent for each decontamination operation.
A33	The outer door is not closed.	Make sure to securely close the outer door. During H_2O_2 decontamination, the outer door is automatically locked with the electric lock.
A34	The amount of H_2O_2 decontamination reagent is insufficient.	Use the entire contents of one bottle of the H_2O_2 decontamination reagent for each decontamination operation.
	The amount of H_2O_2 decontamination reagent is excessive. Therefore, the gas injection did not finish within the allocated time.	Use the dedicated H_2O_2 decontamination reagent and use entire contents of one bottle of the H_2O_2 decontamination reagent for each decontamination operation.
	The duct is not installed correctly.	Check that the 4 pins at the rear side of the chamber securely fit into the 4 holes of the duct (see page 24).
	The fan is not attached properly.	Check that the fan is fully pushed in the projection of the motor shaft (see page 24).
	The H_2O_2 generator is not installed properly.	Check that the positioning fittings of the H_2O_2 generator is set into the slits of the duct (page 80).
	The H_2O_2 generator has exceeded its cumulative usage limit.	Replace the H_2O_2 generator after 200 uses.
A35	A power failure occurred during H_2O_2 decontamination.	Decontamination may not have been completed correctly in steps other than step 6. In such a case, restart the decontamination process from the beginning.
A38	The outer door was opened during H_2O_2 decontamination.	For safety, decontamination is forcibly terminated. Restart decontamination from the beginning.
A81	No measuring instrument is connected to the analog output terminal.	▪Check that a measuring instrument is connected. ▪When you are not using the analog output, set it to OFF.
A89 or A90	The chamber temperature does not match the set value immediately after turning on the power switch.	Run the incubator for about 4 hours until it stabilizes.
	The temperature set point has been changed.	The difference between the current chamber temperature and the changed chamber set point exceeded the value set for the temperature alarm. In this case, the alarm status will clear automatically when the chamber temperature reaches the set value.
	The door has been open for a long time, or high or low temperature contents were placed in the chamber.	The alarm is caused by a temporary change in chamber temperature. In this case, the alarm status will clear automatically if left as is.
	The outer door is open.	Close the outer door and check that "Door: Closed" is displayed on the control panel.
	The incubator has resumed operation after a power failure.	Due to temperature/CO_2 concentration changes during the power failure, alarms may be triggered when resuming due to differences from the set values. In this case, the alarm status will clear automatically when the actual value reaches the set value.
A91 or A92	The CO_2 gas density in the chamber is not at the set point.	Run the incubator for about 4 hours until it stabilizes.
	The CO_2 gas density set point has been changed.	Check the changed CO_2 gas density set point and the CO_2 gas density alarm value that is set.
	The outer door is open.	Close the outer door and check that "Door: Closed" is displayed on the control panel.
	The incubator has resumed operation after a power failure.	Due to temperature/CO_2 density changes during the power failure, alarms may be triggered when resuming due to differences from the set values. In this case, the alarm status will clear automatically if left as is.
A97	The overheat alarm temperature is not set at least 1°C higher than the chamber temperature set point.	Set the overheat alarm temperature at least 1°C higher than the chamber temperature set point.

DOMINIQUE DUTSCHER SAS

TROUBLESHOOTING

Problems related to Notification messages

Code number	Cause	Solution
N1	CO_2 gas will not be injected for about one hour immediately after turning on the power switch of this product or after recovering from a power failure, until the temperature of the CO_2 sensor stabilizes sufficiently.	Run the incubator for more than 1 hour until it stabilizes.
N3	Overcurrent was applied to the USB connector.	Check that no devices, such as a computer, are connected via the USB cable.
N4	The outer door is open.	Close the outer door and check that "Door: Closed" is displayed on the control panel.

Other problems

Problem	Cause	Solution
The incubator does not operate at all.	The incubator is not connected to the power supply properly.	Securely insert the power plug into the outlet.
	The main breaker has tripped.	Turn on the main breaker.
	A power failure has occurred.	Check the cause of the power failure and wait for the power to be restored. If you are unsure when the power will be restored, unplug the incubator and wait.
	The plug of the removable power supply cord is not connected to the power supply port on the lower right side of the incubator.	Securely insert the plug into the power inlet on the incubator.
An alarm is activated.*	An electrical device that generates electromagnetic waves is being used nearby.	Turn off the device or change its location to avoid electromagnetic interference.
The chamber temperature does not match the set value.	The ambient temperature is not at least 5°C lower than the temperature set for the chamber.	Set the chamber temperature at least 5°C above ambient temperature.
	The outer door is closed with the inner door left open.	Close the inner door before closing the outer door.
	An electrical device that generates electromagnetic waves is being used nearby.	Turn off the device or change its location to avoid electromagnetic interference.
The CO_2 gas density does not match the set value.	The secondary pressure for the gas regulator of CO_2 gas cylinder is not at the specified value.	Set the secondary pressure for the gas regulator of CO_2 gas cylinder to 0.03 MPa(G) - 0.1 MPa(G) (0.3 kgf/cm^2(G) to 1 kgf/cm^2(G), 4.4 psi(G) to 14.5 psi(G)).
	The gas tube is blocked.	Remove any foreign objects causing the blockage or replace the gas tube.
	The duct is not attached securely.	Check that the 4 pins at the rear side of the chamber securely fit into the 4 holes of the duct (page 24).
	The fan is not attached properly.	Check that the fan is fully pushed in the projection of the motor shaft (see page 24).
	The output from the CO_2 sensor has fluctuated relative to its initial value at the start of the incubator operation.	Contact our sales representative or agent.
	An electrical device that generates electromagnetic waves is being used nearby.	Turn off the device or change its location to avoid electromagnetic interference.
A large quantity of CO_2 gas is being consumed.	The outer and inner doors are frequently opened and closed.	Open and close the doors swiftly, and do not open the doors frequently. Opening and closing the doors let the outer air into the incubator, causing the CO_2 density change in the chamber.
	The gas tube may be damaged or deteriorated.	Check for gas leaks due to deterioration at the connection points and for any leaks from pinholes. The gas tube is a consumable part, and replacing it once a year is recommended.
	The door gasket on the inner door is coming off.	There is a risk of air leaking from the chamber. Ensure that the door gasket is properly in place and not slack.
	The door gasket for the inner door is defective.	Replace the door gasket. A damaged door gasket may cause air leakage from the chamber.
	The access port is open.	Cover the access port with the silicone caps on both the rear side and the chamber side.
	The output from the CO_2 sensor has fluctuated relative to its initial value at the start of the incubator operation.	Contact our sales representative or agent.

* If the control panel is unresponsive while an alarm is active, try turning the incubator's power off and back on. If the problem persists, contact our sales representative or agent.

DOMINIQUE DUTSCHER SAS

Problem	Cause	Solution
Cannot culture properly and the CO_2 gas density in the chamber may be the cause.	The air environment around the incubator is abnormal. There is a source of polluted gas in the vicinity.	Use the incubator in an appropriate air environment.
	The output from the CO_2 sensor has fluctuated relative to its initial value at the start of the incubator operation.	Contact our sales representative or agent.
CO_2 gas is not being injected.	This incubator adopts ON-OFF method for controlling CO_2. The CO_2 gas is intermittently injected as the gas density in the chamber approaches the set value. Injections may stop for approximately 15 seconds, but that is not an error.	-
	-	When the code number N1 is displayed on the message display area, see the description in “Problems related to Notification messages”
The CO_2 gas density is slow to recover.	A dust filter is used for the incubator gas piping. If CO_2 gas density is slow to recover when the gas pressure is normal, the dust filter may be clogged.	Contact our sales representative or agent.
	There is little gas remaining in the CO_2 gas cylinder.	A primary pressure reading of 3.8 MPa(G) (38 kgf/cm^2(G), 551 psi(G)) or lower is a sign that there is little gas remaining. Replace the cylinder soon.
	The gas tube is bent.	Install the gas tube to ensure it does not bend.
	The secondary pressure for the gas regulator of CO_2 gas cylinder is not at the specified value.	Set the secondary pressure for the gas regulator of CO_2 gas cylinder to the specified value of 0.03 MPa(G) to 0.1 MPa(G) (0.3 kgf/cm^2(G) to 1 kgf/cm^2(G), 4.4 psi(G) to 14.5 psi(G)).
	The gas tube is blocked.	Remove any foreign objects causing the blockage or replace the gas tube.
	The duct is not securely attached.	Check that the 4 pins at the rear side of the chamber securely fit into the 4 holes of the duct (page 24).
	The fan is not attached properly.	Check that the fan is fully pushed in the projection of the motor shaft (see page 24).
The outer door will not open.	The outer door is automatically locked by the electric lock when the power is not supplied to the incubator.	Turn on the incubator or use the key that came with this product to unlock the outer door.
	The outer door is automatically locked by the electric lock during H_2O_2 decontamination.	Wait until the process is complete.
Data cannot be exported to the USB flash drive.	The USB flash drive is not inserted properly. (Error message window showing “USB flash drive is not recognized.” is displayed.)	▪Check that a USB flash drive is inserted into the USB port. ▪Check that the USB flash drive is fully inserted. ▪The USB flash drive should be formatted in FAT16 or FAT32 format. Please check the format.
	Data during the specified time period does not exist. (Information message window showing “No data was found.” is displayed.)	Check the exporting setting and try again.
	The USB flash drive has no free space or very little available space. (Error message window showing “There is not enough space on the USB flash drive.” is displayed.)	Check that the USB has enough free space for saving the data.
	The export failed. (Error message window showing “An Error has occurred. Please retry.” is displayed.)	The export failed for some reason. Please check that the USB flash drive is fully inserted and try exporting again.
	The USB flash drive that requires password is used.	Use a USB flash drive that has no password function.
	The USB flash drive with capacity of more than 32 GB is used.	Use a USB flash drive with capacity of 32 GB or less.
H_2O_2 decontamination cannot be performed.	The outer door was opened during the decontamination process.	For safety reasons, the decontamination process will be forcibly terminated. If it happens at any step other than STEP 6, restart the decontamination process from the beginning. Do not open the outer door during decontamination. If the outer door is accidentally opened during STEP 4 or 5, promptly close it as H_2O_2 gas is being sprayed inside the chamber.
	-	When one of the following code numbers is displayed in the message display area, see the relevant description in the “Problems related to Alarm messages” section. A31, A32, A33, A34, A35, A38

Note:
If the problem is still not solved after trying the above solutions, or for any problems not covered here, contact our sales representative or agent.

DOMINIQUE DUTSCHER SAS

SPECIFICATIONS

Product name	CO_2 Incubator MCO-173AIC	CO_2 Incubator MCO-173AICUV	CO_2 Incubator MCO-173AICUVH
External dimensions	W620 mm x D730 mm x H905 mm (W24.4 in. x D28.7 in. x H35.6 in.)		
Internal dimensions	W490 mm x D523 mm x H665 mm (W19.3 in. x D20.6 in. x H26.2 in.)		
Interior volume	165 L (5.8 cu.ft.)		
Exterior	Painted steel (Rear cover has no paint)		
Interior	Stainless steel containing copper		
Outer door	Painted steel		
Inner door	Tempered glass		
Trays	4 trays made of stainless steel containing copper W470 mm x D450 mm x H12 mm (W18.5 in. x D17.7 in. x H0.47 in.) Maximum load: 7 kg (15 lbs.)/tray		
Access port	Inner diameter: 30 mm (1.18 in.), On the back side		
Insulation	Expandable styrene-acrylonitrile beads		
Heating system	DHA system (heater jacket + air jacket system)		
Heater	300 W		
Humidifying system	Natural evaporation with humidifying pan		
Temperature controller	PID control system		
Temperature display	Digital display		
CO_2 controller	PID control system		
CO_2 density display	Digital display		
Air circulation	Fan assisted		
Gas line filter	0.3 µm, Efficiency: 99.97 % or higher		
UV LED	-	UVC-LED x 1 (ozone-free emission)	
Alarms	Temperature alarm, CO_2 density alarm, Over heat alarm, CO_2 gas, various sensor/heater alarms		
Remote alarm contacts	Allowable contact capacity: DC 30 V, 2 A *		
CO_2 inlet connection	Soft Polyurethane tube can be connected (ID 4 mm, OD 6 mm (ID 0.157 in. ,OD 0.236 in.))		
CO_2 inlet pressure	0.03 MPa(G) - 0.1 MPa(G) (0.3 kgf/cm^2(G) - 1 kgf/cm^2(G), 4.4 psi(G) - 14.5 psi(G))		
Weight	69 kg (152 lbs.)		
Accessories	1 removable power supply cord for UK, 1 removable power supply cord for EU countries, 1 power cord cover plate, 4 trays 1 gas tube, 1 humidifying pan, 1 conversion joint		2 keys to outer door

*A maximum length of 30 meters is recommended for the cable used for the connection.

DOMINIQUE DUTSCHER SAS

SPECIFICATIONS

Product name	CO_2 Incubator		
	MCO-173AIC	MCO-173AICUV	MCO-173AICUVH
Optional accessories required for H_2O_2 decontamination			
H_2O_2 generator (MCO-50HP)	✓	✓	✓
H_2O_2 reagent (MCO-20H2O2)	✓	✓	✓
H_2O_2 decon kit (MCO-HK)	✓	✓	Standard
Electric lock (MCO-170EL)*1	✓	✓	Standard
Optional accessories required for stacking incubator (see Table 1 below)			
Double stacking bracket (MCO-170PS)	✓	✓	✓
Stacking plate (MCO-170SB/MCO-230SB)	✓	✓	✓
Optional accessories			
UV system kit (MCO-LUVS)	✓	Standard	Standard
Gas regulator (MCO-010R)	✓	✓	✓
Gas auto changer (MCO-21GCP)	✓	✓	✓
Reinforced tray (MCO-170RT)	✓	✓	✓
Tray (MCO-170ST)	✓	✓	✓
Half tray (MCO-25ST)	✓	✓	✓
Small door kit (MCO-173ID)	✓	✓	✓
Roller base (MCO-170RB)	✓	✓	✓

*1: Electric lock can also be used for the auto lock function.

Note:
Refer to the updated catalog when ordering an optional component.

Table 1 Combination of incubators and required stacking bracket/plate (other combinations are not recommended).

		Upper incubator		
		MCO-50 series	MCO-170/171/173 series	MCO-230/233 series
Lower incubator	MCO-50 series	MCO-170PS	—	—
	MCO-170/171/173 series	—	MCO-170PS or MCO-170SB	—
	MCO-230/233 series	—	MCO-230SB	MCO-170PS or MCO-230SB

DOMINIQUE DUTSCHER SAS

PERFORMANCE

Product name	CO_2 Incubator MCO-173AIC	CO_2 Incubator MCO-173AICUV	CO_2 Incubator MCO-173AICUVH
Model number	MCO-173AIC-PE	MCO-173AICUV-PE	MCO-173AICUVH-PE
Temperature control range	Ambient temperature plus 5°C to max. 50°C* (ambient temperature: 5°C to 35°C)		
Temperature distribution	±0.25°C (ambient temperature: 23°C, setting: 37°C, CO_2: 5%, no load)		
Temperature variation	±0.1°C (ambient temperature: 23°C, setting: 37°C, CO_2: 5%, no load)		
CO_2 setting range	0% to 20%		
CO_2 variation	±0.15% (ambient temperature: 23°C, setting: 37°C, CO_2: 5%, no load)		
Chamber humidity	95% RH ± 5% RH		
Environmental conditions	Ambient temperature: 5°C to 35°C, Humidity: max. 80% RH (The designed performance may not be obtained if the ambient temperature is 15°C or lower)		
Noise level	25 dB (A scale)		
Power consumption	Max. 350 W		
Heat emission	Max. 1,180 kJ/h		
Rated voltage, frequency	AC 220-240 V, 50/60 Hz		
Amperage	Max. 1.5 A		

*When set temperature is 37°C, ambient temperature must be 32°C or lower. Regardless of ambient temperature, the maximum of temperature control range is always 50°C.

Notes:

- Based on our measuring method.
- Default calibration conditions: 37°C, CO_2: 5%
 When using under other conditions, we recommend performing calibration under the conditions of use.
- We recommend performing calibration every year.

DOMINIQUE DUTSCHER SAS

SAFETY ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):

■ Indoor use;

■ Altitude up to 2,000 m;

■ Temperature 5˚C to 40˚C

■ Maximum relative humidity 80% for temperature up to 31˚C decreasing linearly to 50% relative humidity at 40˚C;

■ Mains supply voltage fluctuations up to ±10% of the nominal voltage;

■ Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;

■ Temporary OVERVOLTAGES occurring on the mains supply;

■ Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases);

* Above conditions do not indicate the performance of this product. For the performance of this product, refer to the "PERFORMANCE" section.

DOMINIQUE DUTSCHER SAS

EMC SPECIFICATION

This product is suitable for use in specialized healthcare facility environment such as hospitals and medical research institutions. This product has been confirmed as complying with IEC 60601-1-2 in such environment.

Electromagnetic Emission	IEC 60601-1-2
Electromagnetic Immunity	IEC 60601-1-2

1. The compliance with the requirements for the Electromagnetic Emissions

Emissions test	Basic EMI standard	Compliance LEVELS
Conducted emissions	CISPR 11	Group 1, Class A
Radiated emissions		
Harmonic current emissions	IEC 61000-3-2	Class A
Voltage changes, voltage fluctuations and flicker emissions	IEC 61000-3-3	—

Note:
The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

DOMINIQUE DUTSCHER SAS

EMC SPECIFICATION

2. The compliance with the general requirements for Electromagnetic Immunity

Phenomenon	Basic EMS standard	Compliance LEVELS	
Electrostatic Discharge	IEC 61000-4-2	Enclosure SIP/SOP	±8 kV contact ±15 kV air
Radiated RF EM field	IEC 61000-4-3	Enclosure	3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz
Electrical fast transients / bursts	IEC 61000-4-4	Input AC	± 2 kV 5/50 ns, 100 kHz
		SIP/SOP	± 1 kV 5/50 ns, 100 kHz
Surges	IEC 61000-4-5	Input AC	± 1 kV Line-to-line ± 2 kV Line-to-ground
Conducted disturbances induced by RF fields	IEC 61000-4-6	Input AC SIP/SOP	3 V (0.15 MHz to 80 MHz) 6 V (in ISM bands between 0.15 MHz to 80 MHz) 80 % AM at 1 kHz
Rated power frequency magnetic fields	IEC 61000-4-8	Enclosure	30 A/m, 50 Hz/60 Hz
Voltage dips	IEC 61000-4-11	Input AC	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1cycle and 70 % UT; 25/30 cycle Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	Input AC	0 % UT; 250/300 cycle
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Enclosure	See next page

Note:
Keep this product at least 15 cm away from the AC power cables of other devices. Otherwise, the performance of this product may degrade.

Do not use this product adjacent to the following systems. Otherwise, the performance of this product may degrade.

– electronic article surveillance (EAS) systems,
– access control systems (radio-frequency identification – RFID),
– magnetic levitated rail systems,
– radio-frequency (RF) medical procedures such as high frequency surgery and ablation therapy,
– metal detectors,
– use of transponders for traffic control,
– wireless charging/inductive power transfer

DOMINIQUE DUTSCHER SAS

3. The compliance with the requirements for IMMUNITY from RF wireless communications equipment

Test frequency [MHz]	Band [MHz]	Service	Modulation	Maximum power [W]	Distance [m]	Immunity test level [V/m]
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704 - 787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth; WLAN 802.11 b/g/n; RFID 2450; LTE band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

Note:
Keep portable RF communication devices (including peripheral devices such as antenna cables and external antennas) at least 30 cm (12 inches) away from this product. Otherwise, the performance of this product may degrade.

4. Additional information

This product also complies with IEC 61326-1 in basic electromagnetic environment (including the environment of specialized healthcare facilities) and intended to be used in such environment.

DOMINIQUE DUTSCHER SAS

SAFETY CHECK SHEET

⚠ CAUTION	**Please copy and fill out this form before servicing. Hand over the form to the service engineer for their and your safety.**

Safety check sheet

1. Stored material

Risk of infection:	☐Yes	☐No	☐Maybe
Risk of toxicity:	☐Yes	☐No	☐Maybe
Risk from radioactive sources:	☐Yes	☐No	☐Maybe

List all potentially hazardous materials that have been stored in this unit:

2. Contamination in the unit

a) Contamination ☐Yes ☐No ☐Maybe

Types of contamination (if any):_____

b) Decontaminated ☐Yes ☐No

Methods used for the decontamination work:

3. Status of the unit

a) The unit is now safe to work on ☐Yes ☐No

b) If the answer is "No,"

Details on the danger:_____

Measures we should take to reduce the danger:

Date:

Signature:

Address, Division:

Telephone:

Product name: CO_2 incubator	Model No. MCO-	Serial number:	Date of Installation:

Please decontaminate the unit yourself before calling the service engineer.

DOMINIQUE DUTSCHER SAS

MEMO

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DISPOSAL OF UNIT

Decontamination of unit

Before disposal of this CO_2 incubator, contact our sales representative or agent for further information.
Improper handling of biohazardous waste can result in accidental exposure to infectious agents. If there is a danger of biohazard, decontaminate the incubator as thoroughly as possible before disposal.

Disposal of Old Equipment and Batteries
Only for European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

Note for the battery symbol (bottom symbol):
This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

Recycle of battery

■ This label is a recycle mark complying with the Taiwanese battery regulation.

Original Operating Instructions

■ Medical Device Regulation ((EU) 2017/745) CE 0123

■ Low Voltage Directive (2014/35/EU)
EMC Directive (2014/30/EU)
RoHS Directive (2011/65/EU) CE

■ Electrical Equipment (Safety) Regulations 2016 (LV)
The Electromagnetic Compatibility Regulations 2016 (EMC) UKCA

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„Time Zone“ (Zeitzone)	
UTC -12:00: Internationale Datumslinie West	UTC +04:00: Golf, Samara
UTC -11:00: Samoa	UTC +04:30: Kabul
UTC -10:00: Hawaii-Aleuten	UTC +05:00: Pakistan, Yekaterinburg
UTC -09:30: Marquesas Inseln	UTC +05:30: Neu-Delhi, Mumbai, Colombo
UTC -09:00: Alaska	UTC +05:45: Katmandu
UTC -08:00: Pazifik	UTC +06:00: Bangladesch, Omsk
UTC -07:00: Mountain	UTC +06:30: Yangon
UTC -06:00: Central	UTC +07:00: Indochina, Novosibirsk
UTC -05:00: Eastern	UTC +08:00: China, Australian Western
UTC -04:30: Venezuela	UTC +08:45: Eucla
UTC -04:00: Atlantik	UTC +09:00: Japan, Korea
UTC -03:00: Brasilia, Grönland	UTC +09:30: Adelaide
UTC -02:00: Koordinierte Weltzeit-02	UTC +10:00: Australian Eastern, Vladivostok
UTC -01:00: Azores, Kap Verde	UTC +10:30: Lord Howe Insel
UTC 00:00: Koordinierte Weltzeit, GMT	UTC +11:00: Salomonen, Magadan
UTC +01:00: Mitteleuropäisch	UTC +12:00: Neuseeland, Fidschi
UTC +02:00: Osteuropäisch	UTC +12:45: Chatham Inseln
UTC +03:00: Moskau, Ostafrika	UTC +13:00: Koordinierte Weltzeit +13
UTC +03:30: Tehran	UTC +14:00: Weihnachtsinsel

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*UTC: Tempo coordinato universale

Elenco dei fusi orari

Fuso orario	
UTC -12:00: Linea internazionale del cambio di data provenendo da ovest	UTC +04:00: Golfo, Samara
UTC -11:00: Samoa	UTC +04:30: Kabul
UTC -10:00: Hawaii-Aleutine	UTC +05:00: Pakistan, Ekaterinburg
UTC -09:30: Isole Marquesas	UTC +05:30: Nuova Delhi, Mumbai, Colombo
UTC -09:00: Alaska	UTC +05:45: Kathmandu
UTC -08:00: Pacifico	UTC +06:00: Bangladesh, Omsk
UTC -07:00: Mountain	UTC +06:30: Yangon
UTC -06:00: Centrale	UTC +07:00: Indocina, Novosibirsk
UTC -05:00: Orientale	UTC +08:00: Cina, Australia occidentale
UTC -04:30: Venezuela	UTC +08:45: Eucla
UTC -04:00: Atlantico	UTC +09:00: Giappone, Corea
UTC -03:00: Brasilia, Groenlandia	UTC +09:30: Adelaide
UTC -02:00: Tempo coordinato universale-02	UTC +10:00: Australia orientale, Vladivostok
UTC -01:00: Azzorre, Capo Verde	UTC +10:30: Isola di Lord Howe
UTC 00:00: Tempo coordinato universale, GMT	UTC +11:00: Isole Salomone, Magadan
UTC +01:00: Europa centrale	UTC +12:00: Nuova Zelanda, Fiji
UTC +02:00: Europa orientale	UTC +12:45: Isola Chatam
UTC +03:00: Mosca, Africa orientale	UTC +13:00: Tempo coordinato universale +13
UTC +03:30: Teheran	UTC +14:00: Isola di Natale

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© PHC Corporation 2025