Sigma-Aldrich®

www.sigmaaldrich.com

1.01424.0100
1.01424.0103
1.01424.0500
1.01424.0503
1.01424.1000
1.01424.1022
1.01424.2500
1.01424.9025

Microscopy

May-Grünwald's eosin methylene blue solution modified

for microscopy

For professional use only

IVD *In Vitro* Diagnostic Medical Device

CE

Intended purpose

This "May-Grünwald's eosine-methylene blue solution modified - for microscopy" is used for human-medical cell diagnosis and serves the purpose of the hematological and clinico-cytological investigation of sample material of human origin. It is a staining solution that when used together with other *in vitro* diagnostic products from our portfolio makes target structures evaluable for diagnostic purposes (by fixing, embedding, staining, where necessary counterstaining, mounting) in hematological and clinico-cytological specimen materials, for example whole blood or bone-marrow smears.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further tests must be carried out according to recognized, valid methods to reach a definitive diagnosis.

Principle

When used in hematological applications, May-Grünwald's stain is frequently used in combination with other staining solutions, e. g. with Giemsa's solution for Pappenheim (MGG) overview staining. This staining solution generally stains the nuclei red to violett, based on the molecular interaction between the eosin Y dye and an Azure B-DNA complex. Both dyes assemble to an Eosin Y - Azure B-DNA complex and the intensity of the resulting stain depends on the content of Azure B and the ratio of Azure B : Eosin Y. Furthermore, the resulting stain can vary depending on the influence of fixation, staining times, pH-value of the solutions or buffer substances.

Sample material

air-dried blood or bone-marrow smears just as clinical cytological material like urine sediment, sputum, smears from fine needle aspiration biopsies (FNAB), rinses, imprints

Reagents

Cat. No. 1.01424
May-Grünwald's eosine-methylene blue solution modified for microscopy — 100 ml, 500 ml, 1 l, 2.5 l, 25 l

Also required:

for the staining of air-dried smears:

Cat. No. 1.09468 Buffer tablets pH 7.2 for preparing buffer solution acc. to WEISE for staining of blood smears — 100 tabs

or

Cat. No. 1.11373 Buffer tablets pH 6.4 for preparing buffer solution acc. to WEISE for the staining of blood smears — 100 tabs

or

Cat. No. 1.11374 Buffer tablets pH 6.8 for preparing buffer solution acc. to WEISE for the staining of blood smears — 100 tabs

for the Pappenheim staining:

Cat. No. 1.09204 Giemsa's azur eosin methylene blue solution for microscopy — 100 ml, 500 ml, 1 l, 2.5 l

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.
All samples must be clearly labeled.
Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Reagent preparation

May-Grünwald's eosine-methylene blue solution modified

The solution is supplied as a concentrated staining solution and before use must be diluted with a buffer solution as described below. The diluted staining solution should be filtered before use.

Buffer solution

For preparation of approx. 1000 ml of solution, add and dissolve:

Buffer tablet, Cat. No. 1.11373 (pH 6.4), Cat. No. 1.11374 (pH 6.8), or Cat. No. 1.09468 (pH 7.2) depending on the required reaction color	1 tablet
Distilled water	1000 ml

Dilute May-Grünwald's staining solution for manual staining

For preparation of approx. 200 ml solution mix:

May-Grünwald's eosine-methylene blue solution modified	30 ml
Buffer solution	20 ml
Distilled water	150 ml
mix and filter	

Dilute May-Grünwald's staining solution for May-Grünwald's stain with automatic stainer (pH 7.2)

For preparation of approx. 300 ml solution mix:

May-Grünwald's eosine-methylene blue solution modified	50 ml
Buffer solution pH 7.2	30 ml
Distilled water	220 ml
Mix, leave to stand for 10 min, and filter	

Dilute Giemsa's staining solution for manual staining

For preparation of approx. 200 ml solution mix:

Giemsa's azur eosin methylene blue solution	10 ml
Buffer solution	190 ml
Mix well, leave to stand for 10 min, and filter if necessary	

Dilute Giemsa's staining solution for Pappenheim's stain with automatic stainer

For preparation of approx. 300 ml solution mix:

Giemsa's azur eosin methylene blue solution	15 ml
Buffer solution pH 6.8 or 7.2	285 ml
Mix well, leave to stand for 10 min, and filter if necessary	

DOMINIQUE DUTSCHER SAS

May-Grünwald's stain

Procedure

Air-dried smears

Staining in the staining cell

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear	
May-Grünwald's eosine-methylene blue solution modified	3 min
Dilute May-Grünwald's staining solution for manual staining	6 min
Buffer solution	1 min
Buffer solution	1 min
Air-dry (e.g. over night or at 50 °C in the drying cabinet)	

Staining on the staining rack

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear			
May-Grünwald's eosine-methylene blue solution modified	1 ml	cover completely	3 min
Buffer solution	1 ml	mix	6 min
Buffer solution		rinse	
Air-dry (e.g. over night or at 50 °C in the drying cabinet)			

Staining in the automatic stainer (pH 7.2)

The stated times should be adhered to in order to guarantee an optimal staining result.

	Time	Station	Dip
Slide with air-dried smear			
May-Grünwald's eosine-methylene blue solution modified	3 min	2	on
Dilute May-Grünwald's staining solution for May-Grünwald's stain with automatic stainer (pH 7.2)	20 min	3	on
Buffer solution pH 7.2	1 min	4	on
Dry	3 min	6	-

All diluted solutions should be replaced after each working day. Only in the case of May-Grünwald's eosin methylene blue solution modified should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated May-Grünwald's eosin methylene blue solution modified (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

Note: The best staining result in the automatic stainer is achieved at pH 7.2.

Covering with non-aqueous mounting media (e.g. Neo-Mount™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

	Buffer solution pH 6.4	Buffer solution pH 6.8	Buffer solution pH 7.2
Nuclei	red-violet	red-violet	violet
Cytoplasm of lymphocytes	blue	blue	blue
Cytoplasm of monocytes	grey-blue	grey-blue	grey-blue
Neutrophilic granules	light violet	light violet	light violet
Eosinophilic granules	brick-red to red-brown	brick-red	red-brown
Basophilic granules	dark violet to black	dark violet to black	dark violet to black
Thrombocytes	violet	violet	violet
Erythrocytes	reddish	reddish	pink to brownish

Pappenheim's stain

with May-Grünwald's solution and Giemsa's solution

Procedure

Air-dried smears

Staining in the staining cell

The slides must be immersed and moved briefly in the solutions, simple immersion alone yields inadequate staining results.

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear	
May-Grünwald's eosin-methylene blue solution modified	3 min
Dilute Giemsa's staining solution for manual staining	20 min
Buffer solution	1 min
Buffer solution	1 min
Air-dry (e. g. over night or at 50 °C in the drying cabinet)	

Staining on the staining rack

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear			
May-Grünwald's eosin-methylene blue solution modified	1 ml	cover completely	3 min
Buffer solution	1 ml	mix	
Dilute Giemsa's staining solution for manual staining		cover completely	20 min
Buffer solution		rinse	
Air-dry (e. g. over night or at 50 °C in the drying cabinet)			

Staining in the automatic stainer

The stated times should be adhered to in order to guarantee an optimal staining result.

	Time	Station	Dip
Slide with air-dried smear			
May-Grünwald's eosine-methylene blue solution modified	3 min	1	on
Dilute Giemsa's staining solution for Pappenheim's stain with automatic stainer	20 min	2	on
Buffer solution	2 min	3	on
Buffer solution	2 min	4	on
Dry	3 min	6	-

All diluted solutions should be replaced after each working day. Only in the case of May-Grünwald's eosin methylene blue solution modified should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated May-Grünwald's eosin methylene blue solution modified (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

DOMINIQUE DUTSCHER SAS

Covering with non-aqueous mounting media (e.g. Neo-Mount™, Entellan™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

	Buffer solution pH 6.4	**Buffer solution pH 6.8**	**Buffer solution pH 7.2**
Nuclei	red-violet	purple to violet	violet
Cytoplasm of lymphocytes	blue	blue	blue
Cytoplasm of monocytes	grey-blue	grey-blue	grey-blue
Neutrophilic granule	light violet	light violet	violet
Eosinophilic granule	brick-red	brick-red	red-brown
Basophilic granule	dark violet	dark violet to black	dark violet to black
Thrombocytes	violet	violet	violet
Erythrocytes	reddish	reddish	reddish-grey

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.
When using automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software.
The freshly prepared staining solutions should be filtered before use.
Remove surplus immersion oil before filing.

Analytical performance characteristics

"May-Grünwald's eosine-methylene blue solution modified" stains and thereby visualizes biological structures, as described in the "Result" chapters of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100 %:

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
May-Grünwald's staining				
Nuclei	20/20	20/20	6/6	6/6
Cytoplasm of lymphocytes	20/20	20/20	6/6	6/6
Neutrophilic granule	20/20	20/20	6/6	6/6
Eosinophilic granule	20/20	20/20	6/6	6/6
Erythrocytes	20/20	20/20	6/6	6/6

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel.
Valid nomenclatures must be used.
This method can be supplementarily used in human diagnostics.
Further tests must be selected and implemented according to recognized methods.
Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Store the May-Grünwald's eosine-methylene blue solution modified - for microscopy at +15 °C to +25 °C.

Shelf-life

The May-Grünwald's eosine-methylene blue solution modified - for microscopy can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

Capacity

approx. 1500 stainings / 500 ml

Additional instructions

For professional use only.
In order to avoid errors, the application must be carried out by qualified personnel only.
National guidelines for work safety and quality assurance must be followed.
Microscopes equipped according to the standard must be used.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.
Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No.	1.00579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No.	1.00974	Ethanol denatured with about 1 % methyl ethyl ketone for analysis EMSURE®	1 l, 2.5 l
Cat. No.	1.03699	Immersion oil Type N acc. to ISO 8036 for microscopy	100-ml dropping bottle
Cat. No.	1.04699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No.	1.06009	Methanol for analysis EMSURE® ACS,ISO,Reag. Ph Eur	1 l, 2.5 l, 5 l
Cat. No.	1.07960	Entellan™ rapid mounting medium for microscopy	500 ml
Cat. No.	1.07961	Entellan™ new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No.	1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No.	1.09016	Neo-Mount™ anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml
Cat. No.	1.09204	Giemsa's azur eosin methylene blue solution for microscopy	100 ml, 500 ml, 1 l, 2.5 l
Cat. No.	1.09468	Buffer tablets pH 7.2 for preparing buffer solution acc. to WEISE for staining of blood smears	100 tabs
Cat. No.	1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l
Cat. No.	1.11373	Buffer tablets pH 6.4 for preparing buffer solution acc. to WEISE for the staining of blood smears	100 tabs
Cat. No.	1.11374	Buffer tablets pH 6.8 for preparing buffer solution acc. to WEISE for the staining of blood smears	100 tabs

Hazard classification

Cat. No. 1.01424

Please observe the hazard classification printed on the label and the information given in the safety data sheet.
The safety data sheet is available on the website and on request.

Main components of the product

Cat. No. 1.01424

C.I. 52015	0.7 g/l
C.I. 45380	0.5 g/l

contains CH_3OH
1 l = 0.79 kg

DOMINIQUE DUTSCHER SAS

Other IVD products

Cat. No.	1.00869	Entellan™ new for cover slipper for microscopy	500 ml
Cat. No.	1.01383	Wright's eosin methylene blue solution for microscopy	100 ml, 500 ml, 2.5 l
Cat. No.	1.02439	Eosin Y-solution 0.5%, alcoholic for microscopy	500 ml, 2.5 l
Cat. No.	1.05174	Hematoxylin solution modified acc. to Gill III for microscopy	500 ml, 1 l, 2.5 l
Cat. No.	1.05387	Leishman's eosin methylene blue solution modified for microscopy	500 ml
Cat. No.	1.09844	Eosin Y-solution 0.5% aqueous for microscopy	1 l, 2.5 l
Cat. No.	1.11661	Hemacolor® Rapid staining of blood smear staining set for microscopy	1 set

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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H225: Highly flammable liquid and vapor.

H301 + H311 + H331: Toxic if swallowed, in contact with skin or if inhaled.

H370: Causes damage to organs (Eyes, Central nervous system).

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

DOMINIQUE DUTSCHER SAS

Consult instructions for use

Manufacturer

Catalog number

Batch code

Caution, consult accompanying documents

Use by YYYY-MM-DD

Temperature limitation

Status: 2023-Mar-06

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

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