Silica Test Kit

SI-5 (1455400)

DOC326.98.00041

Test preparation


CAUTION: ⚠ Review the Safety Data Sheets (MSDS/SDS) for the chemicals that are used. Use the recommended personal protective equipment.

- Put the color disc on the center pin in the color comparator box (numbers to the front).
- Use sunlight or a lamp as a light source to find the color match with the color comparator box.
- Rinse the tubes with sample before the test. Rinse the tubes with deionized water after the test.
- If the color match is between two segments, use the value that is in the middle of the two segments.
- If the color disc becomes wet internally, pull apart the flat plastic sides to open the color disc. Remove the thin inner disc. Dry all parts with a soft cloth. Assemble when fully dry.
- Undissolved reagent does not have an effect on test accuracy.
- To verify the test accuracy, use a standard solution as the sample.
- The sample temperature must be 15–25 °C (59–77 °F) for accurate results.
- Phosphate levels of more than 50 mg/L PO_4^{3-} interfere with the 0–30 mg/L test procedure. If 60 mg/L PO_4^{3-} is in the sample, the error is –2% silica . If 75 mg/L PO_4^{3-} is in the sample, the error is –11% silica . Use the 0–600 mg/L test procedure if high levels of phosphate are in the sample.

Replacement items

Description	Unit	Item no.
Acid Reagent Powder Pillows, 5 mL	100/pkg	1454599
Citric Acid Powder Pillows, 5 mL	100/pkg	1454999
Molybdate Reagent Powder Pillows, 5 mL	100/pkg	1454699
Silica 3 Reagent Powder Pillows	100/pkg	27169
Bottle, square, 29 mL, with 10, 15, 20 and 23-mL marks	6/pkg	232706
Color disc, silica, 0–30 mg/L	each	9262000
Color comparator box	each	173200
Dropper, glass, 0.5- and 1.0-mL marks	5/pkg	1419705
Glass viewing tubes, glass, 18 mm	6/pkg	173006
Stoppers for 18-mm glass tubes and AccuVac Ampuls	6/pkg	173106
Water, deionized	100 mL	27242

Test procedure—Silica (0–30 mg/L SiO_2)


1. Fill two tubes to the first line (5 mL) with sample.

2. Add one Acid Reagent Powder Pillow and one Molybdate Reagent Powder Pillow to one tube. Swirl to mix.


3. Wait 10 minutes. A yellow color develops if silica or phosphate is in the sample.

4. Add one Citric Acid Powder Pillow to the same tube. Swirl to mix.


5. Wait 2 minutes. Interference from phosphate is removed.


6. Add one Silica 3 Reagent Powder Pillow to the same tube. Swirl to mix.


7. Wait 5 minutes. A blue color develops if silica is in the sample.

8. Put the blank tube into the left opening of the color comparator box. Put the prepared sample tube into the color comparator box.

9. Hold the color comparator box in front of a light source. Turn the color disc to find the color match.

10. Read the result in mg/L in the scale window.

DOWNLOADED FROM CHEMSAVE


Test procedure—Silica (0–600 mg/L SiO_2)


1. Use the dropper to add 1 mL of sample to the bottle.


2. Fill the bottle to the 20-mL mark with deionized water. Swirl to mix.

3. Fill two tubes to the first line (5 mL) with the diluted sample.


4. Put one tube into the left opening of the color comparator box.

5. Add one Acid Reagent Powder Pillow and one Molybdate Reagent Powder Pillow to the second tube. Swirl to mix.


6. Wait 10 minutes. A yellow color develops if silica or phosphate is in the sample.

7. Add one Citric Acid Powder Pillow to the same tube. Swirl to mix.


8. Wait 2 minutes. Interference from phosphate is removed.


9. Add one Silica 3 Reagent Powder Pillow to the same tube. Swirl to mix.


10. Wait 5 minutes. A blue color develops.

11. Put the second tube into the color comparator box.

12. Hold the color comparator box in front of a light source. Turn the color disc to find the color match.

13. Read the value in the scale window.

14. Multiply the value by 20 to get the result in mg/L.

DOMINIQUE DUTSCHER SA

© Hach Company/Hach Lange GmbH, 2016. All rights reserved. Printed in U.S.A. www.hach.com 800-227-4224 techhelp@hach.com

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