TPP Techno Plastic Products AG
Zollstrasse 7, CH-8219 Trasadingen, Switzerland
Tel. +41 (0)52 687 01 87, Fax +41 (0)52 687 01 77
info@tpp.ch, www.tpp.ch

Gebrauchsanleitung TPP - Zellkulturtestplatten / *Instruction for Use TPP Tissue Culture Test Plate*

Anleitung

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Instruction

TPP tissue culture plates are designed for cell culture and for manual handling.

- *Take product packaging and look for intact bag and product.*
- *Open sleeve in sterile area remove a complete system (lid and bottom) for usage.*
- *Take off lid.*
- *Introduction of medium: Use a pipette to dispense the medium over the sidewalls into the dish, or directly onto the cell growth surface (base).*
- *Cover wells with specially designed lid separates each well, thus preventing liquid condensation or cross-contamination.*
- *Identification: clear alpha-numerical labeling between the wells for easy orientation during microscopy, easy to see black alpha-numeric on the side, yellow writing area*
- *Incubation: growth area on spherical zone only, lid with air-venting system guarantees controlled gas exchange with low evaporation.*
- *Cell removal:*
 - *Adherent cells can be harvested using trypsinisation or a cell scraper.*
 - *Cells in suspension can be removed using normal medium change procedures such as aspiration.*

Advice and Recommendation

- *Test plates with 96 wells:*
 - *Manufactured according to ANSI/SBS-3d-standards.*
 - *Standardized with an interleaved barcode 2/5.*
 - *2 well geometries available: F-base and U-base.*
- *A controlled level base surface and regular use of feeding medium are essential to uniform cell growth.*
- *The necessary gas-exchange is by means of the lid where spacers on the inside of the dish lid supply a constant and optimal movement of air, which guarantees optimum aeration of the culture.*
- *ELISA: The TPP plates feature a low binding capacity and are therefore not suited for the use in ELISA tests.*

DOWNLOADED BY BIOSCHER DIAGNOSTICS

TPP Techno Plastic Products AG
Zollstrasse 7, CH-8219 Trasadingen, Switzerland
Tel. +41 (0)52 687 01 87, Fax +41 (0)52 687 01 77
info@tpp.ch, www.tpp.ch

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DOMINIQUE DUTSCHER SAS