ITALIANO

CONTAM SWAB - LISTERIA

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ENGLISH

CONTAM SWAB - LISTERIA

Detection of *Listeria* spp. directly from surfaces.

DESCRIPTION

CONTAM SWAB - LISTERIA is a swab containing a culture medium used for the detection of *Listeria* spp. directly from surfaces.

PACKAGE CONTENTS

The product consists of an envelope made of transparent paper/plastic bag containing a culture medium and a sterile swab.
Each package contains:

- 30 sealed envelopes (swab + culture medium).
- 1 instructions sheet.

TEST PROCEDURE

1. Take the sterile swab contained in the envelope.
2. Moisten the tip of the swab by immersing it in a tube containing sterile physiological solution.
3. Press the tip of the swab against the wall of the tube to remove excess liquid.
4. Streak the swab on a surface (10cm x 10cm) horizontally and vertically.
5. Insert the swab in the tube containing the culture medium.
6. Cap the tube and record date and the sampling point.
7. Incubate at 37°C for 18-24 hours.

RESULTS INTERPRETATION

Watch for the color change of the medium and interpret the results as shown in the table.

Amber	Black
	
Absence of *Listeria* spp.	Presence of *Listeria* spp.

QUALITY CONTROL

Each batch of CONTAM SWAB - LISTERIA is submitted to quality control using the following bacterial:
Escherichia coli ATCC 25922, *Streptococcus faecalis* ATCC 19433, *Listeria monocytogenes* CIP 7836, *Listeria innocua* CIP 8011.

STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

CONTAM SWAB - LISTERIA does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

PRESENTATION

Product	REF	Σ
CONTAM SWAB - LISTERIA	86102	30 Tests

REFERENCES

- ISO 18593 - Microbiology of food and animal feeding stuffs - Horizontal method for sampling techniques from surfaces using contact plates and swabs.
- ISO 17604 - Microbiology of food and animal feeding stuffs - Carcass sampling for microbiological analysis.

TABLE OF SYMBOLS

LOT	Batch code		Do not reuse		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation	Σ	Contains sufficient for <n> tests		Caution,consult accompanying documents		

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DOMINIQUE DUTSCHER SAS