ITALIANO

CONTAM SWAB - MRSA

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ENGLISH

CONTAM SWAB - MRSA

Detection of Methicillin-resistant *Staphylococcus aureus* directly from surfaces.

DESCRIPTION

CONTAM SWAB - MRSA is a swab containing a culture medium used for the detection of methicillin/oxacillin resistant *Staphylococcus aureus* directly from surfaces.

PACKAGE CONTENTS

The product consists of an envelope made of transparent paper/plastic bag containing a culture medium and a sterile swab.
Each package contains:

- 30 sealed envelopes (swab + culture medium).
- 1 instructions sheet.

TEST PROCEDURE

1. Take the sterile swab contained in the envelope.
2. Moisten the tip of the swab by immersing it in a tube containing sterile physiological solution.
3. Press the tip of the swab against the wall of the tube to remove excess liquid.
4. Streak the swab on a surface horizontally and vertically. A Sampling Template 10x10 (ref. 96762) can be used to delimit and sample an area of 100 cm^2.
5. Insert the swab in the tube containing the culture medium.
6. Cap the tube and record date and the sampling point.
7. Incubate at 30°C for 24-48 hours.

RESULTS INTERPRETATION

Watch for the color change of the medium and interpret the results as shown in the table.

Yellow Medium	Red Medium
	
Absence of MRSA	Presence of MRSA

QUALITY CONTROL

Each batch of CONTAM SWAB - MRSA is submitted to quality control using the following bacterial:
Escherichia coli ATCC® 25922, *Proteus mirabilis* ATCC® 25933, *Pseudomonas aeruginosa* ATCC® 27853, *Staphylococcus aureus* (MRSA) ATCC® 43300, *Staphylococcus aureus* (MSSA) ATCC® 6538.

STORAGE

2-8°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

CONTAM SWAB - MRSA does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

PRESENTATION

Product	**Ref.**	Σ
CONTAM SWAB - MRSA	86104	30 Tests

REFERENCES

- ISO 18593 - Microbiology of food and animal feeding stuffs - Horizontal method for sampling techniques from surfaces using contact plates and swabs.
- ISO 17604 - Microbiology of food and animal feeding stuffs - Carcass sampling for microbiological analysis.

TABLE OF SYMBOLS

LOT	Batch code		Do not reuse		Manufacturer		Use by		Fragile, handle with care
REF	Catalogue number		Temperature limitation		Contains sufficient for <n> tests		Caution,consult accompanying documents		

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DOMINIQUE DUTSCHER SAS