Making Testing Easier

PROTEON GLUTEN EXPRESS

ZE/PR/GL10

Immunochromatographic test
for gluten detection

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DOMINIQUE DU

TEST PRINCIPLE

PROTEON GLUTEN EXPRESS is a qualitative and immunochromatographic rapid test for detection of gluten traces in food and working surfaces. The assay is based on the reaction of the prolamins with a specific antibody coated with red particles. The appearance of two lines, red and blue, on the strip test is a positive result.

Limit of detection in food: 3 ppm of gluten.

Limit of detection in surfaces: 1 µg of gluten.

Consult the technical sheet or contact ZEULAB for more information.

SAFETY

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. This kit should be used following good laboratory practices.

STABILITY AND STORAGE

Store the test under dry conditions and temperatures between 15 and 25ºC. Check the expiry date on the package.


NOTES

- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed.
- The test should be performed between 15 and 25 ºC.
- Use a filter bag, a pipette, tube and swab for each sample.
- To handle the strip, hold it from the opposite end of the arrow area.
- The test has been validated in a large number of representative matrices. It is recommended to verify that the test works with your specific matrix before using it. To do this, analyze at least one positive sample of your matrix and one negative to ensure that the test works correctly. It is also recommended to re-verify it if there are changes in the ingredients or in the production processes of the matrix. Please, contact ZEULAB for further information.
- WARNING: The kit detects allergen traces. The pure allergen should not be tested. Concentrations above 10 g gluten proteins/ Kg food product can lead to an invalid or a negative result, due to "prozone or Hook effect". For this reason, gluten cereal samples can not be directly analysed. Dilutions higher than 1/1000 are recommended for the analysis.

DOMINIQUE DU

NOTES for food analysis and rinse water

- Additional material:
 - Beater/grinder
 - Balance

- For using the filter bag, add sample and buffer inside the inner bag. Close it with the help of the wire, folding it upon itself.
- Samples must be representative (at least 50 g) and homogeneous.
- The provided syringe should be used to add the analysis buffer to each sample. Extreme care is needed not to touch the bag content to avoid contamination.
- Make sure that the sample does not reach the area of the test strip arrows.
- For **samples containing polyphenols, tannins or antioxidants** special procedure must be followed:

 0.4 g of Polyvinylpyrrolidone (Sigma PVP-10) and 1 g of gelatine (fish gelatine) to the extraction buffer and continue the procedure in step 2. These samples must be incubated at 50°C for 40 minutes. The minimum level of detection for these samples is 10 ppm.

NOTES for surfaces analysis

- The area should be visually clean.
- Swab a representative area of at least 100 cm^2.
- The result obtained from the test is the presence or absence of the allergen on the analysed surface, and it cannot be extrapolated to ppm values.
- The test limit of detection on surfaces is subjected to the use of the swabs provided by ZEULAB.

PROCEDURE

(*Please, read NOTES before running the assay*)

Analysis of food and rinse water

For liquid samples such us rinsing water or juices the analysis can start in step 4.

1. Sampling
 - For solid samples, take a representative sample (50g) and mash or crush to homogenise. Weigh 1 g inside of the filter bag and add 10 ml of EXTRACTION BUFFER using the 10 ml syringe.

DOMINO DE DU

- For liquid samples, place 1 ml of sample (or weigh 1 g) into de filter bag and add 9 ml of EXTRACTION BUFFER using the 10 ml syringe.

2. Rub the bag with your hands to obtain a homogeneous mixture (1-2 min approximately).
3. The filtered sample (placed between the filter and the outer bag) is ready to be analysed.
4. Dilute the filtered sample with ANALYSIS BUFFER. The dilution will depend on the desired level of detection (3, 10 or 20 ppm). See table below.

 For detecting 20 ppm of gluten, take 100 µl of the filtered sample and add 10 ml of ANALYSIS BUFFER.
5. Pour the diluted sample into the test tube cap.
6. Introduce the strip test into the tube cap following the arrows' direction.
7. Wait for 10 min and read the test result following indications on "Interpretation of the results".

Dilutions for levels of detection of 3, 10 and 20 ppm

Level of detection	**3 ppm**	**10 ppm**	**20 ppm**
Dilution	1/15	1/50	1/100
Sample Volume	100 µl	100 µl	100 µl
ANALYSIS BUFFER Volume	1.5 ml	5 ml	10 ml

Analysis of surfaces

1. Add 0.5 ml of ANALYSIS BUFFER into a tube, use the 1 ml syringe.
2. Introduce a swab into the tube and soak it in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
3. Rub the swab along in all possible directions on the analysed surface. Rotate from time to time the swab's point in contact with the surface.
4. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
5. Dispose of the swab.
6. Introduce the strip test into the tube following the arrows' direction.
7. Wait for 10 min and read the test result following indications on "Interpretation of the results".

DOMINIQUE DU

INTERPRETATION OF RESULTS

Visual reading

NEGATIVE: Only a blue line (control line) is shown. This is the control line and it indicates that the test works properly.

POSITIVE: A blue line (control line) and a red line (test line) are shown.

INVALID: The blue line (control line) is not shown. This may occur when the assay has not been performed properly, when the matrix interferes or when the reagents are damaged. Repeat the assay with a new strip and contact ZEULAB for further information.

Objective reading – Reader IRIS

Place the strip in the proper tray (named at number 5) and introduce it in the reader IRIS. The NEGATIVE/POSITIVE result will be displayed on the screen. Please, check the IRIS instructions for further information.

Note: There is a specific reading method for each batch. The method is transferred IRIS using a QR code that is available in the kit, or using IRIS Lite. When different batches are used, pay special attention to the method uploaded on IRIS.

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DOMINIQUE DU

PROTEON is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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