Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

AREX 5 Heating Magnetic Stirrer
F20500570, F20510570, F20580580

ARE 5 Heating Magnetic Stirrer
F20500560, F20510560, F20580560

RC 5 Heating Plate
F20700610, F20710610, F20780610

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
使用前请仔细阅读以下使用说明书。

Caution, hot surface!
Attenzione, superficie calda!
Attention, surface chaude!
Prudencia, superficie caliente!
Vorsicht, heiße Oberfläche!
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.
Il prodotto può essere utilizzato con liquidi infiammabili.
Le produit peut être utilisé avec des liquides inflammables.
El producto puede utilizarse con líquidos inflamables.
Das Produkt kann mit brennbaren Flüssigkeiten verwendet warden.
本产品可与易燃液体混合使用。

DOMINIQUE DUTSCHER SAS

This unit must be used for indoor laboratory applications only.
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with non specified accessories, product's safety may be compromised.

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
NON UNIQUE DUTSCHER SAS

Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids (< 50ml) or flammable liquids that have a fire point higher than 825°C and a flash point higher than 800°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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DOMINIQUE DUTSCHER SAS

1. Introduction

The AREX / ARE 5 heating magnetic stirrer is used to heat and mix liquids inside a suitable container placed on the aluminium plate. The stirring is possible thanks to a magnetic drive stir bar placed into the vessel.

The RC 5 heating plate is used to heat liquids inside a suitable container placed on the aluminium plate.

The instrument is used in the laboratory for general use and for all those applications that require precise regulation of the heating plate temperature (and, for AREX / ARE 5, also of the stirring speed).

The Aluminum alloy heating plate ensures:
- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock

The AREX 5 heating plate has a special ceramic coating that ensures high resistance to chemical agents and surface abrasion. Using this heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

For AREX / ARE 5, the powerful brushless motor is able to stir from 100 up to 1500 rpm.
SpeedServo technology maintains constant speed as viscosity changes.

***Note**: it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).*

The AREX 5 can work in combination with thermoregulator VTF for precise and accurate control of the sample temperature.

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking

- The box includes
 - AREX / ARE 5 heating magnetic stirrer or RC 5 heating plate
 - Instruction manual
 - Power cable
 - Magnetic stir bar Ø 8 x 40 mm (A00000356)

- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and of the power supply line are the same
 - Make sure that the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

 ***Note**: keep the power cable far away from the hot plate.*

 - Place the container with the liquid and the stir bar (only for AREX / ARE 5) on the stirring plate

3. Operation

Power-up
- Turn the unit on using the main switch
- The right led between the two knobs lights up

Stirring (only for AREX / ARE 5)
- Rotate the right knob to set the stirring speed, the led over the speed knob lights up
- The real speed increases until setpoint achievement
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo)
- Turn the knob to the 0 positions to stop stirring

DOMINIQUE DUTSCHER SAS

Heating

- Rotate the left knob to set the heating plate temperature
- If the set temperature is higher than the actual one, the led over the temperature knob lights up
- When the heating plate temperature exceeds 50 °C, the left led between the two knobs lights up until the temperature drops below 50 °C

 ***Note**: This warning is active even if the main switch is off but not if the instrument is not powered.*

- Turn the knob to the 0 positions to stop heating

Operation with thermoregulator VTF (only for AREX 5)

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed into the sample.
- Plug the VTF into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch. When the VTF is correctly connected, the display of the VTF lights up showing the temperature of the sample
- Select the temperature set point on the VTF
- Turn the temperature knob to the max value to start heating

 ***Note:** The AREX 5 always has primary control of the heating plate temperature. The temperature control function of the heating plate can also be used as a safety thermostat. In this case, the maximum temperature of the heating plate will not exceed the temperature setting on the magnetic stirrer, meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature set, thus reducing temperature oscillation at the setpoint value.*

4. Error and warning messages

When the control panel shows one of the following error messages, the instrument functions automatically stop.

Display indicator	**Malfunction:**
Flashing Heating led and Hot plate led ON	Heating plate overtemperature (T > 330°C)
Flashing Heating led and Hot plate led OFF	Heating plate safety probe overtemperature (T >370°C)
Flashing Stirring led (only for AREX / ARE 5)	The stirring system doesn't run correctly

To remove the error message, restart the instrument from the main switch. For the overtemperature of the safety probe is necessary to disconnect the instrument from the power supply.
If the alarm persists on the control panel, please contact VELP Scientifica's technical assistance service.

5. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

Cleaning
Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair
Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g. palletizing).
It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

DOMINIQUE DUTSCHER SAS

6. Technical data

		AREX 5			ARE 5			RC 5		
General features	Model	F20500570	F20510570	F20580570	F20500560	F20510560	F20580560	F20700610	F20710610	F20780610
	Voltage	230 V 50/60 Hz	115 V 60 Hz	100 V 50/60 Hz	230 V 50/60 Hz	115 V 60 Hz	100 V 50/60 Hz	230 V 50/60 Hz	115 V 60 Hz	100 V 50/60 Hz
	Power input	630 W			630 W			600 W		
	Current consumption	2,8 A	5,5 A	6,3 A	2,8 A	5,5 A	6,3 A	2,6 A	5,2 A	6 A
	Dimensions (WxHxD)	160x85x270 mm (6,3x3,3x10,6 in)								
	Weight	1,7 kg (3,7 lb)						1,4kg (3,1 lb)		
	Construction material	Aluminum body –Technopolymer enclosure								
	Working in continuous	Admitted								
	Maximum load on the plate	25 kg								
	Noisiness	<< 80 dBa								
	Environmental temperature admitted	+5…+40 °C								
	Storage temperature admitted	-10…+60 °C								
	Level of electrical protection CEI EN60529	IP 42								
	Max humidity	80%								
	Overvoltage category	II								
	Pollution degree CEI EN61010-1	2								
	Max altitude	2000 m								
Heating	Heating plate power output	600 W								
	Heating plate dimensions	Ø 135 mm								
	Temperature range	50 ÷ 310 °C								
	Heat control type	Analogic								
	Heating plate construction material	CerAlTop™			Aluminum					
	Safety circuit	Separated with dedicated probe								
	Hot plate alarm	Over 50°C								
	Overtemperature alarm	Over 330°C								
VTF	Temperature range	-10 ÷ 300 °C (1 °C step)								
	Temperature resolution	0,2 °C								
	Accuracy	+/-0,5°C*								
	Electrical data	12Vdc – 1.2W (max)								
Stir	Stirring capacity	20 l H_2O			15 l H_2O					
	Speed range	100 ÷ 1500 rpm								
	Stir control type	Analogic								
	Motor type	BLDC								
	Motor rating output	10 W								

* in the following conditions: 800 ml of water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50°C.

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
EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:
e soddisfa i requisiti essenziali delle direttive:
et qu'il satisfait les exigences essentielles des directives:
y cumple con los requisitos esenciales de las directivas:
und den Anforderungen folgender Richtlinien entspricht:
并满足以下指令的基本要求:

2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung / 加修改.

DOMINIQUE DUTSCHER SAS

11. Declaration of conformity UK CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

DOMINIQUE DUTSCHER SAS

Thank you for having chosen Velp!

Established in 1983, Velp is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. Velp works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.
Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments
Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment
Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortexers and Shakers
Dispersers
Dry Block Heaters and COD
BOD and Respirometers
Cooled Incubators
Flocculators
Turbidimeter
Open Circulating Baths
Pumps

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Velp Scientifica Srl