Medicom

Medicom SAS
Boulevard de la chanterie | 49124 Saint Barthélemy d'Anjou
France
Tel +33 0 241 963 434 | sales.support.pro.eu@medicom.com | www.medicom-eu.com

CH REP
Finagest SA
Rue Pierre-Fatio 12 | 1204 Genève
Switzerland

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Medicom

ENGLISH – EN

Description and indications: The shoe covers are category I personal protective equipment (Regulation EU 2016/425) and class I medical device (Regulation EU 2017/745 and Swiss Ordinance RS 812.213). They are designed:

- As personal protective equipment to provide protection against minor risks only (small splashes of cleaning materials of weak action or prolonged contact with water and superficial mechanical injuries), and to limit contamination in sensitive areas.
- As a medical device to cover the shoes of the healthcare professional and/or the feet of the patient during surgical procedures, medical care or examinations, in order to prevent the risk of cross-contamination.

Reference	Size	Colour
10.001D-10	L	White
H01038A	L	White
10.040R	XL	White with blue sole
10.041R	XL	White with white sole

Limitations on use: The user will be the sole judge of the suitability of the protective accessories for the situation and the protection required. The manufacturer cannot be held responsible for the improper or inappropriate use of the accessories. The shoe covers are for single use only and must be disposed of after any job. The shoe covers must be used on intact skin only. They are for professional use and do not protect against biological and chemical hazards. Do not use them if any defects are noticed (e.g.: holes, damaged seams, ...). The shoe covers are flammable, keep away from fire.
Biocompatibility: The shoe covers are non-cytotoxic, non-irritant, and non-sensitizing according to the requirements of the standards EN ISO 10993.
Performances: The shoe covers meet the requirements according to the EN ISO 13688:2013+A1:2021 and the following standards:

Tests	Standards	Results
Tests on the nonwoven		
Resistance to water penetration – increasing hydrostatic pressure	EN ISO 811 :2018	13.4 cm H_20
Abrasion resistance	EN ISO 12947-2 : 2017	6 000 cycles
Tensile strength	EN ISO 9073-3 :2023	Machine direction : 70 N Cross direction : 55 N
Elongation at maximum force	EN ISO 9073-3 :2023	Machine direction : 47 % Cross direction : 55 %
pH	EN ISO 3071 :2020	6.0
Tests on the sole		
Resistance to water penetration – increasing hydrostatic pressure	EN ISO 811 :2018	> 2 000 cm H_20
Abrasion resistance	EN ISO 12947-2 : 2017	45 000 cycles
Tensile strength	EN ISO 9073-3 :2023	Machine direction : 134 N Cross direction : 93 N
Elongation at maximum force	EN ISO 9073-3 :2023	Machine direction : 63 % Cross direction : 66 %

Shelf life: Use the product within 5 years of its date of manufacture.
Storage and disposal: Shoe covers should be stored in their original packaging in a dry place and away from sunlight. They must be disposed off in accordance with the regulations in force.
Incident notification: Any serious incident arising from the use of the product must be reported to the manufacturer and to the competent authority of the Member State in which the user is established.
The internet address at which the declaration of conformity can be accessed is: www.medicom-eu.com

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DOMINIQUE DUTSCHER SAS