Medicom

INSTRUCTIONS FOR USE

SINGLE USE
NITRILE GLOVES

1142/1144

FR EN DE ES IT

PT SV NL PL AR

CH REP

Medicom SA
Boulevard de la Chanterie
49124 Saint-Barthélemy-d'Anjou
France

Finagest SA
Rue Pierre-Fatio 12
1204 Genève
Switzerland

medicom-eu.com

Rev 02 - 2024/07

SYMBOLS & PICTOGRAMS

A

B

C

D

E

F

G

H

(01)

EN ISO 374-5:2016
VIRUS

(02)

EN ISO 374-1/Type B
JKT

(03)

EN 421:2010

(04)

CE
2927

(06)

MD

(07)

(08)

LATEX

(09)

(10)

(11)

(12)

(13)

(14)

(15)

(16)

(17)

(18)

(19)

(20)

(21)

FR

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EN

1142_MEDICOM KOLMI COVERFEEL CARE & 1144_MEDICOM KOLMI COVERFEEL CARE LONG - INSTRUCTIONS FOR USE FOR CATEGORY III PERSONAL PROTECTIVE GLOVES, LETHAL OR IRREVERSIBLE RISKS - CLASS I MEDICAL EXAMINATION GLOVES - SINGLE USE - POWDER-FREE - NON-STERILE. These gloves are recommended for protection against micro-organisms (bacteria, fungi and viruses), tested chemicals and chemotherapy drugs. These gloves are also intended for medical use (examination gloves) by healthcare professionals. These gloves are suitable for food contact. These gloves are available in the following sizes: 6/XS - 7/S - 8/M - 9/L - 10/XL - 11/XXL. The length of glove 1142 = 240mm and glove 1144 = 290 mm. **INSTRUCTIONS:** Choose the appropriate size. Put the glove(s) on as tightly as possible. To achieve the desired protection, these gloves can be used with other suitable PPE (overalls, footwear, respiratory protection, etc.) - Carry out a risk analysis according to the intended use. Ensure a perfect seal at the junction points of the various accessories (compatible adhesives can be used). Special precautions for donning and doffing may be required depending on the area of application (e.g. double gloving). Donning procedure: follow figures A-B-C-D and doffing procedure: follow figures E-F-G-H. These gloves comply with the standards illustrated by the relevant pictograms. **EXPLANATION OF SYMBOLS, PICTOGRAMS AND PERFORMANCE LEVELS:** (01) EN ISO 21420:2020 - Read these instructions before using the products. (02) EN ISO 374-5: 2016 - Protection against bacteria, fungi and viruses. Resistance to penetration has been assessed under laboratory conditions and relates only to the test specimen of the test. (03) EN ISO 374-1:2016+A1:2018 TYPE B - Protection against chemicals: (Type B = breakthrough time > 30 minutes for at least 3 chemicals) - Chemical permeation data, obtained using the test method in EN 16523-1:2015+A1:2018 and degradation data, obtained using the test method in EN ISO 374-4:2019 ; Breakthrough time/performance level: >10min/level 1; >30min /level 2; >60min /level 3; >120min /level 4; >240min /level 5; >480min /level 6 – Results: J = n-heptane, Permeation level 6 and degradation 46,7%, K = Sodium hydroxyde - Permeation level 6 and degradation -16.8%, and T = formaldehyde - Permeation level 6 and degradation 33.9%. ASTM D6978-05: Resistance to permeation by a chemotherapy drug (Breakthrough time = maximum duration of use of the cytostatic): Carmustine = 43,8min, Doxorubicin Hydrochloride >240min, Cyclophosphamide >240min, Paclitaxel >240min, Etoposide >240min, Fluorouracil >240min, Thiotepa = 26,7min, Busulfan >240min, Irinotecan >240min, Vidaza >240min, Fentanyl Citrate >240min, Cytarabine >240min, Ifosfamide >240min, Methotrexate >240min, Vincristine Sulfate >240min. **WARNING!** Make sure the gloves are suitable for their intended use. Before use, we recommend that you inspect your gloves to ensure that they are free from defects or imperfections. When worn, protective gloves may offer less resistance to hazardous chemicals and micro-organisms, due to changes in their physical properties. Movement, snags, rubbing or degradation caused by contact with chemicals, etc. can considerably reduce the actual duration of use of a chemical protective gloves. Degradation may be the most important factor to consider when choosing chemical-resistant gloves. Chemical resistance has been assessed under laboratory conditions using samples taken from the palm only and relates only to the chemical being tested. It may be different if used in a mixture. This information does not reflect the actual duration of protection in the workplace, nor the differentiation between mixtures and pure chemicals. This is because conditions in the workplace may differ from those in the type test, depending on temperature, abrasion and degradation. As a result, and insofar as Medicom SAS has no control over the conditions of use, the responsibility for controlling the actual conditions of use of its products, these recommendations are purely indicative and cannot engage its responsibility. (04) EN 421:2010 - Protection against radioactive contamination. These gloves do not protect against microbial risk. According to EN 388:2016+A1:2018 dexterity (level 5) is the highest parameter for our special purpose glove (protection against micro-organisms and viruses, radioactive contamination and chemical risks). Does not protect against ionising radiation. Do not use in confinement enclosures. (05) CE2927= Product certified to comply with the provisions of European Regulation (EU) 2016/425 on personal protective equipment. The CE marking is followed by a four-digit code that corresponds to the identification number of the notified body in charge of conformity assessment of category III PPE, for products intended to protect against serious risks. Module B carried out by the notified body CTC (0075), Parc Scientifique Tony Garnier - 4 rue Hermann Frenjel - 69367 Lyon Cedex 07 - France and Module D (Conformité en vue de l'assurance de la qualité du procédé de production) carried out by the notified body QUINTIN CERTIFICATIONS (2927), 825 route de Romans - 38160 Saint-Antoine-l'Abbaye - France. (06) Medical Device (MD) marking, class I MD according to European regulation (EU) 2017/745. (07) Pictogramme Food contact = Apte au contact alimentaire. Conforme aux règlements européens (CE) 1935/2004 et (CE) 2023/2006. (8) DO NOT RE-USE (9) DOES NOT CONTAIN NATURAL RUBBER LATEX (10) KEEP DRY (11) KEEP AWAY FROM SUNLIGHT (12) DO NOT USE IF PACKAGING IS DAMAGED AND CONSULT INSTRUCTION FOR USE (13) EXPIRY DATE (14) DATE OF MANUFACTURE (15) MANUFACTURER (16) BATCH NUMBER (17) BATCH NUMBER (18) REFERENCE (19) UNIQUE DEVICE IDENTIFIER (20) DISTRIBUTOR (21) CH RESPONSIBLE PERSON. **PRECAUTION FOR USE:** The glove must ONLY be used for the risks and performances described, by trained personnel. These gloves do not protect against mechanical risks. The combined use of this product with other PPE of a higher class does not extend the scope of application of this product. The user will be solely responsible for choosing the correct combination and adaptation of the glove with other equipment (overalls, footwear, PPE, etc.) and the period during which the glove may be worn, depending on the task to be performed, the level of protection required and comfort of use. Keep product and packaging out of reach of children. For use on intact skin only. The gloves are for single use only and are not designed to be washed out or reused. Any reuse may lead to post-contamination and post-operative complications. The manufacturer has not validated the cleaning and re-sterilization of this product. Any damaging process likely to compromise the integrity of the product. The materials and components making up the glove do not contain substances at levels known or suspected to have adverse effects on the hygiene or health of the user under the foreseeable conditions of use. These gloves may contain residues of chemicals used in their manufacture which may cause allergic skin reactions in some people. In the event of a skin reaction, please discontinue use of the glove. Any serious incident arising from the use of the product must be reported to the manufacturer and to the competent authority of the Member State in which the user is established. **STORAGE:** Store the product in its original packaging, away from humidity, sunlight and light. In this way, their performance will remain intact and their characteristics will not change significantly. Use the product within 3 years of its date of manufacture. **DISPOSAL:** They should be discarded as soon as the slightest sign of deterioration during use appears (tear, hole, discolouration or embrittlement, for example). Depending on the area of application, the user must comply with the national regulations in force for the product and its packaging.
The internet address at which the declaration of conformity can be accessed is: www.medicom-eu.com

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