Making Testing Easier

RC BOVINO

ZE/RCB96
ZE/RCB960

Immunoenzymatic test for detection of cow's milk in sheep or goat's milk

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ENGLISH

SCOPE

RC-BOVINO is a sandwich ELISA test based on the specific detection of bovine IgG. It is a test for the determination of cow´s milk in sheep´s milk, sheep-cheese, goat´s milk and goat-cheese. Working range of RC Bovino is 0-20% of cow's milk when samples are diluted 1:200.

KIT COMPONENTS

	ZE/RCB96	ZE/RCB960
Microtiter plate for 96 tests	1	10
SET STANDARDS (0%, 0.01%, 0.05% and 0.10%)	1	3
CONJUGATE x20	1	3
CONJUGATE DILUENT	1	3
SUBSTRATE	1	1
STOP SOLUTION	1	2

STABILITY AND STORAGE

All kit components must be stored refrigerated (4-12ºC) and protected from light. The kit is provided with a minimum shelf life of 6 months. See kit expiring date on the package.

ADDITIONAL MATERIAL (NOT PROVIDED)

- ✓ Micropipettes.
- ✓ Microtiter spectrophotometer with 450 nm filter.
- ✓ Pestle and mortar (for cheese samples).

SAFETY

Safety clothing should be worn and skin contact with the reagents avoided. Do not ingest. A SAFETY DATA SHEET is available from your local distributor on request.

*Warning: *STOP SOLUTION* contains sulfuric acid. Gloves and other protective clothing must be worn when handling this solution.

NOTES

- The microtiter plate is provided dry in a sealed bag. Take out only the strips to be used. Store the remaining strips in the provided bag at 4ºC.
- The SUBSTRATE solution has a pink colour under normal circumstances. A blue colouring of the SUBSTRATE indicates contamination of the reagent and should be discarded. Protect the reagent from exposure to light and air. Avoid dispensing this solution pipetting directly from the tube because it could be contaminated.
- Do not allow wells to dry out while running the assay.
- Let reagents to reach room temperature before using the kit (at least 30 min). Return it to refrigerator after use.

DOMINIQUE DUTSCHER SAS

- Whole, skimmed, fresh and pasteurised milk, and whey samples can be used with this kit. However UHT and sterilization treatments can destroy milk immunoglobulins and results can be underestimated. Therefore, it is not recommended to use this test with samples treated at those temperatures.
- Lab temperature should be between 20 and 24ºC while running the test. Higher or lower temperatures could produce non-lineal calibration curves or low absorbance values respectively.
- For cheese analysis, only qualitative assays should be considered. Positive results will indicate the presence of cow´s milk and negative ones cow´s milk absence in cheese. During the cheese making, the lost of soluble proteins is variable, depending on the variety of cheese and the production system of each cheese maker. Therefore, for a quantitative estimation in cheese, it is recommended to develop specific standards for each variety of cheese.

PREPARATION OF SAMPLES AND REAGENTS

SUBSTRATE , STOP SOLUTION and STANDARDS are ready to use.

Preparation of milk samples

- Ranges of mixtures between 0 and 10%. Dilution 1/100 as follows:
 Mix 0.1 mL of milk sample with 9.9 mL of distilled water
- Ranges of mixtures between 10 and 20%. Dilution 1/200 as follows:
 Mix 1 ml of milk sample with 9 ml of distilled water
 Mix 0.5 ml of first dilution with 9.5 ml of distilled water

Preparation of cheese samples:

- Weight 1 ± 0.1 g of homogenised cheese sample and add 10 mL of distilled water.
- Homogenize the mixture with a pestle and a mortar until a cheese paste is obtained.
- Centrifuge at 3000 x g for 10 minutes (alternatively, it can be filtered through a filter paper).
- The supernatant will be used in the assay.

Preparation of conjugate solution:

Prepare the quantity of conjugate solution required by mixing 1 volume of *CONJUGATE x20* with 19 volumes of *CONJUGATE DILUENT*. For instance, for analysing 8 wells it is necessary 8 x 0.05 mL = 0.4 mL of conjugate solution. Prepare 0.5 mL of solution to have a small volume of excess, mixing 0.025 mL of *CONJUGATE x20* with 0.475 mL of *CONJUGATE DILUENT*.

In case of ZE/RCB960, for analysing large number of samples it is recommended to prepare conjugate solution adding the volume of 1 bottle of CONJUGATE x20 to 1 bottle of CONJUGATE DILUENT. This conjugate solution is stable for 5 days when stored under refrigeration (4-12ºC).

DOMINIQUE DUTSCHER SAS

TEST PROCEDURE


1. Leave the kit components at room temperature for at least 30 min.
2. Add 50 µL of sample or standard into each well. Incubate at 20-24ºC for 30 minutes. Analyse samples and standards by duplicate.
3. Empty the strip / plate by turning the plate upside down and wash the plate/wells with distilled water by adding 300 µL/well of distilled water, when using a micropipette, or filling the well up, when using a squeeze bottle. Empty the wells/plate by turning the plate upside down on top of an absorbent paper to remove the excess of water. Repeat the washing step 2 to 3 times
4. Add 50 µL of conjugate solution into each well. Incubate at 20-24ºC for 30 minutes.
5. Repeat 4 times the washing sequence as in step 3.
6. Add 50 µL of SUBSTRATE into each well. Incubate at 20-24ºC for 30 minutes.

 Note: Pour the volume of SUBSTRATE to be used in a clean container to prevent contamination
7. Add 50 µL of *STOP SOLUTION* into each well and mix gently.
8. Read immediately the absorbance of each well at 450 nm with a microtiter spectrophotometer.

DATA ANALYSIS

Note: A worksheet to calculate the results is available from ZEULAB free of charge. Please, contact us for further information.

1. Obtain a standard curve by plotting % of cow milk in standard samples in a X-axis and absorbance values in the Y-axis. Draw standard curve using a lineal or polynomial fitting. Standard curve goes from 0 to 0.10% of cow's milk. An example of standard curve is shown to the right of this text.

2. To calculate the percentage of cow's milk in the sample is needed to interpolate the absorbance values obtained into the calibration curve and multiply by the diluting factor used. If samples have been diluted 1:100 the working range will be up to 10% of cow milk. If dilution used is 1:200, then the working range is up to 20%.

 Any sample that produces optical density above the value of the highest standard should be further diluted with distilled water and tested again. Dilutions higher than 1:200 are not recommended.

For cheese analysis, only qualitative assays should be considered (see Notes).

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PRINTING DUTSCHER SAS

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1. After preparation of samples, add Standards/Samples.

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2. Incubation 30 minutes at room temperature. Wash.

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3. Add conjugate and incubate 30 min at room temperature. Wash.

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4. Add substrate and incubate 30 min at room temperature.

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5. Add stop solution and read the results at 450 nm.

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DOMINIQUE DUTSCHER SAS

RC BOVINO is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be revaluated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility

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