Sonda PHC201 PT-BR

Introdução

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


PHC201



www.hach.com

© Hach Company/Hach Lange GmbH, 2011, 2013, 2021. All rights reserved.
Printed in U.S.A. 09/2021, Edition 3

DOMINIQUE DUTSCHER SAS

PHC201 probe

EN

Introduction

The probe is a non-refillable, semi-liquid gel filled combination pH probe with a built-in temperature sensor for general aqueous applications.

Warning

⚠ *Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.*

User manual

A user manual with all instructions is on the manufacturer's website.

Collect the following items:

Compatible meter; pH standard solutions; beakers; DI water; lint-free cloth.

1 **Connect the probe**

2 **Calibration**

Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

3 **Measurement (direct method)**

4 **Storage**

Technical information

pH range	pH 0 to 14
Slope	-59 mV/pH (85 to 115% at 25 °C (77 °F) per Nernstian theoretical value
Operating temperature	0 to 80 °C (32 to 176 °F)
Temperature accuracy	±0.3 °C (±0.54 °F)
Minimum sample depth	15 mm (0.59 in.)
Storage solution	Hach pH electrode storage solution or 3 M KCl

Warranty

6 months on the probe. This warranty covers manufacturing defects, but not improper use or wear.

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