Sonda MTC101 PT-BR

Introdução

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


MTC101



www.hach.com

© Hach Company/Hach Lange GmbH, 2010-2011, 2013, 2014, 2021. All rights reserved.
Printed in U.S.A. 09/2021, Edition 5

DOMINIQUE DUTSCHER SAS

MTC101 probe

EN

Introduction

The probe is a non-refillable, gel-filled combination oxidation reduction potential (ORP/Redox) probe with a built-in temperature sensor.

Warning

⚠ Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.

User manual

A user manual with all instructions is on the manufacturer's website.

Collect the following items:

Compatible meter; ORP/Redox solution; beakers; DI water; lint-free cloth.

1 Connect the probe

2 Calibration

Make sure that the platinum electrode is clean before sample measurement.

3 Measurement

⚠Make sure that the platinum electrode is clean before sample measurement. Soaking in reducing solution or sample will speed up measurements with negative mV samples, as with waste water samples. Keep the shroud installed on the rugged probe during sample measurements to prevent damage to the sensor.

4 Storage

Technical information

Range	±1200 mV
Reference potential versus Standard Hydrogen Electrode	207 mV at 25 °C
Operating temperature	0 to 80 °C (32 to 176 °F)
Temperature accuracy	±0.3 °C (±0.54 °F)
Minimum sample depth	20 mm (0.79 in.)
Storage solution	Hach pH electrode storage solution or 3 M KCl

Warranty

6 months on the probe. This warranty covers manufacturing defects, but not improper use or wear.

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