Sonda LDO101 PT-BR

Introdução

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

LDO101



www.hach.com


© Hach Company/Hach Lange GmbH, 2010, 2013, 2017, 2021–2023. All rights reserved.

Printed in U.S.A. 01/2023, Edition 6

DOMINIQUE DUTSCHER SAS

LDO101 probe

EN

Introduction

The probe measures the dissolved oxygen concentration in wastewater, drinking water and general applications.

Warning

⚠ *Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.*

User manual

A user manual with all instructions is on the manufacturer's website.

Collect the following items:

Compatible meter; timer; BOD bottle; DI water; lint-free cloth.

1 Connect the probe

2 Calibration

Note:
Allow up to 30 minutes for contents to equilibrate to room temperature.

3 Measurement

⚠ *The protective shroud is recommended during measurement to protect from damage to the sensor cap and light shielding for field applications.*

	Do not submerge the pressure sensor module.

4 Storage

The storage conditions for the LDO probe can be dry, wet, or humid.

Technical information

Dissolved oxygen range	0.1 to 20.0 mg/L (ppm) 1 to 200% saturation
Dissolved oxygen accuracy	±0.1 mg/L for 0 to 8 mg/L ±0.2 mg/L for greater than 8 mg/L
Operating temperature	0 to 50 °C (32 to 122 °F)
Minimum sample depth	25 mm (0.984 in.)

Warranty

3 years on the probe, 1 year on the sensor cap. This warranty covers manufacturing defects, but not improper use or wear.

Note:
An iButton is no longer supplied with the probe. It is not necessary for new probes to have an iButton, because the calibration data is installed on the probe at the factory.

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