Sonda ISENA381 — PT-BR

Introdução

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

ISENA381

www.hach.com

© Hach Company/Hach Lange GmbH, 2010-2011, 2013, 2021. All rights reserved.
Printed in U.S.A. 09/2021, Edition 5

DOMINIQUE DUTSCHER SAS

ISENA381 probe

EN

Introduction

The probe is a digital, combination ion selective probe with a double junction, refillable junction reference and a built-in temperature sensor.

Warning

⚠ *Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current safety data sheets (MSDS/SDS) for safety protocols.*

User manual

A user manual with all instructions is on the manufacturer's website.

Collect the following items:

Compatible meter; beakers; Na^+ standard solutions; rinse bottle; ISA (0.4 g) powder pillows; prepared ISA rinse solution; lint-free cloth.

1 **Connect the probe**

2 **Calibration**

Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

3 **Measurement (direct method)**

4 **Storage**

⚠*Solution: 25 mL of 100 mg/L Sodium standard plus one ISA pillow (0.4 g).*

Technical information

Range	0.023 mg/L ($1x10^{-6}$ M) to 23,000 mg/L (1 M) Na^+
Slope	-59 mV/pNa (90 to 110% at 25 °C per Nernstian theoretical value)
Operating temperature	0 to 50 °C (32 to 122 °F)
Minimum sample volume	25 mL
Filling solution	0.02 M NH_4Cl

Warranty

1 year on the probe. This warranty covers manufacturing defects, but not improper use or wear.

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