
52 03


pH

HACH

DOC022.98.90535
10-2020, Edition 2a

52 03

DOMIN

Puntos clave. *Punti chiave.* Points clé. *Key points.*

1. Tapón protector.
 Tappo protettore
 Bouchon protecteur.
 Protective cap.

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DOMIN

pH electrode 52 03

Introduction
The **52 03** is a pH electrode with screw cap S7, for general use up to pH 14.

Electrode preparation
Remove the storage protector containing electrolyte by unscrewing it ⑤.
Verify that there are not air bubbles inside the membrane ⑧. They can be removed by shaking the electrode as a clinical thermometer. Remove the cap from the electrolyte filling port ③.

Calibration
⚠ *See pH-meter instruction manual.*

Measurements. Recommendations
- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragm ⑦.
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ *Do not leave the electrode immersed in distilled water.*

Maintenance
- Periodically refill the electrode with electrolyte. Maintain the level as high as possible ④.
- Maintain the membrane's surface clean ⑧. Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for HACH regeneration solutions).
- Re-hydration. A "dry" electrode recovers by immersing it in HACH cleaning electrode solution or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting
... the pH-meter always measures the same pH value.
Short circuit in the electrode or connector.
Substitute the electrode.
... unstable reading.
Air bubble in the membrane.
Shake the electrode as a thermometer.
... impossible to perform calibration.
... very slow response.
Dirt on the diaphragm or the membrane.
De-hydrated membrane.
Electrode ageing or wear.
Clean, regenerate or substitute the electrode.

ORIGINAL

Storage
The connector should be protected from humidity by the protective cap ①.
Place the electrode in the storage protector ⑤ with KCl 3M.
Close the refilling port with the stopper ③.

Life of an electrode
The average life expectancy of a pH electrode is approximately one year depending on the working conditions.
An electrode used for more than six months that does not response appropriately after regeneration, must be replaced by a new one.
To avoid accidental breakage of the glass body electrodes, the use of the working protector is recommended.

⚠ *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

Measuring range	0... 14
Operating temp.	0... 80 °C
Reference element	encapsulated Ag/AgCl crystals
Diaphragms	2 ceramic
Electrolyte	KCl 3M
Body material	glass
Min. immersion depth	20 mm

Applications
For aqueous media in general.

Limitations
Low-conductivity or very viscous solutions. Products with colloids or solids in suspension.
Temperatures > 80 °C.

Warranty
HACH guarantees this electrode against manufacturing defects.
Validity: 6 months from shipment date.
Limitations:
- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

DOMINI

QUALITY CERTIFICATE

The attached electrode **52 03** has passed the quality test and meets the specifications.

HACH COMPANY World Headquarters
P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
Tel. (970) 669-3050
(800) 227-4224 (U.S.A. only)
Fax (970) 669-2932
orders@hach.com
www.hach.com

HACH LANGE GMBH
Willstätterstraße 11
D-40549 Düsseldorf, Germany
Tel. +49 (0) 2 11 52 88-320
Fax +49 (0) 2 11 52 88-210
info@de.hach.com
www.de.hach.com

HACH LANGE Sàrl
6, route de Compois
1222 Vésenaz
SWITZERLAND
Tel. +41 22 594 6400
Fax +41 22 594 6499

© Hach Company/Hach Lange GmbH, 2015–2020
All rights reserved. Printed in Germany.