97001 X4610A Anaerogen 28/2/07 12:48 pm Page 1

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OXOID **AnaeroGen** is a convenient and reliable product for the rapid generation of an anaerobic environment which is essential for the isolation and growth of fastidious anaerobes.

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DO NOT DISTRIBUTE

DESCRIPTION

When an AnaeroGen sachet is placed in a sealed jar, the atmospheric oxygen in the jar is rapidly absorbed with the simultaneous generation of carbon dioxide. This novel method differs from those commonly used in that the reaction proceeds with **no** evolution of hydrogen, and therefore does **not** require a catalyst. Furthermore, **no** addition of water is needed to activate the reaction.

When used as directed, the AnaeroGen sachet will reduce the oxygen level in the jar to below 1% within 30 minutes. The resulting carbon dioxide level will be between 9% and 13%.

COMPONENTS

Each box contains:
10 AnaeroGen paper sachets which are individually foil packed.
1 Product Insert.

The active component within each AnaeroGen sachet is ascorbic acid.

PRECAUTIONS

This product is for in-vitro use only.
As soon as the AnaeroGen paper sachet is exposed to air, the reaction will start. It is therefore essential that the paper sachet is placed in the jar and the jar sealed within one minute.
The reaction of the ascorbic acid with oxygen is exothermic. However, the temperature of the AnaeroGen paper sachet will not exceed 65°C.

STORAGE

Store at 2-25°C. Under these conditions, the AnaeroGen sachets will retain their reactivity until the expiry date declared on the outer box and on the foil sachet.

DIRECTIONS

AN 35 is designed for use in 3.5 litre jars. It is therefore suitable for the OXOID Anaerobic jar HP11 and for other jars of similar capacity.

AN 25 is designed for use in 2.5 litre jars. It is therefore suitable for the Oxoid AnaeroGen.

1) Place the inoculated media plates in the appropriate anaerobic jar. Disposable plastic Petri dishes should be of the vented variety to aid gas transfer between the interior and exterior of the plates.
2) Tear open an AnaeroGen foil sachet at the tear-nick indicated, and remove the AnaeroGen paper sachet from within.
3) Immediately place the AnaeroGen paper sachet in the appropriate clip on the plate carrier within the jar.

 N.B. The AnaeroGen paper sachet will become warm to the touch on exposure to air.
4) Close the jar lid immediately.

 N.B. The time taken between opening the foil sachet and sealing the jar should not exceed 1 minute. Extended exposure will result in loss of reactivity, and full anaerobic conditions may not be achieved in the jar.
5) After the appropriate incubation period remove the plates and examine for the presence of anaerobes. If the plates require re-incubation then a fresh AnaeroGen sachet must be used following steps 2)-5) described above.
6) After incubation, the exhausted AnaeroGen sachet should be discarded with the appropriate laboratory waste.

CONTROL TESTING

It is recommended that an OXOID Anaerobic Indicator (BR55) is also used in the jar as a visual check that anaerobic conditions have been achieved and maintained.

The user should check their Anaerobic system periodically for its ability to provide adequate conditions for the growth of appropriate bacteria. The following strains are recommended:

Clostridium novyii	ATCC 9690 growth
Micrococcus luteus	ATCC 9341 no growth

DISPOSAL

On removal from the jar after incubation, the AnaeroGen paper sachet will retain a small amount of reactivity and will warm up. The sachets should therefore be allowed to cool to room temperature prior to disposal alongside the appropriate laboratory waste.

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DOMINIQUE DUTSCHER S.A.S.
DOMINIQUE DUTSCHER SAS