Instruction for Use (GB)

3M™ Comply™
LF Steam Chemical Indicator Label

Class 1: Process indicator for STEAM as classified by ISO 11140-1.
The white indicator stripe changes to a dark brown/black endpoint after exposure to the critical parameters of the steam sterilization process: moisture, time and temperature.
An alpha-numerical printing on the dispenser box label provides the unique lot code/expiry date. The expiry date is based on storage in the sealed package in a dry condition at room temperature.
Description: The 3M™ Comply™ 1322CH LF Steam Chemical Indicator Label consists of a 100mm by 25mm label composed of a 25mm by 25mm piece of 3M™ Comply™ 1322 Lead-Free Steam Indicator Tape, applied to a blue adhesive coated backing delivered on a clear liner. The adhesive is designed to adhere to metal or glass surfaces during steam sterilization. After sterilization, the label is designed to remove easily and cleanly from the goods. The blue backing provides an easy to write on surface for identification of contents. A paper tab is provided to make removal from both the liner and the container following sterilization easy.
Indications For Use: Use the steam chemical indicator label to demonstrate that goods have been exposed to a steam sterilization process and to distinguish between processed and unprocessed goods.
Contraindication:
DO NOT use with dry heat, ethylene oxide, or other low temperature sterilization processes.
DO NOT use the steam chemical indicator label on packs wrapped with any fabric, paper or paper/plastic wraps.
Potential problems may occur with marginal adhesion and/or adhesive residue.

Instructions For Use:
Open the dispenser box and loosen the cut end of the label roll.
Remove one of the 3M™ Comply™ 1322CH LF Steam Chemical Indicator Labels from the clear plastic liner.
Adhere to the outside of a sterilizing container.
If desired, information regarding contents of the load can be written on the blue section of the label with any permanent marker.
After sterilization, allow containers to cool according to established procedures before removing label.
After use the indicator will not change visually within 6 months when stored at <50% RH and 15-30°C.

Explanation of Symbols:

LOT 2000-10AZ • The lot in a box and the hourglass symbols represent lot number and expiration date. The hourglass is followed by a year and month which represent the expiration date (year and month: 2000-10). The entire line after the hourglass represents the lot number. (2000-10AZ).

STEAM • All steam sterilization processes

50% • Humidity limitation, store below 50% RH

Made in Poland
3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss
Germany

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DOMINIQUE DUTSCHER SAS