Wasserdestillierapparate

Water Stills

2001/2 – 2001/4

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Operating Instructions

DOMINIQUE DUTSCHER SAS

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GFL Water Stills models 2001/2 and 2001/4 produce distillate of good quality with a low conductivity (approx 2.3 µS/cm at 25°C), dependent on quality and composition of the feed water.

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Before installation, please check whether contents of package are in good order and complete.

Should you note any damages or have any reasons for complaint, please contact your supplier or directly:

LAUDA-GFL Gesellschaft für Labortechnik mbH
Schulze-Delitzsch-Strasse 4+5
D-30938 Burgwedel
Telefon +49 (0)5139 9958 0
Telefax +49 (0)5139 9958 21
E-Mail info@lauda-gfl.de
Internet www.GFL.de

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Contents –

1	Intended and Improper Use	Page	11
2	Warranty	Page	11
3	Before Installation	Page	12
4	Set-up and Location of the Water Still	Page	12
5	Operating Voltage	Page	12
6	Water Connections	Page	12
7	Initiation and Setting the Cooling Water Temperature	Page	13
8	Functional Description	Page	14
9	Maintenance, Service and Removal of Malfunctions	Page	14
10	Disposal of Old Units	Page	16
11	Technical Data	Page	17
12	Circuit Diagram	Page	18
13	Connection to the Mains Supply	Page	19
14	Spare Parts	Page	21
15	EC Declaration of Conformity	Page	23
16	Accessories	Page	24

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1.Use

1.1 Intended Use

> GFL Water Stills have been tested by TÜV SÜD for compliance to the basic requirements of the laboratory equipment regulations EN 61010-1:2010 and EN 61010-2-010:2014 as well as to the general engineering standards. The manufacturing site is regularly controlled by the testing institution. The units may be marked with the depicted quality symbol of TÜV-SÜD Product Service GmbH.

> In Mono Water Stills 2001/2 and 2001/4, an electric heating element serves to boil water and transform it to steam. The steam is led into the condenser that is placed on the boiler and condenses on the cooling coil. The produced distillate drains through a drain tube on the front of the condenser. Depending on the model, the distillate quantity is 2 litres (2001/2) or 4 litres (2001/4).

> The Water Stills must be operated within the user's field of vision.

> Depending on the tap water quality, the produced distillate has a conductivity of approx. 2.3 µS/cm at 25°C. To feed the Water Still, preferably use tap water of drinking quality. Please also observe the information in chapters 16.1 to 16.3 of these instructions on the possible use of filters and cartridges to pretreat the tap water.

> 

The information contained in these operating instructions must **by all means** be read and observed. Only then can a perfect operation of the **Water Still** be guaranteed. The units may only be installed and operated by persons who have made themselves familiar with these operating instructions.

> 

Both the outer housing of the boiler and the steam condenser get strongly heated up during operation. Both parts may only be touched after they have cooled down, or when using suitable safety gloves.

1.2 Improper Use

> GFL Water Stills, operated in a laboratory, are **no Medical Devices**. They fall neither under national nor international Medical Device Directives and have to be used and applied accordingly.

> The Water Still may not be used in potentially explosive surroundings.

> The Water Still may neither be set up nor operated in laboratory areas with aggressive or corrosive ambient conditions.

2. Warranty

The warranty period for all laboratory apparatus and their accessories by GFL Gesellschaft für Labortechnik mbH is 24 months from the date of the invoice made out by GFL. The warranty period for repairs and modifications, carried out by GFL, as well as spare parts is 6 months from the date of the invoice made out by GFL. In order to identify defective units, we require both model and serial number on the nameplate at the back, left-hand side of the Water Still and, if applicable, a copy of the invoice.
A free-of charge guarantee repair will not be granted for defects that are due to improper installation and handling.

DOMINIQUE DUTSCHER SAS

3. Before Installation

Important information are marked in bold letters in these instructions, safety indications are additionally marked by the following warning symbols and mandatory signs.

Warning of hot liquids and vapour

Warning of hot surfaces

Read and observe the operating instructions

Warning of dangerous electrical voltage

Before maintenance and repair disconnect the unit all-pole from the electrical mains (pull the plug from the socket).

All free openings of the water connections on the right-hand side of the unit are taped close. These transport protections have to be removed before installation.

4. Set-up and Location of the Water Still

Place on solid, even and level surfaces **inside buildings** only. Make sure to place the unit only on a watertight, temperature-resistant, non-flammable surface. The location must provide sufficient space as well as the necessary carrying capacity for the unit (unit weight as per technical data, chapter 11 of this manual, plus weight of the filling.
The unit is **not** suitable for use in **potentially explosive areas**.

5. Operating Voltage



The Water Still must be connected to the mains through a correctly installed shock-proof socket or through an on-site main switch. **The Still is a protection class I electrical appliance, a connection to the earth conductor (PE) must be ensured.** For information on the required mains fuse please view Technical Data, chapter 11 of this manual.

The electrical connection must ensure an all-pole separation of the Water Still from the mains at any time.
The mains connection cable must not touch any hot surfaces of the unit. It may not lead underneath the unit.
The main switch of the unit must be **off** (position O)**. The voltage on the nameplate (at the back of the unit) must be identical to the mains voltage. If they are identical, connect the unit to the mains.**
For further information, please also refer to chapter 13 of these operating instructions "Connection to the Mains Supply".

6. Water Connections

All water connections of the Water Still, with the exception of the distillate withdrawal at the front, are situated at the right-hand side of the unit. Hoses for water inlet and outlet are not included in the standard scope of supply.

DOMINIQUE DUTSCHER SAS

6.1 Tap water inlet (ill. A)

Water is fed to the Water Still through the tap water inlet. Connect a ½" pressure hose between the tap for cooling water inlet and an on-site manual stop valve of the mains water supply.

Make sure to secure both hose connections with hose clips!

6.2 Cooling water outlet (ill. B)

Connect a ½"

Connect a ½" **temperature-resistant hose** with a maximum length of 1.5 m to the cooling water outlet. The hose must lead into a drain on a lower level and must have a slope on its complete length.

Make sure the cooling water can drain without back draughts.

Caution, the cooling water leaves the Water Still with a temperature of up to 70 °C.
Danger of scaldings!

6.3 Distillate outlet


Use a temperature-resistant hose (not depicted) to connect the distillate outlet tube on the front of the condenser with a storage tank (to be provided on-site). When choosing the storage tank, observe the special storage conditions for distilled water according to your requirements as well as the high temperature of freshly produced distillate.

Caution, distilled water leaves the Water Still with a temperature of up to 100 °C.
Danger of scaldings!

7. Initiation and Setting the Cooling Water Temperature

Lift the condenser from the boiler (ill. E, upwards arrow) and fill in water until the heating element in the boiler is completely covered by water.

7.1 Control Thermometer (ill. C)


Insert the supplied thermometer into the tube support on the right-hand of the upper part.

DOMINIQUE DUTSCHER SAS

7.2 Water inlet and Setting the Cooling Water Temperature (ill. D and E)

Open the tap of the cooling water supply on the unit completely (handle in the direction of water flow), and the on-site stop valve of the mains water supply a little.
Check the power supply and switch on the main switch. The pilot lamp in the main switch lights up.
As soon as the water in the boiler reaches boiling point, regulate the cooling water supply on the tap of the mains water supply so as to set the temperature of the draining cooling water to approx. 60°C to 70°C. The temperature can be checked with the inserted thermometer. The Water Still achieves its maximum efficiency with a cooling water temperature of 60°C to 70°C.

7.3 Degassing and Distillate Withdrawal (ill. F)

Carbon dioxide that forms during the distillation process, is released through a vent in the condenser. Produced distillate is to be led through a hose (not shown) from the drain tube into a storage tank.

8. Functional Description

The tap water flows through the cooling coil in the condenser to the water level regulator on the right-hand side of the unit. The mechanical water level regulator determines the water level in the boiler so as to ensure that the heating element is constantly kept below water level. Water that is not used for evaporation flows off through the cooling water outlet.
When the heating element is switched on, it brings the water in the alembic to the boil. A thermostatic **low water cut-off** protects the heating element from running dry. The produced steam is led into the condenser, condenses on the cooling coil and drips off as distillate through the drain tube on the front. By setting the cooling water temperature to 60°C to 70°C, to be checked with the inserted thermometer, the unit's efficiency and the produced quantity of distillate is optimised.

9. Maintenance, Service and Removal of Malfunctions


Caution! Before servicing the unit
let the Water Still cool down!
Danger of burns!

Caution, before opening the Water Still as well as before cleaning procedures, disconnect the unit from the mains!
(Unplug the unit or switch off the main switch).
Danger of electrical shock!

DOMINIQUE OUTSCHIER SAS

9.1 Descaling

Dependent on the hardness of the tap water, the Water Still must be cleaned from scale formations in regular intervals. To do so, close the on-site water supply valve and lift the condenser from the boiler.
A suitable descaling agent is a mixture of 10% formic acid, 10% acetic acid and 80% distilled water. Fill the solvent into the boiler until it covers the topmost scale deposits and **heat up to maximum approx. 70°C**. Make absolutely certain not to bring the mixture to a boil.
Pour the scale-solvent mix from the boiler after approx. 30 minutes and **then rinse the boiler thoroughly with tap water**. Commercial descaling agents, suitable for use with stainless steel, can also be used according to the manufacturer's instructions (e.g. rea-calc® of M/s CHEMOTEC GmbH, 63486 Bruchköbel, Germany).

> **Never use any products containing hydrochloric acid!**
> **These will lead to damages to heating element, boiler, temperature sensor and the ducts of the screw connections.**

Restart the unit as described in **chapter 7**.
After descalings, the first few litres of distillate should not be used as it might contain traces of evaporated descaling agent.

9.2 Cleaning

Depending on the degree of impurities in the feed water and the growing contamination of the water in the boiler, caused by the distilling process itself, the boiling process will cause foam formation. To prevent these impurities from migrating into the pure water area of the condenser, the water in the boiler must be exchanged at least twice a week. To do so, close the on-site water inlet valve and lift the condenser from the boiler. Then empty the contents from the boiler and rinse thoroughly with tap water.
Restart the unit as described in **chapter 7**.
The "Cleaning" procedure serves to exchange contaminated water in the boiler against clean water.

9.3 Re-Initiation after Low Water

In case of low water, the Water Still's temperature limiter will switch off the unit. In order to restart, allow the **unit to cool down** and, as described in chapter 7, fill water into the boiler until the heating element is below water level.
The triggered low water cut-off must be released. To do so, loosen the black cap nut at the right-hand bottom of the Water Still. Inside the thread a small white plastic pin can be seen that has to be gently pressed inside (e.g. with a pen) until a clicking sound can be heard.

DOMINIQUE DUTSCHER SAS

The **GFL Water Still** is produced with first-class materials. Nevertheless, the unit should only be subjected to mechanical strain in sensible limits.

> **Make sure that no fluids come into contact with cable connections or the electrical parts inside the unit!**

The powder-coated surfaces of the unit can be cleaned with mild detergents, if necessary.

> **Servicings, repairs or modifications must be carried out according to the General Engineering Standards by a competent electrician.**

Only original spare parts must be used. Always demand a detailed confirmation by the person in charge (company, date, signature) on the kind and volume of the tasks carried out.

10. Disposal of Old Units

GFL will take responsibility, within the scope of the legal directives, for an environmentally sound handling and disposal of all used GFL units as of the production year 1995 that are returned to us **free of charge** and will have it materially recycled. Before the unit is returned, a legally binding declaration must be provided from the sender confirming that the unit is free from harmful and/or hazardous contaminations as well as from hazardous substances caused by the previous use of the unit.
GFL laboratory apparatus are exclusively designed for industrial use and may not be disposed of through public waste disposal authorities.
EAR Registration Number WEEE-ID.NO.DE 67770231

DOMINIQUE DUTSCHER SAS

11. Technical Data


Exterior dimensions (W x D x H)			270 x 250 x 510 mm
Distillation capacity		Model 2001/2 Model 2001/4	2 l / h distillate 4 l / h distillate
Distillate quality			Mono distillate approx. 2.3 µS/cm at 25°C **The conductivity of the distilled water is directly related to the chemical composition of the feed water. Components of the feed water having the same of a lower evaporation point than water may deteriorate the conductivity.**
Cooling water requirement		Model 2001/2 Model 2001/4	20 litres / hour 40 litres / hour
Low water cut-off (thermostat)			Electromechanical temperature limiter, with capillary tube sensor. Switch-off temperature 135°C / -15K.
Required cooling water pressure			1-2 bar 14.5 psi / 29.0 psi
Electrical connection / Mains connection		Model 2001/2	230V +/-10%, 50/60 Hz, 2.0 kW Schuko plug CEE 7/7
Caution! **Mains voltage deviations**, even within the mentioned limits, **influence the quantity of distillate produced!**		Model 2001/4	230V +/-10%, 50/60 Hz, 3.0 kW Schuko plug CEE 7/7
Mains fuse	On-site Internal On-site Internal	Model 2001/2 Model 2001/2 Model 2001/4 Model 2001/4	10 Amp. (max. 16A) 10 Amp. inert 16 Amp. 15 Amp. inert
Protection class / - type			I / IP20
Surrounding conditions			Use only inside buildings **(not in potentially explosive areas)** Up to 2000 m MSL +10 °C to +40 °C max. 80% rel. humidity to 31°C, decreasing to 40% rel. humidity at 40°C
Weight Net / with water filling		Model 2001/2 Model 2001/4	7.7 kg / 12.2 kg 8,0 kg / 12.5 kg

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
* s. Typenschild

13. Examples for connection to the mains supply.

Examples for connection to the mains supply page 20

Water Stills models 2001/2 and 2001/4 are supplied with a pre-assembled, cast-on shock-proof plug (CEE 7/7). Make sure to connect to a protective conductor terminal.

Colour coding of the mains cable	**Mains supply**
ge/gr – yellow/green	PE (Protective earth)
bl – blue	N
sw – black	L1

Models 2001/2 and 2001/4 **supplied for 230V** (see information on the nameplate) can be connected **to all power supplies of 220V or 230V**.

Electrical fuses

Model	**Power**	**Power consumption at mains voltage***	**Mains fuse** (F4, F5)
2001/2	2,0 kW	8,7 Amp at 230 V	**10 Amp**
2001/4	3,0 kW	13.0 Amp at 230 V	**16 Amp**

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

Models 2001/2 and 2001/4 for 230V with mains supply 230V / N / PE / 50 / 60 Hz connected through 3-pole shock-proof plug system.

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
Models 2001/2 and 2001/4 for 230V with mains supply 400V / 3 / N / PE / 50 / 60 Hz connected through on-site main switch

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Models 2001/2 and 2001/4 for 230V with mains supply 220V / 3 / PE / 50 / 60 Hz connected through 3-pole shock-proof plug system

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declare that the below stated **Water Stills** models:

2001/2 and 2001/4

with the technical data: **230 V, 50/60 Hz**
2,0 kW (2001/2)
3.0 kW (2001/4)

are in conformity with the following EC Directives:

I	**2014/35/EU**	**(Low Voltage Directive)**
II	**2014/30/EU**	**(EMC Directive)**
III	**2011/65/EU**	**(RoHS 2 Directive)**

For conformity **with I** the following standard was applied:

EN 61010–1:2010
EN 61010-1-010:2014

For conformity **with II** the following standard was applied:

EN 61326-1:2013

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Andreas Degmayr
Managing director

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Dechlorite Filter eliminates chlorine particles in tap water added by the local waterworks. Complete with connection for pressure hose ½ inch, with first filling.

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Phosphate cartridge prevents scale formation in the condenser by phosphatising of tap water. Can be used from 4 to 15°dH (German hardness), equivalent to approx. 0.7 to 2.7 mMol/l. With connection for pressure hose ½ inch, with first filling.

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Pre-Filter 1 µm, for pre-cleaning the tap water, and to protect the unit from premature contamination. Complete with connections for pressure hose ½ inch, including filter candle. The candle should be replaced after maximum 6 months.

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Wall bracket, for one filter or for filter combinations of two or three filters, including sleeves to connect the filters as well as screws to fix the filters to the wall bracket.
A data sheet on possible connection variants of articles 2904, 2906 and 2912 can be supplied on request.

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Hose Set, consisting of hoses for water inlet and outlet (length 1.5 m), including hose clips.

Schlauchsatz Bestell Nr. 2940	Hose Set Order no 2940

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

LAUDA-GFL Gesellschaft für Labortechnik mbH
Postfach / P.O. Box 1152 · D-30927 Burgwedel
Schulze-Delitzsch-Straße 4+5 · D-30938 Burgwedel
Telefon +49 (0)5139 9958 0 · Telefax +49 (0)5139 9958 21
E-Mail: info@lauda-gfl.de
Internet: www.GFL.de