MOINEAU Instruments – ZA P. Pierre - 8 Allée des Erables – 79110 CHEF-BOUTONNE – FRANCE
Tel : (33) 05.49.29.81.25 – Fax : (33) 05.49.29.66.21
Site internet : www.moineau-instruments.com - E-Mail : contact@moineau-instruments.com

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Instructions 2143D – Maxi/Mini Thermometer / Hygrometer

Safety advice
- Clean regularly with a damp cloth. Do not use abrasive chemicals
- Never immerse the thermometer in water. Do not allow water to run down the thermometer.
- Close the battery compartment securely to prevent moisture ingress. Store in a sheltered place.
- The device must not be exposed to temperatures above +60°C or used inside an oven

Installation
1. Open the battery compartment on the back of the device.
2. Remove the protective seal and close the battery compartment cover.
3. When the display becomes dim, replace the faulty battery with a 357 - LR44 – 1.5 V button cell battery observing proper polarity.

Note : If the device is to be stored for long periods, remove the battery.

Positioning
Choose one of the following :
- Apply 3M double-sided adhesive tape to the top of the thermometer, then stick the thermometer onto a clean, flat surface.
- Lift the stand at the back of the device to place it on the desired surface.

Units
- Select °C / °F unit on the front case of the device.

Temperature reading
- The thermometer continuously displays the temperature and humidity.
- Allow it to stabilise to obtain an accurate reading.

Note : After changing the battery, the maximum and minimum temperatures are reset.
You can safely leave the thermometer switched on: its power consumption

Maximum/minimum readings
The thermometer has a memory for the minimum and maximum readings recorded since it was last reset:
1. Press the MAX/MIN button. The MAX icon and the maximum temperature reading will appear on the display.
2. Press the MAX/MIN button again. The MIN icon and the minimum temperature reading will appear on the display.

Note : To reset the memory, press and hold the MAX/MIN button. The temperature and humidity readings will return to the current values.

Technical Specifications
- Ambient temperature range : -10 ~ +60°C
- Accuracy : ± 1°C / Resolution: 1°
- Ambient humidity range : 10 ~ 99 %RH
- Accuracy : ± 5 %RH / Resolution: 1 %RH
- Dimensions : 52x39x15 mm / 30 g
- Display : 37x17 mm
- ABS
- IPX2 waterproof
- Battery: AG13 – LR44 – 1.5 V

Compliance
IHM – MOINEAU INSTRUMENTS, an importer of time-keeping and measuring instruments, declares that this equipment complies with EC Directive 2014/30/EU and other French and European directives governing the placing of products on the market. This device complies with EU standards on the management of waste electrical and electronic equipment (WEEE) and the restriction of hazardous substances in electrical and electronic equipment (RoHS).
This product has been manufactured using high-quality materials that can be recycled and reused.

IMPORTANT : Please ensure that used batteries and equipment are disposed of in an environmentally responsible manner.

DOMINIQUE DUTSCHER SAS

RoHS CE