SARTORIUS

1000062289

EN **Operating Instructions for disposable small scale cassette devices used in Biopharmaceutical Applications**

Sartocon® Slice 50

1. Important information

Observe the following when setting up the system and during operation:

- Do not use denatured ethanol for storage.
- Do not exceed pressure at the feed inlet P_{inlet} of 3.5 bar (50 psig). Higher pressures may damage the filter.
- Avoid abrupt changes in pressures (approx. > 1 bar/sec.), like those that can occur when the pump is switched on or the valve settings are changed.

2. Technical data on the filter

Materials	
Membrane Materials	Hydrosart® \| Polyethersulfone
Sealing and Gaskets	Silicone, Epoxide, Polypropylene
Spacer	Polypropylene
Encapsulation	Reinforced ABS
Luer Lock Connectors	ABS
Filtration area	
Sartocon® Slice 50	50 cm²
Operating Data	
Feed pressure, P_{in}	max. 50 psi \| 3.5 bar
Max. backpressure	7 psi \| 0.5 bar for Hydrosart® 14.7 psi \| 1 bar for PESU
Operating temperature	max. 40°C
pH stability	1–14 Pesu \| 2–14 Hydrosart®
Integrity limit at 36 psi \| 2.5 bar	≤ 3 mL/min
Recommended cross flow rate E Screen	50 – 120 mL/min
Recommended cross flow rate ECO Screen	20 – 60 mL/min

3. Possible setup of a Cross Flow System

The following picture shows a Sartoflow® Smart System with a Sartocon® Slice 50 device installed.

Sartocon® Slice 50

4. Filter holder (optional)

No specific filter holder is required to operate the Sartocon® Slice 50. To ensure vertical configuration, a non-compression device holder can be ordered from Sartorius Stedim Biotech.

Ordering Information	
Description	**Order Number**
Filter holder for Sartocon® Slice 50	3ZG---0108

5. Installing the filter in the System

To install the Sartocon® Slice 50 in a filtration system, proceed as follows:

- Remove the filter from the packaging.
- Check whether the filter type is correct.
- If you use an optionally available device holder, insert the filter vertically.
- Unscrew all of the port covers.
- Connect the individual filter connections to the corresponding ports on the filtration system (i.e. feed to pump outlet, retentate to tank recirculation port and permeate to drain).

Sartocon® Slice 50

6. Flushing the Filter

The filters are supplied with glycerol wetted membranes. Before using the filters for the first time, they should be flushed with 500 mL ultrapure water, isotonic saline solution and/or buffer per unit. To Flush the cassette add flushing solution to the reservoir and then turn on pump and adjust the flow to the desired flow rate with retentate and permeate lines being directed to drain.

7. Filtration

To start a filtration process, proceed as follows:

- Make sure that the filtration system is filled with one of the liquids listed below:
 - The final rinse water
 - Isotonic saline solution
- Close the permeate valve.
- Introduce the fluid into the system. Generate the desired crossflow rate. The inlet pressure will increase.
- Open the permeate valve completely and adjust the retentate valve until the desired TMP is reached.
- Make sure that the desired crossflow rate for the process is maintained.
- Make sure that the crossflow rate is sufficient after adjusting the back pressure valve.

For information about the reference values for your solution, ask our applications specialists.

8. Integrity testing

After flushing and wetting, an Integrity test can be conducted. Please proceed as follows:

- Drain the system
- Close the retentate valve
- Pump air into the system until 2.5 bar (37 PSI) is reached
- Measure diffusive air flow after two minutes stabilization time.

9. Short term storage

Flush cassettes with water and with either 0.1N NaOH or 20% EtOH. Leave cassette wetted with storage solution and re-install all luer lock caps.

10. Cleaning

Circulate 1N NaOH for 60 minutes as follows:
Slice 50 ECO (part numbers ending MLLPU); 50 mL/min
Slice 50 E (part numbers ending ELLPU); 80 mL/min
p_{in} = 29 psi | 2.0 bar | 200 kPa
p_{ret} = 7.3 psi | 0.5 bar | 50 kPa
p_{per} = 0 psi | 0 bar | 0 kPa

The cleaning agent must flow from both permeate ports.

11. Part numbers

Part Number	Membrane Material	Nominal molecular weight cut-off
3D91441950ELLPU	Hydrosart®	2 kD
3D91442950ELLPU	Hydrosart®	5 kD
3D91443950ELLPU	Hydrosart®	10 kD
3D91443950MLLPU*	Hydrosart®	10 kD
3D91445950ELLPU	Hydrosart®	30 kD
3D91445950MLLPU*	Hydrosart®	30 kD
3D91446850ELLPU	Hydrosart®	100 kD
3D91446850MLLPU*	Hydrosart®	100 kD
3D91447950ELLPU	Hydrosart®	300 kD
3D91447950MLLPU	Hydrosart®	300 kD
3D91460950ELLPU	PESU	1 kD
3D91462950ELLPU	PESU	5 kD
3D91463450ELLPU	PESU	8 kD
3D91463950ELLPU	PESU	10 kD
3D91463950MLLPU*	PESU	10 kD
3D91465950ELLPU	PESU	30 kD
3D91465950MLLPU*	PESU	30 kD
3D91465050ELLPU	PESU	50 kD
3D91466850ELLPU	PESU	100 kD
3D91466850MLLPU*	PESU	100 kD
3D91467950ELLPU	PESU	300 kD
3D91467950MLLPU*	PESU	300 kD

* ECO Screen

12. Disposal (Status: November 2018)

If non-hazardous residual impurities are present, the EWC waste code 150203 (European Waste Catalogue) for absorbents and filter materials, wiping cloths and protective clothing can be applied.
If the filters are contaminated with hazardous materials, the EWC waste code 150202* for absorbents and filter materials contaminated with hazardous materials should be used.
If you are uncertain, contact your waste disposal authority.

* Hazardous waste that requires special monitoring under the German laws on the recovery and recycling of waste.

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© 2022
DOMINIQUE DUTSCHER SAS
Sartorius Stedim Biotech GmbH