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
General Guidelines for Freezing and Long-Term Storage of Microbial Cultures with CLEARStable

CLEARStable is a system designed for the freezing and long-term storage of microbial cultures, allowing the preservation of microorganism viability without the need for lyophilization or frequent subculturing.
The system is available in two configurations:

CLEARStable with cryoprotectant: 2 mL cryotube containing large or small porous beads, immersed in a glycerol-based cryoprotective solution, suitable for the preservation of bacteria, bacterial spores, and fungi, including molds and yeasts.

CLEARStable Dry: 2 mL cryotube containing large or small porous beads, without cryoprotective solution, allowing the operator to choose the most appropriate cryoprotectant based on the microorganism to be preserved.

Usage Procedure

Preliminary preparation for CLEARStable Dry only:

Before beginning the preservation procedure, add a suitable cryoprotectant (selected by the operator), depending on the microorganism or strain to be preserved.
Ensure that the porous beads are completely immersed in the cryoprotective solution before proceeding.
Once the cryoprotectant has been added, continue by following the standard instructions for use provided below.

1. Preservation of the Microbial Strain:

1.1. Preparation of the microorganisms

- For greater efficiency, ensure that the culture (either liquid suspension or plate) is prepared the day before use, allowing it to reach the logarithmic growth phase (typically 20-24 hours, depending on the microorganism).

1.2. Preparation of the Cryotube

- Remove a sterile cryotube from its packaging, ensuring sterility is maintained.
- Label the cryotube with the microbial strain name, date, and any other relevant information.

1.3. Inoculation


- From liquid suspension: using a sterile inoculation loop or pipette, transfer approximately 100-200 µl of actively growing culture into the cryotube.
- From plate Culture: Take a sterile bead from the cryotube, gently streak it across the inoculated plate containing microorganisms. Ensure proper adherence of the microorganisms to the bead before transferring it back into the cryotube. Alternatively, you can also pick colonies directly from the plate and transfer them into the cryotube.

1.4. Mixing and Distribution

- Cap the cryotube and gently invert it several times to ensure the beads are fully immersed in the cryoprotective solution and that the microorganisms are evenly distributed. **Do NOT vortex**.
- Let the vial stand for 30-60 seconds.
- Using a sterile pipette, remove the excess cryoprotective solution, leaving the beads as free of liquid as possible.
- Close the vial finger tight.

1.5. Freezing the Cryotube

- Immediately store the inoculated cryotube in a *-20°C freezer* for short-term storage or a *-80°C freezer* for long-term storage, or in the *vapor phase of nitrogen* for optimal preservation, ensuring the freezer is well-maintained and monitored to maintain a constant temperature.

2. Retrieval of the Microbial Strain:

2.1. Removal of the cryotube from the nitrogen vapor phase or the freezer


- When necessary, promptly remove the cryotube from storage.

2.2. Removal of a Bead

- Open the cryotube and, using sterile tweezers or a similar tool, (e.g., Microloops ref BSV120 from Biosigma), remove a frozen bead.
- Immediately close the cryotube and return it to the freezer to minimize temperature fluctuations.

2.3. Thawing and Inoculation

- Transfer the frozen bead into a fresh, appropriate culture medium (e.g., nutrient broth or agar plate). If using a liquid medium, immerse the bead into the broth; if using a solid medium, gently streak the bead across the surface.
- Incubate at the optimal temperature for the recovery and growth of the microorganisms.
- Safely discard the bead after use, according to the established guidelines for contaminated material disposal.

Precautions and Warnings:

- **Sterility**: Always work under aseptic conditions to prevent contamination of microbial cultures.
- **Personal Protective Equipment (PPE)**: Wear appropriate PPE throughout all handling procedures.
- **Temperature Control**: Avoid thermal shocks and the thawing of the beads during retrieval. The use of a pre-cooled cryoblock is recommended to limit temperature fluctuations in the cryotube..
- **Vial Integrity**: Do not use the product if there are leaks, damages or abnormal turbidity of the cryoprotective solution.
- **Storage Conditions**: Select the most appropriate temperature depending on the micro-organism and the expected shelf life.
- **Product Suitability**: CLEARStable was tested with reference strains. It is recommended to check their suitability for the specific strains to be stored.

Additional product Information:

- **Quality Control**: Each batch is subjected to sterility and performance tests, documented in the Certificate of Analysis (CoA).
- **Color Variants**: CLEARStable is available with porous beads and matching color caps in white, red, blue, yellow, green and mixed variants. The cryoprotective formulation is identical for all stains.
- **Available Configurations**: CLEARStable is available with large or small porous beads, with or without cryoprotective solution, in rack or bulk packs and with 1D or 2D barcode identification systems.
- **Storage before use**: Store unused CLEARStable at temperatures between +4 °C and +30 °C.