

Now Available from **BIOLGIX**  
**NEW --Solution Basins!**



**Features and Benefits**

- **Volume**-- Available in 25ml and 100ml
- **Materials**-- Made from polystyrene
- **Graduated**-- Marks inside wall
- **Multi-purpose**-- Can be used with single or multi-channel pipettors
- **Stable**-- A wide base and reinforced sidewalls avoids spills and provides stability
- **Radiation sterilization, Individually Wrapped**

**New Design**

Graduated marks on inside wall



Pour-off spouts on all four corners



Angled sides and V-shaped bottom permits easy access to solution

**Do not forget our other Solution Basins**

Catalog #	Material	Capacity	Color	Sterile	Unit
25-0052	PVC	55ml	Clear	No	100 Pieces/Pack, 8 Packs/Case
25-0051	PS	55ml	White	Yes	50 Pieces/Case
<b>25-1025</b>	<b>PS</b>	<b>25ml</b>	<b>White</b>	<b>Yes</b>	<b>50 Pieces/Case</b>
<b>25-1100</b>	<b>PS</b>	<b>100ml</b>	<b>White</b>	<b>Yes</b>	<b>50 Pieces/Case</b>
25-1202	PP	50ml	Assorted Colors	No	25 Pieces/Pack, 4 Packs/Case

