

# Cell Strainers



**Sterile**



**Sturdy**




15-1040: 40µm




15-1070: 70µm



15-1100: 100µm




**40µm, 70µm, and 100µm**



**Individual package**

**Features:**

1. Made of sturdy polypropylene and strong nylon mesh
2. Strainers are color-coded according to pore size for convenience
3. Fits perfectly into Biologix 50ml Centrifuge Tubes
4. Radiation sterilized, non-pyrogenic



DOMINIQUE DUTSCHER SAS

**Biologix Group Limited**

Phone: +86-531-67802668

[Inquiry@BiologixGroup.com](mailto:Inquiry@BiologixGroup.com)

[www.BiologixGroup.com](http://www.BiologixGroup.com)