

CATALOG NO.	NOMINAL CAPACITY. (mL)				MATERIALS	
	10	28	30	50 ¹	85 ²	TUBE CLOSURE
3118						PC PP

SIZE CODE.....	-0010		-0028		-0030	
NOM. CAP. (mL).....	10		28		30	
CLOSURE SIZE (mm)...	13		20		20	
	INCH		MM		INCH	
A (NECK I.D.).....	0.32	8.1	0.54	13.7	0.54	13.7
B (HT. W/CLOSURE)...	3.21 ± 0.03	81.5 ± 8	3.70 ± 0.03	94.0 ± 8	3.70 ± 0.03	94.0 ± 8
C (TUBE O.D.).....	0.635 ⁺ _{-0.012}	16.13 ⁺ _{-0.30}	1.016 ⁺ _{-0.040}	25.81 ⁺ _{-1.02}	1.016 ⁺ _{-0.040}	25.81 ⁺ _{-1.02}


SIZE CODE.....	-0050 ¹		-0085 ²	
NOM. CAP. (mL).....	50		85	
CLOSURE SIZE (mm)...	24		33	
	INCH		MM	
A (NECK I.D.).....	0.64	16.3	1.00	25.4
B (HT. W/CLOSURE)...	4.19 ± 0.03	106.4 ± 8	4.14 ± 0.03	105.2 ± 8
C (TUBE O.D.).....	1.135 ⁺ _{-0.023}	28.83 ⁺ _{-0.58}	1.503 ⁺ _{-0.023}	38.18 ⁺ _{-0.58}

¹ NOTE: ACTUAL BRIM CAPACITY OF CATALOG NO. 3118-0050 IS 43mL.

² NOTE: ACTUAL BRIM CAPACITY OF CATALOG NO. 3118-0085 IS 81mL.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponoreno Creek Drive
Rochester, NY 14625-2385, USA

NAME	TL	REV
TL	REV	1009
CR 9	DV	
DATE	REV	
10/23/08		
PSS 00-3118-XXXX-1		

OAK RIDGE CENTRIFUGE TUBES