

Corning® Gosselin™ Square HDPE Storage Bottles

High density polyethylene bottles with wide neck

CORNING

- ▶ HDPE storage bottles are easier to handle, require less space (13% to 20%) on the shelf, and are ideal for sampling and storage
- ▶ HDPE resin is lightweight and break-resistant
- ▶ Bottles are leak-proof

Key Features

- ▶ Good for storage of liquid, semi-liquid, and solid forms
- ▶ Available in six sizes: 150 mL, 250 mL, 500 mL, 1L, 2L, and 2.5L
- ▶ Large neck opening provides easier access for liquid transfer
- ▶ Graduated (except 150 mL)
- ▶ Sterile
- ▶ Raw materials are food contact compliant



Ordering Information

Corning® Gosselin™ Square HDPE Storage Bottles

Cat. No.	Description	Qty/Cs
LC150A-05	150 mL, 37 mm red screw cap, sterile	300
LC250-06	250 mL, graduated, 37 mm red screw cap, sterile	210
LC500-11	500 mL, graduated, 58 mm red screw cap, sterile	175
LC1000-13	1L, graduated, 58 mm red screw cap, sterile	90
LC2000-03	2L, graduated, 58 mm red screw cap, sterile	50
LC2500A-07	2.5L, graduated, 77 mm white screw cap, sterile	35

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400