

Address Im Hackacker 15, 8902 Urdorf, Switzerland
Mail address P.O. Box, CH-8902 Urdorf
Phone +41-44-729 62 11
Fax +41-44-729 66 36
Bank Credit Suisse, 8070 Zurich, BC 4835 / SWIFT CRESCHZZ80A
Account no. 370501-21-4 CHF/IBAN CH65 0483 5037 0501 2100 4

www.mt.com/pro

To whom it may concern

Certificate of Quality

We

Mettler-Toledo GmbH, Process Analytics
Im Hackacker 15
8902 Urdorf
Switzerland

confirm that we hold the ISO 9001 quality certificate and the ISO 14001 environmental certificate, and that the following product of our company has been manufactured in accordance with the material specifications mentioned below:

Description of Goods

Redox Buffer Solution 220 mV/pH 7 (9881)

Confirmation

Shelf life, unopened from date of production: two years

After opening use immediately (Oxygen sensitive: store at room temperature)

Traceable to certified reference materials

T°C	U H(mV) ± 5 mV	U Ag/AgCl ± 5 mV	pH ± 0.05
10	462	245	7.05
20	439	228	7.02
25	427	220	7.00
30	415	212	6.98
40	391	195	6.95
50	366	178	6.95
60	340	160	6.96
70	314	142	6.97

Mettler-Toledo GmbH, Process Analytics



Alessandro Deambrosi
Head of Customer Support and Service



Peter Rowing
Head of Quality Management

Place and date of issue

Urdorf, 10.07.2023

This original must not be copied as it may be subject to technical changes.

QC_Redox_buffer_220mV_pH7_Feb21.docx