

Nolato Treff Certificate of Conformity

CoC Certificate according to DIN EN 10204 2.1



Article: Tip clear 0.5 – 5 ml (160mm)

Article Number: 96.08755.9.01

Purchase Number:

Lot- Number:

Quantity:

Statement

Nolato Treff AG confirms that TreffLab consumables and products marked with the CleanRoom® Pure quality symbol have been manufactured in microbiologically controlled clean rooms class 8, ISO 14644-1. To avoid direct human contact, employees are fully protective dressed and masked.

Quality Parameters

Quality standard for original sealed packages

- Free of RNase (< 5x10E-9 Kunitz units)
- Free of DNase (< 5x10E-6 Kunitz units)
- Free of DNA (< 2 pg DNA)
- Free of PCR inhibitors (< 10 copies)
- Free of endotoxins (< 0.001 EU/ml)
- Free of ATP (< 30 RLU)

Analyses are performed by an independent Swiss laboratory for molecular biological & microbiological analytics, accredited ISO/IEC17025 for real-time PCR.

The above products are free of mold release agents, Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE). Nolato Treff is manufacturing under ISO 13485, and therefore uses raw materials which are in compliance with the actual valid RoHS and REACH regulations strictly (heavy metals). It also includes Art. 33 of REACH concerning SVHC substances.

Date: 14.04.2025

Place: CH-9113 Degersheim
Nolato Treff AG

Customer Service
Sandra Boesch