


## Technical Data Sheet

Eppendorf Deepwell Plate 96/1000 µL

English (EN)



### Plate dimensions [mm]

| Length |       |       | Width |      |      | Height |     | Stack height     |
|--------|-------|-------|-------|------|------|--------|-----|------------------|
| l      | l1    | l2    | w     | w1   | w2   | h      | fh  | sh (without lid) |
| 127.8  | 124.5 | 123.2 | 85.5  | 82.2 | 81.0 | 44.1   | 6.1 | 86.8             |

### Well dimensions [mm]

| Well dimensions [mm] |     |          |     | Bottom [°] | Positions of well center [mm] |     |      |     | Bottom space [mm] |     |
|----------------------|-----|----------|-----|------------|-------------------------------|-----|------|-----|-------------------|-----|
| Depth                |     | Diameter |     | Angle      |                               |     |      |     |                   |     |
| h1                   | h2  | d        | d1  | a          | P1                            | P2  | P3   | P4  | b                 | b1  |
| 38.7                 | 2.4 | 7.0      | 6.0 | 90         | 14.4                          | 9.0 | 11.2 | 9.0 | 3.2               | 1.4 |

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)

Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)