

# KIMBLE® Screw Thread Cap, PTFE-Faced Rubber, 20-400

- Specially formulated phenolic resin and liner adhesive to withstand the effects of repeated autoclaving
- PTFE faced/general purpose white rubber liners are highly resistant to chemical effects
- Excellent for general laboratory use
- Identified in ASTM Specification E982, Class A requirements

Catalog Number: **73802-20400**

## Specifications

Material	Phenolic
Liner Material	PTFE-Faced Rubber
Color	Black
Thread	20-400
Cap Style	Solid Top
Quantity	144

Learn about the physical and chemical properties of our glass, plastic and rubber materials.

[View Technical Information](#)

*This information is accurate as of 14/08/2025 and is subject to change.*