





A-A (3 : 1)

Undefined radius: R 0.2 mm

Tolerances unless otherwise specified:

Without decimal: ± 0.5 mm

1 place decimal: ± 0.1 mm

2 place decimal: ± 0.05 mm

Thermo Fisher
S C I E N T I F I C

Nunc A/S · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.thermoscientific.com

SUBJECT:

Cryotube 1,8 ml, round
bottom
int. thread

DRAWN BY/DATE:

PS/24.11.93

CHANGE ORDER NO.:

CO-5368

VALID FROM:

31-03-2014

APPROVED BY/DATE:

ROSPBE 31-03-2014

FORMAT:

A4

ITEM NUMBER:

-

MATERIAL:

PP

DRAWING NO.:

2152

VERSION:

12-b

PAGE TYPE:



-C

SCALE:

3:1

METRIC

First Angle
Projection



Undefined radius: R 0.2 mm
 Tolerances unless otherwise specified:
 Without decimal: ± 0.5 mm
 1 place decimal: ± 0.1 mm
 2 place decimal: ± 0.05 mm

Thermo Fisher
 S C I E N T I F I C

Nunc A/S · Kamstrupvej 90 · P.O.Box 280
 DK-4000 Roskilde · Denmark
 Tel +45 4631 2000 · Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

www.thermoscientific.com

SUBJECT:

Stopper with gasket
 for cryotube, int. thread

DRAWN BY/DATE:

PBE/2012-10-12

CHANGE ORDER NO.:

CO-4444

VALID FROM:

22-07-2014

APPROVED BY/DATE:

ROSPBE 22-07-2014

FORMAT:

A3

ITEM NUMBER:

MATERIAL:

PP/TPE

DRAWING NO.:

5285

VERSION:

10-a

PAGE TYPE:

-C

SCALE:

5:1

METRIC

First Angle
 Projection