






A (4 : 1)




Cavity no.
Text font: Arial
Text size: 1 mm
0,05 positive



123,52 ± 0,25

18,90



85,48 ± 0,25

127,76 ± 0,25

NUNC logo to be polished
Text font: Arial
Text size: 3 mm
0,05 positive

 = Surface frosting, charmillles 25

Thermo Fisher
S C I E N T I F I C

Nunc A/S · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.thermoscientific.com

Undefined radius: R 0.2 mm
Tolerances unless otherwise specified:
Without decimal: ±0.5 mm
1 place decimal: ±0.1 mm
2 place decimal: ±0.05 mm

SUBJECT:

MD24 2nd generation

DRAWN BY/DATE:

AVG/12.03.07

CHANGE ORDER NO.:

CO-2859

VALID FROM:

2011-06-16

APPROVED BY/DATE:

ROSPSI 2011-05-30 A3

FORMAT:

ROSPSI 2011-05-30 A3

PAGE TYPE:

-C

SCALE:

1:1

METRIC:

METRIC