





Tolerances unless otherwise specified:
Without decimal: ± 0.5 mm
1 place decimal: ± 0.1 mm
2 place decimal: ± 0.05 mm

SUBJECT: Upper part Easy Flask 25cm2				www.nunchrand.com			
Resale Site - Kimsrupvej 90 - P.O.Box 280 DK-4000 Roskilde - Denmark Tel: +45 4631 2042 Fax: +45 4631 2042 Email: kimsrup@nunchrand.com							
DRAWING DATE: TSS/22.09.99		DRAWING NO.: CO-101		DATE OF ISSUE: 2008.03.07		APPROVED BY: JR 2008.03.04	
TYPE NUMBER: PS		MATERIAL: PS		QUANTITY: 2521		SCALE: 10	
						-C	
						2:1	
						1/2	