

KIMBLE® Screw Thread Cap, PTFE-Faced Rubber, 22-400

- Specially formulated phenolic resin and liner adhesive to withstand the effects of repeated autoclaving
- PTFE faced/general purpose white rubber liners are highly resistant to chemical effects
- Excellent for general laboratory use
- Identified in ASTM Specification E982, Class A requirements

Catalog Number: **73802-22400**

Specifications

| | |
|----------------|-------------------|
| Material | Phenolic |
| Liner Material | PTFE-Faced Rubber |
| Color | Black |
| Thread | 22-400 |
| Cap Style | Solid Top |
| Quantity | 144 |

Learn about the physical and chemical properties of our glass, plastic and rubber materials.

[View Technical Information](#)

This information is accurate as of 14/08/2025 and is subject to change.