



Chemical resistance guide

LEVEL 0	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	LEVEL 6
< 10 min	10 > 29 min	30 > 59 min	60 > 119 min	120 > 239 min	240 > 479 min	> 480 min

SMARTLine™ Powder-Free Blue Nitrile Gloves, 240



- Length: 240 mm/ 9.4"
- Palm thickness: 0.08 mm/ 3.1 mil
- Chemical performance: Type C
- Biological risk: AQL 1.5 / Level 2
- Particles level: N/A / N/A
- Virus resistant
- Allergies: Latex protein < 50 µg/g / Latex proteins < 200 µg/g / Latex-free / Free of Thiazoles and Thiurams / Free of Thiurams / Free of Thiazoles / Accelerator-free / Single chlorination / Double chlorination / Silicone free
- Design: Ambidextrous / Powder-free
- Colour: Blue
- Mechanical risk: N/A
- Applications: Laboratory / General workplace / Healthcare

1310-73-2 Sodium Hydroxide 40%	LEVEL 6 480 min DR -5%
7722-84-1 Hydrogen peroxide 30%	LEVEL 6 480 min DR 17%
50-00-0 Formaldehyde 37%	LEVEL 2 34 min DR 34%

DISCLAIMER: The data provided was based on gloves tested under laboratory conditions, in accordance with EN 16523-1:2015 (formerly EN 374-3:2003) and EN 374-4:2013. The information is for guidance only and may not reflect the user's application. A risk assessment should always be made by purchaser to assess the suitability of gloves for a specific application.